# Oxygen Saturation Monitor

# **PULSOX™-2**

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## 1.DESPECTION OF THE ANALYZER / PARTS

Please open and unfold the appendix page with diagrams of the analyser, at the end of this instruction manual and follow the different points with the pictures:

### A1 Name of parts

- |                        |                                                                                                                             |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| 1. Battery cover lock: | Locks the battery cover.                                                                                                    |
| 2. Battery cover:      | Open or close this cover when replacing the batteries with new ones.                                                        |
| 3. Strap eyelet:       | Attach the neck strap to this hook.                                                                                         |
| 4. Finger holder:      | Holds finger when inserted for measurement.<br>A light source is located in the upper side, and a sensor in the lower side. |
| 5. Display:            | Displays measured value, pulse level meter and error messages.                                                              |
| 6. Thumb rest          | Place thumb here during measurement.                                                                                        |

### A2 Measurement Values

- |                                                |                                                                           |
|------------------------------------------------|---------------------------------------------------------------------------|
| 1. Battery indication:                         | The battery mark will begin to blink when the battery power level is low. |
| 2. Pulse rate (P.R.)                           |                                                                           |
| 3. Oxygen saturation value (SpO <sub>2</sub> ) |                                                                           |
| 4. Pulse level meter                           |                                                                           |

Calculation method for displayed values

Oxygen saturation value (SpO<sub>2</sub>):

Value obtained by performing moving averaging for the last 5 seconds is updated and displayed every second.

Pulse rate (P.R.):

Value obtained by performing moving averaging for the last 8 pulse rates is updated and displayed every second.

## 2. OPERATING METHOD

### A4 Measurements

This instrument measures with the forefinger inserted in the finger cuff. To ensure stable measurement, hold the instrument properly by placing your thumb on the thumb rest. Most important is to hold the instrument body securely and to not place your thumb on the bottom of the instrument (which would apply a force to the forefinger).

1. Open the finger holder.
2. Insert the index finger and release the finger holder.
  - Opening the finger holder will start measurement automatically.

### Notes on Measurements

- Remove nail polish before taking measurements.
- If no measurement values are displayed, attach the finger holder to a thinner finger to take measurements.
- It may not be possible to take measurements in case of blood circulation or blood stream problems. In this case, remove the finger holder, rub or warm the finger to improve blood circulation, then re-attach the finger holder.
- Do not use the instrument in areas where it is exposed to strong light such as direct sunlight.
- Check that the pulse level meter changes in synchronization with the pulse rate during measurements. In the following cases, the level meter may not operate properly, hindering display of correct values.
  - Rapid body movement
  - Finger Holder is not attached to the patient properly.
  - Arm or finger is pressed, causing bad blood circulation.
- When using the instrument, take care not to expose the finger holder to strong light.

### Ending Measurements

1. Open the finger holder and remove the finger.
  - Measurement will stop automatically, and the display will go off after 10 seconds.

### Notes on Use

- If the instrument will not be used for a long period of time, remove the batteries from the battery chamber.

### 3.ERROR MESSAGES

Error Messages	Cause	Solution
	<ul style="list-style-type: none"> <li>• Insufficient light for measurement.</li> <li>• Strong light enters the finger holder's sensor directory.</li> </ul>	<ul style="list-style-type: none"> <li>• Check that the finger holder is properly attached to the patient.</li> <li>• Clean the light-source and sensor sections.</li> <li>• Take measures to prevent direct entry of strong light to the finger holder.</li> </ul>
	<ul style="list-style-type: none"> <li>• Pulse signal required for measurement is not being received.</li> </ul>	<ul style="list-style-type: none"> <li>• Check that the finger holder is properly attached to the patient. If the finger holder is properly attached but this message still appears, the circulation in the area being measured may not be good. If this occurs, warm up the measurement location.</li> </ul>
<p>The Measurement Value and</p>  <p>flash alternately</p>	<ul style="list-style-type: none"> <li>• The pulse is weak and the pulse signals required for the measurement are not obtained. The measurement value is displayed, however, the accuracy of the measurement values cannot be guaranteed. Thus, the measured value must be considered as guidance only.</li> </ul>	<ul style="list-style-type: none"> <li>• Attach the finger holder to the body correctly. If this message is displayed when the unit is attached to the patient correctly, the measurement conditions are regarded as blood circulation or other circulatory problems. If this occurs, warm up the measurement location.</li> <li>• In particular, when the pulse is weak, this error tends to occur. For improved measurement accuracy, if the finger holder is exposed to strong light, wrap it with a piece of black cloth, etc. Also make sure that the patient rests calmly to avoid influence of body movement.</li> </ul>
<p>The Measurement Value and</p> 	<ul style="list-style-type: none"> <li>• Motion artifact</li> </ul>	<ul style="list-style-type: none"> <li>• Attach the finger holder to the patient properly. If this message still reappears even though the finger holder is attached properly, keep the measuring point as stationary as possible.</li> </ul>
	<ul style="list-style-type: none"> <li>• The pulse rate exceeds the upper limit (250 bpm) of the measurement range.</li> </ul>	<ul style="list-style-type: none"> <li>• Measurement cannot be performed when the upper limit of the measurement range is exceeded.</li> </ul>

## 4.SAFETY PRECAUTIONS

To ensure correct use of this instrument, read the following points carefully and adhere to them. After you have read this manual, keep it in a safe place where it can be referred to anytime a question arises.

### WARNING

(Failure to adhere to the following points may result in death or serious injury.)

- Do not use the instrument in places where flammable or combustible gases (anesthetic gas etc.) are present. Doing so may cause a fire.
- Do not put the batteries on a fire, short-circuit them, heat them or disassemble them. Doing so may cause explosion or heat generation, resulting in fire or injury.
- This device is designed to measure the oxygen saturation level (SpO<sub>2</sub>) and the pulse rate. It should not be used as a warning device to detect cessation of breathing or irregular breathing during sleep.

### CAUTION

(Failure to adhere to the following points may result in injury or damage to the instrument or other property.)

- Do not use batteries other than those specified by MINOLTA. When installing batteries in the instrument, make sure that they are correctly oriented according to the (+) and (-) marks. For a model which use more than one battery, do not mix new and old batteries, or mix batteries of different types. Failure to adhere to these instructions may cause explosion of the batteries or leakage of electrolyte, resulting in fire, injury or corrosion.
- Do not place the instrument on an unstable or sloping surface. Doing so may result in its dropping or overturning, causing injury. Take care not to drop the instrument when carrying it.
- Do not operate the instrument for long periods of time with a finger holder attached to a patient. Low-temperature burn, redness or rash may result. If you feel pain or itchiness, stop use of the instrument immediately and consult a doctor. A doctor should also be consulted before using the instrument on infants, those with peripheral blood circulation problems or those with sensitive skins.
- If alkali fluid from the battery comes in contact with eyes, skin, or clothing, immediately wash the affected area and see a physician for treatment.
- Do not use wet batteries. do not use instrument if the battery compartment contains water. Doing so may cause explosion or heat generation, resulting in fire or injury.
- Do not disassemble or modify the instrument. Doing so may cause a fire or electric shock.
- Do not disassemble or modify this instrument: Otherwise, it may cause a fire or electric shock.

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© 2002 Minolta Co., Ltd.
**Minolta Co.,Ltd.**