
- MANUALE ISTRUZIONI D'USO Aerosol per terapia generica
- INSTRUCTIONS FOR USE MANUAL Aerosol therapy unit for seasonal use
- MODE D'EMPLOI Aérosol pour thérapie générique
- **GEBRUIKSAANWIJZING** Aerosoltherapie voor seizoensgebonden gebruik
- **D** BEDIENUNGSANLEITUNG Aerosol-Inhalationsgerät für den allgemeinen Gebrauch
- MANUAL DE INSTRUCCIONES Aerosol inhalador para terapia general
- РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ Терапевтический аппарат для аэрозольтерапии

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Neb-aid

Mod. P0603EM F400 Aerosol Therapy Device

Congratulations on your purchase and thank you for choosing our unit.

Our goal is to fully satisfy our customers by offering them cutting-edge products for treating ailments of the respiratory system.

Please read these instructions carefully and keep them handy for future reference. Use this device only as described in this manual. This is a home medical device for the nebulization and administration of medication prescribed or suggested by your doctor.

Please visit www.flaemnuova.it to see our entire range of Flaem products www.flaem.it

THIS DEVICE COMES EQUIPPED WITH:

(A) - Aerosol Therapy Device

A1 - Switch

A2 - Air outlet port

A3 - Air filter

A4 - Nebulizer port

A5 - Socket for power cord

A6 - Power cord

A7 - Carrying handle

B - Connection tube (main unit / nebulizer)

C - Accessories

C1 - Pediatric nose piece

C2 - Adult mask

C3 - Mouthpiece

C4 - RF7 Dual Speed Nebulizer

C4.1 - Top piece

C4.2 - Nozzle

C4.3 - Bottom piece

C4.4 - Speed selector

C5 - Manual nebulization control

D - Convenient, roomy carrying bag

NIMPORTANT SAFEGUARDS

- Before using the unit for the first time, and periodically throughout its life, check the power supply cord for any damage; if you detect any damage, do not plug it in. Take the appliance immediately to an authorised service centre or your dealer.
- The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect.
- Children and disabled persons should always use the unit under strict supervision of an adult who
 has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
- · Keep the cord away from heated surfaces.
- Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Only use the device in dust-free environments; otherwise, therapy may be compromised.
- The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to high temperatures.
- Do not put the device near heat sources, direct sunlight, or in hot environments.
- Do not obstruct or insert any objects in the filter or in its housing.

- Do not obstruct the air vents placed on both sides of the unit.
- During use, always place the unit on a firm surface clear of obstacles.
- Check that there is no material obstructing the air vents before each use.
- Do not place any objects inside the air vents.
- The average duration provided for compressor ranges is: F400: 400 hours, F700: 700 hours, F1000: 1000 hours. F2000: 2000 hours.
- Repairs must be performed by authorized personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.
 - WARNING: Do not modify this equipment without authorization of the manufacturer.
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- Interactions: The materials used that come into contact with medications have been tested with a
 wide range of medications. However, due to the variety and continuous evolution of medication,
 the possibility of a chemical interaction cannot be excluded. We recommend that once the
 medication has been opened you should use it as soon as possible and avoid prolonged contact of
 the medication with the nebulizer.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

OPERATING INSTRUCTIONS

Before each use, clean hands thoroughly and clean the device as described in the section on "CLEANING, SANITIZATION AND DISINFECTION". During use, it is advisable to protect yourself from any dripping. It is recommended that each person use their own nebulizer cup and accessories to prevent risk of infection due to contamination.

- 1. Plug the power cord (A6) into the socket (A5) and plug the power cord into a power socket corresponding to the voltage of the unit. This must be positioned so that it is not difficult to make the disconnection from the mains.
- 2. Open the nebulizer by turning the top piece (C4.1) anticlockwise.
- 3. Pour the medicine prescribed by your physician in the bottom piece (C4.3). Close the nebulizer by turning the top piece (C4.1) clockwise.
- 4. Connect accessories as indicated in the "Connection diagram" on the cover. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nose piece or mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).

- 5. Turn on the device using the switch (A1) and then breathe in and exhale deeply.
- 6. When the treatment is finished, shut off the device and unplug the power cord.
- 7. If after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the device's compressor fan; this operation will prevent mold from forming inside the tube.

How to use the "RF7 Dual Speed" nebulizer with the speed selector

Professional and fast, this device is suitable for administering all types of medications, including more expensive ones, even in patients with chronic diseases. The geometry of the internal lines of the RF7 Dual Speed nebulizer ensures the ideal granulometry for active treatment all the way down to the lower respiratory tract.

To speed up the inhalation therapy, move the speed selector button (C4.4) by pressing on MAX with your finger.

To increase the effectiveness of the inhalation therapy, press with your finger on the wording MIN of the speed selector (C4.4). In this position, the speed selector acts as a valve and allows nebulizing the optimal amount of drug to the lower respiratory tract, reducing its dispersion in the environment.

SoftTouch Mask

SoftTouch masks have an outer edge made of **soft biocompatible material** that ensures excellent adherence to the face, and is also equipped with an **innovative Dispersion Limiting Device**. These distinctive elements that distinguish it allow greater sedimentation of medication in the patient and also **limit dispersion**.

During the inspiratory phase, the **tab**, acting as a **Dispersion Limiting Device**, bends inwards towards the mask.

During the expiratory phase, the **tab**, acting as a **Dispersion Limiting Device**, bends outwards from the mask.

USE OF NEBULIZATION MANUAL CONTROL

To achieve continuous nebulisation action you should not use the manual nebulisation control (C5), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (C5) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited,

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

CLEANING, SANITIZATION AND DISINFECTION

Turn off the unit before any cleaning and disconnect the power cord from the mains.

We recommend a personal use of both the nebulizer and accessories to avoid any risk of infection or contagion.

CLEANING THE DEVICE AND THE TUBE'S OUTER SURFACE

Only use a cloth dampened with antibacterial detergent (non-abrasive and free of solvents of any kind).

SANITIZATION AND DISINFECTION OF THE ACCESSORIES

Open the nebulizer by rotating the upper part (C4.1) counter-clockwise, detach the nozzle by pressing as indicated by the 2 arrows in the "Assembly diagram" section C4.

SANITIZATION

Before and after each use, sanitize the nebulizer and accessories according to one of the methods descrive below.

(Method A): The accessories that can be sanitized with the method A are (C1-C2-C3-C4.1-C4.2-C4.3-C4.4). Sanitize the accessories in warm (approximately 40° C) drinking water with mild dishwashing liquid (non abrasive) or in the dishwasher with hot cycle.

(Method B): The accessories that can be sanitized with the method B are (C1-C2-C3-C4.1-C4.2-C4.3-C4.4). Sanitize the accessories by immersion in a solution consisting of 60% water and 40% white vinegar. Finally rinse thoroughly with warm (approximately 40° C) drinking water.

After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

DISINFFCTION

The accessories that can be disinfected are (C1-C2-C3-C4.1-C4.2-C4.3-C4.4)

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized. The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).
- Dispose of the solution as recommended by the disinfectant manufacturer.


REPLACING THE FILTER

The device has a suction filter (A3) that should be replaced when it is dirty or changes color. Do not wash or reuse the same filter. Regular replacement of the filter is required to help and ensure correct performance of the compressor. Check the filter on a regular basis. Contact your local dealer or authorized service center for spare filters.

To replace the filter:

Insert a flathead screwdriver between the edge of the filter and the body. Lift the filter and remove it by turning and pulling it. The filter has been designed so that it is always fixed in its seat.

Use original FLAEM accessories only.

SYMBOLS

Class II device

Important: check the operating instructions

Switch off

Switch on

Complies with: European standard EN 10993-1 "Biological Evaluation of Medical Devices" and European Directive 93/42/ EEC "Medical Devices"

IP21

Enclosure protection rating: IP21. (Protected against solid bodies over 12 mm. Protected against access with a finger; Protected against vertically falling water drops.)

Manufacturer

Type BF applied part

Alternating current

Hazard: electrocution. Consequence: Death. Do not use device while taking a bath or shower

CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments

TÜV Certification ref. EN 60601-1 3° Ed. EN 60601-1-11

SN

Serial number of device

AVAILABLE SPARE PARTS

Description Code RF7 Dual Speed Nebulizer Set composed of: nebulizer and mouthpiece ACO461P Mask Set composed of: adult mask, pediatric mask and elastic band ACO462P 1 m connection tube ACO35 Spare Air fi Iter kit per unit (2 pcs) ACO439P Nose piece Set composed of: adult nose piece and pediatric nose piece ACO463P

TECHNICAL DATA

Mod. P0603FM F400

Voltage:	230V~ 50Hz 130VA	115V ~ 60Hz	220V ~ 60Hz	
In conformity with:	C€ ₀₀₅₁			
Safety certifications:	- account dependent total			

Max Pressure: $1.8 \pm 0.3 \, \text{bar}$ Compressor air output: 9 I/min approx Noise level (at 1 m): 54 dB (A) approx Continuous use Operation: Dimensions (L)x(P)x(H) 18x14x13 cm Weight 1.3 Kg

Operating conditions: Temperature: min 10°C: max 40°C Air humidity: min 10%: max 95% Atmospheric pressure: min 69KPa: max 106KPa

Storage conditions:

Temperature: min -25°C: max 70°C Air humidity: min 10%: max 95% Atmospheric pressure: min 69KPa: max 106KPa

APPLIED PARTS

Type BF applied parts include: (C1,C2,C3)

RF7 Dual Speed nebulizer

Medication minimum capacity 2 ml Medication maximum capacity: 8 ml Operating pressure (with neb.) 0.65 bar

	Speed selector C4.4 in pos Max	Speed selector C4.4 in pos Min
(1) Delivery:	0.55 ml/min approx.	0.25 ml/min approx.
(2) MMAD:	4.58 μm	3.78 μm
(2) Breathable fractions < 5 μm (FPF):	54.4%	63%

DEVICE DISPOSAL

In conformity with Directive 2002/96/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore deliver (or have delivered) the above waste to an appropriate waste collection centre provided by the local authorities, or deliver it to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects to the environment and public health due to incorrect waste management. Unlawful disposal of the product by the user could result in administrative fines pursuant to Article 50 and subsequent amendments of Legislative Decree no. 22/1997.

ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60 601-1-2:2007). Electrical medical devices require special care. During installation and use with respect to EMC requirements, it therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit www.flaemnuova.it. Flaem reserves the right to make technical and functional changes to the product without notice.

Odata detected according to Flaem 129-P07.5 internal procedure
In vitro characterization certified by TÜV Rheinland Product Safety GmbH - Germany in compliance with the new European Standard for aerosol therapy units, Standard EN 13844-1, ANNEX CC. Further details are available on request.

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© 2013 FLAFM NUOVA® All rights riserve cod. 14086C0-0 Rev. 11/2013 (TÜV 4)