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

For full details, please read and take note of the instructions for use.

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*) TouchSpray® Technology made under license from the Technology Partnership PLC.

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
2 INHALE

Sit in an upright position.
Hold nebuliser handset horizontally.

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


Disinfect

Disinfect in the steam disinfector for 6 minutes, or boil thoroughly in distilled water for 5 minutes.

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Dry

Allow to air dry (for about 4 hours).
Store in a clean, dry place.

de	Gebrauchsanweisung	3
en	Instructions for use	29
da	Brugsanvisning.....	53
es	Instrucciones de uso.....	77
fr	Notice d'utilisation.....	103
it	Istruzioni per l'uso	129
nl	Gebruiksaanwijzing	155
no	Bruksanvisning.....	181
pt	Manual de instruções	205
sv	Bruksanvisning.....	231

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For eFlow[®]rapid Type 178G1005

Table of contents

1	IMPORTANT INFORMATION	31
	Organisation of warnings	31
	Handling the nebuliser system	31
	Treatment of babies, children and individuals requiring assistance	31
2	PRODUCT DESCRIPTION	32
	Intended use	33
	Proper use	33
	Intended patient groups	33
	Medications	33
	Life cycle	33
3	PREPARING FOR INHALATION	34
	Information for the first application	34
	Preparing the controller	34
	Operation with the supplied batteries or standard commercial rechargeable batteries	34
	Operation with the power adapter	34
	Preparing the nebuliser	35
4	PERFORMING INHALATION	38
	Pausing during inhalation	38
	Display indicators during inhalation	39
	End of the inhalation session	39
5	CLEANING AND DISINFECTION	40
	Preparation	40
	Control unit and connection cord	41
	Nebuliser	41
	Cleaning	41
	Disinfection	41
	Storage	43
6	TROUBLESHOOTING	44
7	DISPOSAL	48
8	SPARE PARTS AND ACCESSORIES	48

9 TECHNICAL DATA 49
 General 49
 Classification according to DIN EN 60601-1 49
 Materials used 50
 Ambient conditions for operation 50
 Storage and transportation conditions 50
10 EXPLANATION OF SYMBOLS 51
11 TERMS AND CONDITIONS OF WARRANTY 52
 Service and manufacturer 52

1 IMPORTANT INFORMATION

Read these instructions for use and the instructions for use of the accompanying accessories completely before using for the first time. Keep them safely so that you can refer to them again later. If these instructions for use are not followed, injuries or damage to the product cannot be ruled out.

You should have a medical examination before treating illnesses.

Organisation of warnings

In these instructions for use, safety-related warnings are organised in hazard levels:

- The signal word CAUTION is used to indicate dangers which may result in slight to moderate injuries or impair treatment if precautionary measures are not taken.
- The signal word NOTE is used to indicate general precautionary measures for working with the product which should be followed to avoid damaging the product.

Handling the nebuliser system

Check the nebuliser system before every application. Replace all broken, misshapen and heavily discoloured parts. Damaged parts can impair the function of the nebuliser system and thus also affect treatment.

Do not operate the nebuliser system if the power adapter is visibly damaged, since otherwise there is a danger from contact with live parts (e.g., electric shock).


In this case, follow the instructions below:

- Always remove the power adapter from the socket to ensure that the power is completely cut off.
- Never remove the power adapter from the socket with wet hands. There may be a danger of electric shock.
- Do not allow pets (e.g., rodents) near the cables.
- Do not use the nebuliser system in areas where there is a risk of explosion or in the presence of gases that promote combustion (e.g., oxygen, nitrous oxide, flammable anaesthetics).
- Do not carry out an inhalation session while you are driving a motor vehicle (risk of accident).

Treatment of babies, children and individuals requiring assistance

- Babies, children and anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. This is the only way to ensure safe and effective treatment. These individuals often underestimate the hazards associated with the equipment (e.g. risk of strangulation with the cable), and this can result in risk of injury.
- The product contains small parts. Small parts can block the airways and pose a choking hazard. Therefore, be sure to keep them out of the reach of children at all times.
- The nebuliser handsets are only suitable for use by patients who can breathe unassisted and are conscious. This is the only way to ensure effective treatment and avoid a choking hazard.

2 PRODUCT DESCRIPTION

- 1 Control unit
- 2 2 complete nebuliser handsets (incl. aerosol heads):
- 3 Aerosol head*
- 4 Medication cap (preassembled), consisting of:
 - 4a Cap seal
 - 4b Cap
- 5 Medication reservoir and
- 5a Aerosol chamber (preassembled)
- 6 Inspiratory valve
- 7 Mouthpiece with
- 7a Expiratory valve (preassembled)

- 8 Connection cord (connection between control unit and nebuliser handset)
- 9 AC adapter
 - 9a Interchangeable adapter (3 items)
 - Transport case and nebuliser bag
 - easycare cleaning aid for aerosol heads

Check that all components are contained in your package. If they are not, contact your local service partner (see “Service and manufacturer”, page 52).

*) TouchSpray® Technology made under licence from the Technology Partnership PLC

Intended use

The eFlow[®]*rapid* inhalation system is intended for the treatment of diseases of the airways and lungs.

Proper use

The nebuliser system must only be operated in accordance with its intended purpose. The frequency of application may vary depending on the clinical picture and the medication used. The nebuliser handset is designed for use by one patient only.

Info:

The eBase Controller control unit can also be used for a medication-specific nebuliser handset. Follow the instructions for use of the nebuliser handset as well.

CAUTION:

Medications that are approved for inhalation with a medication-specific nebuliser handset must not be used with an eFlow[®]*rapid* nebuliser handset (risk of incorrect dosage).

Only original PARI accessories must be used to operate the system.

NOTE:

Using accessories (components or cables) other than those supplied as spare parts by the manufacturer may cause increased electromagnetic emissions or reduced electromagnetic immunity, leading to a malfunction.

PARI Pharma GmbH is not responsible for damage if the system is used improperly or contrary to the instructions.

Intended patient groups

The nebuliser system is only suitable for use by patients who can breathe unassisted and are conscious.

Inhalation with the mouthpiece is the most effective form of inhalation, since it results in minimum medication loss on the way to the lung. For babies and infants who cannot use the mouthpiece to inhale, PARI offers appropriate masks as accessories.

Medications

The eFlow[®]*rapid* inhalation system must only be used to nebulise permitted inhalation solutions and suspensions. These are usually specified by the doctor. Take note of any limitations in the package leaflet accompanying each medication.

CAUTION:

Never attempt to inhale essential oils, as this can cause irritation of the airways (intense coughing).

Life cycle

- Control unit: 3 - 5 years
- Nebuliser handset (without aerosol head): 1 year
- Aerosol head: 3 - 6 months

The expected life cycle of the aerosol head depends on how often and for how long it is used. With two inhalation treatments lasting a total of 30 minutes' inhalation time and one disinfection per day, a life cycle of up to 6 months can be expected.

3 PREPARING FOR INHALATION

Information for the first application

- To avoid irritation due to inhaling aerosol when it is too cold, do not use medications straight from the refrigerator.
- Clean the nebuliser handset and aerosol head before the first application (see CLEANING AND DISINFECTION, page 40).
- The controller can be operated with non-rechargeable batteries or rechargeable batteries, or it can be connected to the mains power supply via the power adapter supplied.


Preparing the controller

Operation with the supplied batteries or standard commercial rechargeable batteries


i Info:

i Info: Batteries cannot be recharged in the control unit!

- Open the battery compartment by lifting the clip:

- Insert the batteries as indicated by the polarity symbols:

- Close the battery compartment.


Operation with the power adapter

Check that the power adapter is intact before every application. Do not use a power adapter that is damaged or faulty.


i Info:

The batteries or rechargeable batteries can be left in the control unit while the device is operated on mains power with the power adapter. This does not deplete the energy reserves of the batteries or rechargeable batteries.

- Replace the country-specific adapters for the power adapter as required. Make sure that the adapter (9a) engages in the power adapter housing.

- Connect the power adapter to the power cord socket on the controller:

When the power adapter is plugged into a socket, the controller is ready for use.

Preparing the nebuliser


CAUTION:

Check the parts of your nebuliser handset before every application. Replace all broken, misshapen and heavily discoloured parts. Damaged parts and incorrectly assembled nebuliser handsets can impair the function of the nebuliser and thus also affect the treatment.

Info:

The medication reservoir and the aerosol chamber are already preassembled.


- Open the aerosol chamber (5a):


- Insert the aerosol head (3).
It must snap noticeably into place.

NOTE:


Do not touch the membrane (metal part in the middle of the aerosol head), it can be damaged very easily.

- Insert the inspiratory valve (6).
When doing so, check that the two valve flaps are seated correctly. They must lie perfectly flush.


- Close the aerosol chamber:

NOTE:


If the button closure does not close, ensure that the inspiratory valve is seated correctly. It must lie flat against the aerosol chamber.

- Attach the mouthpiece (7) to the aerosol chamber.

i *Info:*
For inhalation with a mask, follow the instructions for use of the mask you are using.


- Connect the connection cord (8) to the controller and the nebuliser.

- Add the quantity of medication prescribed by your doctor to the medication reservoir:


CAUTION:

Ensure that the medication is filled no higher than the upper scale graduation mark (6 ml). Do not overfill under any circumstances!


i *Info:*
If you intend to inhale several medications one after the other, rinse all nebuliser handset parts thoroughly with tap water before refilling.

- Close the medication reservoir:

- Check that all parts are connected tightly and that the medication reservoir is sealed:

4 PERFORMING INHALATION

- Sit in an upright position and relax. This makes inhaling easier and improves uptake of the medication in the airways.
- Hold the nebuliser handset in your hand.


CAUTION:

Hold the nebuliser handset horizontally throughout the inhalation session. If held in a tilted position, an incorrect dosage may occur.

- Hold the mouthpiece between your teeth and enclose it with your lips without blocking the blue expiratory valve.

- Press the ON/OFF button on the control unit to start the inhalation session:

A green LED lights up and an audible signal (1 beep) is emitted to indicate proper functioning.

- Inhale and exhale slowly and deeply through the mouthpiece. The mouthpiece can also remain in the mouth while exhaling. Do not breathe through your nose. You may use a nose clip if your doctor agrees.

i Info:
It is normal for aerosol to escape through the expiratory valve in the mouthpiece during exhalation; this is not a malfunction.

Pausing during inhalation

If it is necessary to interrupt the inhalation session briefly, use the Pause function:

- Press the ON/OFF button for approx. 1 second. Pause mode is activated if the LED flashes green.
- Press the ON/OFF button for approx. 1 second again to resume the inhalation session.

i Info:
The Pause function is available for use after an operating time of at least 5 seconds.

If pausing for a longer time, switch the controller off by pressing the ON/OFF button for approx. 2 seconds.

As long as there is enough liquid in the medication reservoir, the nebuliser can be restarted by pressing the ON/OFF button again.


Display indicators during inhalation

The following indicators appear on the display during inhalation:


Start

During inhalation


Nebulising with battery power


Nebulising with mains power

Activated Pause function

Pause function when operating with batteries

Pause function with mains power

End of the inhalation session


The inhalation session is finished as soon as no more aerosol is generated and the following message appears in the display:

The controller switches off automatically.

i Info:
The nebulising time may become longer with increasing length of use. This does not constitute a malfunction of the nebuliser handset.

i Info:
Note that about 1 ml of medication cannot be nebulised and remains unused in the medication reservoir; this residue must be discarded. This is necessary to ensure a safe quantity of medication is delivered and does not constitute a malfunction:

- When the inhalation session has ended, disconnect the power adapter plug from the socket.

5 CLEANING AND DISINFECTION


The nebuliser handset (and the aerosol head) must be cleaned immediately after each use, and disinfected once a day.

Preparation

Disassemble the nebuliser system into its constituent parts:


- Pull the connection cord out of the nebuliser handset:
- Remove the mouthpiece from the nebuliser.
- Carefully pull the blue expiratory valve (7a) on the mouthpiece (7) out of the slot.

i Info: *The expiratory valve is held in place by a thin fixing clip so that it cannot get lost. This clip must not be detached from the mouthpiece.*


- Open the medication reservoir and pour out the medication residue.

- Disassemble the medication cap. To do this, press the cap seal (4a) downwards and out of the cap (4b).

- Open the aerosol chamber.
- Remove the aerosol head. To do this, gently press the engaging hooks on the aerosol head together from the side:

- Remove the inspiratory valve from the aerosol chamber.

Control unit and connection cord


Clean the housing surface of the control unit and the cable as necessary with a clean, damp cloth.

NOTE:

Do not hold the control unit under running water and do not use liquid cleaning agents. If liquid gets into the control unit, the electronics may be damaged and the control unit may no longer function properly.

If liquid does get into the control unit, contact the local service partner immediately.

- Shake all parts to remove excess water more quickly.
- Place the individual parts on a dry, clean surface and allow to dry completely:

CAUTION:

Moisture promotes the growth of bacteria. Please ensure adequate drying after each cleaning.

Nebuliser

Cleaning

Clean the nebuliser handset (incl. aerosol head) immediately after each application.

- Place all nebuliser parts as well as the aerosol head in warm tap water with a little washing up liquid for approx. 5 minutes. If the parts are heavily soiled, move the nebuliser parts gently back and forth in the washing water.

Disinfection

Disinfect the nebuliser handset parts and the aerosol head at least once a day, e.g., immediately after cleaning at the end of the day (a nebuliser handset cannot be disinfected effectively unless it has been cleaned first).

NOTE:

Do not attempt to disinfect the aerosol head in a microwave oven as it cannot be ensured that the product components will not be damaged.


- Rinse all nebuliser parts and the front and back of the aerosol head thoroughly under running tap water.

With a commercially available disinfectant

- To ensure effective disinfection, use a thermal disinfectant with an operating time of at least 6 minutes. With regard to performing the disinfection, the duration of the disinfection process and the quantity of water required, follow the instructions for use of the disinfectant you are using.

CAUTION:

Inadequate disinfection promotes the growth of bacteria and thus increases the risk of infection. If disinfection is not carried out for the time specified for the device concerned, it will not be effective. Therefore, do not switch the device off before the disinfection cycle is complete. Also, make sure that the device is clean and check it regularly for proper operation.

- Take the parts out of the disinfectant immediately after disinfection is completed. Place them on a dry, clean surface and allow them to dry completely:
- Check the nebuliser parts after every disinfection. Replace all broken, misshapen and heavily discoloured parts.


Possible alternative: Boiling in distilled water

- Put all parts in boiling water for at least 5 minutes. Use a clean saucepan and distilled water for this.

NOTE:

Plastic will melt if it comes into contact with the hot bottom of the saucepan. Therefore, make sure that there is sufficient water in the pan. This will prevent the parts from being damaged.

- Take the parts out of the saucepan.
- Shake all parts to remove excess water more quickly.
- Place the individual parts on a dry, clean surface and allow to dry completely:

CAUTION:

Moisture promotes the growth of bacteria. Please ensure adequate drying after every disinfection.

- Check the nebuliser parts after every disinfection. Replace all broken, misshapen and heavily discoloured parts.

Possible alternative: Chemical disinfection

In general, disinfectants based on quaternary ammonium compounds are suitable for use in disinfecting the nebuliser, including the aerosol head.

When choosing the chemical, please ensure that it is suitable for disinfecting medical devices made from the specified materials (see "Materials used", page 50). Follow the instructions for use of the disinfectant, particularly the dosage information and instructions for safe use.

The effectiveness of this method has been tested with a 2% solution of Bomix® plus (Bode) for an application time of 5 minutes.

Storage

Between applications, particularly for prolonged breaks in therapy, keep the nebuliser handset in a clean, dust-free place (e.g., nebuliser bag).




Pack the control unit, the power adapter and the connection cord in the case provided for this purpose.

NOTE:

Leaking batteries may cause damage to the control unit. Therefore, remove the batteries or rechargeable batteries if the nebuliser system is not expected to be used for a relatively long period of time.

6 TROUBLESHOOTING

You can correct malfunctions that occur during daily operation yourself. The following tables will help you to locate and correct the causes of malfunction.

Display indicator (flashes alternately)	Possible cause	Solution
	Cable connection to the control unit or on the nebuliser handset interrupted.	Check cable connections.
	No medication added.	Add medication.
	Inhalation interrupted.	Resume inhalation by restarting the device (press ON/OFF button).
	The maximum operating time of 20 minutes for one application has been exceeded.	Resume inhalation by restarting the device (press ON/OFF button).
	Battery charge level 50%.	Have new batteries ready.
	Battery flat.	Insert new batteries or operate with the power adapter connected.

Fault	Possible cause/Remedy
The controller cannot be switched on (no beep, no red or green light signal).	<ul style="list-style-type: none"> - To switch on, press the ON/OFF button for approx. 2 seconds. - Is power plug plugged into mains and the device (for mains operation)? - Are batteries inserted correctly? - Check the charge level of the batteries or rechargeable batteries! - Check the connection between the control unit and the nebuliser.
After switching the controller on, no aerosol is generated, and/or the controller switches off again after a few seconds.	Has the medication reservoir been filled with medication?
<p>The controller switches off even though there is still medication left in the medication reservoir.</p> <p>i <i>Info:</i> <i>Note that about 1 ml of medication cannot be nebulised and remains unused in the medication reservoir; this residue must be discarded. This is intentional and does not represent a malfunction.</i></p>	<p>More than 1 ml is left over in the medication reservoir:</p> <ul style="list-style-type: none"> - Check the charge level of the batteries or rechargeable batteries! - The nebuliser handset was not being held horizontally. Hold the nebuliser handset horizontally. - The controller switches off automatically after 20 minutes. Press the ON/OFF button again to resume treatment, so that you can inhale the prescribed quantity of medication.
The device does not switch off automatically when no more aerosol is being produced and there is no more than about 1 ml medication left in the medication reservoir.	The treatment session has been completed successfully and the device can be switched off by pressing the ON/OFF button.

Fault	Possible cause/Remedy
<p>The nebulisation period has become noticeably longer.</p>	<p>If the nebuliser system is used regularly, the nebulisation period will become longer. This is not a malfunction. However, if you notice that the nebulisation period has become considerably longer, this may be due to one of the following causes:</p> <ul style="list-style-type: none"> - You are inhaling a different medication. - The fill volume has been increased. - The medication reservoir cap was not closed tightly. - The aerosol head was not cleaned immediately after the last application, or it is visibly damaged. - You are using the aerosol head for longer or more often than is described in section Life cycle. <p>If none of the points listed above apply, carry out the following function test for the aerosol head:</p> <ul style="list-style-type: none"> • Pour 2.5 ml isotonic sodium chloride solution (0.9%) into the medication reservoir. • Switch the controller on and note how long nebulisation takes without inhaling. If the nebulisation period is longer than 4 minutes, use the <i>easycare</i> cleaning aid. Take note of the information regarding frequency of use of the <i>easycare</i> cleaning aid in the accompanying instructions for use. • If the nebulisation period is still longer than 5 minutes after cleaning with the <i>easycare</i> cleaning aid, you must replace the aerosol head. <p>i <i>Info:</i> The time specified above only applies for 2.5 ml isotonic sodium chloride solution (0.9%).</p>
<p>Aerosol continuously escapes in large volumes from the opening slots in the medication reservoir.</p>	<ul style="list-style-type: none"> • Check whether the nebuliser handset has been assembled correctly. When doing so, check that the two valve flaps of the inspiratory valve are seated correctly. They must lie perfectly flush.

Fault	Possible cause/Remedy
Little or no aerosol is escaping.	The ambient temperature should be at least 10 °C. <ul style="list-style-type: none">• Wait until the device has stabilised at the lowest ambient temperature.
The medication cap cannot be removed from the medication reservoir after the inhalation treatment.	The negative pressure in the aerosol chamber is too great. <ul style="list-style-type: none">• Open the aerosol chamber and gently pull the aerosol head so that air is able to flow in. The cap can now be removed easily.

If the malfunction persists, contact your local service partner.

7 DISPOSAL

The live parts (such as the control unit, the cable, the power adapter and the batteries) as well as the aerosol head must not be disposed of with normal domestic waste.

The device components must be disposed of in accordance with local regulations.

8 SPARE PARTS AND ACCESSORIES

Description	Part no.
Nebuliser handset (incl. aerosol head)	678G8222
Aerosol head	678B2620
Mouthpiece	078B3600
PARI universal mouthpiece (alternative)	022E3050
AC adapter (100-240 V~, 50/60 Hz)	078B7106
Connection cord	178G6009
Carrying case	078E8005
PARI Filter/Valve Set	041G0500
PARI SMARTMASK® (adults)	041G0730
SMARTMASK® Kids (children aged 2 and older)	078G5000
SMARTMASK® Baby (babies weighing 2.5 kg and more)	078G5026
<i>easycare</i> cleaning aid for aerosol heads	078G6100

9 TECHNICAL DATA

General

Electrical connection	Power adapter (REF 078B7106, Type FW7555M/12): Input: 100 V - 240 V~, 50/60 Hz Output: 12 V $\overline{=}$
Batteries	4 x 1.5 V Mignon (AA LR6/Alkaline)
Rechargeable batteries	4 x 1.2 V Mignon (AA rechargeable, at least 2100 mAh)
Nebuliser handset weight	approx. 55 g
Weight eFlow [®] rapid (incl. batteries)	approx. 300 g
Nebuliser handset dimensions (W x H x D)	5.0 cm x 6.3 cm x 14.5 cm
Control unit housing dimensions	H 4.0 cm, Ø 11.6 cm
Min. fill volume	2.0 ml
Max. fill volume	6.0 ml
TOR (Total Output Rate) ^{a)}	610 mg/min ^{b)}
MMD (Mass Median Diameter) ^{a)}	4.1 $\mu\text{m}^{\text{b)}$
RF (Mass fractions) < 5 $\mu\text{m}^{\text{a)}$	69% ^{b)}

a) Measurement using laser diffraction (Helos SympaTec) at 23 °C and 50% relative humidity. Nebulised medium: 0.9% NaCl (4 ml). Inspiratory flow: 20 l/min.

b) May vary depending on medication and the aerosol head used.

Technical electromagnetic compatibility data is available in table form upon request from PARI Pharma GmbH or on the Internet at www.eflowrapid.info.

Classification according to DIN EN 60601-1

Type of electric shock protection (power adapter)	Protection class II
Degree of protection from electric shock of part used (nebuliser)	Type BF
Degree of protection against ingress of water as per DIN EN 60529 (IP rating)	IP 21
Degree of protection when used in the presence of flammable mixtures of anaesthetics with air, with oxygen or with nitrous oxide	No protection
Operating mode	Continuous operation

Materials used

- Nebuliser: Polypropylene, thermoplastic elastomers, polyoxymethylene
- Mouthpiece: Polypropylene, silicone rubber

The eFlow[®]rapid contains no parts made from natural rubber (latex).

Ambient conditions for operation

- Ambient temperature: +5 °C to +40 °C
- Relative humidity of ambient air: 15% to 93% (non-condensing)
- Air pressure: 700 hPa to 1060 hPa

Within professional care settings, the device must only be used in in-patient wards and intensive care units. The device must never be used in areas exposed to elevated magnetic or electrical radiation (e.g. in the vicinity of an MRI scanner).

Users can operate the device both at home and when out and about. Its use on trains and aeroplanes is restricted to passenger areas only. If using the device in a car, it must be operated with battery power.

NOTE:

When using portable wireless communication devices, you must ensure that a distance of at least 30 cm is maintained between all parts of the nebuliser system and these devices (including their accessories, e.g. antenna cables or external antennas). Otherwise, the performance of your nebuliser system may deteriorate.

Storage and transportation conditions

- Temperature: -25 °C to +70 °C
- Relative humidity of ambient air: 0% to 93% (non-condensing)
- Air pressure: 500 hPa to 1060 hPa

NOTE:

Moisture from condensation can impair the controller's functional capabilities. Therefore, you should avoid extreme temperature fluctuations. Wait until the controller has reached room temperature before beginning an inhalation session.

If the device is not expected to be used for a prolonged period, the batteries or rechargeable batteries must be removed.

Info:

i Batteries vary considerably in terms of storage life and output depending on battery quality. The life cycle of about 90 minutes can only be achieved with the recommended batteries (4x).

10 EXPLANATION OF SYMBOLS

	Please follow the instructions for use (background = blue, symbol = white)
	Alternating current
	Direct current
	Item number
	Serial number
	Protection class of the part used: Type BF
	Protection class II appliance
	Minimum and maximum ambient temperature
	Minimum and maximum humidity
	Minimum and maximum air pressure
	CE marking: This product is compliant with the requirements set forth in Directives 93/42/EEC (Medical Devices) and 2011/65/EU (RoHS).
	The medical device was distributed commercially after 13 August 2005. This product must not be disposed of with normal domestic waste. The symbol of a crossed-through refuse bin indicates that separate collection is required.
	Manufacturer

11 TERMS AND CONDITIONS OF WARRANTY

PARI guarantees to you that if used for its intended purpose your nebuliser system will be free from defects in material or workmanship caused by the manufacturing process for the warranty period of two years.

However this warranty does not cover wearing parts, that is to say device parts that are exposed to normal wear, such as the nebuliser handset and aerosol head. The warranty shall no longer be valid if:

- the device has been operated or used in a manner inconsistent with the guidelines in the instructions for use,
- damage is evident that is attributable to external factors such as water, fire, lightning strike or similar,
- the damage has been caused by transporting the device incorrectly, or resulted from a fall,
- the device has been handled or serviced incorrectly,
- the serial number on the device has been changed, removed or rendered illegible;

- repairs, adaptations or modifications have been made to the device by persons not authorised by PARI.

If, exceptionally, a defect does occur, PARI will replace the device. The device may be replaced with an identically or at least comparably equipped model.

Replacement of the device shall not serve as grounds for new claims under the warranty. All exchanged used devices or parts will become the property of PARI. All other claims under the warranty are excluded to the extent permitted by law.

In the event of a claim under the warranty, please contact your local service partner.

The original purchase receipt from the dealer must be presented as proof of warranty and ownership. This warranty applies solely for the initial purchaser of the device.

Service and manufacturer

Manufacturer	PARI Pharma GmbH Moosstr. 3 82319 Starnberg Germany
Contact for technical problems or questions regarding the device	Local service partner Contact: www.eflowrapid.info

The device must only be serviced by PARI Pharma GmbH or a service location expressly authorised by PARI Pharma GmbH.

We reserve the right to make technical changes.

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Reparation av apparaten får endast utföras av PARI Pharma GmbH eller av ett behörigt serviceställe.