Oxy-One Plus

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Mode d'emploi

Shenzhen Creative Industry Co., Ltd.

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Chers clients,

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CE 0123

Shenzhen Creative Industry Co., Ltd.
Floor 5, BLD 9, Baiwangxin High-Tech
Industrial Park, Songbai Road, Xili Street,
Nanshan District, 518110 Shenzhen, P. R.
China

Tel: +86-755-2643 3514

Fax: +86-755-2643 0930

E-mail: info@creative-sz.com

Website: www.creative-sz.com

Shanghai International Holding Corp.
GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg Germany

Oxy-One Plus
Fingertip Oximeter
User Manual

Shenzhen Creative Industry Co., Ltd.

Instructions to User

Dear Customers,

Thank you for purchasing this quality product. Please read the manual very carefully before using this device. Failure to follow these instructions can cause measuring abnormality or damage to the Oximeter.

The manual is published in English and we have the ultimate right to explain the Manual. No part of this manual may be photocopied, reproduced or translated into another language without the prior written consent. We reserve the right to improve and amend it at any time without prior notice.

Version of the Manual: Ver 1.1

Revised date: January 21, 2021

All rights reserved.

Notes

- The contents contained in this manual are subject to change without notice.
- Information furnished by our company is believed to be accurate and reliable. However, no responsibility is assumed by us for its use, or any infringements of patents or other rights of third parties that may result from its use.

Instructions for Safe Operation

- Check the device to make sure that there is no visible damage that may affect user's safety or measurement performance with regard to sensors and clips. It is recommended that the device should be inspected minimally before each use. If there is obvious damage, stop using the device.
- Special attention should be paid while the Oximeter is used constantly at ambient temperature above 37°C, burns may occur because of over-heating of the sensor at this situation.

-
- ⚠ Necessary maintenance must be performed only by qualified service technicians. Users are not permitted to service this device.
 - ⚠ The Oximeter must not be used with devices and accessories not specified in User Manual.

Cautions

- ⚠ Explosive hazard—**DO NOT** use the Oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- ⚠ **DO NOT** use the Oximeter while the patient is under MRI or CT scanning. This device is NOT MRI Compatible.

Warnings

- ⚠ Discomfort or pain may appear if using the Oximeter continuously on the same location for a long time, especially for patient with poor microcirculation. It is recommended that the Oximeter should not be applied to the same location for longer than 2 hours. If any abnormal condition is found, please change the position of Oximeter.

-
- ⚠ DO NOT clip this device on edema or tender tissue.
 - ⚠ The light (the infrared light is invisible) emitted from the device is harmful to the eyes. Do not stare at the light.
 - ⚠ The Oximeter is not a treatment device.
 - ⚠ Local laws and Regulations must be followed when disposing of the device.

Attentions

- 🔔 Keep the Oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- 🔔 The device should be kept out of the reach of children.
- 🔔 If the Oximeter gets wet, please stop using it and do not resume operation until it is dry and checked for correct operation. When it is carried from a cold environment to a warm and humid environment, please do not use it immediately. Allow at least 15 minutes for Oximeter to reach ambient temperature.
- 🔔 DO **NOT** operate the button on the front panel

with sharp materials or sharp point.

-  DO **NOT** use high temperature or high pressure steam disinfection on the Oximeter. Refer to Chapter 7 for instructions regarding cleaning and disinfection.
-  The equipment is IP22 with protection against harmful solid foreign objects and ingress of liquid. So that means the equipment is protected against solid foreign objects of 12.5mm and greater, and protected against vertically falling water drops when enclosure tilted up to 15°.
-  Please pay attention to the effects of lint, dust, light (including sunlight), etc.

Declaration of Conformity

The manufacturer hereby declares that this device complies with the following standards:

IEC 60601-1: 2005 Medical electrical equipment-Part 1: General requirements for basic safety and essential performance;

ISO 80601-2-61:2017 - Medical electrical equipment --

Part 2-61: Particular requirements for basic safety and essential performance of pulse Oximeter equipment. And it also follows the provisions of the council directive MDD 93/42/EEC.

Caution: U.S. federal law restricts this device to sale or use by or on the order of a physician.

Table of Contents

1 OVERVIEW	1
1.1 APPEARANCE	1
1.2 NAME AND MODEL	2
1.3 INTENDED USE	2
1.4 FEATURE LIST	2
2 BATTERY INSTALLATION	5
3 OPERATION	6
4 TECHNICAL SPECIFICATIONS	11
5 PACKING LIST	13
6 REPAIR AND MAINTENANCE	14
6.1 MAINTENANCE	14
6.2 CLEANING AND DISINFECTING INSTRUCTION	16
7 TROUBLESHOOTING	16
8 KEY OF SYMBOLS	17
APPENDIX I EMC	19

1 Overview

1.1 Appearance


Figure 1 Front View


Figure 2 Rear View

1.2 Name and Model

Name: Fingertip Oximeter

Model: Oxy-One Plus

1.3 Intended Use

This Fingertip Oximeter is intended for measuring the pulse rate and functional oxygen saturation (SpO_2) through a patient's finger. It is applicable for spot-checking SpO_2 and pulse rate of adult and pediatric patients in homes and medical clinics.

1.4 Feature List

Explanation of abbreviations:

<i>Model</i>	<i>Oxy-One plus</i>
<i>Function</i>	
<i>Display type</i>	<i>OLED</i>
<i>SpO₂, PR, PI</i>	✓
<i>Plethysmogram</i>	✓
<i>Auto on/off</i>	✓
<i>Pulse bar graph</i>	✓

<i>PR and PI shifts</i>	√
<i>Over-limits indication</i>	√
<i>Over-limits indication sound</i>	√
<i>Indication sound mute</i>	<i>Manual</i>
<i>Four directions display</i>	√
<i>Setting menu</i>	√
<i>Over-limits setting</i>	√
<i>Pulse beep</i>	√
<i>Measuring mode</i>	<i>Continuous, Spot check</i>
<i>Record list</i>	√

Description:

➤ ***Indication sound mute***

For the model with both over-limits indication function and pulse beep function, when beep is on and over-limits indication sound is activated, then Display key will work as the Mute key, and short time pressing it can mute the over-limits indication sound and pulse beep for 90 seconds.

➤ **Measuring mode**

① *Spot check mode: the measurement starts automatically when the finger is inserted into the finger clip properly, the measuring time lasts 30 seconds with a counting-down indication. The SpO₂ and PR readings will freeze at the end of 30 seconds, the analysis result for the pulse rhythm will be displayed on the screen as well. Once the finger is out, the display will be cleared and the Oximeter shuts down automatically.*

② *Continuous mode: measurement will start automatically when finger is inserted into the oximeter, SpO₂ and PR readings will be displayed until the finger is removed from the oximeter. Once finger is removed the Oximeter will automatically turn off.*

➤ **Record list**

① *A single group of stable readings will be recorded in the record list each time when the Oximeter shuts down regardless of spot-check or continuous mode. However, if the time from displaying valid readings to the end of measurement is less than 5 seconds, then*

no recording will be done.

② Up to 12 groups of records can be stored in the record list, the newest record is marked as M1, and the oldest record is marked as M12. The new record will override the previous record.

③ When batteries are removed from the device all readings will be deleted.

④ On power off status, long pressing the Display key brings up the record recall screen. On record recall screen, short time pressing the Display key can shift the records display, and if there is no key operation for 6 seconds, then the Oximeter will power off automatically again.

2 Battery Installation

1. Refer to Figure 3, insert two AAA size batteries into the battery compartment properly, and note the polarity markings.


Figure 3 Battery Installation

2. Replace the cover.

- Please make sure that the batteries are correctly installed. Incorrect installation may cause the device not to work.
- Please remove batteries if the device is not being used for more than 7 days to prevent and avoid potential damage from the battery leaking. Any such damage is not covered under the product warranty.

3 Operation

1. Open the clip and put finger inside the

rubber cushions of the clip (make sure the finger is in the correct position), and then clip the finger, as shown in figure 4.

Figure 4

2. Next enter into data display screen:

- ✧ *A short time press of the Display Key can change display direction, the four display directions are as shown in figure 5A, 5B, 5C, & 5D. For display screens of figure 5B and 5D, the PI% display value*

will be replaced with PR display value after 20 seconds if no key operation.

Figure 5A


Figure 5B

Figure 5C

Figure 5D

✧ Icon “ 25 ” on display screen means the counting-down time if the Oximeter works at Spot check mode. The total measuring time for Spot check mode is 30 seconds.

No irregularity found
%SpO ₂ :99 PR:78

S:	98	99	98	97
P:	68	77	82	75
	M1	M2	M3	M4

Figure 5E

- ✧ When the measurement ends up for Spot check mode, the measured SpO₂, PR value and the analysis result of pulse rhythm will be displayed on the screen, as shown in figure 5E.

Figure 6

Other result descriptions see Appendix I.

Menu

Long time pressing display key can enter the setup menu screen.

SpO ₂ alm Lo	89
PR alm Hi	100
PR alm Lo	30
Setting menu >>	

Mode	Continuous
Beep	On
Exit	
<< Setting menu	

“Beep” : Pulse beep option. If it is set to on, every pulse beat makes a beep.

“Mode”: to set the measuring mode. “Continuous” and “Spot check” for optional, the default is “Spot check”.

Menu setup

Short time press Display Key to choose the setting item; Longtime press Display Key to active the setting item, then short time press it to modify the setting parameter; Next, longtime press Display Key to confirm the modification and exit from this setting item. At last, move the setting item to “Save, exit menu”, and long time pressing Display Key to store the modification and exit from the setup menu.

Attention to the operation

- The finger should be put into the sensor correctly.
- Do not shake the finger and relax during measurement.
- Do not put wet finger directly into sensor.
- Avoid placing the device on the same limb which is wrapped with a cuff for blood pressure measurement or during venous infusion.
- Do not let anything block the emitting light from device, i.e. do not use finger nail polish/paints.
- Vigorous exercise and electrosurgical device interference may affect the measuring accuracy.
- The orientation-sensor works on the basis of the gravity. A small movable metal ball is built in the orientation-sensor for detecting the orientation of

the Oximeter. When you want to change the Oximeter's display direction, if you move the Oximeter too slowly, the movable metal ball will also move slowly because of not enough acceleration. Consequently the response of orientation detection would be delayed. Acceleration needs to be provided to the orientation-sensor for quick sensing the orientation change.

- Nail polish may affect the measuring accuracy, and too long fingernail may cause failure of measurement or inaccurate result.
- Existence of high intensive light sources, such as fluorescence light, ruby lamp, infrared heater or strong sunshine, etc. may cause inaccuracy of measurement result. Please put an opaque cover on the sensor or change the measuring site if necessary.
- If the first reading appears with poor waveform (irregular or not smooth), then the reading is unlikely true, the more stable value is expected by waiting for a while, or a restart is needed when necessary.

4 Technical Specifications

A. SpO₂ Measurement

Transducer: dual-wavelength LED sensor with wavelength:

Red light: 663 nm, Infrared light: 890 nm.

Maximal average optical output power: $\leq 2\text{mW}$

SpO₂ display range: 35%~100%

SpO₂ measuring accuracy:

$\leq 2\%$ for SpO₂ range from 70% to 100%

B. Pulse Rate measurement

PR display range: 30bpm~240bpm

PR measuring accuracy: $\pm 2\text{bpm}$ or $\pm 2\%$
(whichever is greater)

F. Perfusion Index(PI) Display range

0%~20%

G. Preset over-limits

SpO₂ low limit: 90%

Pulse Rate: high limit: 120bpm
low limit: 50bpm

H. SpO₂ Over-limit settings

Low limit setting range: 85%~99%, step: 1%

Default setting: 90%

Pulse Rate:

Low limit setting range: 30~60bpm, step: 1bpm;

High limit setting range: 100~240bpm, step: 5bpm;

Default setting: high: 120bpm; low: 50bpm

I. Audible & visual alert function

When measuring, if SpO₂ value or pulse rate value exceeds the preset limit, the device will alert with beep automatically and the value which exceeds limit will flash on the screen.

H. Power supply requirement:

2 x LR03 (AAA) alkaline batteries

Supply voltage: 3.0VDC

Operating current: ≤40mA

N. Environmental Conditions:

Operating Temperature: 5°C ~40°C

Operating Humidity: 30%~80%

Atmospheric pressure: 70kPa~106kPa

O. Low Perfusion Performance:

The accuracy of SpO₂ and PR measurement still meet the precision described above when the modulation amplitude is as low as 0.6%.

P. Ambient Light Interference:

The difference between the SpO₂ value measured in the condition of indoor natural light and that of darkroom is less than ±1%.

Q. Dimensions: 56 mm (L) × 34 mm (W) × 30 mm (H)

Net Weight: approx. 60g (including batteries)

R. Classification

The type of protection against electric shock:
Internally powered equipment.

The degree of protection against electric shock:
Type BF applied parts.

The degree of protection against harmful solid foreign objects and ingress of liquid:

The equipment is IP22 with protection against harmful solid foreign objects and ingress of liquid.

Electro-Magnetic Compatibility: Group I, Class B

5 Packing List

- 1) Fingertip Oximeter

- 2) User Manual
- 3) Batteries
- 4) Pouch (optional)
- 5) Lanyard (optional)

Note: the items and its quantity are subject to change, please refer to your subject in hand.

6 Repair and Maintenance

6.1 Maintenance

The expected service life (not a warranty) of this device is 5 years. In order to ensure its long service life, please pay attention to the maintenance.

- Please change the batteries when the low-voltage indicator lightens.
- Please clean the surface of the device before using, with 75% alcohol wipes, then let it air dry or wipe it dry. Do not allow liquid to enter the device.
- Please take out the batteries if the Oximeter will not be used any more than 7 days.
- The recommended storage environment of the device:

ambient temperature: -20°C ~60°C, relative

humidity 10%~95%, atmospheric pressure:
50kPa~107.4kPa.

- The Oximeter is calibrated in the factory before sale, so there is no need to calibrate it during its life cycle. Any SpO₂ simulators should not be used to validate the accuracy of the Oximeter, they can only be used as functional testers to verify its precision. The SpO₂ accuracy claimed in this manual is supported by the clinical study conducted by inducing hypoxia on healthy, non-smoking, light-to-dark skinned subjects in an independent research laboratory.
- If it is necessary to verify the precision of the Oximeter routinely, the user can do the verification by means of SpO₂ simulator, or it can be done by the local third party test house. Please note that the specific calibration curve (so called R-curve) should be selected when use of SpO₂ simulator, e.g. for Index 2 series SpO₂ simulator from Fluke Biomedical Corporation, please set "Make" to "DownLoadMake: KRK", then the user can use this particular R-curve to test the Oximeter. If the SpO₂ simulator does not contain "KRK" R-curve, please

ask the manufacturer for helping to download the given R-curve into the SpO₂ simulator.

- ⚠ **High-pressure sterilization cannot be used on the device.**
- ⚠ **Do not immerse the device in liquid.**
- ⚠ **It is recommended that the device should be kept in a dry environment. Humidity may reduce the life of the device, or even damage it.**

6.2 Cleaning and Disinfecting Instruction

- Surface-clean sensor with a soft cloth damped with a solution such as 75% isopropyl alcohol, if low-level disinfection is required, use a mild bleach solution.
- Then surface-clean with a cloth damped **ONLY** with clean water and dry with a clean, soft cloth.

Caution: Do not sterilize by irradiation steam, or ethylene oxide.

Do not use the Oximeter if it is damaged.

7 Troubleshooting

Problem:

1. The SpO₂ and Pulse Rate display instable
2. Cannot turn on the device

3. No display
4. Display direction doesn't change or changes insensitively.

Solution

7. Place the finger correctly inside and try again.
8. Changing batteries.
9. Let the patient keep calm.
10. Please shake the Oximeter with a certain force to make the movable metal ball move freely. If the problem still exists, maybe the orientation-sensor is not working properly.
11. Hardware failure of wireless transmission function.
12. If the above problem still exists please contact the local service center.

8 Key of Symbols

Symbol	Description
%SpO₂	Pulse oxygen saturation
♥BPM/PR	Pulse rate (beats per minute)
PI%	Perfusion Index (%)

Pulse Strength Bar Graph

Low battery voltage

CE mark

SN

Serial number

Date of manufacture

EC REP

Authorised representative in the European community

Manufacturer (including address)

BF type applied part

Attention – refer to User Manual

Follow WEEE regulations for disposal

Appendix I EMC

The equipment meets the requirements of IEC 60601-1-2:2014.

Table 1

Guidance and manufacturer's declaration-electromagnetic emission		
The Fingertip Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Oximeter should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The Fingertip Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The Fingertip Oximeter suitable for use in all establishments, including domestic establishments and those directly network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC61000-3-3	N/A	

Table 2

Guidance and manufacturer's declaration-electromagnetic emission			
The Fingertip Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Oximeter should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance

<p>Electrostatic discharge(ESD) IEC61000-4-2</p>	<p>± 8 kV contact ± 15kV air</p>	<p>± 8 kV contact ± 15kV air</p>	<p>Floors should be wood, concrete or ceramic tile. if floors are covered with synthetic material, the relative humidity should be at least 30%</p>
<p>Electrical fast transient/ burst IEC61000-4-4</p>	<p>± 2kV for power Supply lines ± 1 kV for input/output lines</p>	<p>N/A</p>	<p>N/A</p>
<p>Surge IEC 61000-4-5</p>	<p>± 1kV line (s) to line(s) ± 2kV line(s) to earth</p>	<p>N/A</p>	<p>N/A</p>

<p>Voltage dips, short interruptions and voltage variations on power supply input lines</p> <p>IEC61000-4-1 1</p>	<p><5% U_T (>95% dip in U_T) for 0.5 cycle</p> <p><40% U_T (60% dip in U_T) for 5 cycles</p> <p><70% U_T (30% dip in U_T) for 25 cycles</p> <p><5% U_T (>95% dip in U_T) for 5 s</p>	<p>N/A</p>	<p>N/A</p>
<p>Power frequency(50Hz/60Hz) magnetic field</p> <p>IEC61000-4-8</p>	<p>3A/m</p>	<p>3A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Table 3

**Guidance and manufacturer's declaration –
electromagnetic immunity**

The Fingertip Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of The Fingertip Oximeter should assure that it is used in such an electromagnetic environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment -guidance
Conducted	3 Vrms		Portable and mobile RF communications equipment should be used no closer to any part of The Fingertip Oximeter,

<p>RF IEC61000-4-6</p>	<p>150 kHz to 80 MHz</p>	<p>N/A</p>	<p>including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p>
<p>Radiated RF IEC61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	<p>Recommended separation distance $d=1.2 \sqrt{P}$ $d=1.2 \sqrt{P}$ 80MHz to 800MHz $d=2.3 \sqrt{P}$ 800MHz to 2.5GHz Where P is the maximum output power rating of</p>

			<p>the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p>
--	--	--	--

			Interference may occur in the vicinity of equipment marked with the following symbol.
--	--	--	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a: Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should

be considered. If the measured field strength in the location in which The Fingertip Oximeter is used exceeds the applicable RF compliance level above, The Fingertip Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating The Fingertip Oximeter.

b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 4

Recommended separation distances between portable and mobile RF communication the equipment

The Fingertip Oximeter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of The Fingertip Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Fingertip Oximeter as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W(Watts)	Separation distance according to frequency of transmitter M(Meters)		
	150kHz to 80MHz	80MHz to 800MHz	80MHz to 2,5GHz
	$d=1.2 \sqrt{P}$	$d=1.2 \sqrt{P}$	$d=2.3 \sqrt{P}$
0,01	N/A	0.12	0.23
0,1	N/A	0.38	0.73
1	N/A	1.2	2.3

10	N/A	3.8	7.3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 : At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Quality Certificate

Name: Fingertip Oximeter

Model: _____

Date: _____

QA: _____

This product has been inspected in accordance with the standards specified in the User Manual.

CE 0123

**Shenzhen Creative Industry Co., Ltd.
Floor 5, BLD 9, Baiwangxin High-Tech
Industrial Park, Songbai Road, Xili Street,
Nanshan District, 518110 Shenzhen, P. R.
China**

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