


STUDY SUMMARY

Use of the *AeroEclipse*® Brand of
Breath Actuated Nebulizer

CONFIDENCE IN AEROSOL DELIVERY



Trudell Medical International*

AeroEclipse* Breath Actuated Nebulizer (BAN) and AeroEclipse* II BAN

Study Summary

Table of Contents

| | |
|---|-----------|
| AeroEclipse* BAN Equivalence to AeroEclipse* II BAN | 3 |
| Study Summary by Active Pharmaceutical Ingredient | 6 |
| Ventolin™ (salbutamol sulfate/albuterol), GSK™ Inc. | 6 |
| Pulmicort™ (budesonide), AstraZeneca™ | 16 |
| Combivent™ (ipratropium bromide & albuterol sulfate), Boehringer Ingelheim™ Pharmaceuticals Inc. | 18 |
| Atrovent™ (ipratropium bromide), Boehringer Ingelheim™ Pharmaceuticals Inc. | 19 |
| Alupent™ (metaproterenol sulphate), Boehringer Ingelheim™ Pharmaceuticals Inc. | 20 |
| Intal™ (cromolyn sodium), Fisons™ Pharmaceuticals Ltd. | 20 |
| Methacholine chloride..... | 20 |
| Ablecet™ (Amphotericin), Enzon™ Pharmaceuticals | 21 |
| Measles Vaccine (Placebo) | 23 |
| Interferon Gamma..... | 23 |
| Saline..... | 23 |
| Xopenex™ (levalbuterol), Sepracor™ Inc. | 24 |
| Actiq® (fentanyl citrate), Abbott Laboratories..... | 26 |
| AeroLEF™ (Liposome-Encapsulated Fentanyl), YM Biosciences, Inc. | 26 |
| Study Summary by Nebulizer | 28 |
| Ventstream™, Medic-Aid™ | 28 |
| UpDraft Neb-U-Mist™, Hudson Oxygen Therapy Sales Co..... | 28 |
| Updraft-II Nebu-U-Mist™, Hudson Oxygen Therapy Sales Co. | 30 |
| Micromist™, Hudson RCI™ | 31 |
| VixOne™, Westmed™ Corp..... | 32 |
| Airlife™ Misty-Neb™, Baxter™ Healthcare Corp..... | 33 |
| Misty Max 10™, Cardinal Health | 38 |
| Pari LC-JET™, PARI™ Respiratory Products Inc. | 38 |
| Pari LC-Star™, PARI™ Respiratory Equipment..... | 39 |
| Pari LC Plus™, PARI™ Respiratory Equipment..... | 40 |
| Pari LC-D™, PARI™ Respiratory Equipment | 40 |
| Sidestream™, MedicAid™ | 43 |
| Circulaire™, Westmed™ Corp..... | 44 |
| Salter HDN™, Salter Labs™ | 48 |
| Salter Model 8900™, Salter Labs™ | 50 |
| PediNeb™, Westmed™ Corp. | 50 |
| AeroChamber* Valved Holding Chamber and MDI | 51 |
| Large Volume Nebulizers | 51 |
| General Information | 53 |

AeroEclipse* BAN Equivalence to AeroEclipse* II BAN

1. **TRANSFER FROM THE MALVERN MASTERSIZER-X TO MALVERN SPRAYTEC LASER DIFFRACTOMETERS: EXPERIENCE WITH TWO BREATH-ACTUATED NEBULIZERS (BAN).** J Mitchell, K Wiersema, C Doyle, M Nagel, P Kippax, H Krarup. Presented at Respiratory Drug Delivery (RDD), Boca Raton, FL, 2006.

Introduction: Laser diffraction is widely used for the measurement of droplet sizes of aqueous solution aerosols from nebulizers on account of its rapidity and size resolution capability (1), and is indicated in an Informative Annex of a European standard for the evaluation of this class of inhalers (2). The second generation Malvern Spraytec laser diffractometer (LD) (Malvern Instruments Ltd., Malvern, UK) has recently been introduced for the purpose of size-characterizing aerosols and droplet sprays, replacing earlier instruments. We describe our recent experience transferring from a Mastersizer-X LD to the Spraytec LD at the same time as bringing a second-generation breath-actuated nebulizer (**AeroEclipse* II BAN**, Trudell Medical International, London, Ontario, Canada) to market.

TRANSFER FROM MASTERSIZER-X TO SPRAYTEC LD SYSTEMS

In the first part of the study, we compared droplet size distributions of normal saline (0.9% w/v NaCl, 5 mL fill) determined by Mastersizer-X and Spraytec LDs, using first generation **AeroEclipse*** BANs (n=3 devices, 2 measurement per device) operated at 7 to 8 L/min by compressed air supplied at 345 kPa (50 psi). The complex refractive index (RI) for saline was defined as $1.33 + 0i$, with air (RI = 1.00) as support medium. Measurements were made with the Mastersizer LD in the open bench configuration with a 100-mm focal length range lens, delivering an additional flow of 20 L/min through the cap of the nebulizer containing the air entrainment entry passages to move the droplets through the measurement zone without risk of recirculation. In contrast, the aerosol from the nebulizer was drawn via the inhalation cell of the Spraytec (300-mm range lens) at 28 L/min using an external vacuum source. This arrangement is more representative of the process of inhalation.




Figure 1. LD-measured size distributions from the **AeroEclipse*** BAN.

The cumulative volume (mass)-weighted size distributions (Figure 1) were comparable for droplets > 3 µm, so that the Mastersizer-X-determined fine droplet fraction < 4.8 µm ($84.0 \pm 1.2\%$ (mean \pm SD)) compared with $83.5 \pm 1.9\%$ < 4.6 µm for the Spraytec system. The cause of the 'tail' of fine droplets present in the Mastersizer data requires further investigation. Preliminary studies suggest that the cause was not multiple scattering, even though obscurations in excess of 25% were obtained. It may, however, be associated with the way the aerosol was transported to the measurement zone and the working range of the optical system. Here the Spraytec offers advantages over the Mastersizer-X in that the working range is 150-mm compared with 2.4-mm. The angular range of the scattering measurements made using the Spraytec is also greater than for the Mastersizer-X so that the former would be expected to provide a more accurate measure of the fine particle fraction.

FIRST AND SECOND GENERATION BAN COMPARISON

In the second part of the study we compared saline droplet size distributions from the original **AeroEclipse*** BAN with those produced by a second generation BAN (**AeroEclipse* II**) designed to improve actuation capability for low inhalation flow rate patients. 5 nebulizers of each type were evaluated, with the Spraytec system configured as described in the first part of the investigation. The entire size distribution profiles from the two nebulizer types were substantially similar (Figure 2), so that the fine droplet fraction < 4.6 µm from the

AeroEclipse[®] BAN ($85.2 \pm 1.5\%$) compared with $80.7 \pm 2.7\%$ for the second generation nebulizer. In both cases, the volume (mass) median diameter was 2.5 to 2.7 μm .




Figure 2. Spraytec LD-measured size distributions for BANs.

These measurements were made with only one solution (saline), and further work with other solution formulations is therefore merited.

2. **ARE FIRST AND SECOND GENERATION, MECHANICALLY-OPERATED BREATH-ACTUATED NEBULIZERS (BAN) COMPARABLE BASED ON *IN-VITRO* PERFORMANCE?** James Schmidt, Jennifer Pevler, Cathy Doyle, Kimberly Wiersema, Mark Nagel, and Jolyon Mitchell. Presented at Respiratory Drug Delivery (RDD), Boca Raton, FL, 2006.

Introduction: The original **AeroEclipse**[®] nebulizer (Monaghan Medical Corp., Plattsburgh, NY) introduced a few years ago was the first mechanically-operated BAN with dosimetric capability, providing a near constant delivery rate of medication from a variety of solution formulations and volume fills (1). This nebulizer required an inhalation flow rate close to 25 L/min to operate the breath-actuation mechanism. The second generation **AeroEclipse**[®] II BAN now actuates at flow rates as low as 15 L/min, making it potentially more suitable for younger patients. At the same time, a control located on the nebulizer cap enables a smooth transition to be made from breath-actuated to continuous operation. We report a study in which the delivery of albuterol sulfate solution from the new BAN was evaluated with a 3 mL fill, corresponding to a single unit dose ampoule (0.83 mg/mL albuterol sulfate) in widespread use within the US (1). Previously published data for the original BAN (1) were used as a benchmark for demonstrating *in vitro* equivalence. The study was extended to examine comparative behavior with a low volume (1 mL) fill, used to reduce treatment time. **Materials and Methods:** In the first part, we evaluated 5 **AeroEclipse**[®] II nebulizers (n=3 replicates/device) using a piston-driven breathing simulator (Compas[®], PARI GmbH, Starnberg, Germany) set at tidal volume of 600-mL, inspiratory/expiratory ratio of 1:2, rate of 10 breaths/minute, based on a previous study simulating adult use (2). Each nebulizer was operated at 8.0 ± 0.2 L/min with compressed air supplied at 50 ± 0.5 psig. 3 mL albuterol solution obtained by diluting respirator solution (5 mg/mL albuterol base equivalent, Hi-Tech Pharmacal, Amityville, NY) with normal saline to the desired concentration (0.83 mg/mL) was placed in the reservoir of the nebulizer prior to test. The measurement protocol to determine the total mass of drug delivered on a minute-by-minute basis was as described previously (1). Fine droplet fraction < 4.8 μm diameter (FDF<4.8 μm) was also determined by laser diffraction (Mastersizer-X, Malvern Instruments plc, UK) as described previously (1). At each minute, the mass of drug delivered as fine particles was calculated as the product of total mass and the mean (FDF<4.8 μm). Measurements were made at comparable conditions ($22 \pm 2^\circ\text{C}$, $30 \pm 5\%$ RH) to those of the original study. In the second part, we followed the same protocol, except that the fill volume was decreased to 1 mL, diluting respirator solution with normal saline to achieve an albuterol concentration of 2.5 mg/mL. The delivery rate of fine droplets from the BAN was compared with that produced by the LC PLUS[®] (PARI Respiratory Equipment Inc.), chosen as a benchmark high output, continuous breath-enhanced nebulizer. **Results:** Comparable fine droplet delivery with both the original and new BAN was achieved throughout the 10 min. delivery period (Figure 1).




Figure 1. Comparative delivery of albuterol solution (0.83 mg/mL) with 3-mL fill in reservoir.

Mean FDF<4.8 μ m for both nebulizers was within $80 \pm 2\%$. The rate of delivery of albuterol was constant, as might be expected for a solution formulation. The cumulative mass of fine droplets from the new BAN by the time that audible sputtering occurred was 842 ± 50 μ g compared with 810 ± 34 μ g for the original BAN. In the case of the measurements made with the 1 mL fill (2.5 mg/mL albuterol), the new BAN operated for about 3 minutes before sputtering, delivering 544 ± 54 μ g albuterol as fine droplets, in comparison with 576 ± 49 μ g in a similar time from the original BAN. In contrast, only 67 ± 10 μ g of albuterol was obtained as fine droplets from the LC PLUS[®] (mean FDF<4.8 μ m also ~80%), which operated for just over 1 minute before sputtering. The LC-Plus[®] operated throughout each breathing cycle, reducing delivery time, but medication emitted during exhalation was not collected since it would be wasted in normal use. **Conclusions:** The **AeroEclipse® II BAN** has similar *in vitro* performance with albuterol as the original version, and treatment time can be significantly shortened by reducing the volume fill to 1 mL. The breath-actuation feature avoids the escape and therefore waste of medication during patient exhalation, with attendant concerns concerning possible exposure of the care-giver to medication. These considerations could be important when used with more expensive medications.

Study Summary by Active Pharmaceutical Ingredient


Ventolin™ (salbutamol sulfate/albuterol), GSK™ Inc.

1. **IN VITRO PERFORMANCE COMPARISON OF A BREATH-ACTUATED NEBULIZER (BAN) FOR THE DELIVERY OF ALBUTEROL OPERATED WITH COMPRESSED HELIOX OR AIR.** D Coppolo, J Mitchell, V Avvakoumova, M Nagel. Presented at the American Council of Clinical Pharmacy (ACCP) Annual Meeting, Philadelphia, PA, 2008.

Purpose: The NAEPP Guidelines for the Diagnosis and Management of Asthma were revised in 2007 to include the use of Heliox (21%v/v oxygen/79%v/v helium) for treatment of severe exacerbations that are unresponsive to initial treatments. We report data for delivery of a beta-2 adrenergic agonist by BAN as guidance to clinicians. **Methods:** **AeroEclipse® II** BANs (n=5 devices, Monaghan Medical Corp., Plattsburgh, NY) were operated simulating adult tidal breathing (tidal volume = 600-ml, 10 bpm, 33% duty cycle) and delivering 3-ml albuterol (0.83 mg/ml). Each nebulizer was powered at 50 psig by compressed air at 8 L/min (condition A, maximum achievable); Heliox at 8 L/min (condition B); Heliox at 16 L/min (condition C, maximum achievable). Emitted droplets were collected on separate filters at the mouthpiece of the BAN at 1-min intervals and recovered albuterol assayed by HPLC-UV spectrophotometry. The nebulizers were operated until onset of sputtering to determine total emitted mass (TEM). In a parallel study the emitted fine droplet fraction < 4.7 µm diameter obtained at each condition (FDF<4.7µm) was determined by laser diffractometry (n=3 replicates with 1 device). Total fine droplet delivery (FDM<4.7µm) was calculated as the product of TEM and FDF<4.7µm. **Results:** FDF<4.7 µm (mean ± SD) was 78.4 ± 1.8% (condition A); 68.7 ± 2.9% (condition B) and 84.8 ± 3.2% (condition C). The BANs operated for 10-min, 19-min and 11-min with corresponding values of FDM<4.7 µm (mean ± SD) of 90.2 ± 3.3, 28.8 ± 2.0 and 80.3 ± 4.5 µg/min at conditions A, B and C respectively. **Conclusion:** Fine droplet delivery from the BAN can be maintained at a near equivalent delivery rate with Heliox if the flow rate is set to maximum. The reduction in aerosol output if flow rate is unchanged between air and Heliox reflects the lower density of the latter driving gas. **Clinical Implication:** Clinicians should be mindful of the need to set the flow rate of Heliox to the BAN at maximum to maintain aerosol delivery characteristics established for air. **Disclosures:** The authors are employees of Trudell Medical Group, manufacturer of the **AeroEclipse® II** BAN.

2. **RAPID DELIVERY OF BRONCHODILATOR MEDICATION IS POSSIBLE USING A BREATH-ACTUATED SMALL VOLUME NEBULIZER AS AN ALTERNATIVE TO EXTENDED DELIVERY OF MEDICATION BY LARGE VOLUME NEBULIZER.** J Mitchell, DCoppolo, C Doyle, MW Nagel, KJ Wiersema. Presented at the American Association for Respiratory Care (AARC) Open Forum, Orlando, FL, 2007.

Background: Inhaled beta-2 adrenergic agonist bronchodilators are often given to patients with severe reversible airways disease by continuous nebulization in extended treatments. However data are limited as to whether or not shorter, but higher concentration delivery is as an effective treatment modality. The development of a new breath-actuated nebulizer (**AeroEclipse® II**, Monaghan Medical Corp., Plattsburgh, NY (AEIIBAN)) provided an opportunity to compare the two treatment methods in a laboratory study before undertaking a clinical comparison. We investigated the delivery of diluted generic respirator solution albuterol by a widely used continuous jet nebulizer (MiniHeart® Hi-Flo, Westmed Corp., Tucson, AZ (CONT)) with that from the AEIIBAN. **Method:** The continuous nebulizers (n=5) were operated with 8 L/min air supplied at 50 psig with a 20-ml fill (albuterol concentration of 0.5 mg/mL). A similar number of AEIIBANs were operated with ca. 8.0 L/min air at 50 psi with a 1-ml fill (albuterol concentration of 5 mg/mL). Aerosol from both nebulizers was sampled onto electret filters using a breathing simulator mimicking small child use (250-ml tidal volume, inspiratory/expiratory ratio 1:2, rate 12 cycles/min) until onset of sputtering. Assay for albuterol was undertaken by UV spectrophotometry. In a parallel study, droplet size distributions were determined by laser diffractometry, so that the fine droplet fraction (mass % < 4.7 µm diameter) likely to penetrate to the airways of the lungs (FDF) could be determined. **Results:** Values of FDF for the AEIIBAN and CONT were 78.4% and 62.0% respectively. The AEIIBAN delivered 758 ± 36 µg as fine droplets after 4-min (delivery rate of 190 ± 9 µg/min), compared to 180 ± 76 µg in the same period by CONT (delivery rate of 45 ± 19 µg/min). **Conclusions:** The faster delivery rate from the AEIIBAN/high albuterol concentration modality (un-paired t-test, p < 0.001) may offer an important clinical alternative to CONT/low concentration treatment modality.



3. **DELIVERY OF ALBUTEROL VIA A NEW BREATH ACTUATED NEBULIZER: COMPARISON WITH CONTINUOUS JET NEBULIZERS.** DP Coppolo, MW Nagel, CC Doyle, VA Avvakoumova, and JP Mitchell. Accepted for presentation at the American Thoracic Society International Conference, San Francisco, CA, 2007.

A new breath actuated nebulizer (**AeroEclipse*** II BAN, Monaghan Medical Corp., Plattsburgh, NY) has been developed to deliver medication only when the patient inhales. This study sought to determine the delivery of albuterol (3-ml fill of diluted solution (0.83 mg/ml)) as fine droplets < 4.7 μ m aerodynamic diameter, and compare this fine droplet mass (FDM) with equivalent data from 4 widely available continuous jet nebulizers as benchmark devices. Each nebulizer (n=5; 3 replicates/device) was operated with compressed air at 50 psig at ca. 8 L/min to simulate hospital wall outlet conditions. The nebulizer on test was coupled to a breathing simulator set to mimic adult use (tidal volume = 600 ml, rate = 10 breaths/min; duty cycle = 0.33), and the emitted droplets were collected on an electret filter at the mouthpiece. The total mass of albuterol (TM) was assayed subsequently by HPLC-UV spectrophotometry. In a separate study, the droplet size distribution was determined by laser diffractometry so that the fine droplet fraction (FDF) could be obtained. FDM was determined as the product of TM and FDF. FDM (mean SD) from the BAN was 791 84 g, delivered in 8 minutes. Corresponding values (FDM in time from start to sputter) for the VixOne™ (Westmed, Tucson, AZ), MicroMist™ (Hudson RCI, Temecula CA), Misty Max 10™ (Cardinal Health, McGaw Park (IL) and model 8900™ (Salter Labs, Arvin, CA) were 267 11 g in 4 min, 133 8 g in 4 min, 249 10 g in 6 min and 161 10 g in 5 min. Aside from dosage assurance imparted by breath-actuation, the **AeroEclipse*** II BAN delivered substantially more FDM/min than the other devices. The clinician is now able to treat either for extended high dose delivery (potentially eliminating the need for additional therapy), or titrate to a shorter interval based on response.

4. **A BREATH-ACTUATED SMALL VOLUME NEBULIZER (BAN) OFFERS A RAPID ALTERNATIVE TREATMENT MODALITY FOR THE DELIVERY OF BRONCHODILATORS FOR ASTHMATIC PATIENTS IN A SEVERE EXACERBATION.** DP Coppolo, CC Doyle, JP Mitchell, MW Nagel, KJ Wiersema. American Association for Respiratory Care (AARC), Las Vegas, NV, 2006.

Large volume continuous nebulizers (LVNs) are often used for the delivery of beta-2 adrenergic agonist bronchodilators in the emergency department to treat severe, reversible airways disease, in particular asthma 1. Treatment time, however, can be lengthy for delivery of the typical LVN fill volume from 20- to 120-ml. Quick delivery of a bronchodilator with an efficient nebulizer may help relieve symptoms from bronchospasm in a shorter period of time. We report a study in which the delivery of diluted generic respirator solution albuterol by LVN (Hope, B&B Medical Technologies Inc., Loomis, CA) was compared with that from a small volume breath-actuated nebulizer (BAN) (**AeroEclipse***, Monaghan Medical Corp., Plattsburgh, NY). The LVNs (n=5) were operated with 10 L/min air supplied at 50 psig with a 20-ml fill (albuterol concentration of 0.167 mg/ml). A similar number of BANs were operated with 8.0 L/min air at 50 psi with a 3-ml fill (albuterol concentration of 0.833 mg/ml). The aerosol from the LVNs was sampled continuously until onset of sputtering at 12 L/min via a Dreschel filter/bottle where the albuterol was captured quantitatively. Aerosol from the BANs was sampled onto electret filters using a breathing simulator (600-ml tidal volume, inspiratory/expiratory ratio 1:2, rate 10 cycles/min) until onset of sputtering, so that operation of the breath actuation mechanism was effected. Assay for albuterol was undertaken by UV spectrophotometry. In a parallel study droplet size distributions were determined by laser diffractometry, so that the fine droplet fraction < 4.8 μ m diameter likely to penetrate to the airways of the lungs could be determined. Fine droplet albuterol delivery rates were constant as a function of time for all nebulizers. After 15-min, the LVNs had supplied 127.3 \pm 37.4 μ g as fine droplets at a rate of 8.5 \pm 2.5 μ g/min. In contrast, the BANs delivered 810.0 \pm 20.4 μ g in a 10-min period, equivalent to a rate of 81.0 \pm 2.0 μ g/min. The significantly higher delivery rate from the BAN group (un-paired t-test, p < 0.001) offers an important clinical alternative to the LVN in the emergency department where rapid delivery of a bronchodilator is critical.

Reference: M McPeck, R Tandon, K Hughes, GC Smaldone. *Aerosol delivery during continuous nebulization*. CHEST 1997;111:1200-1205.

5. **A RANDOMIZED CONTROLLED TRIAL COMPARING A BREATH ACTIVATED NEBULIZER TO STANDARD INTERMITTENT AND ONE-HOUR CONTINUOUS ALBUTEROL IN THE TREATMENT OF EMERGENCY ROOM PEDIATRIC ASTHMA.** Katie Sabato MS RRT, Priscilla Ward RRT, William Hawk MD, Jeanette Asselin RRT MS, Children's Hospital and Research Center at Oakland. Respiratory Care Journal. November 2005; 50(11):1489.

Background: Bronchodilator treatments for asthma can be provided by a various number of aerosol generating devices and methods. To date, there are few large randomized, controlled trials comparing the efficacy, effectiveness and safety of undiluted and continuous diluted administration of albuterol in the treatment of pediatric asthma. Data are also limited on whether certain nebulizers and their masks are more effective than others and whether blow-by treatments area at all effective. Children's Hospital and Research Center at Oakland (CHRCO) Respiratory Care Department is currently conducting a large randomized controlled study comparing the efficacy of a one-time treatment with the **AeroEclipse*** breath actuated small volume nebulizer (BA SVN) used with mask or mouthpiece, to a one-time treatment with a standard small volume nebulizer (SSVN) or a one-hour continuous treatment (CONT) for asthmatics presenting to the emergency room (ER). **Methods:** Patients were eligible for inclusion if they were admitted to the ER for respiratory distress, were between 0 months to 18 years of age, and had wheezing or status asthmaticus. Patients were objectively assessed utilizing a CHRCO designed clinical asthma score (CAS) and peak flows when possible. The CAS scores clinical wheezing on a scale from 0 to 11, with 11 representing the most severe distress. Patients were stratified by CAS score (CAS < 4 and > 4) and weight (< 20 kg and > 20kg). Patients were randomized to receive their first bronchodilator treatment in the ER via the BA SVN or standard therapy (CONT or SSVN). Bronchodilator doses for the BA SVN and SSVN were: 0.5cc (2.5 mg) Albuterol in 0.5cc normal saline for patients < 20 kg, and 1cc (5.0 mg) undiluted Albuterol for patients > 20 kg. Bronchodilators given via the CONT method used 2.0cc (10 mg) Albuterol in 18cc normal saline. Patients were evaluated at baseline and again 10 minutes after completion of the assigned treatment. Primary endpoints include change in CAS pre/post treatment, need for additional bronchodilator treatments, and time spent in the emergency room. Secondly, we evaluated the ability of infants to breath activate the BA SVN, the effectiveness of different aerosol

interface adapters (patients utilizing the mouthpiece, vented and non-vented aerosolized masks versus blow-by administration), and incidence of side effects documented with each of the approaches. **Results:** Between 10/14/04 and 11/11/05, we enrolled 151 patients into the study. 2 patients were dropped due to consent issues. The remaining 149 represented 90 male and 59 female patients with an average age of 5.5 years. 84 patients were randomized to the BAN and 65 were randomized to CONT/SSVN (57 CONT and 8 SSVN). There were no differences in demographics between the groups. Initial CAS scores were 5.3 and 5.2 for the BAN and CONT/SSVN groups respectively. After treatment, the BAN group showed significant improvement in their CAS (38% vs 24%, $p < 0.003$), and the number of patient requiring admission (31 vs 37, $p = 0.03$). Other than a significant decrease in respiratory rate in the BAN group (-3.9 vs 0.5 , $p = 0.002$), there were no differences in side effects. **Conclusions:** Use of the Monaghan breath-actuated **AeroEclipse*** nebulizer resulted in significant improvements in CAS ($p < 0.003$), need for admission ($p = 0.03$), and decrease in respiratory rate ($p = 0.002$) as compared to our standard treatments (CONT/SSVN). 66% of the BAN patients were able to breath-activate their treatment. We contend that the Monaghan **AeroEclipse*** is a safe and effective nebulizer for the administration of bronchodilator aerosols in pediatrics and may be more effective than continuous aerosols in the treatment of Emergency Room pediatric asthma.

6. **PERFORMANCE COMPARISON OF NEBULIZER DESIGNS: CONSTANT-OUTPUT, BREATH-ENHANCED, AND DOSIMETRIC.** Joseph L Rau PhD RRT FAARC, Arzu Ari MSc CRT CPFT, and Ruben D Restrepo MD RRT. Respiratory Care 2004; 49(2): 174-179.

Introduction: Design differences among pneumatically powered, small-volume nebulizers affect drug disposition (percentage of the dose delivered to the patient, lost to deposition in the equipment, and lost via exhalation to ambient air) and thus affect drug availability and efficacy. **Objective:** Evaluate *in vitro* the dose disposition with 5 nebulizer models, of 3 types (constant-output, breath-enhanced, and dosimetric), using simulated normal, adult breathing. **Methods:** We compared 5 nebulizer models: 2 constant-output (Misty-Neb and SideStream), 1 breath-enhanced (Pari LCD), and 2 dosimetric (Circulaire and **AeroEclipse***). Each nebulizer was filled with a 3-mL unit-dose of albuterol sulfate and powered by oxygen at 8 L/min. The nebulizers were connected to an induction throat, connected to a breathing simulator. We measured (1) inhaled drug (subdivided into mass deposited in the induction throat and mass deposited in the filter at the distal end of the induction throat), (2) exhaled drug (lost to ambient air), (3) drug lost to deposition in the apparatus, and (4) drug left in the unit-dose bottle. The duration of nebulization (until sputter) was measured with a stopwatch. All drug amounts were analyzed via spectrophotometry and expressed as a percentage of the total dose. **Results:** The mean \pm SD inhaled drug percentages were: Misty-Neb $17.2 \pm 0.4\%$, SideStream $15.8 \pm 2.8\%$, Pari LCD $15.2 \pm 4.2\%$, Circulaire $8.7 \pm 1.0\%$, and **AeroEclipse*** $38.7 \pm 1.3\%$. The mean \pm SD percentages of drug lost to ambient air were: Misty-Neb $26.8 \pm 0.7\%$, SideStream $17.3 \pm 0.4\%$, Pari LCD $18.3 \pm 0.8\%$, Circulaire $12.3 \pm 0.8\%$, and **AeroEclipse*** $6.6 \pm 3.3\%$. The mean \pm SD percentages of drug lost to deposition in the apparatus were: Misty-Neb $52.3 \pm 0.6\%$, SideStream $63.4 \pm 3.0\%$, Pari LCD $62.5 \pm 4.0\%$, Circulaire $75.8 \pm 0.5\%$, and **AeroEclipse*** $51.0 \pm 2.1\%$. Duration of nebulization was shortest with the Circulaire and longest with the **AeroEclipse*** ($p < 0.05$ via 1-way analysis of variance). **Conclusions:** The nebulizers we tested differ significantly in overall drug disposition. The dosimetric **AeroEclipse*** provided the largest inhaled drug mass and the lowest loss to ambient air, with the test conditions we used.

7. **COMPARISON OF BREATH-ACTUATED JET NEBULIZER (BAN) IN 'CONTINUOUS DELIVERY' MODE WITH OTHER CONTINUOUS DELIVERY NEBULIZERS.** JP Mitchell, KJ Wiersema, CC Doyle and MW Nagel. Trudell Medical Aerosol Laboratory, London, Canada. Respiratory Care 2003; 48 (11): S1077.

The **AeroEclipse*** BAN (Monaghan Medical Corp., Plattsburgh, N.Y.) has been equipped with an optional blue cap whose purpose is to retain the actuator piston in the position it would occupy during inhalation in breath-actuated mode, so that the nebulizer operates continuously. The present study compared the delivery of a bronchodilator from diluted albuterol sulfate respirator solutions (3-ml of 0.83 and 1-ml of 2.5 mg/ml albuterol in physiologically normal saline (0.9% w/v NaCl)), via this nebulizer, the Micromist® (Hudson RCI, Temecula, CA), Misty-Neb™ (Allegiance Healthcare Corp., McGaw Park, IL) and the LCD™ (PARI Respiratory Equipment, Monterey, CA). Each nebulizer was tested using a breathing simulator set to the following parameters representative of adult use: tidal volume = 600-ml, rate = 10 breaths/min, inspiratory/expiratory ratio 1:2. The total mass of albuterol (TM) delivered to the first sputter was determined by filter collection at the mouthpiece of the nebulizer operated with compressed air supplied at 50 psig at 8 L/min ($n = 5$ devices/group, 3 replicates/device). The fraction of the aerosol contained in droplets finer than $4.8 \mu\text{m}$ aerodynamic diameter (FPF) was determined by laser diffractometry in a parallel study, so that the fine droplet mass (FM) could be calculated as the product of TM and FPF. Values of FM (mean \pm SD) and time to deliver medication (T_{med}) were as follows:

| Solution(mg/ml) | AeroEclipse* | | LCD™ | | Micromist® | | MistyNeb™ | |
|------------------------|---------------------|-----------------|-----------------|----------------|-----------------|-------------------|---------------|---------------|
| | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 |
| FM (μg) | 360 ± 22 | 263 ± 26 | 149 ± 16 | 108 ± 4 | 209 ± 12 | 15.4 ± 5.9 | 82 ± 9 | 31 ± 5 |
| T_{med} (min) | 3 | <1 | 2 | <1 | 7 | <1 | 4 | <1 |

The **AeroEclipse*** nebulizer delivered significantly more FM in continuous delivery mode than the other nebulizers when operated in continuous mode with either solution strength (1-way repeated measures ANOVA, $p < 0.05$). T_{med} from the **AeroEclipse*** nebulizer was comparable with the best performing continuous nebulizer (LCD™).

8. **BREATH-ACTUATED NEBULIZER DELIVERS BRONCHO-DILATOR MORE EFFICIENTLY THAN CONVENTIONAL JET NEBULIZER IN A SIMULATION OF AN ADULT TIDAL-BREATHING PATIENT.** MW Nagel and JP Mitchell. Am. J. Resp. Crit. Care Med., 2002;165(8):A189.

Rationale: To compare delivery of albuterol sulfate inhalation solution (2.5 mg/3 ml vial equivalent to 0.083% w/v albuterol, Zenith Goldline Pharmaceuticals, Miami, FL) by conventional and breath-actuated nebulizer (BAN), simulating adult use. **Methods:** Each SVN (n = 5/group, 3 replicates/nebulizer) was operated with 8 l/min air at 50 psig and simulating breathing at tidal volume, I:E ratio and rate of 600-ml, 1:2 and 10/min respectively. Total emitted dose (TED) was determined for 5-**AeroEclipse*** BANs (Monaghan Medical Corp., N.Y., 1.5 ml solution) and 5 Micromist® nebulizers (Hudson RCI, Temecula, CA, 3.0 ml solution) by filter collection, and droplet size distributions were measured in a parallel study by laser diffractometer. Fine particle dose (FPD) was calculated as the product of TED and the percentage by mass of droplets finer than 4.8 mm aerodynamic diameter. **Results:** After 3 minutes, the **AeroEclipse*** BAN delivered 282 ± 10 mg FPD (mean ± SD) and the Micromist® delivered 209 ± 12 mg albuterol after 7 minutes. **Conclusion:** Dose delivery and patient compliance are assured by virtue of the breath actuation feature of the **AeroEclipse*** nebulizer and the reduced time to deliver a specific equivalent dose of medication compared with a conventional nebulizer will improve cost effectiveness of treatment.

9. **SAFETY AND EFFICACY OF FIVE-MINUTE TIMED AEROSOL ADMINISTRATION WITH THE AEROECLIPSE* BREATH ACTUATED NEBULIZER: COMPARISON OF LEVALBUTEROL WITH RACEMIC ALBUTEROL.** RS Pikarsky, R Acevedo, C Roman, W Fascia, T Farrell. Respiratory Care, September 2002; 47(9):1075.

Purpose: Beta₂-agonist Racemic Albuterol has been used extensively in the performance of pre & post bronchodilator studies in the pulmonary function laboratory. This study evaluated the safety and efficacy of timed nebulization of the two dosages of Levalbuterol (Sepracor Inc., Marlborough, MA) as compared to Racemic Albuterol (Dey, Napa, CA) with the use of the **AeroEclipse*** Breath Actuated Nebulizer (BAN) (Monaghan Medical Corp. Plattsburgh, N.Y.). **Methods:** A consecutive, non-randomized, mostly COPD population (n = 93) receiving pre & post bronchodilator testing in our Pulmonary Function Lab were studied. Two different Levalbuterol medication dosages were administered: 0.63mg Levalbuterol UD or 1.25mg UD Levalbuterol. The Racemic Albuterol dosage was 2.5mg UD. All 5 minute timed aerosol treatments were administered using the BAN with an oxygen flow rate of 8L/min. The Sensormedics Vmax 22 Pulmonary Function System was utilized to measure both FEV1 and PEFR. A standardized subjective questionnaire to determine side effects was completed. **Results:** The table shows the Levalbuterol and Racemic Albuterol dosages, mean % change of FEV1 and PEFR from pre-treatment to 10-minute post treatment, administration time, tremulousness and increase in heart rate. There was no significant difference in % change in FEV1 or PEFR. There was a significant increase in heart rate with the 1.25mg Levalbuterol UD group (7.2 vs. 3.4, p<.05*; 7.2 vs. 2.2, p<.01**). There was no difference in respiratory rate, tremulousness, or nausea.

| Nebulizer (n) | Dose | % Change FEV1 | % Change PEFR | Time (min) | Trem. | HR (Inc.) |
|------------------------|-----------|---------------|---------------|------------|-------|-----------|
| Levalbuterol (38) | 0.63mg UD | 7.8 | 6.2 | 5 | 4 | 3.4* |
| Levalbuterol (29) | 1.25mg UD | 7.7 | 16.6 | 5 | 2 | 7.2 |
| Racemic Albuterol (26) | 2.25mg UD | 12.2 | 10.5 | 5 | 0 | 2.2** |

Conclusion: Five minute timed administration of Levalbuterol and Racemic Albuterol using the BAN was equally efficacious and had similar safety profiles. The change in FEV1 and PEFR are consistent with our mostly COPD population. The increase in heart rate was greatest with the Levalbuterol 1.25 mg dosage. **Clinical Implications:** Five minute timed administration of Levalbuterol and Racemic Albuterol using the BAN is a safe and efficient alternative to the use of small volume nebulizers. Additional caution should be taken when administering Levalbuterol at the 1.25 mg dosage utilizing the BAN in cardiac patients. The efficiency of timed aerosol administration could have significant impact on resource utilization while maintaining the quality of aerosol delivery. This may be one of several strategies to address the problems of Respiratory Care staff shortages or high seasonal effect in the acute care facility.

10. **COMPARISON IN RATES OF BREAKTHROUGH TREATMENTS DURING A CONVERSION FROM RACEMIC ALBUTEROL TO LEVALBUTEROL.** RS Pikarsky, RA Acevedo, C Roman and T Farrell. CHEST 2002; 122(4):146S.

Purpose: in order to meet our patient care demands, Crouse Hospital approved an automatic conversion from Racemic Albuterol to Levalbuterol. This study compares the breakthrough rates of Racemic Albuterol and Levalbuterol, with and without Ipratropium. **Methods:** Racemic Albuterol (Alb) 2.5 mg Q4h was converted to either Levalbuterol (Lev) 0.63 mg Q6h or Levalbuterol 1.25 mg Q8h. If ordered, Ipratropium (Ipra) 0.5 mg was administered at the same frequency as the Levalbuterol. Patients with acute coronary syndromes, need for cardiac monitoring, or requiring more frequent aerosol administration received the lower Levalbuterol dose Q6h. A majority of aerosol therapy was provided with the use of the **AeroEclipse*** Breath Actuated Nebulizer (BAN). All aerosol treatments, including breakthrough treatments, delivered between July 1, 2001 and February 28, 2002 were recorded. **Results:** Tx/Pt/day represents the number of treatments delivered per patient per day. Rate/100 Pt/days = (Breakthrough) / (Total Tx / Tx/Pt/day) x 100. Rate/100 Pt/days corrects for the differences in daily administration frequency, and may better reflect the daily impact of the breakthrough rate. The breakthrough rate of the combined Albuterol group was significantly greater than both Levalbuterol groups (5.29 vs. 2.29, 5.29 vs. 2.43, p<.001)*. The breakthrough rate with Albuterol was significantly reduced with the addition of Ipratropium (p<.001)**. Ipratropium did not significantly change the breakthrough rate when added to Levalbuterol groups.

| Medication | Total Tx | Breakthrough | Rate/1000 | Tx/Pt/day | Rate/100 Pt/day | |
|----------------------|----------|--------------|-----------|-----------|-----------------|-------|
| Alb Q4h | 3832 | 47 | 12.27 | 6 | 7.36** | 5.29* |
| Alb/lpra Q4h | 3767 | 20 | 5.31 | 6 | 3.19** | |
| Lev 0.63mg Q6h | 3592 | 24 | 6.68 | 4 | 2.67 | 2.29* |
| Lev 0.63 mg/lpra Q6h | 1821 | 7 | 3.84 | 4 | 1.54 | |
| Lev 1.25mg Q8h | 1791 | 17 | 9.49 | 3 | 2.85 | 2.43* |
| Lev 1.25mg/lpra Q8h | 678 | 3 | 4.42 | 3 | 1.33 | |

Conclusions: The conversion from Racemic Albuterol to Levalbuterol allowed for a decreased frequency of daily medication administrations and a significant decrease in breakthrough requirements. Ipratropium showed a significant benefit in breakthrough reduction for the Racemic Albuterol group. **Clinical Implications:** The efficiencies gained by decreasing the daily frequency of aerosol administration can have a significant impact on resource utilization. The conversion to Levalbuterol allows for decreased respiratory therapy time or the re-allocating of workforce needs while maintaining, or improving, quality of aerosol administration, as evidenced by the decrease in breakthrough requirements.

11. **THE DELIVERY TIME, EFFICACY, AND SAFETY OF BETA AGONIST BRONCHODILATOR ADMINISTRATION WITH THE AEROECLIPSE* BREATH - ACTUATED NEBULIZER ("BAN").** RS Pikarsky, T Farrell, R Acevedo, W Fascia and C Roman. CHEST 2001; 120(4): 218S.

Purpose: Aerosol delivery consumes the highest level of Respiratory Care resources. This study evaluated the delivery time, efficacy, and safety of rapidly nebulized Albuterol with the use of the **AeroEclipse*** Breath Actuated Nebulizer as compared to both an MDI with **AeroChamber*** VHC (both from Monaghan Medical Corp. Plattsburgh, N.Y.) and the Airlife Misty-Neb Nebulizer (SVN) (Allegiance Healthcare Corporation). **Methods:** A consecutive, non-randomized, mostly COPD population receiving pre & post bronchodilator testing in our Pulmonary Function Lab were studied. Three different Albuterol medication dosages were administered with the BAN: 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline, 1.0 ml (5 mg) of undiluted Albuterol, and 0.75 ml Albuterol (3.75 mg) using an oxygen flow rate of 8 L/min. Two puffs of Albuterol were administered by MDI with **AeroChamber*** VHC. Treatments with the SVN consisted of nebulizing 2.5 mg of Albuterol diluted with 3 ml of Normal Saline Unit Dose (UD) using an oxygen flow rate of 8 L/min. The Sensormedics Vmax 22 Pulmonary Function System was utilized to measure FEV1. A standardized subjective questionnaire to determine side effects was completed.

| Nebulizer (n) | Dose | % Change FEV1 | Time(min) | Tremulousness |
|---------------------------------|-----------------------|---------------|-----------|---------------|
| AeroEclipse* BAN (12) | 0.5 ml + 0.5 ml NS | 8.2% | 2.67* | 0 |
| AeroEclipse* BAN (64) | 1.0 ml undil. | 10.9% | 3.29* | 17 |
| AeroEclipse* BAN (23) | 0.75 ml undil. | 5.6% | 1.30* | 5 |
| MDI (21) | 2 puffs | 8.5% | 2.86** | 1 |
| Misty-Neb (52) | 2.5 mg UD | 9.1% | 8.33 | 2 |


Results: The table shows the Albuterol dosages, mean % change of FEV1 from pre-treatment and 10 minute post treatment, mean administration time and tremulousness. The mean treatment time with all BAN patients was 2.78 minutes as compared to 8.33 minutes with the SVN ($p < .001$) *. The mean treatment time with the MDI was 2.86 minutes as compared to 8.33 minutes with the SVN ($p < .001$) **. The changes in FEV1 were not significant. There was no difference in heart rate, respiratory rate or nausea. Seventeen patients receiving the 1.0 l undiluted Albuterol indicated an increase in tremulousness. **Conclusion:** The rapid administration of Albuterol in the 0.5 ml + 0.5 ml NS and 1.0 ml undiluted doses using the BAN was equally efficacious as the MDI with **AeroChamber*** VHC and SVN UD. The 1.0 ml Albuterol dosage has the highest incidence of tremulousness. The 0.75 ml Albuterol dosage under-performed. Delivering 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline using the BAN offered the best delivery time, efficacy and safety profile of the nebulizer trials. The BAN performance was comparable to the MDI with **AeroChamber*** VHC. **Clinical Implications:** In a health care facility that delivers large volumes of aerosol treatments, the decrease in delivery time could have a significant impact on resource utilization. The results supported changes in the Respiratory Care practice throughout Crouse Hospital. Further studies evaluating additional medication dosing regimens measuring safety, efficacy and resource utilization are needed.

12. **THE CLINICAL EFFICACY OF USING THE AEROECLIPSE* BREATH ACTUATED NEBULIZER ("BAN") IN PULMONARY LAB TESTING AND IMPLICATIONS FOR GENERAL USE.** YM Christensen, CJ Flanigan, SA Ravenscraft. Respiratory Care 2001;46(10):1084.

Purpose: To compare the clinical efficacy and delivery time of nebulization of beta agonist bronchodilator with the use of the **AeroEclipse*** Breath Actuated Nebulizer ("BAN") (Monaghan Medical Corp.) as compared to the Airlife Misty-Neb Nebulizer (SVN) (Allegiance Healthcare Corporation). **Methods:** Adult patients (n=40) presenting with Asthma (50%), COPD (10%) and other pulmonary disorders (40%); receiving pre and post bronchodilator spirometry testing in our Pulmonary Function Lab were included in the study. Each patient received both nebulizers on two separate visits (less than 24 hours apart). Patient received a nebulizer treatment with the BAN (n=40) 2.5mg Albuterol (0.5ml) in 0.5cc saline run to sputter, or the SVN (n=40) 2.5mg Albuterol in 2.5cc saline (3ml unit dose) run to sputter. FVC, FEV1, FEV1% ratio and FEF 25-75% spirometry was conducted using the Medical Graphics 1085DX pre and 5 minutes post treatment with the BAN and 10 minutes post treatment with the SVN. **Results:** The results demonstrated that FVC, FEV1 and FEF 25-

75% for patients using the BAN were substantially higher while FEV1% ratio favored the SVN (Table and Chart). Importantly, total nebulization time was reduced from 22 minutes (SVN) to 7 minutes (BAN), and total test time was reduced from 30 minutes (SVN) to 15 minutes (BAN).

| SPIROMETRY RESULTS | | | | |
|-----------------------------|------|------|------------------|-------|
| Absolute % Change by Device | | | % Difference BAN | |
| | SVN | BAN | | BAN |
| FVC | 5.3 | 10.2 | FVC | 91.3 |
| FEV1 | 7.3 | 13.1 | FEV1 | 79.8 |
| FEV1%ratio | 3.0 | 2.3 | FEV1% | -25.1 |
| FEF 25-75% | 29.8 | 57.7 | FEF 25-75% | 93.3 |



Conclusion: The administration of 2.5mg of albuterol with the BAN produced improved results in FVC, FEV1 and FEF 25-75%. Substantially shorter test times delivered by the BAN would allow for more tests and associated revenue. These data support the thesis that the BAN can reduce costs of care by delivering clinically acceptable outcomes in significantly less time.

- BREATH-ACTUATED VS RESERVOIR NEBULIZERS FOR UNDILUTED ALBUTEROL.** D Geller, B Kesser. The Nemours Children's Clinic Aerosol Research Lab, Orlando FL, USA. Presented at the 13th International Congress on Aerosols in Medicine, Interlaken, Switzerland, September 17-21, 2001.

Aim: Some Emergency Departments use undiluted albuterol in nebulizers designed to conserve drug during exhalation. We compared the *in vitro* performance of 4 devices to estimate which would be most effective clinically: **AeroEclipse*** Breath-Actuated Nebulizer ("BAN"); Circulaire® (C) and AeroTee™ (AT) which use a 750 ml reservoir bag to conserve drug during exhalation; and Salter HDN™ (S) with a 50 ml tower reservoir. **Method:** We studied 4 units of each nebulizer type in duplicate, using a Pari Proneb Turbo compressor. Nebulizers were filled with undiluted 0.5% albuterol, 1 ml (5 mg) or 2 ml (10 mg). Particle size distributions were measured by laser diffraction (Malvern SprayTec). Drug output (1 minute after "sputter") was captured on a filter between the device mouthpiece and a Pari breath-simulator, which used a recorded waveform from a 9 yr old male. Albuterol was measured by spectrophotometry, and fine particle dose (FPD) (mg of drug < 5 mm in size) was calculated.

Results:

| Neb | MMAD | FPD (1cc) | Minutes | FPD (2cc) | Minutes |
|-----|------|-----------|---------|-----------|---------|
| AE | 3.9 | 0.60 | 3.8 | 2.41 | 11.0 |
| AT | 4.8 | 0.03 | 2.0 | 0.62 | 3.2 |
| C | 2.5 | 0.09 | 2.0 | 0.65 | 3.7 |
| S | 8.5 | 0.08 | 2.0 | 0.57 | 3.7 |

Conclusions: The AE was superior to the reservoir-type nebulizers in fine-particle output for each fill volume. The AT and C had large dead volumes, and the S produced larger particles. These shortcomings were overcome with larger nominal doses. Each nebulizer produced 0.6-mg FPD of albuterol over 3½ minutes, but the AE required only half the starting dose. Albuterol 0.6 mg is a reasonable clinical respirable dose in a child with acute asthma. These findings must be taken into account when designing clinical treatment protocols for acute asthma. **Background:** Many nebulizers are designed to decrease the amount of drug that is lost during exhalation. The Circulaire® (Westmed) and AeroTee™ (Hudson) incorporate a 750 ml bag on the expiratory side of the nebulizer that collects aerosol while the patient exhales, making it available for inhalation on the next breathing cycle. The Salter HDN™ (Salter) has a 50 ml tower that acts as a reservoir. The **AeroEclipse*** BAN (Trudell/Monaghan) has a spring mechanism that allows generation of aerosol during inhalation

only, so no drug waste occurs during exhalation. We recently reported the aerosol characteristics with these devices nebulizing unit-dose albuterol sulfate (2.5 mg/3 ml).¹ Delivery time with unit-dose (0.083%) albuterol can be long, which may increase personnel costs. To maintain lung-dose delivery and minimize the treatment time, some hospitals use drug-conserving nebulizers with small fill-volumes of undiluted (0.5%) albuterol for patients presenting with acute bronchospasm. We measured the particle size distributions and used a child's breathing pattern to compare albuterol output of these 4 drug-conserving nebulizers, using unit-dose albuterol 2.5 mg (3ml), 0.5% albuterol 5 mg (1ml) and 10 mg (2ml) nominal doses. We calculated the fine particle dose and measured the dose of drug remaining within the nebulizer and all attachments to determine the residual dose. For reference, we compared these results to those of a T-piece (Hudson Micromist) nebulizer using unit-dose albuterol to simulate conventional dosing. **Materials and Methods:** Drug: Albuterol Sulfate 0.083% unit-dose (2.5mg/3ml); Albuterol Sulfate 0.5% (5mg/ml) 1 & 2cc fill volumes. Nebulizers: Circulaire[®] (Model 0260), AeroTee[™] with Micromist Nebulizer (Model 1002), Salter HDN[™] (Model 8960), and **AeroEclipse[®] BAN** (Figure 1). Compressor: PARI PRONEB TURBO. 4 nebulizers of each type studied in duplicate; Particle size by laser diffraction (Malvern Insotec); Breathing pattern from 9 year old male volunteer, using the PARI breath simulator (RR 19 bpm, Vt 421 cc, Ti 1.3 seconds). **Definitions:** Inspired dose = drug on inspiratory filter; Residual dose = drug collected from nebulizer and accessory components after completion of nebulization; Fine particle dose (FPD) = (Inspired dose) x (% of particles <5 µm) Figure 1; Duration = time (minutes) from the beginning of nebulization to 1 minute past the onset of sputter; Samples assayed with spectrophotometer at 228 λ.

Results:

| | | AeroEclipse[®] BAN | AeroTee[™] | Circulaire[®] | Salter HDN[™] |
|-------------------------|-------------------------|------------------------------------|----------------------------|-------------------------------|-------------------------------|
| Particle Sizing | MMD | 3.87 | 4.80 | 2.47 | 8.46 |
| | GSD | 2.3 | 2.0 | 2.1 | 2.0 |
| | % < 5 µm | 61.7% | 52.9% | 83.6% | 30.0% |
| 2.5 mg Unit Dose* | Duration (minutes) | 14.7 | 7.2 | 7.0 | 3.6 |
| | Inspired Dose (mg) | 0.77 | 0.37 | 0.14 | 0.30 |
| | Residual Dose (mg) | 1.5 | 1.8 | 2.1 | 1.9 |
| | Fine Particle Dose (mg) | 0.52 | 0.19 | 0.12 | 0.10 |
| | | | | | |
| 5 mg (1 ml) Dose | Duration (minutes) | 3.8 | 2.0 | 2.0 | 2.0 |
| | Inspired Dose (mg) | 0.97 | 0.06 | 0.11 | 0.28 |
| | Residual Dose (mg) | 3.5 | 4.9 | 4.6 | 4.4 |
| | Fine Particle Dose (mg) | 0.60 | 0.03 | 0.09 | 0.08 |
| | | | | | |
| 10 mg (2 ml) Dose | Duration (minutes) | 11.0 | 3.2 | 3.7 | 3.7 |
| | Inspired Dose (mg) | 3.9 | 1.2 | 0.8 | 1.9 |
| | Residual Dose (mg) | 5.8 | 8.7 | 8.6 | 6.9 |
| | Fine Particle Dose (mg) | 2.40 | 0.60 | 0.60 | 0.60 |
| | | | | | |

* Unit-dose data presented at ATS 2001¹

For comparison, the Hudson Micromist conventional T-Piece Nebulizer (with Unit-Dose 2.5 mg Albuterol) produced a fine-particle dose of 0.14 mg in 7.0 minutes.

Discussion:

- **AeroEclipse[®] BAN** had highest FPD with all nominal doses:
 - FPD was 2.7 to 5.2 times higher with unit-dose; 6.7 to 20 times higher with 5 mg dose; 4 times higher with 10 mg dose.
 - Lowest residual dose
 - Higher fine particle fraction except for Circulaire[®]
- **Nebulizer Inefficiencies:**
 - AeroTee[™] and Circulaire[®] had high residual doses in part due to valves and collection bags.
 - Salter HDN[™] produces larger particles
 - These inefficiencies were partially compensated for by increasing nominal dose to 10 mg (2 ml)
- **Duration of Nebulization:**
 - **AeroEclipse[®] BAN** had longer delivery time because it is breath actuated; no waste during exhalation
 - Using 0.5% albuterol, all nebulizers produced 0.6 mg fine-particle dose in < 4 minutes, but the **AeroEclipse[®] BAN** only required half the nominal dose to accomplish this.
- **Comparison to Unit-Dose 2.5 mg:**
 - **AeroEclipse[®] BAN** produced comparable FPD with unit-dose and 5 mg (1 ml) nominal dose, but delivery time was less than a third with undiluted drug.
- **Comparison to conventional nebulizers:**
 - The FPD with the Hudson and unit-dose drug was 0.14 mg, similar to the reservoir-type nebulizers with unit-dose.
 - The higher FPD with **AeroEclipse[®] BAN** (all doses) and the reservoir nebs (10 mg dose) may result in better and longer lasting bronchodilation than the Hudson with conventional dosing, thus reducing number of treatments, therapist time, and total costs.

Conclusion:

- **AeroEclipse*** BAN was superior to the reservoir-type nebulizers at all nominal doses.
- **AeroEclipse*** BAN has the additional advantage of being a dosimetric device, i.e. it will not operate or waste drug while the patient is coughing or resting. The patient and health care providers get visual feedback of adequate inspiratory effort necessary to actuate the nebulizer.
- Use of undiluted 0.5% albuterol may result in higher lung doses in a shorter amount of time. These results can be used as a guide when developing bronchodilator protocols for the hospital or E.D. setting.

Funded by the Nemours Foundation

¹ Geller D, Kesser B. Am J Respir Crit Care Med 2001; 163:A444 Journal of Aerosol Medicine 14 (3) 2001; 395:1-41.

14. **THE DELIVERY TIME, EFFICACY, AND SAFETY OF BETA AGONIST BRONCHODILATOR ADMINISTRATION WITH THE AEROECLIPSE* BREATH ACTUATED NEBULIZER ("BAN") VERSUS A CONVENTIONAL T-TYPE SMALL VOLUME NEBULIZER.** RS Pikarsky, T Farrell, R Acevedo, W Fascia and C Roman. Respiratory Care 2001; 46(10):1085.

Purpose: Aerosol delivery consumes the highest level of Respiratory Care resources. This study evaluated the delivery time, efficacy, and safety of rapidly nebulized albuterol with the use of a novel breath actuated nebulizer compared to a standard small volume nebulizer. **Methods:** A consecutive, non-randomized, mostly COPD population receiving pre & post bronchodilator testing in our Pulmonary Function Lab were studied. 0.5 ml albuterol (2.5 mg) with 0.5 ml Normal Saline (NS) was administered with the **AeroEclipse*** Breath Actuated Nebulizer ("BAN") (Monaghan Medical Corp. Plattsburgh, N.Y.) using an oxygen flow rate of 8 L/min. Treatments with the AirLife™ brand Misty-Neb™ small volume nebulizer (SVN) (Allegiance Healthcare Corporation) consisted of nebulizing 2.5 mg of albuterol diluted with 3 ml of Normal Saline Unit Dose (UD) using an oxygen flow rate of 8 L/min. The Sensormedics Vmax 22 Pulmonary Function System was utilized to measure FEV1. A standardized subjective questionnaire to determine side effects was completed. **Results:** The table shows the albuterol dosages, mean % change of FEV1 from pre-treatment and 10 minute post treatment, mean administration time and tremulousness. The mean treatment time with all BAN patients was 2.67 minutes as compared to 8.33 minutes with the SVN ($p < .001$)*. The changes in FEV1 were not significant. There was no difference in heart rate, respiratory rate or nausea. **Conclusion:** The rapid administration of albuterol in the 0.5 ml + 0.5 ml NS dose using the BAN was equally efficacious as the SVN UD. Delivering 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline using the BAN offered the best delivery time, efficacy and safety profile between the two devices. **Clinical Implications:** In a health care facility that delivers large volumes of aerosol treatments, the decrease in delivery time achieved with the BAN could have a significant impact on resource utilization. The results supported changes in the Respiratory Care practice throughout Crouse Hospital. Further studies evaluating additional medication dosing regimens measuring safety, efficacy and resource utilization are needed.

| Nebulizer (n) | Dose | % Change FEV1 | Time (min) | Tremulousness |
|---------------------------------|---|---------------|------------|---------------|
| AeroEclipse* BAN (12) | 2.5mg (0.5 ml albuterol + 0.5 ml NS) | 8.2% | 2.67* | 0 |
| Misty-Neb™ (52) | 2.5mg (3 ml unit dose) | 9.1% | 8.33 | 2 |

15. **COMPARISON OF DRUG OUTPUT FROM 4 DIFFERENT RESERVOIR TYPE NEBULIZERS.** DE Geller, B Kesser. American Journal of Respiratory Care & Critical Care Medicine, 2001; 163(5): A444.

Rationale: Many nebulizers currently being marketed utilize different techniques to conserve drug that would normally be lost during exhalation. The Circulaire™ and Aero Tee™ nebulizers use a 750 cc reservoir bag to accumulate nebulized drug, while the Salter HDN™ uses a 50ml tower to serve as a reservoir. The **AeroEclipse*** nebulizer uses breath actuated nebulization to deliver drug only during inspiration. We evaluated all 4 nebulizers using a recorded pediatric breathing pattern to measure total drug output. We additionally measured the particle size characteristics of each type with the laser diffraction technique. **Methods:** 4 nebulizers of each type were studied in duplicate for sizing and total output characteristics. Each nebulizer was charged with a unit dose of 2.5 mgs albuterol sulfate in 3cc's. Sizing studies were averaged values preformed over 5 minute runs on each nebulizer with a Malvern Spray Tec laser. Drug output was as calculated as the assayed amount of albuterol collected on a filter distal to the mouthpiece of the nebulizer. Simulated breathing was performed through the nebulizer by a Pari breath simulator from waveforms originally recorded from a healthy 9-year-old male.

Results:

| | Inspired dose | %>1 & <5M | Respirable Dose | Residual Dose |
|---------------------|----------------|------------|-----------------|---------------|
| AeroEclipse* | 0.64 ± 0.06 mg | 52.7 ± 2.0 | 0.34 ± 0.03 mg | 1.27 ± 0.09 |
| Aero Tee | 0.31 ± 0.09 mg | 41.2 ± 7 | 0.13 ± 0.04 mg | 1.51 ± 0.11 |
| Circulaire | 0.12 ± 0.03 mg | 61.9 ± 1 | 0.07 ± 0.02 mg | 1.72 ± 0.13 |
| Salter HDN | 0.25 ± 0.05 mg | 24.7 ± 5 | 0.06 ± 0.02 mg | 1.59 ± 0.10 |

Conclusion: The **AeroEclipse*** delivers a greater total dose of drug as well as a greater amount of drug in the fine particle range, most likely to deposit in the lower airways.

16. **CLINICAL EVALUATION OF A BREATH ACTUATED SMALL VOLUME NEBULIZER (BA-SVN).** S Klopff, N Scneiderman, H Payne, C Schramm, MW Nagel, and JP Mitchell. Respiratory Care, August 2000; 45(8).

Background: In prior *in-vitro* studies using laser diffractometry, the aerosol produced by a novel breath-actuated nebulizer (BAN), the **AeroEclipse*** (Monaghan Medical Corp. Plattsburgh, NY) has been shown to contain a high proportion of droplets < 4.8 µm diameter (80.9% ± 2.4%). Such droplets are more likely to penetrate beyond the oro-pharyngeal region where bronchodilation is achieved. These *in-vitro* results should therefore be predictive of improved *in-vivo* delivery of nebulized medications to the respiratory tract. This study explored the clinical performance of the **AeroEclipse*** BAN in the delivery of a beta2-agonist (albuterol 2.5 mg/ml) accompanied by anticholinergic (ipratropium bromide 250 µg/ml) bronchodilator in some cases. **Methods:** Patients (n=48) with a previous diagnosis for asthma presenting to the Emergency Department for acute exacerbation of asthma were included in this study. Upon presentation, an asthma care path, an assessment driven, algorithm-based tool was used to place patients in one of three stages of severity as recommended by the NIH-NAEPP Guidelines for the Diagnosis of Asthma. Each patient was assigned to receive inhaled aerosol treatment using the **AeroEclipse*** BAN. Stage 1 asthmatics were given 0.5-ml of albuterol with 0.5-ml normal saline delivered until sputter. Patients categorized in stage two and three were given 0.5-ml albuterol with the addition of 1.5-ml of ipratropium bromide unit dose. Treatments repeated every 20 minutes times three if necessary by protocol.

Results:

| Asthma Severity | Stage 1 | Stage 2 | Stage 3 |
|-----------------------------------|-----------|----------------|---------------|
| Number | 10 | 30 | 8 |
| Treatments Given | 2.4 | 2.03 | 2.25 |
| Treatment Duration (min) | 3.7 | 3.78 | 5 |
| Increase in PEF (mean, range (%)) | 44(0-120) | 67.7(-2.7-580) | 120.7(28-420) |

Four patients had greater than 20% increase in heart rate, three patients noted tremor following treatment. Twenty four patients had positive comments about the device focused on shorter treatment time and improved relief from dyspnea. Two imminent intubations were avoided with the use of the BA-SVN. **Conclusions:** Use of the **AeroEclipse*** BAN appears to result in good clinical outcomes. Minimum number of treatments, shorter treatment duration and minimal side effects were noticed with this device. Further outcome studies are needed to assess this impact on other groups of patients.

17. **EVALUATION OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) FOR THE DELIVERY OF ALBUTEROL SULFATE AND CROMOLYN SODIUM.** JP Mitchell, MW Nagel, A Archer, and D. Presented at ALA/ATS International Conference, San Diego, CA, 1999. (Am. J. Resp. Crit. Care Med., 1999; 159(3): A120)

Purpose: To evaluate the delivery of Ventolin® (0.2% v/v, albuterol sulfate, GlaxoSmithKline, Canada) and Intal® (1.0% v/v cromolyn sodium, Fisons Pharmaceuticals Ltd., Canada) by a prototype AE-SVN (Trudell Medical International) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 AE-SVNs were tested using an Andersen Mark II Cascade Impactor operated at 28.3±0.5 l/min to determine the size distribution of droplets emitted at the mouthpiece during the first 10 seconds following nebulization. The mass of drug emitted was determined directly by HPLC-UV spectrophotometry. **Results:** Total (TM) and fine particle ((FPM), droplets finer than 4.7 µm diameter) mass output rates and droplet mass median diameter (MMD) were as follows:


| Drug | TM (µg/s) | FPM (µg/s) | MMD (µm) |
|-----------|--------------|-------------|-----------|
| Ventolin® | 32.4 ± 3.1 | 27.6 ± 1.3 | 3.0 ± 0.1 |
| Intal® | 138.6 ± 10.2 | 109.7 ± 8.3 | 3.2 ± 0.1 |

Conclusion: The fine MMD produced from the AE-SVN resulted in an improved FPM output rate, which is likely to produce increased lung deposition.

18. **EFFECT OF NEBULIZER DESIGN ON FINE PARTICLE MASS.** D Hess, JP Mitchell, D Coppola, MW Nagel, AD Archer, R Blacker. Presented at Open Forum, Annual Meeting of the American Association For Respiratory Care (AARC), Las Vegas, NV, 1999. (Resp Care, Oct 1999; 44(10): 1289)

Background: Nebulizer design is known to affect performance. In this study, we compared fine particle mass from nebulizers of four designs. **Methods:** We tested traditional disposable nebulizers (Baxter Misty-Neb™, Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™), nebulizers with collection bags (Westmed Circulaire™), and a Trudell **AeroEclipse*** (with breath actuation disabled). Five of each device with three replicates (n = 15) were tested using an *in-vitro* model of spontaneous breathing. A rigid bar was placed between the two compartments of a test lung (Michigan Instruments TTL). The drive lung was attached to a ventilator (Infrasonics Infant Star™) to simulate spontaneous breathing (tidal volume 0.6 L, rate 10/min, T_i 2 s). A bacterial/viral filter (Trudell MT3000) was placed between the nebulizer and slave lung. Flow was monitored between the test lung and filter (Novamatrix Ventcheck™). Albuterol solution (0.625 mg/mL) was placed into the nebulizers (4 mL), which were powered with air (8 L/min). Filters were replaced at one minute intervals (flow to the nebulizer was discontinued during filter replacement) until sputtering occurred. The filter was washed with methanol (20 mL) and albuterol concentration was measured with HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer. Fine particle mass was calculated as the product of mass % <4.7 µm and total nebulizer output.

Results: Fine particle mass from the **AeroEclipse*** nebulizer was greater than that from the other nebulizers ($P < 0.001$) (see figure). **Conclusions:** Fine particle mass was affected by nebulizer design. The clinical relevance of this finding awaits further investigation. Further evaluation of the breath-actuated feature of the **AeroEclipse*** is warranted.



19. **PERFORMANCE OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER.** A Archer, JP Mitchell, MW Nagel, AMW Verdun. Eur. Resp. J.1998; 12(28): 68.

We report an *in vitro* investigation in which the performance of a new disposable AE-SVN ($n = 3$ devices) has been assessed with salbutamol sulphate (Ventolin®: 5 µg/2.5 ml, GlaxoSmithKline Inc.), metaproterenol sulphate (Alupent®: 10 µg/2.5 ml, Boehringer Ingelheim Pharmaceuticals Inc.) and cromolyn sodium (Intal®: 20 µg/2 ml, Fisons Pharmaceuticals) nebulizers. Each AE-SVN was filled with 2 nebulizers and operated continuously with oxygen supplied at 50 psig and 8 l/min. The AE-SVN was coupled directly to an Andersen cascade impactor, sampling at 28.3 l/min. Total and fine particle (< 4.7 µm aerodynamic diameter) delivery rates were 33.5 ± 1.8 µg/s and 27.6 ± 1.3 µg/s (Ventolin®); 54.2 ± 10.6 µg/s and 45.0 ± 7.8 µg/s (Alupent®); 138.6 ± 10.2 µg/s and 109.7 ± 8.3 µg/s (Intal®) over a 10 s period following the start of nebulization. The mass median aerodynamic diameter (MMAD) and mass % contained in fine droplets were 3.0 ± 0.1 µm and $82.4 \pm 1.2\%$ (Ventolin®); 2.9 ± 0.2 µm and $83.3 \pm 2.6\%$ (Alupent®); 3.1 ± 0.1 µm and $79.2 \pm 1.9\%$ (Intal®). This new nebulizer appears to perform well with all three formulations.

20. **THE EFFECT OF SMALL VOLUME NEBULIZER (SVN) DESIGN ON FINE PARTICLE MASS DELIVERY OF A BRONCHODILATOR.** R Blacker, RW Morton, JP Mitchell, MW Nagel and DR Hess. Drug Delivery to the Lungs-X, London, UK, 1998. (J. Aerosol Med., 1998;13(1):65)

Fine particle mass delivery was compared from six different SVNs, including continuous un-enhanced flow designs (Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™, Medic-Aid Sidestream®), nebulizers with aerosol collection bag (Westmed Circulaire™), and an **AeroEclipse*** with breath actuation disabled (Trudell Medical International). Five of each type of SVN were tested operating with air (8 l/min, 50 psig), using an *in-vitro* model that simulated spontaneous breathing by an adult (tidal volume 0.6 l, rate 10/min, TI = 2 s). A bacterial/viral filter was placed between the nebulizer and breathing simulator. In each case, salbutamol sulphate (Ventolin®) respirator solution (0.625 mg/ml, 4 ml) was placed into the reservoir of the SVN. The filters were replaced at one-minute intervals until sputtering occurred. The salbutamol collected on the filter was assayed by HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer laser diffractometer. Fine particle mass delivery rates varied significantly from each of the SVNs from more than 110 µg/min (**AeroEclipse***) to ca. 20 µg/min (Circulaire™).

21. **PERFORMANCE OF A NEW BREATH-ACTUATED SMALL VOLUME NEBULIZER UNDER CONDITIONS THAT SIMULATE USE BY AN ADULT PATIENT.** R Blacker, JP Mitchell, MW Nagel, and AMW Verdun. Eur. Resp. J. 1997; 10(25): 235.

The development of pneumatic small volume nebulizers (SVNs) in which atomization is enabled during the inhalation portion of a patient's breathing cycle has important ramifications in terms of the efficiency at which medication can be delivered. We report an investigation in which the effectiveness for the delivery of salbutamol (Ventolin® nebulizer: 5 mg/2.5 ml, GlaxoSmithKline, Canada) via a prototype breath-actuated SVN (Trudell Medical, Canada (TRU)) was compared with that of a high performance closed-system SVN (Ventstream™, Medic-Aid, Pagham, U.K. (VEN)). Each device was connected in turn to a ventilator-test lung apparatus in such a way that aerosol delivered on inhalation (800 ml tidal volume, I/E of 1/1, 15 breaths/min) was collected on a filter (Filtrete™, 3M Corp., St Paul, MN) located at the mouthpiece. Oxygen (440 kPa, 8 l/min) was supplied to operate each SVN, and the contents of a single nebule (2.5 ml) were added to the reservoir at the start of each test. Over a 5 minute period of use, the TRU SVN provided 1.74 ± 0.04 mg salbutamol to the filter ($n=5$ replicates). In comparison, the VEN delivered 1.28 ± 0.01 mg in 3.5 min after which the device sputtered dry ($n = 5$ replicates). These data indicate that the new breath-actuated device has important benefits in reducing wastage of medication by operating more efficiently, as well as an optimal impact on the environment.

22. **A NOVEL BREATH-ACTUATED SMALL VOLUME NEBULIZER UNDER SIMULATED ADULT USE CONDITIONS.** R Blacker, JP Mitchell, MW Nagel, AMW Verdun. Presented at Open Forum, Ann. Meet. of the American Association For Respiratory Care (AARC), New Orleans, LA, 1997.

Pneumatic small volume nebulizers (SVNs) in which atomization only occurs during the inhalation phase of the breathing cycle have important ramifications in terms of the efficiency of medication delivery. We report an investigation in which the effectiveness for the delivery of salbutamol (Ventolin[®] nebulizer: 5 mg/2.5 ml, GlaxoSmithKline, Canada) via a prototype breath-actuated SVN (Trudell Medical, Canada (TRU)) was compared with that of a high performance closed-system SVN (Ventstream[™], Medic-Aid, Pagham, U.K. (VEN)). Each nebulizer was connected in turn to a dual-chambered test lung with one chamber driven by a ventilator and the other connected to the SVN mouthpiece. Aerosolized salbutamol delivered on inhalation (800 ml tidal volume, I/E of 1/1, 15 breaths/min) was collected on a filter (Filtrete[™], 3M Corp., St Paul, MN) located at the mouthpiece. Oxygen (440 kPa, 8 l/min) was used to operate each SVN, and the contents of a single nebule (2.5 ml) were added to the reservoir at the start of each test. Over a 5 minute period of use, the TRU SVN provided 1.74 ± 0.04 mg salbutamol to the filter (n=5 replicates), significantly more than the VEN which delivered 1.28 ± 0.01 mg in 3.5 min (Mann Whitney Rank Sum Test, p = 0.008), after which the device sputtered dry (n = 5 replicates). These data indicate that the new breath-actuated device may have important benefits in reducing wastage of medication by operating more efficiently, as well as reducing exposure to the care-giver.

Pulmicort[™] (budesonide), AstraZeneca[™]

1. **DELIVERY OF BUDESONIDE INHALATION SOLUTION (BIS) THROUGH AN INFANT UPPER AIRWAY MODEL.** DE Geller, KC Kesser, HM Janssens, HAWM Tiddens. Am J Respir Crit Care Med, 2003; 167(7):A508.

We investigated variables that may be important in the delivery of BIS to the lungs of infants, a challenging population for aerosol delivery. **Methods:** The Sophia Anatomical Infant Nose Throat (SAINT) airway model mounted on a breath simulator mimicked the breathing pattern of a 9-mo old infant (RR=30, Vt=100 ml, I:E ratio=1:1.3). Nebulizers were charged with BIS 0.25 mg and run continuously until dry. Drug captured on a filter distal to the SAINT model was the lung dose. Compressor: PARI PRONEB TURBO. Nebulizer/mask systems studied: VIX1/aerosol mask (AM), PediNeb pacifier device (PN) or blow-by (BB); **AeroEclipse^{*}** neb and mask (AE); PARI LC+ and PARI LC*/ PARI Baby mask (PB), Fish mask (FM), and AE masks. The AE neb/mask was also studied with an ill breathing pattern (RR=50, Vt=100, I:E=1:2). **Results:** Lung dose ranged from 2.0 to 7.6% of the neb charge. Lung dose was AE (5.0%) > VIX1 (3.5%), LC+/FM (3.2%), LC*/PB (2.9%), and LC+/PB (2.8%). Also, VIX1/AM (3.5%)>VIX1/PN (2.5%)>VIX1/BB (2.0%). The lung dose of the LC+ and LC* more than doubled (6.8 and 6.3%) when used with the AE mask. Lung dose increased with the ill breath pattern in proportion to increased minute ventilation (7.6%). **Conclusion:** 1) The AE system provided higher lung dose than other nebulizers with standard masks. 2) Mask design and fit can substantially impact nebulizer performance. 3) PN performed better than BB, but not as good as a mask. If crying decreases lung dose by 75%, we speculate that the PN and BB (non-crying) may improve lung dose vs mask with a crying infant. 4) An increase in lung dose may occur in ill infants if minute ventilation is elevated. *This abstract is funded by: AstraZeneca*

2. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA A BREATH-ACTUATED SMALL VOLUME NEBULIZER (SVN): A COMPARATIVE IN-VITRO ASSESSMENT.** MW Nagel, KJ Wiersema, SL Bates and JP Mitchell. American Journal of Respiratory Care & Critical Care Medicine, 2001;163(5):A442.

Rationale: To compare the delivery of budesonide suspension in terms of fine particle dose (< 4.7 μ m aerodynamic diameter (FPD)) from a breath-actuated (BA) SVN with that from a continuous flow air entrainment (AE) SVN. **Methods:** FPD values were determined for 5-**AeroEclipse^{*}** BA SVNs (Monaghan Medical Corp., Plattsburgh, N.Y.) and 5-LC-D[™] AE SVNs (PARI Respiratory Equipment, Inc., Monterey, CA), nebulizing 4ml of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). Each SVN was operated with air at 50 psig, 8 l/min until sputtering occurred. Breathing parameters were: tidal volume= 600 ml, I:E=1:2 rate= 10/min. FPD was determined by cascade impactor at 28.3 ± 0.5 l/min. **Results:** From the beginning of nebulization until sputtering, the **AeroEclipse^{*}** and the LC-D[™] SVNs produced 164 ± 3 and 71 ± 4 μ g FPD of budesonide respectively. During the first 5 minutes (after which time the LC-D[™]s sputtered), values of FPD for the **AeroEclipse^{*}** and the LC-D[™] SVNs were 76 ± 4 and 71 ± 4 μ g budesonide respectively.

Conclusion: The **AeroEclipse*** was more efficient than the LCD™ SVN for this suspension formulation [Mann-Whitney rank sum test, $p < 0.001$]. Almost no medication delivery took place from the **AeroEclipse*** SVN during the exhalation portion of the breathing cycle, thereby providing important benefits to both patient and care giver.

| RESULTS | | |
|-------------------------|-----------|----------|
| Nebulizer | FILT (µg) | ENV (µg) |
| AeroEclipse* BAN | 283 ± 33 | 80 ± 11 |
| LCD™ | 97 ± 7 | 305 ± 2 |

mean ± SD

- DELIVERY OF A SUSPENSION CORTICOSTEROID FORMULATION BY SMALL VOLUME NEBULIZERS: A COMPARATIVE BENCH STUDY.** JP Mitchell, MW Nagel, KJ Wiersema and SL Bates. Accepted for presentation ERS Annual Congress, Berlin, Germany, September 2001, Abstract #290.

We report a study of the delivery of 0.25% mg/ml budesonide suspension (Pulmicort®, Nebuamp® (2 x 2-ml), Astra-Zeneca, Canada) by two types of small volume nebulizer (SVN), simulating adult breathing conditions ((tidal volume = 600-ml, duty cycle = 1:2 (2-s inspiration), PIFR = 31 l/min). Each SVN was operated by compressed air (8 l/min at 50 psig). Budesonide mass delivery was determined by filter collection (n = 5 SVN/group, 3-replicates/device). The **AeroEclipse*** BANs (Trudell Medical International, London Canada) delivered 283 ± 32 mg prior to sputtering, and 80 ± 11 mg were lost to the environment. Corresponding data for the LCD™ SVN (Pari Respiratory Equipment Inc., Richmond, VA, USA) were 97 ± 7 mg and 305 ± 2 mg respectively. The breath-actuation feature of the **AeroEclipse*** SVN minimizes aerosol release to the environment during exhalation, which may cause adverse effects to both patient and health care provider.

- ENHANCED *IN VITRO* DELIVERY OF BUDESONIDE VIA CONTINUOUS AND BREATH-ACTIVATED NEBULIZATION.** Saldone GC. European Respiratory Journal, August 2000; 16(31):540s.

In vitro bench testing designed to mimic clinical aerosol delivery is predictive of *in vivo* delivery of nebulized medications to the respiratory tract. This study tested a new nebulizer designed for either continuous or breath-actuated use (**AeroEclipse*** BAN, Monaghan/Trudell International). Using a piston pump and Pari Master compressor, a range of breathing patterns were utilized to estimate drug delivery [Inhaled mass (IM)] to pediatric patients over a wide range of breathing patterns. 500mg of budesonide comprised the nebulizer charge (0.25mg/ml in 2ml) delivered via three patterns of breathing (Vt f: 50ml, 40; 200ml, 25; 440ml, 19; duty cycle 0.50). The 50 and 200ml Vt patterns were delivered using continuous nebulization, while 440 was breath-actuated. IM was measured at 1 min intervals using a low deadspace filter with drug activity analyzed by HPLC. Low flow cascade impaction measured aerodynamic diameters (MMAD) and fine particle fraction (FPF, cutpoint 6.0µm). For the three breathing patterns IM averaged (mean ±SD), 11.1±0.74%, 22.9±2.74%, and 36.3±1.22% respectively. These values exceed by 35% those previously reported for the most efficient devices (J. Aerosol Med. 1998, 11:113-125). MMAD averaged 3.55±0.07µm, GSD 2.55 FPF 0.72. When corrected for FPF, pulmonary delivery is estimated to be 60% higher than that reported for conventional and air-entrained nebulization.

- THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS: A COMPARATIVE *IN-VITRO* ASSESSMENT.** JP Mitchell, MW Nagel, and AD Archer. Presented at Proc. Ann. Meet. Amer. College of Chest Physicians (ACCP), Toronto, November 1998, in Chest, 114 (4S), 295, (1998); and at the World Asthma Meeting, Barcelona, December 1998, in Eur. Resp. J., 12S29, 7, (1998).

Purpose: To compare the performance of a new air entrainment small volume nebulizer (AE-SVN, Trudell Medical International with other widely used SVN (LC-Star™ (PARI Respiratory Equipment), Updraft™ Neb-U-Mist™ (Hudson Oxygen Therapy Sales Co.), Circulaire™ (Westmed), Sidestream™ (Medic-Aid), AirLife™ Misty-Neb™ (Baxter Healthcare Corp.)) for the delivery of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). **Methods:** Each SVN (n = 5 devices for each group, 3 replicates per device) was operated with compressed air at 50 psig at a flow rate of 8 l/min. The total mass of budesonide nebulized from 2 x 2 ml ampoules was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The SVN was operated until it spluttered, was then tapped gently to dislodge droplets back to the reservoir. Nebulization was deemed complete 20 seconds later. The mass of budesonide collected was determined by HPLC-UV spectrophotometry. **Results:** The delivery rate ((mean ± 1 S.D) µg budesonide/min) from the AE-SVN (102 ± 9) was significantly greater than with the other groups: (LC-Star™ (91 ± 6), Misty-Neb™ (49 ± 2), Sidestream™ (46 ± 4), Circulaire™ (26 ± 4) and Neb-U-Mist™ (25 ± 6)), (1-way ANOVA, $p < 0.02$). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with LC-Star™ (229 ± 10 s), Sidestream™ (365 ± 19 s), Circulaire™ (420 ± 84 s), Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). **Conclusions:** The new AE-SVN is highly efficient at entraining the budesonide particles into the liquid droplets at these conditions. **Clinical Implications:** The good delivery rate combined with comparatively short duration of delivery offers the potential for rapid treatment and patient convenience.

6. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS - THE RELATIONSHIP BETWEEN NEBULIZED DROPLET SIZE AND THE PARTICLE SIZE OF THE SUSPENSION.** JP Mitchell, MW Nagel, and AD Archer. Presented at Drug Delivery to the Lungs-IX, London, UK, 1998, in J. Aerosol Med. 12(3), 208, (1999).

A new air entrainment small volume nebulizer (AE-SVN) has been compared with two other SVN's (Neb-U-Mist™ and Misty-Neb™) for the delivery of a suspension of 0.25 µg/ml budesonide. Each SVN was operated at 8 l/min with compressed oxygen (50 psig). The total mass of budesonide was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The time-averaged delivery rate over the period of nebulization ((mean ± 1 S.D.) µg budesonide/min) from the AE-SVN (102 ± 9) was greater than with the Misty-Neb™ (49 ± 2), or Neb-U-Mist™ (25 ± 6). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with the Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). The mass median diameter (MMD) of the droplets from the AE-SVN measured using a laser diffractometer (2.9 ± 0.1 µm), was significantly finer compared with those from the Misty-Neb™ (4.5 ± 0.9 µm) and Neb-U-Mist™ (5.6 ± 0.6 µm) and closest to the size of the micronized budesonide particles in the original suspension. The efficient delivery of medication formulated as micronized powder in aqueous suspension necessitates that the droplets produced upon nebulization are large enough so that single particles are efficiently entrained during atomization, but not so coarse that they cannot leave the nebulizer, extending nebulization time.

Combivent™ (ipratropium bromide & albuterol sulfate), Boehringer Ingelheim™ Pharmaceuticals Inc.

1. **COMPARISON IN RATES OF BREAKTHROUGH TREATMENTS DURING A CONVERSION FROM RACEMIC ALBUTEROL TO LEVALBUTEROL.** RS Pikarsky, RA Acevedo, T Farrell, R Bear and W Fascia. Respiratory Care, Nov. 2003;48(11):1080.

Purpose: In order to meet our adult patient care demands, Crouse Hospital approved an automatic conversion from Racemic Albuterol to Levalbuterol. This study compares the breakthrough rates of Racemic Albuterol and Levalbuterol, with and without Ipratropium. Different dosing schedules for Levalbuterol were evaluated. **Methods:** Racemic Albuterol (Alb) 2.5 mg Q4h was converted to either Levalbuterol (Lev) 0.63 mg Q6h or Levalbuterol 0.63 mg Q8h. Patients dosed Q8h who required more frequent aerosol administration received Levalbuterol 0.63 mg Q6h (cardiac patients) or Levalbuterol 1.25 mg Q8h (all others). If ordered, Ipratropium (Ipra) 0.5 mg was administered at the same frequency as the Levalbuterol. A majority of aerosol therapy was provided with the use of the **AeroEclipse*** Breath Actuated Nebulizer (BAN). All aerosol treatments, including breakthrough treatments, delivered between June 1, 2002 and September 30, 2002 were recorded. **Results:** Tx/Pt/day represents the number of treatments delivered per patient per day. Rate/100 Pt/days = (Breakthrough) / (Total Tx / Tx/Pt/day) x 100. Rate/100 Pt/days corrects for the differences in daily administration frequency, and may better reflect the daily impact of the breakthrough rate. The breakthrough rate of the combined Albuterol group was significantly greater than all Levalbuterol groups (25.8 vs. 18.43, 25.8 vs. 18.43, 25.8 vs. 5.96 p<.001)*. The breakthrough rate with Albuterol was significantly reduced with the addition of Ipratropium (40.76 vs. 13.35 p<.001)**. The 1.25 mg dose of Levalbuterol outperformed both 0.63 mg dosage groups (3.78 vs. 13.48 p<.02, 3.78 vs. 21.36 p<.001)***. Ipratropium did not significantly change the breakthrough rate when added to Levalbuterol groups.

| Medication | Total Tx | Breakthrough | Rate/1000 | Tx/Pt/day | Rate/100 Pt/day |
|----------------------|----------|--------------|-----------|-----------|-----------------|
| Alb Q4h | 898 | 61 | 67.93 | 6 | 40.76** |
| Alb/Ipra Q4h | 1079 | 24 | 22.24 | 6 | 13.35** |
| Lev 0.63mg Q6h | 2047 | 69 | 33.71 | 4 | 13.48*** |
| Lev 0.63 mg/Ipra Q6h | 2728 | 151 | 55.35 | 4 | 22.14 |
| Lev 0.63mg Q8h | 660 | 47 | 71.21 | 3 | 21.36*** |
| Lev 0.63 mg/Ipra Q8h | 707 | 37 | 52.33 | 3 | 15.70 |
| Lev 1.25mg Q8h | 238 | 3 | 12.61 | 3 | 3.78*** |
| Lev 1.25mg/Ipra Q8h | 215 | 6 | 27.91 | 3 | 8.37 |

Conclusions: The conversion from Racemic Albuterol to Levalbuterol allowed for a decreased frequency of daily medication administrations and a significant decrease in breakthrough requirements. Levalbuterol at the 1.25 mg dose performed better than the 0.63 dose for Q8h administration. Ipratropium showed a significant benefit in breakthrough reduction for the Racemic Albuterol group. **Clinical Implications:** The efficiencies gained by decreasing the daily frequency of aerosol administration can have a significant impact on resource utilization. The conversion to Levalbuterol allows for decreased respiratory therapy time or the re-allocating of workforce needs while maintaining, or improving, quality of aerosol administration, as evidenced by the decrease in breakthrough requirements.

2. **PERFORMANCE OF A NEW BREATH-ACTUATED SMALL VOLUME NEBULIZER (SVN) FOR THE DELIVERY OF A COMBINATION ANTICHOLINERGIC/BRONCHODILATOR.** MW Nagel, KJ Wiersema, SL Bates and JP Mitchell. American Journal of Respiratory and Critical Care Medicine, 2001;163(5):A443.

Purpose: To compare the delivery of ipratropium bromide (IPR) and albuterol sulfate (ALB) as fine droplets (< 4.8 µm diameter (FPD)) and as total emitted dose (ED) from a breath-actuated (BA- SVN) with that from a continuous flow air entrainment (AE-SVN) after 5-minutes of operation. **Methods:** FPD and ED were determined for 5-**AeroEclipse*** BAN (Monaghan Medical Corp., N.Y.) and 5-PARI LCD™ SVN's (PARI Respiratory Equipment, Inc., CA) nebulizing Combivent® (2.5-ml, 0.2 mg/ml IPR and 1.0 mg/ml ALB; Boehringer-Ingelheim (Canada) Inc.). Each SVN was operated

with 8 l/min air at 50 psig, simulating breathing at tidal volume, I:E ratio and rate of 750-ml, 1:2 and 10/min respectively. Droplet size distributions were measured by laser diffractometer. **Results:** (ED) and (FPD) were as follows:

| | | | |
|-----|-------------------------|------------------|-------------------|
| IPR | AeroEclipse* BAN | ED = 102 ± 7 µg | FPD = 82 ± 6 µg |
| IPR | PARI LCD™ SVNs | ED = 55 ± 7 µg | FPD = 45 ± 5 µg |
| ALB | AeroEclipse* BAN | ED = 581 ± 17 µg | FPD = 471 ± 14 µg |
| ALB | PARI LCD™ SVNs | ED = 279 ± 33 µg | FPD = 226 ± 26 µg |

Differences in ED and FPD between SVN for IPR and ALB components were statistically significant (unpaired t-test for each variable, $p < 0.001$). Mass median aerodynamic diameters were close to 2.8 µm for both SVN groups. **Conclusion:** The **AeroEclipse*** BAN is significantly more efficient for the delivery of this combination anticholinergic/bronchodilator than a conventional AE-SVN.

3. **CLINICAL EVALUATION OF A BREATH ACTUATED SMALL VOLUME NEBULIZER (BA-SVN).** S Klopff, N Schneiderman, H Payne, C Schramm, MW Nagel and JP Mitchell. Respiratory Care, August 2000, 45(8).

Background: In prior *in-vitro* studies using laser diffractometry, the aerosol produced by a novel breath-actuated nebulizer (BAN), the **AeroEclipse*** (Monaghan Medical Corp. Plattsburgh, NY) has been shown to contain a high proportion of droplets < 4.8 µm diameter (80.9% ± 2.4%). Such droplets are more likely to penetrate beyond the oropharyngeal region where bronchodilation is achieved. These *in-vitro* results should therefore be predictive of improved *in-vivo* delivery of nebulized medications to the respiratory tract. This study explored the clinical performance of the **AeroEclipse*** BAN in the delivery of a beta2-agonist (albuterol 2.5 mg/ml) accompanied by anticholinergic (ipratropium bromide 250 µg/ml) bronchodilator in some cases. **Methods:** Patients (n=48) with a previous diagnosis for asthma presenting to the Emergency Department for acute exacerbation of asthma were included in this study. Upon presentation, an asthma care path, an assessment driven, algorithm-based tool was used to place patients in one of three stages of severity as recommended by the NIH-NAEPP Guidelines for the Diagnosis of Asthma. Each patient was assigned to receive inhaled aerosol treatment using the **AeroEclipse*** BAN. Stage 1 asthmatics were given 0.5-ml of albuterol with 0.5-ml normal saline delivered until sputter. Patients categorized in stage two and three were given 0.5-ml albuterol with the addition of 1.5-ml of ipratropium bromide unit dose. Treatments repeated every 20 minutes times three if necessary by protocol.

Results:

| Asthma Severity | Stage 1 | Stage 2 | Stage 3 |
|-----------------------------------|-----------|----------------|---------------|
| Number | 10 | 30 | 8 |
| Treatments Given | 2.4 | 2.03 | 2.25 |
| Treatment Duration (min) | 3.7 | 3.78 | 5 |
| Increase in PEF (mean, range (%)) | 44(0-120) | 67.7(-2.7-580) | 120.7(28-420) |

Four patients had greater than 20% increase in heart rate, three patients noted tremor following treatment. Twenty four patients had positive comments about the device focused on shorter treatment time and improved relief from dyspnea. Two imminent intubations were avoided with the use of the BA-SVN. **Conclusions:** Use of the **AeroEclipse*** BAN appears to result in good clinical outcomes. Minimum number of treatments, shorter treatment duration and minimal side effects were noticed with this device. Further outcome studies are needed to assess this impact on other groups of patients.

Atrovent™ (ipratropium bromide), Boehringer Ingelheim™ Pharmaceuticals Inc.

1. **A MECHANICALLY OPERATED BREATH-ACTUATED NEBULIZER ENABLES BOTH IMPROVED CONTROL OF DOSING AND DELIVERY EFFICIENCY.** JP Mitchell, MW Nagel. Presented at Drug Delivery to the Lungs Conference, December 2005.

A mechanically operated, breath-actuated nebulizer (BAN) offers the clinician the prospect of being able to control the rate and duration of medication delivery dosimetrically, providing greater precision when titrating patients to establish an appropriate treatment regimen. We describe an *in vitro* study obtained with two formulations that are representative of formulations available for nebulization (amphotericin-B and ipratropium bromide), in which a BAN (**AeroEclipse***) delivered slightly more medication as fine droplets < 4.8 µm aerodynamic diameter with approximately one-half of the dose in the reservoir compared with a continuously operating nebulizer (VixOne™). These measurements were made simulating use by an adult (500-ml tidal volume, inspiratory/expiratory ratio 1:2, 20 breaths/minute). Significant cost savings are therefore possible with the BAN with expensive medications, such as antibiotics, if less volume fill is required per treatment.

Alupent™ (metaproterenol sulphate), Boehringer Ingelheim™ Pharmaceuticals Inc.

1. **PERFORMANCE OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER.** A Archer, JP Mitchell, MW Nagel and AMW Verdun. Eur. Resp. J., 1998; 12(28): 68.

We report an *in vitro* investigation in which the performance of a new disposable AE-SVN (n = 3 devices) has been assessed with salbutamol sulphate (Ventolin®: 5 µg/2.5 ml, GlaxoSmithKline Inc.), metaproterenol sulphate (Alupent®: 10 µg/2.5 ml, Boehringer Ingelheim Pharmaceuticals Inc.) and cromolyn sodium (Intal®: 20 µg/2 ml, Fisons Pharmaceuticals) nebulizers. Each AE-SVN was filled with 2 nebulizers and operated continuously with oxygen supplied at 50 psig and 8 l/min. The AE-SVN was coupled directly to an Andersen cascade impactor, sampling at 28.3 l/min. Total and fine particle (< 4.7 µm aerodynamic diameter) delivery rates were 33.5 ± 1.8 µg/s and 27.6 ± 1.3 µg/s (Ventolin®); 54.2 ± 10.6 µg/s and 45.0 ± 7.8 µg/s (Alupent®); 138.6 ± 10.2 µg/s and 109.7 ± 8.3 µg/s (Intal®) over a 10 s period following the start of nebulization. The mass median aerodynamic diameter (MMAD) and mass % contained in fine droplets were 3.0 ± 0.1 µm and 82.4 ± 1.2% (Ventolin®); 2.9 ± 0.2 µm and 83.3 ± 2.6% (Alupent®); 3.1 ± 0.1 µm and 79.2 ± 1.9 % (Intal®). This new nebulizer appears to perform well with all three formulations.

Intal™ (cromolyn sodium), Fisons™ Pharmaceuticals Ltd.

1. **EVALUATION OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) FOR THE DELIVERY OF ALBUTEROL SULFATE AND CROMOLYN SODIUM.** JP Mitchell, MW Nagel, A Archer and D Coppola. Presented at 1999 ALA/ATS International Conference, San Diego, CA, 1999, in Am. J. Resp. Crit. Care Med, 1999;159(3):A120.

Purpose: To evaluate the delivery of Ventolin® (0.2% v/v, albuterol sulfate, GlaxoSmithKline, Canada) and Intal® (1.0% v/v cromolyn sodium, Fisons Pharmaceuticals Ltd., Canada) by a prototype AE-SVN (Trudell Medical International) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 AE-SVNs were tested using an Andersen Mark II Cascade Impactor operated at 28.3±0.5 l/min to determine the size distribution of droplets emitted at the mouthpiece during the first 10 seconds following nebulization. The mass of drug emitted was determined directly by HPLC-UV spectrophotometry. **Results:** Total (TM) and fine particle ((FPM), droplets finer than 4.7 µm diameter) mass output rates and droplet mass median diameter (MMD) were as follows:

| Drug | TM (µg/s) | FPM (µg/s) | MMD (µm) |
|-----------|--------------|-------------|-----------|
| Ventolin® | 32.4 ± 3.1 | 27.6 ± 1.3 | 3.0 ± 0.1 |
| Intal® | 138.6 ± 10.2 | 109.7 ± 8.3 | 3.2 ± 0.1 |

Conclusion: The fine MMD produced from the AE-SVN resulted in an improved FPM output rate, which is likely to produce increased lung deposition.

2. **PERFORMANCE OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER.** A Archer, JP Mitchell, MW Nagel, and AMW Verdun. Eur. Resp. J., 1998; 12(28): 68.

We report an *in vitro* investigation in which the performance of a new disposable AE-SVN (n = 3 devices) has been assessed with salbutamol sulphate (Ventolin®: 5 µg/2.5 ml, GlaxoSmithKline Inc.), metaproterenol sulphate (Alupent®: 10 µg/2.5 ml, Boehringer Ingelheim Pharmaceuticals Inc.) and cromolyn sodium (Intal®: 20 µg/2 ml, Fisons Pharmaceuticals) nebulizers. Each AE-SVN was filled with 2 nebulizers and operated continuously with oxygen supplied at 50 psig and 8 l/min. The AE-SVN was coupled directly to an Andersen cascade impactor, sampling at 28.3 l/min. Total and fine particle (< 4.7 µm aerodynamic diameter) delivery rates were 33.5 ± 1.8 µg/s and 27.6 ± 1.3 µg/s (Ventolin®); 54.2 ± 10.6 µg/s and 45.0 ± 7.8 µg/s (Alupent®); 138.6 ± 10.2 µg/s and 109.7 ± 8.3 µg/s (Intal®) over a 10 s period following the start of nebulization. The mass median aerodynamic diameter (MMAD) and mass % contained in fine droplets were 3.0 ± 0.1 µm and 82.4 ± 1.2% (Ventolin®); 2.9 ± 0.2 µm and 83.3 ± 2.6% (Alupent®); 3.1 ± 0.1 µm and 79.2 ± 1.9 % (Intal®). This new nebulizer appears to perform well with all three formulations.

Methacholine chloride

1. **AN IN VITRO STUDY TO INVESTIGATE THE USE OF A BREATH-ACTUATED, SMALL-VOLUME, PNEUMATIC NEBULIZER FOR THE DELIVERY OF METHACHOLINE CHLORIDE BRONCHOPROVOCATION AGENT.** JP Mitchell, MW Nagel, SL Bates and CC Doyle. Respir Care 2003;48(1):46–51.

Background: Current American Thoracic Society and American Association for Respiratory Care guidelines for the delivery of aerosol agents such as methacholine chloride (MC) for bronchoprovocation testing require the use of pneumatic jet nebulizers that have well-defined droplet size and mass output. A recently developed disposable, breath-actuated nebulizer (**AeroEclipse***) may offer bronchoprovocation testers an alternative to existing devices. **Methods:** We studied the performance of 5 **AeroEclipse*** nebulizers with regard to mass of MC delivered with various MC solution concentrations and numbers of inhalations, using a model of adult tidal breathing. Each nebulizer was operated with compressed air (8 L/min at 50 psig) and an initial fill of 2 mL. MC solutions with mass concentrations of 0.25, 0.98, 3.85, and 15.70 mg/mL were tested. The total mass of MC delivered was determined after 5, 10, and 15 complete breathing cycles, by assaying the MC collected on a filter placed at the nebulizer mouthpiece. The aerosol droplet size distribution, fine droplet fraction (FDF) (percentage of droplets < 4.8 µm diameter), and fine droplet mass (FDM) (mass of droplets <

4.8 µm diameter) were determined by laser diffractometry, using physiologically normal saline as a surrogate for MC solution. **Results:** The mean ± SD FDM collected in 5 breathing cycles was 654 ± 29 µg with the 15.70 mg/mL solution, 158 ± 9 µg with the 3.85 mg/mL solution, 37 ± 3 µg with the 0.98 mg/mL solution, and 7 ± 2 µg with the 0.25 mg/mL solution. FDM showed a linear correlation ($r^2 = 0.9999$) with MC concentration, within the range studied. FDM also showed a linear correlation ($r^2 = 0.999$) with the number of breathing cycles. For instance, with the 15.70 mg/mL solution, FDM was 654 ± 29 µg with 5 breathing cycles, 1,228 ± 92 µg with 10 breathing cycles, and 1,876 ± 132 µg with 15 breathing cycles. **Conclusions:** Although the bronchoprovocation test procedure had to be slightly modified from the guidelines to accommodate the operation of the **AeroEclipse***'s breath-actuation feature, our measurements indicate that a predictable dose of MC, within the useful range for bronchoprovocation testing, can be delivered to an adult patient breathing tidally. The green indicator on the **AeroEclipse*** could be used to coach the patient to inhale for a specific period, thereby controlling MC delivery per breathing cycle.

2. **PREDICTING LUNG DEPOSITION WITH A CASCADE IMPACTOR.** S Sangwan, F Hull, R Condos and GC Smaldone. Journal of Aerosol Medicine, 2001; 14(3):421. Presented at the 13th International Congress on Aerosols in Medicine, Interlaken, Switzerland, September 17-21, 2001.

Introduction: In recent deposition studies of interferon-β, we failed to predict the deposition pattern from bench studies of aerosols using multistage cascade impaction (MCI). Recent mass balance studies have identified impaction in connecting tubing and effects of breathing on interpretation of cascade data (Gurses BK et al AJRCC 163; 5(A166). 2001). In the present study we related MCI data using our new bench test protocol directly to lung scans in humans. This protocol emphasizes deposition of large particles in connecting tubing and influence of conditions internal to the nebulizer during breathing. **Methods:** Two devices (Misty-Neb and **AeroEclipse*** Breath-Actuated Nebulizer ("BAN")) were studied. Mass median aerodynamic diameter (MMAD) and mass balance were measured under standing cloud and ventilation using a piston pump. Deposition images were obtained using gamma camera.

Results:

| Nebulizer & method of assessment | | Respirable Mass* (<6µm) | Regional Deposition | |
|----------------------------------|----------------|-------------------------|---------------------|---------------------|
| | | | Lung deposition** | Throat deposition** |
| Misty-Neb | Standing Cloud | 46.2% | 32% | 68% |
| | Ventilated | 24.6% | | |
| AeroEclipse* BAN | Standing Cloud | 48.3% | 72% | 28% |
| | Ventilated | 71.2% | | |

*Calculated by adding T connector deposition to the first stage (>8µm) of cascade ** Expressed as Percent of total deposition in the body

Conclusion: Regional deposition (upper airway vs. lung) was predicted by analysis only when effects of both connecting tubing and breathing were considered in the bench protocol.

Ablecet™ (Amphotericin), Enzon™ Pharmaceuticals

1. **AEROSOLIZED AMPHOTERICIN B LIPID COMPLEX (aABLC) DISTRIBUTION IN LUNG TRANSPLANT RECIPIENTS: A COMPARISON OF CONTINUOUS VERSUS BREATH ACTUATED NEBULIZERS.** ES Dodds, NA Petry, JD Davies, DW Zaas, SM Palmer, SW Shipes, RH Drew, BD Alexander, RE Coleman, JR Perfect. Presented at the American Association for Respiratory Care Congress, Orlando, FL., December 2007.

Background: Aerosolized amphotericin B has become an attractive option for antifungal prophylaxis following solid organ and stem cell transplantation.^{1,2} This therapeutic strategy facilitates localized delivery of antifungal agent, thereby minimizing toxicities and drug-drug interactions associated with currently-available systemic antifungal agents. Determining drug delivery characteristics, including dose and nebulizer system, for aerosol drug administration is important to ensure optimal drug delivery. Newer, breath-actuated nebulizers (BAN's) are available and, in theory, provide the ability to limit environmental exposure and also deliver a higher percentage of the prepared dose to the patient. **Objective:** To characterize the distribution of aerosolized ABLC immediately following nebulization in bilateral lung transplant recipients via 2 different nebulizer systems – continuous nebulizer (CN): Up-Draft, Model 1724 (Hudson RCI, Temecula, CA) and breath actuated nebulizer (BAN): **AeroEclipse*** II (Monaghan Medical, Plattsburgh, NY). ABLC 20 mg/4mL was mixed with prepared 99mTc-ABLC (Ablecet®-Enzon Pharmaceuticals) prior to loading into the radioaerosol delivery system. **Methods:** Nebulizer assignment was performed sequentially with the first 5 subjects receiving treatment via the continuous flow nebulizer and the subsequent 5 subjects receiving study drug treatment via the BAN. Immediately following inhalation, drug product distribution image were obtained with patients in the supine position. Subjects were then placed on the table of a dual-head gamma camera system (General Electric Healthcare, Milwaukee, WI). Total delivered dose (TDD) was calculated by determining the difference in the known starting counts for the medication vial and counts of the nebulizer apparatus, including filter, subject waste materials and empty medication vials, obtained after study medication administration. Gastric activity of 99mTc-ABLC was also measured. Drug exposure was reported as: TDD: total delivered dose; Drug delivery to each of the following lung regions was reported as a

percentage of TDD: right lung (RL), left lung (LL) and GI tract; the two nebulizer groups were compared for differences in mean TDD and regional distribution using student's t-test.

Results: Total drug delivery (as percent of prepared dose) was significantly higher for the BAN (20.7% versus 3.5%, $p=0.01$). Mean regional distribution (as percent of total delivered dose) did not differ between the two nebulizer devices for the left lung, right lung, or GI tract.

| Subject | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 |
|---------------------------|--------------------------------|-----|-----|-----|-----|----------------------------------|------|------|------|------|
| | Continuous Nebulizer | | | | | Breath Actuated Nebulizer | | | | |
| Drug Delivery* | % of total dose in vial | | | | | % of total dose in vial | | | | |
| RL | NR | 1.6 | 1.2 | 0.4 | 1.2 | 7.4 | 9.6 | 5.2 | 5.8 | 11.3 |
| LL | NR | 1.4 | 0.9 | 0.3 | 0.7 | 6.4 | 5.5 | 5.4 | 6.0 | 8.9 |
| GI | NR | 3.6 | 1.3 | 0.6 | 0.5 | 5.1 | 5.1 | 7.2 | 11.2 | 3.5 |
| Total Drug Delivery (TDD) | NR | 6.6 | 3.4 | 1.3 | 2.4 | 18.9 | 20.2 | 17.8 | 23.0 | 23.7 |
| Regional Delivery** | | | | | | | | | | |
| Right | 50 | 24 | 35 | 31 | 49 | 39 | 47 | 29 | 25 | 48 |
| Left | 17 | 21 | 27 | 23 | 29 | 34 | 27 | 30 | 26 | 37 |
| Esophagus and Stomach | 32 | 55 | 39 | 46 | 22 | 27 | 25 | 40 | 49 | 15 |

* As percent of prepared dose

** As a percentage of the total delivered dose

Conclusion: Use of the BAN resulted in a larger portion of the drug being deposited into the lungs. Since GI distribution was similar between the nebulizers, it appeared that more drug was vented to the surrounding atmosphere with the continuous system.

References: ¹Hussain S, Zaldonis D, Kusne S et al. Variation in antifungal prophylaxis strategies in lung transplantation. *Transpl Infect Dis* 2006;213-8. ²Drummer JS. A survey of fungal management in lung transplantation. *Journal of Heart and Lung Transplantation* 2004;23:1376-81.

2. **SIMILAR DELIVERY OF AMPHOTERICIN LIPID COMPLEX IS POSSIBLE AT ONE-HALF DOSE VIA A BREATH-ACTUATED NEBULIZER COMPARED WITH A CONTINUOUSLY OPERATING NEBULIZER.** JP Mitchell, NR MacIntyre, MW Nagel, DP Coppola. American Thoracic Society 101st International Congress, May, 2005.

Delivery of aerosolized antibiotics via continuous nebulizers wastes these expensive medications during patient exhalation. Breath-actuated nebulizers (BAN) can minimize waste with significant cost savings in medication, since they only operate when the patient inhales. Furthermore, medication is not emitted into the environment during exhalation. We describe a study in which dose delivery from a BAN (**AeroEclipse***, Monaghan Medical Corp., Plattsburgh, NY) was compared with that from a continuously operating nebulizer (VixOne, Westmed Corp., Englewood, CO (VIX)) (n=3/group) for the delivery of amphotericin lipid complex ((AMP) Ablecet, Enzon Pharmaceuticals, Piscataway, NY, 5-mg/ml). Each device was operated with air at 50 psig at 7 L/min (BAN) or 8 L/min (VIX), with the mouthpiece connected to a breathing simulator (Compass, PARI, Germany) set to replicate adult use (500-ml tidal volume, 1:2 inspiratory/expiratory ratio, 20-breaths/min). 5-ml AMP was placed in the BAN and 10-ml in the VIX (5-ml initially, followed by a further 5-ml after 4-min). Each nebulizer was operated for 1-min past first sputter. The mass of AMP collected on a filter at the mouthpiece was determined by HPLC-UV spectrophotometry (3-replicates/nebulizer). Droplet size distributions were determined by laser diffractometer in a separate study. Total emitted mass from the BAN was 7274 123 g, delivered in 10-min, of which 5892 100 g was in fine droplets 4.8 m diameter. The VIX delivered a total mass of 5276 557 g in 10-14 min, of which 4326 457 g was contained in fine droplets. The BAN was therefore capable of delivering 36% more medication as fine droplets with only one-half of the dose inserted in the reservoir.

3. **A MECHANICALLY OPERATED BREATH-ACTUATED NEBULIZER ENABLES BOTH IMPROVED CONTROL OF DOSING AND DELIVERY EFFICIENCY.** JP Mitchell, MW Nagel. Presented at Drug Delivery to the Lungs Conference, December 2005.

A mechanically operated, breath-actuated nebulizer (BAN) offers the clinician the prospect of being able to control the rate and duration of medication delivery dosimetrically, providing greater precision when titrating patients to establish an appropriate treatment regimen. We describe an *in vitro* study obtained with two formulations that are representative of formulations available for nebulization (amphotericin-B and ipratropium bromide), in which a BAN (**AeroEclipse***) delivered slightly more medication as fine droplets < 4.8 μ m aerodynamic diameter with approximately one-half of the dose in the reservoir compared with a continuously operating nebulizer (VixOne[™]). These measurements were made simulating use by an adult (500-ml tidal volume, inspiratory/expiratory ratio 1:2, 20 breaths/minute). Significant cost savings are therefore possible with the BAN with expensive medications, such as antibiotics, if less volume fill is required per treatment.

Measles Vaccine (Placebo)

1. **THE DELIVERY OF PLACEBO MEASLES VACCINE BY A MECHANICALLY-OPERATED BREATH-ACTUATED NEBULIZER (BAN).** JP Mitchell, MW Nagel. Presented at European Respiratory Society, Munich, Germany, September 2006.

Nebulizer-delivered vaccination offers the potential for the mass immunization of children. We report the outcome of a study in which the delivery of a placebo measles vaccine by a novel BAN (**AeroEclipse**[®], Trudell Medical International) was evaluated in comparison with a continuously operating jet nebulizer (Aeromist[™], IPI Medical Products Inc., Chicago, USA), used successfully to deliver aerosol in the so-called Classic Mexican Device (CMD) in previous World Health Organization (WHO) - sponsored studies. Each nebulizer (n=5 devices/group) was operated by portable compressor (Pulmomate[™], De Vilbiss Corp.), with a 3-ml fill of reconstituted placebo vaccine in sterile water. The emitted droplets were drawn at 30 L/min \pm 5% through an electret filter located at the distal end of either a 15-cm length of corrugated tubing forming the outlet of the CMD, or a 5-cm tube with inhalation valve attached to the BAN. Mass output rate was quantified gravimetrically, and a laser diffractometer was used to determine droplet size distributions. The aerosol produced by the BAN (mass median diameter (MMD) = $4.3 \pm 0.23 \mu\text{m}$) was finer than the mass output rate of the BAN ($0.40 \pm 0.01 \text{ ml/min}$) significantly exceeded that from the CMD ($0.15 \pm 0.03 \text{ ml/min}$) ($p < 0.001$). The BAN is dosimetric, so that an estimated mass output/breath close to that from the CMD can be anticipated when used by a tidally breathing patient with duty cycle of 33%. Furthermore, the breath actuation feature avoids the risk of exposing the health care giver to medication when the patient is not inhaling.

Interferon Gamma

1. **IMMUNOMODULATION WITH PHARMACOLOGIC IFN-GAMMA AND ITS EFFECT ON THE LUNG-SPECIFIC IMMUNE RESPONSE IN PULMONARY TB.** R Condos, ML Huie, R Dawson, S Ress, C Brauns, CH Tseng, M Weiden, E Bateman, RN Rom. Presented at American Thoracic Society, San Francisco, CA, September 2007.

Background: In a randomized clinical trial of TB patients treated with interferon gamma (IFN- γ), we have shown safety and efficacy (faster culture conversion). We hypothesize that pharmacological IFN- γ stimulates a TH1 environment in situ in the lung. **Methods:** 24 patients with cavitary TB randomized to DOTS alone or DOTS plus IFN- γ (either by aerosol or by sc injection). Bronchoscopy done at baseline and 16 weeks of treatment. BAL cell differential and 24 hour supernatants were prepared and spontaneous expression of cytokines/chemokines were assayed by Beadlyte multiplex assay on the Luminex 200 platform. Results were reported as averages SEM. **Results:** 12 patients were randomized to DOTS plus aerosol IFN- γ ; 5 patients were randomized to DOTS plus sc IFN- γ ; and 5 were randomized to DOTS alone. BAL cell differentials showed an increase in % lymphocytes in all groups (10.3% pre, 22.3% post). Several cytokines/chemokines were differentially expressed between groups. Eotaxin increases with IFN treatment (47.1 pg/ml to 92.4 pg/ml) but not with DOTS alone (64.3 pg/ml to 61.1 pg/ml). IL-4 was low in all patients (pre- 5.1 to post- 9.2 pg/ml). IL-1 decreased with IFN- treatment (186.132 to 21.7 pg/ml), but increased on DOTS alone (20.7 to 163.156 pg/ml) as did TNF- (IFN- group: 119.85 to 43.28 pg/ml ; DOTS alone 13.5 to 202.200 pg/ml) and MIP1. IFN- increased in the aerosol group (148.111 pg/ml to 229.85 pg/ml) and the DOTS only group (38.16 pg/ml to 111.76 pg/ml), but not in the sc group (217.96 pg/ml to 99.40 pg/ml). IP-10 levels increased in all groups (117.55 to 401.93 pg/ml). **Conclusion:** Immunomodulation with IFN- leads to a decrease in pro-inflammatory chemokines/cytokines independent of changes in cell differential or IFN- levels.

Saline

1. **EVALUATION OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) WHEN USED WITH OXYGEN AS A DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** JP MITCHELL and MW NAGEL. Presented at ALA/ATS International Conference, Chicago, 1998.

Purpose: To compare the delivery of saline (0.9% w/v NaCl) by a new AE-SVN (Trudell Medical Int.) with that from two other representative SVN's (UpDraft Neb-U-Mist[®] (Hudson Oxygen Therapy Sales Co.) and Airlife[™] Misty-Neb[™] (Baxter Healthcare Corp.)) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 AE-SVN's were tested using a laser diffractometer (Malvern Mastersizer-X) to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. 5, Neb-U-Mist[®] and a similar number of MistyNeb[™] SVN's were also evaluated. **Results:** Total (TM) and respirable ((RM), droplets finer than $4.8 \mu\text{m}$ diameter) mass output rates and droplet mass median diameter (MMD) were as follows: AE-SVN: TM = $671 \pm 26 \mu\text{g/min}$, RM = $542 \pm 23 \mu\text{g/min}$ ($80.8 \pm 1.3\%$ respirable), MMD = $2.88 \pm 0.09 \mu\text{m}$; Neb-U-Mist[™]: TM = $266 \pm 13 \mu\text{g/min}$, RM = $119 \pm 16 \mu\text{g/min}$ ($42.1 \pm 5.2\%$ respirable), MMD = $5.6 \pm 0.6 \mu\text{m}$; Misty-Neb[™]: TM = $336 \pm 60 \mu\text{g/min}$, RM = $178 \pm 43 \mu\text{g/min}$ ($53.1 \pm 8.5\%$ respirable), MMD = $4.5 \pm 0.9 \mu\text{m}$. **Conclusion:** TM from the new AE-SVN was substantially greater than those from either the Neb-U-Mist[®] or Misty-Neb[™] (1-way ANOVA, $p < 0.001$). The finer MMD produced from the AE-SVN resulted in a significantly greater RM compared with either of the other SVN's ($p < 0.001$).

2. **PERFORMANCE OF A NEW BREATH-ACTUATED SMALL VOLUME NEBULIZER (BA-SVN) WHEN USED WITH OXYGEN AS A DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** AM Verdun, JP Mitchell and MW Nagel. Presented at ALA/ATS International Conference, Chicago, 1998.

Purpose: To compare the delivery of saline (0.9% w/v NaCl) by a new BA-SVN (Trudell Medical Int.) with that from two other representative SVN's (LC-JET™ (PARI Respiratory Products Inc., Canada) and reusable Sidestream™ (MedicAid, UK)) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 BA-SVN's were tested using a laser diffractometer (Malvern Mastersizer-X) to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. 5, LC-JET™ and 5, Sidestream™ SVN's were also tested similarly. The BA-SVN was operated with manual over-ride engaged (continuous delivery of aerosol).

Results: Total (TM) and respirable (RM), droplets finer than 4.8 µm diameter) mass output rates and droplet mass median diameter (MMD) were as follows: BA-SVN: TM = 672 ± 23 µg/min, RM = 545 ± 31 µg (80.9 ± 2.4% respirable), MMD = 2.79 ± 0.15 µm; LC-JET™: TM = 675 ± 69 µg/min, RM = 449 ± 41 µg/min (66.7 ± 1.8% respirable), MMD = 3.39 ± 0.08 µm; Sidestream™: TM = 442 ± 26 µg/min, RM = 358 ± 38 µg/min (80.8 ± 4.2 %respirable), MMD = 2.94 ± 0.03 µm.

Conclusion: Although TM from the new BA-SVN was comparable with that from the LC-JET™ (Mann-Whitney rank sum test, p = 0.84), the finer MMAD produced from the BA-SVN resulted in a significantly greater RM (p < 0.001). Both TM and RM from the BA-SVN were greater than those from the Sidestream™ SVN (p < 0.001).

3. **COMPARISON OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) WITH OTHER SVN'S WHEN USED WITH OXYGEN AS DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** JP Mitchell and MW Nagel. Presented at Ann. Meet. Amer. Assoc. of Asthma, Allergy and Immunology (AAAAI), Washington D.C., 1998.

The performance of a prototype novel AE-SVN (Trudell Medical International (n = 5)) with normal saline (0.9% w/v NaCl) operating at 20 ± 2°C, 50 ± 10% RH, has been evaluated with oxygen (50 psig, 8 l/min) as driving gas to simulate hospital use. Comparison testing was also undertaken with two other representative AE-SVN's, (a) LC-JET™ (Pari Respiratory Equipment Inc.), without inspiratory valve cap which would otherwise restrict aerosol output, (b) SideStream™ (MedicAid, UK). A laser diffractometer (Malvern Mastersizer-X) was used to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. Total (T) and respirable ((R), droplets finer than 4.8 µm aerodynamic diameter) mass output rates and droplet mass median aerodynamic diameter (MMAD) for the new AE-SVN (5 replicate measurements/device) were: 671 ± 26 µg/min (T), 542 ± 23 µg/min (R) and 2.88 ± 0.09 µm (MMAD). Corresponding data for the LC-JET™ were: 675 ± 65 µg/min (T), 450 ± 45 µg/min (R) and 3.39 ± 0.14 µm (MMAD), and for the SideStream™ were: 442 ± 27 µg/min (T), 357 ± 28 µg/min (R) and 2.95 ± 0.13 µm (MMAD). The total aerosol delivery rate from the new AE-SVN matched that of the LC-JET™ (un-paired t-test, p = 0.79) and exceeded that from the SideStream™ (p < 0.001). The finer MMAD of the aerosol provided by the new AE-SVN resulted in a significantly greater respirable mass fraction, increasing the respirable mass delivery rate compared with the other SVN's (p < 0.001).

4. **COMPARISON OF A NEW BREATH ACTUATED SMALL VOLUME NEBULIZER (BA-SVN) WITH AN SVN SUPPLIED WITH COMPRESSOR INTENDED FOR HOME CARE USE.** AM Verdun, JP Mitchell and MW NAGEL. Presented at Ann. Meet. Amer. Assoc. of Asthma, Allergy and Immunology (AAAAI), Washington D.C., 1998.

The performance of a prototype novel BA-SVN (Trudell Medical International (n = 5 devices)) with normal saline (0.9% w/v NaCl) operating at 20 ± 2°C, 50 ± 10% RH, has been evaluated with an air compressor widely used in home care (Proneb™, Pari Respiratory Equipment Inc.). A laser diffractometer (Malvern Mastersizer-X) was used to determine the size distribution of droplets emitted at the mouthpiece (5 replicates per device). The total mass output was determined gravimetrically in a parallel series of tests. The BA-SVN was operated with manual over-ride engaged (continuous delivery of aerosol). Total (T) and respirable ((R), droplets finer than 4.8 µm aerodynamic diameter) mass output rates, respirable mass fraction (RM) and droplet mass median aerodynamic diameter (MMAD) were 167 ± 6 µg/min (T), 96 ± 5 µg/min (R), 57.5 ± 2.1% (RM) and 4.40 ± 0.11 µm (MMAD). In comparison, under similar conditions, a Pari LC-JET™ SVN with Proneb™ (n = 5 replicate measurements) provided 211 ± 3 µg/min (T), 65 ± 4 µg/min (R), 30.9 ± 1.5% (RM) and 6.94 ± 0.20 µm (MMAD). The new BA-SVN provided aerosol having a finer MMAD and greater RM (un-paired t-test, p < 0.001 for each variable) which resulted in an improved respirable mass output rate compared with the LC-JET™ SVN. The BA-SVN also has the advantage that no aerosol is produced to waste during the exhalation portion of each breathing cycle.

Xopenex™ (levalbuterol), Sepracor™ Inc.

1. **SAFETY AND EFFICACY OF FIVE-MINUTE TIMED AEROSOL ADMINISTRATION WITH THE AEROECLIPSE* BREATH ACTUATED NEBULIZER: COMPARISON OF LEVALBUTEROL WITH RACEMIC ALBUTEROL.** RS Pikarsky, R Acevedo, C Roman, W Fascia, T Farrell. Respiratory Care, September 2002; 47(9): 1075.

Purpose: Beta₂-agonist Racemic Albuterol has been used extensively in the performance of pre & post bronchodilator studies in the pulmonary function laboratory. This study evaluated the safety and efficacy of timed nebulization of the two dosages of Levalbuterol (Sepracor Inc., Marlborough, MA) as compared to Racemic Albuterol (Dey, Napa, CA) with the use of the **AeroEclipse*** Breath Actuated Nebulizer (BAN) (Monaghan Medical Corp. Plattsburgh, N.Y.). **Methods:** A consecutive, non-randomized, mostly COPD population (n = 93) receiving pre & post bronchodilator testing in our Pulmonary Function Lab were studied. Two different Levalbuterol medication dosages were administered: 0.63mg Levalbuterol UD or 1.25mg UD Levalbuterol. The Racemic Albuterol dosage was 2.5mg UD. All 5 minute timed aerosol treatments were administered using the BAN with an oxygen flow rate of 8L/min. The Sensormedics Vmax 22 Pulmonary Function System was utilized to measure both FEV1 and PEFR. A standardized subjective questionnaire to determine

side effects was completed. **Results:** The table shows the Levalbuterol and Racemic Albuterol dosages, mean % change of FEV1 and PEFR from pre-treatment to 10-minute post treatment, administration time, tremulousness and increase in heart rate. There was no significant difference in % change in FEV1 or PEFR. There was a significant increase in heart rate with the 1.25mg Levalbuterol UD group (7.2 vs. 3.4, $p<.05^*$; 7.2 vs. 2.2, $p<.01^{**}$). There was no difference in respiratory rate, tremulousness, or nausea.

| Nebulizer (n) | Dose | % Change FEV1 | % Change PEFR | Time (min) | Trem. | HR (Inc.) |
|------------------------|-----------|---------------|---------------|------------|-------|-----------|
| Levalbuterol (38) | 0.63mg UD | 7.8 | 6.2 | 5 | 4 | 3.4* |
| Levalbuterol (29) | 1.25mg UD | 7.7 | 16.6 | 5 | 2 | 7.2 |
| Racemic Albuterol (26) | 2.25mg UD | 12.2 | 10.5 | 5 | 0 | 2.2** |

Conclusion: Five minute timed administration of Levalbuterol and Racemic Albuterol using the BAN was equally efficacious and had similar safety profiles. The change in FEV1 and PEFR are consistent with our mostly COPD population. The increase in heart rate was greatest with the Levalbuterol 1.25 mg dosage. **Clinical Implications:** Five minute timed administration of Levalbuterol and Racemic Albuterol using the BAN is a safe and efficient alternative to the use of small volume nebulizers. Additional caution should be taken when administering Levalbuterol at the 1.25 mg dosage utilizing the BAN in cardiac patients. The efficiency of timed aerosol administration could have significant impact on resource utilization while maintaining the quality of aerosol delivery. This may be one of several strategies to address the problems of Respiratory Care staff shortages or high seasonal effect in the acute care facility.

2. COMPARISON IN RATES OF BREAKTHROUGH TREATMENTS DURING A CONVERSION FROM RACEMIC ALBUTEROL TO LEVALBUTEROL. RS Pikarsky, RA Acevedo, C Roman and T Farrell. CHEST 2002; 122(4):146S.

Purpose: in order to meet our patient care demands, Crouse Hospital approved an automatic conversion from Racemic Albuterol to Levalbuterol. This study compares the breakthrough rates of Racemic Albuterol and Levalbuterol, with and without Ipratropium. **Methods:** Racemic Albuterol (Alb) 2.5 mg Q4h was converted to either Levalbuterol (Lev) 0.63 mg Q6h or Levalbuterol 1.25 mg Q8h. If ordered, Ipratropium (Ipra) 0.5 mg was administered at the same frequency as the Levalbuterol. Patients with acute coronary syndromes, need for cardiac monitoring, or requiring more frequent aerosol administration received the lower Levalbuterol dose Q6h. A majority of aerosol therapy was provided with the use of the **AeroEclipse*** Breath Actuated Nebulizer (BAN). All aerosol treatments, including breakthrough treatments, delivered between July 1, 2001 and February 28, 2002 were recorded. **Results:** Tx/Pt/day represents the number of treatments delivered per patient per day. Rate/100 Pt/days = (Breakthrough) / (Total Tx / Tx/Pt/day) x 100. Rate/100 Pt/days corrects for the differences in daily administration frequency, and may better reflect the daily impact of the breakthrough rate. The breakthrough rate of the combined Albuterol group was significantly greater than both Levalbuterol groups (5.29 vs. 2.29, 5.29 vs. 2.43, $p<.001^*$). The breakthrough rate with Albuterol was significantly reduced with the addition of Ipratropium ($p<.001^{**}$). Ipratropium did not significantly change the breakthrough rate when added to Levalbuterol groups.

| Medication | Total Tx | Breakthrough | Rate/1000 | Tx/Pt/day | Rate/100 Pt/day | |
|----------------------|----------|--------------|-----------|-----------|-----------------|-------|
| Alb Q4h | 3832 | 47 | 12.27 | 6 | 7.36** | 5.29* |
| Alb/Ipra Q4h | 3767 | 20 | 5.31 | 6 | 3.19** | |
| Lev 0.63mg Q6h | 3592 | 24 | 6.68 | 4 | 2.67 | 2.29* |
| Lev 0.63 mg/Ipra Q6h | 1821 | 7 | 3.84 | 4 | 1.54 | |
| Lev 1.25mg Q8h | 1791 | 17 | 9.49 | 3 | 2.85 | 2.43* |
| Lev 1.25mg/Ipra Q8h | 678 | 3 | 4.42 | 3 | 1.33 | |

Conclusions: The conversion from Racemic Albuterol to Levalbuterol allowed for a decreased frequency of daily medication administrations and a significant decrease in breakthrough requirements. Ipratropium showed a significant benefit in breakthrough reduction for the Racemic Albuterol group. **Clinical Implications:** The efficiencies gained by decreasing the daily frequency of aerosol administration can have a significant impact on resource utilization. The conversion to Levalbuterol allows for decreased respiratory therapy time or the re-allocating of workforce needs while maintaining, or improving, quality of aerosol administration, as evidenced by the decrease in breakthrough requirements.

3. IMPROVING RESOURCE UTILIZATION WITH NEW TECHNOLOGIES. MA Lewis, SS Harris, SL Campbell, AL Hodges, DM Clark. Respiratory Care, August 2000; 45(8).

Background: To meet patient care needs during the peak respiratory season using levalbuterol (LEV) (Sepracor Inc., Marlboro, MA) and **AeroEclipse*** Breath Actuated Nebulizer ("BAN") (Monaghan Medical Corp., Plattsburgh, NY). Both pilot projects were approved by the Respiratory Care Advisory Committee. **Methods:** LEV 1.25mg delivered via nebulization q6h was substituted for albuterol 2.5mg ordered q4h in October 1999. Patients could also receive LEV as needed. A standardized subjective questionnaire to determine side effects of LEV was completed. BANs were utilized on patients meeting specified criteria during November 1999. Standard nebulizers were used for all other patients who required nebulized treatments. Treatment times were extracted from the CliniVision Information Management System database. **Results:** LEV was substituted for albuterol in 25 patients. Indications for nebulizer therapy included asthma (8%), COPD (32%), community acquired pneumonia (20%), and other (40%). The average number of LEV treatments per day was 3.7. This compared favorably to albuterol, which historically required = 6 treatments per day. No patients requested breakthrough treatments or noted side effects due to LEV. A total of 298 treatments were delivered using BANs

versus 322 delivered using a standard nebulizer. The average time per treatment using BANs was 9.9 minutes versus 14.76 minutes with the standard nebulizer. The results of these pilot programs prompted changes in respiratory therapy practice throughout the hospital. LEV was added to the Patient Driven Protocols and BANs are now used for nebulizer treatments in patients meeting criteria. Hospital census data indicate a 13.5% increase for 2000 versus 1999. Thus, total treatments for January and February 1999 and 2000 were 30,089 and 32,923, respectively. During this period, 16,000 LEV vials were dispensed from an automated dispensing unit vs 8,900 vials of albuterol. Concurrently, overtime (OT) hours utilized in 2000 were decreased by 693 hours, resulting in a savings of \$16,632, despite the increased number of treatments. Therefore, treatments were delivered to more patients with less OT utilized in 2000. **Conclusions:** These data illustrate the cost-effectiveness of two technologies utilized in our hospital, while patient care and satisfaction were maintained. OT hours decreased by 25% while treatments were delivered to more patients throughout the hospital. The use of LEV has resulted in a 33% decrease in the number of treatments per day with few "prn" treatments, while BAN has decreased the time to deliver therapy by 33%.

Actiq® (fentanyl citrate), Abbott Laboratories

1. **RANDOMIZED CLINICAL TRIAL OF NEBULIZED FENTANYL CITRATE VERSUS IV FENTANYL CITRATE IN CHILDREN PRESENTING TO THE EMERGENCY DEPARTMENT WITH ACUTE PAIN.** JR Miner, C Kletti, M Herold, D Hubbard and MH Biros. Acad Emerg Med 2007 Oct; 14(10):895-8.

Objectives: To compare the pain relief achieved with nebulized fentanyl citrate with intravenous (IV) fentanyl citrate in children presenting to the emergency department (ED) with painful conditions to determine if nebulized fentanyl is a feasible alternative to IV fentanyl for the treatment of acute pain in children. **Methods:** This was a randomized controlled trial in an urban county medical center ED with an annual census of 99,000 visits. ED patients, aged 6 months to 17 years, presenting with acute pain who were going to be treated with IV pain medications, were eligible for enrollment. After the parents had provided informed consent, and children older than 6 years had provided assent, patients were randomized (1:2) to receive either fentanyl citrate IV (1.5 µg/kg) or fentanyl citrate by breath-actuated nebulizer (3.0 µg/kg). Patients aged 6 years and older completed a 100-mm visual analog scale (VAS) describing their pain, and patients younger than 6 years had their pain assessed by the treating physician using the Children's Hospital of Eastern Ontario Pain Scale. Additionally, treating physicians used a 100-mm VAS to describe their perception of the patients' pain. These pain measurements were taken before treatment and every 10 minutes thereafter for 30 minutes. Baseline blood pressure, heart rate, and oxygen saturation were also measured before treatment and every 10 minutes for 30 minutes. After 30 minutes, physicians were asked whether or not they believed the medication provided adequate pain relief for the patient. Parents were asked to rate their satisfaction with the treatment using a five-point scale. Patients who received additional pain medications by any method before the 30-minute measurement period was completed were considered treatment failures. Data were compared using descriptive statistics and 95% confidence intervals; the rates of adequate pain relief between the groups were compared using Fisher exact tests. **Results:** Forty-one patients were enrolled in the study; 14 were randomized to IV fentanyl (ten actually received it), and 27 patients were randomized to nebulized fentanyl (31 actually received it). In the four patients who were randomized to IV fentanyl but received nebulized fentanyl, the parents requested the nebulized medication after being told their child had been randomized to IV fentanyl. Baseline pain VAS scores were 82.8 mm (SD ±14.3, 69–100) in the IV group and 76.2 mm (SD ±20.5, 34–100) in the nebulized group. There were five treatment failures: one who received IV fentanyl and four who received nebulized fentanyl. The four patients who were considered treatment failures in the nebulized fentanyl group were all younger than 3 years and had difficulty triggering the breath-actuated nebulizer. The mean decrease in pain for patients remaining in the study was 55.1 mm (95% CI = 40.3 to 70.0) for the IV group and 77.8 mm (95% CI = 67.4 to 88.4) for the nebulized group. The pain treatment was described as adequate by the treating physician in eight of 14 patients in the IV group and 20 of 27 patients in the nebulized group (p = 0.42). No adverse events were detected. **Conclusions:** Nebulized fentanyl citrate 3 µg/kg through a breath-actuated nebulizer appears to be a feasible alternative to IV fentanyl citrate for a variety of painful conditions in patients older than 3 years.

AeroLEF™ (Liposome-Encapsulated Fentanyl), YM Biosciences, Inc.

1. **A RANDOMIZED CONTROLLED TRIAL DEMONSTRATES THE EFFICACY, SAFETY AND TOLERABILITY OF AEROSOLIZED FREE AND LIPOSOME-ENCAPSULATED FENTANYL (AeroLEF™) VIA PULMONARY ADMINISTRATION.** R Brull, V Chan. Presented at the American Pain Society's 27th Annual Scientific Meeting (APS), Tampa, FL., 2008.

Pain following orthopedic surgery can be severe, requiring rapid onset and prolonged analgesia. The ideal analgesic has rapid onset of action, sustained effect, self titratable dosing and minimal adverse effects (AEs). Inhalation of opioids is conceptually appealing as the alveolar surface permits rapid absorption. We report a prospective randomized, blinded, placebo-controlled study of AeroLEF™ administered via breathactuated nebulizer. Ninety-nine ASA PS I-II patients aged 18-81 years undergoing elective orthopedic surgery under GA were randomized to AeroLEF™ or placebo (2:1 stratification). Nebulizers contained 1500 µg AeroLEF™ (≤1000 µg available for nebulization) or placebo; during each treatment session, a second nebulizer was provided if requested. Treatment was initiated when patients reported ≥ moderate pain. Up to three treatment sessions were permitted over 8-12 hours. Rescue medication was IV morphine. The

primary efficacy endpoint, SPRID4, was better with AeroLEF™ (mean scores of 7.02 vs. 3.35, $P < 0.02$). There was no difference between groups in clinically-significant respiratory depression (< 8 breaths/min or $SpO_2 < 90\%$ for > 20 sec). No patient received opioid antagonists or ventilatory support. Nausea (11% vs. 3%) and vomiting (31% vs. 21%) were more common with AeroLEF™ than with placebo. Following the first dose of study drug, more patients given AeroLEF™ reported mild or no pain (59% vs. 27%; $P < 0.01$). Time to effective pain relief after the first dose of study drug was shorter with AeroLEF™ group ($P < 0.005$). More patients given AeroLEF™ reported moderate-to-complete pain relief (60% vs. 32%, $P < 0.02$). This study suggests that patient-controlled inhalational analgesia with free and liposome encapsulated fentanyl can provide safe and effective pain relief following orthopedic surgery. Industry support provided by YM Biosciences Inc.

2. **AEROSOLIZED LIPOSOME-ENCAPSULATED FENTANYL (AeroLEF™) VIA PULMONARY ADMINISTRATION ALLOWS PATIENTS WITH MODERATE TO SEVERE POST-SURGICAL ACUTE PAIN TO SELF-TITRATE TO EFFECTIVE ANALGESIA.** A Clark, M Rossiter-Rooney, F Valle-Leuttri. Presented at the American Pain Society's 27th Annual Scientific Meeting (APS), Tampa, FL, 2008.

Acute pain is characterized by rapid onset, unpredictable and variable intensity confounded by highly variable patient responses to analgesics. Consequently, a successful dose is difficult to predict and maintain. AeroLEF™, a proprietary combination of free and liposome-encapsulated fentanyl for inhalation provides micro-doses of fentanyl per breath designed to allow real-time patient-controlled dose selection. In this study, nineteen post-surgical patients with moderate to severe pain following ACL surgery, were instructed to self-administer AeroLEF™ via breath actuated nebulizer until they had achieved analgesia, experienced dose-limiting side effects, or completed the maximum available dose (1000µg emitted per nebulizer, ≤ 2 nebulizers allowed). Eighteen (95%) of the patients achieved analgesia following self-administration of AeroLEF™. The median time to first perceptible analgesia was 2.7min. Mean plasma fentanyl concentration at first perceptible analgesia was 0.801ng.mL⁻¹. Median time to effective analgesia was 17min. At analgesia, the mean plasma fentanyl level was 1.30ng.mL⁻¹ but varied widely among patients, covering a 6.5-fold concentration range (0.39 to 2.5 ng.mL⁻¹). The mean duration of analgesia was 3.7h and the request for additional analgesics was associated with a decrease in mean plasma fentanyl levels to 0.887ng.mL⁻¹ (ranging from 0.36ngmL⁻¹ to 1.584ngmL⁻¹), comparable to the concentrations at first perceptible analgesia and consistent with reported ranges for minimal effective plasma fentanyl in post-surgical patients (0.34 to 1.58ng.mL⁻¹). A 9-fold dosing range was selected by patients in order to obtain analgesia with AeroLEF™, emphasizing the inter-patient variability associated with opioid use. AeroLEF™, at doses sufficient to establish analgesia, was well tolerated with no serious adverse events were reported. Adverse events were generally mild and commonly associated with opioid use in the post-operative period. These data suggest that self-titration to analgesia with AeroLEF™ offers a novel and effective approach to address the variability inherent in pain. Industry support provided by YM BioSciences Inc.

3. **COMPARATIVE PHASE I PK STUDY OF AEROSOLIZED FREE AND LIPOSOME-ENCAPSULATED FENTANYL (AeroLEF™) DEMONSTRATES RAPID AND EXTENDED PLASMA FENTANYL CONCENTRATIONS FOLLOWING INHALATION.** O Hung, D Pliura. Presented at the American Pain Society's 27th Annual Scientific Meeting (APS), Tampa, FL, 2008.

AeroLEF is a proprietary combination of free and liposome-encapsulated fentanyl for inhalation via breath-actuated nebulizers. We report the pharmacokinetics, safety, and tolerability of 1500µg AeroLEF vs. 200µg bolus IV fentanyl; values are mean (\pm SD). Healthy, opiate-naïve volunteers inhaled microdoses of AeroLEF ($\leq 5\mu\text{g}/\text{breath}$; total emitted fentanyl dose $\leq 1000\mu\text{g}$) over 7-15 min. Within 4 min of initiating AeroLEF inhalation, subjects attained plasma fentanyl concentrations (Cp) of 0.734 ng.mL⁻¹. Maximum Cp was similar with AeroLEF and IV fentanyl (2.53 vs. 2.80 ng.mL⁻¹). Cmax (mean of 15 min) occurred shortly after completion of AeroLEF™ inhalation (mean of 12 min), indicating rapid absorption from the lung. Cp values in the effective range persisted for several hours with AeroLEF (at 4 hr, Cp was 0.525 ± 0.180 ng.mL⁻¹) but not with IV administration (at 1 hr, Cp was 0.559 ± 0.209 ng.mL⁻¹). Similar inter-subject variability in exposure was observed in both treatment arms: coefficient in variation of AUC was 24% with IV administration vs. 29% with AeroLEF. Subjects were monitored continuously for adverse respiratory events. No severe adverse events were observed. Mild hypoxia was observed in both treatment groups. Mild bradycardia was observed in one subject receiving IV fentanyl. Spirometry measurements (FVC, FEV1 and FEF25%-75%) before and after AeroLEF indicated no significant changes in lung function. In summary, AeroLEF achieves rapid and persistent fentanyl concentrations in the therapeutic range and appears to be well tolerated. Industry support provided by YM BioSciences Inc.

Study Summary by Nebulizer

Ventstream™, Medic-Aid™

1. **PERFORMANCE OF A NEW BREATH-ACTUATED SMALL VOLUME NEBULIZER UNDER CONDITIONS THAT SIMULATE USE BY AN ADULT PATIENT.** R Blacker, JP Mitchell, MW Nagel and AMW Verdun. Eur. Resp. J., 1997; 10(25): 235.

The development of pneumatic small volume nebulizers (SVNs) in which atomization is enabled during the inhalation portion of a patient's breathing cycle has important ramifications in terms of the efficiency at which medication can be delivered. We report an investigation in which the effectiveness for the delivery of salbutamol (Ventolin® nebulizer: 5 mg/2.5 ml, GlaxoSmithKline, Canada) via a prototype breath-actuated SVN (Trudell Medical, Canada (TRU)) was compared with that of a high performance closed-system SVN (Ventstream™, Medic-Aid, Pagham, U.K. (VEN)). Each device was connected in turn to a ventilator-test lung apparatus in such a way that aerosol delivered on inhalation (800 ml tidal volume, I/E of 1/1, 15 breaths/min) was collected on a filter (Filtrete™, 3M Corp., St Paul, MN) located at the mouthpiece. Oxygen (440 kPa, 8 l/min) was supplied to operate each SVN, and the contents of a single nebulizer (2.5 ml) were added to the reservoir at the start of each test. Over a 5 minute period of use, the TRU SVN provided 1.74 ± 0.04 mg salbutamol to the filter (n=5 replicates). In comparison, the VEN delivered 1.28 ± 0.01 mg in 3.5 min after which the device sputtered dry (n = 5 replicates). These data indicate that the new breath-actuated device has important benefits in reducing wastage of medication by operating more efficiently, as well as an optimal impact on the environment.

2. **A NOVEL BREATH-ACTUATED SMALL VOLUME NEBULIZER UNDER SIMULATED ADULT USE CONDITIONS.** R Blacker, JP Mitchell, MW Nagel, AMW Verdun. Presented at Open Forum, Ann. Meet. of the American Association For Respiratory Care (AARC), New Orleans, 1997.

Pneumatic small volume nebulizers (SVNs) in which atomization only occurs during the inhalation phase of the breathing cycle have important ramifications in terms of the efficiency of medication delivery. We report an investigation in which the effectiveness for the delivery of salbutamol (Ventolin® nebulizer: 5 mg/2.5 ml, GlaxoSmithKline, Canada) via a prototype breath-actuated SVN (Trudell Medical, Canada (TRU)) was compared with that of a high performance closed-system SVN (Ventstream™, Medic-Aid, Pagham, U.K. (VEN)). Each nebulizer was connected in turn to a dual-chambered test lung with one chamber driven by a ventilator and the other connected to the SVN mouthpiece. Aerosolized salbutamol delivered on inhalation (800 ml tidal volume, I/E of 1/1, 15 breaths/min) was collected on a filter (Filtrete™, 3M Corp., St Paul, MN) located at the mouthpiece. Oxygen (440 kPa, 8 l/min) was used to operate each SVN, and the contents of a single nebulizer (2.5 ml) were added to the reservoir at the start of each test. Over a 5 minute period of use, the TRU SVN provided 1.74 ± 0.04 mg salbutamol to the filter (n=5 replicates), significantly more than the VEN which delivered 1.28 ± 0.01 mg in 3.5 min (Mann Whitney Rank Sum Test, p = 0.008), after which the device sputtered dry (n = 5 replicates). These data indicate that the new breath-actuated device may have important benefits in reducing wastage of medication by operating more efficiently, as well as reducing exposure to the care-giver.

UpDraft Neb-U-Mist™, Hudson Oxygen Therapy Sales Co.

1. **AEROSOLIZED AMPHOTERICIN B LIPID COMPLEX (aABLC) DISTRIBUTION IN LUNG TRANSPLANT RECIPIENTS: A COMPARISON OF CONTINUOUS VERSUS BREATH ACTUATED NEBULIZERS.** ES Dodds, NA Petry, JD Davies, DW Zaas, SM Palmer, SW Shipes, RH Drew, BD Alexander, RE Coleman, JR Perfect. Presented at the American Association for Respiratory Care Congress, Orlando, FL, December 2007.

Background: Aerosolized amphotericin B has become an attractive option for antifungal prophylaxis following solid organ and stem cell transplantation.^{1,2} This therapeutic strategy facilitates localized delivery of antifungal agent, thereby minimizing toxicities and drug-drug interactions associated with currently-available systemic antifungal agents. Determining drug delivery characteristics, including dose and nebulizer system, for aerosol drug administration is important to ensure optimal drug delivery. Newer, breath-actuated nebulizers (BAN's) are available and, in theory, provide the ability to limit environmental exposure and also deliver a higher percentage of the prepared dose to the patient. **Objective:** To characterize the distribution of aerosolized ABLC immediately following nebulization in bilateral lung transplant recipients via 2 different nebulizer systems – continuous nebulizer (CN): Up-Draft, Model 1724 (Hudson RCI, Temecula, CA) and breath actuated nebulizer (BAN): **AeroEclipse® II** (Monaghan Medical, Plattsburgh, NY). ABLC 20 mg/4mL was mixed with prepared 99mTc-ABLC (Abelcet®-Enzon Pharmaceuticals) prior to loading into the radioaerosol delivery system. **Methods:** Nebulizer assignment was performed sequentially with the first 5 subjects receiving treatment via the continuous flow nebulizer and the subsequent 5 subjects receiving study drug treatment via the BAN. Immediately following inhalation, drug product distribution images were obtained with patients in the supine position. Subjects were then placed on the table of a dual-head gamma camera system (General Electric Healthcare, Milwaukee, WI). Total delivered dose (TDD) was calculated by determining the difference in the known starting counts for the medication vial and counts of the nebulizer apparatus, including filter, subject waste materials and empty medication vials, obtained after study medication administration. Gastric activity of 99mTc-ABLC was also measured. Drug exposure was reported as: TDD: total delivered dose; Drug delivery to each of the following lung regions was reported as a percentage of TDD: right lung (RL), left lung (LL) and GI tract; the two nebulizer groups were compared for differences in mean TDD and regional distribution using student's t-test.

Results: Total drug delivery (as percent of prepared dose) was significantly higher for the BAN (20.7% versus 3.5%, $p=0.01$). Mean regional distribution (as percent of total delivered dose) did not differ between the two nebulizer devices for the left lung, right lung, or GI tract.

| Subject | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 |
|---------------------------|-------------------------|-----|-----|-----|-----|---------------------------|------|------|------|------|
| | Continuous Nebulizer | | | | | Breath Actuated Nebulizer | | | | |
| Drug Delivery* | % of total dose in vial | | | | | % of total dose in vial | | | | |
| RL | NR | 1.6 | 1.2 | 0.4 | 1.2 | 7.4 | 9.6 | 5.2 | 5.8 | 11.3 |
| LL | NR | 1.4 | 0.9 | 0.3 | 0.7 | 6.4 | 5.5 | 5.4 | 6.0 | 8.9 |
| GI | NR | 3.6 | 1.3 | 0.6 | 0.5 | 5.1 | 5.1 | 7.2 | 11.2 | 3.5 |
| Total Drug Delivery (TDD) | NR | 6.6 | 3.4 | 1.3 | 2.4 | 18.9 | 20.2 | 17.8 | 23.0 | 23.7 |
| Regional Delivery** | | | | | | | | | | |
| Right | 50 | 24 | 35 | 31 | 49 | 39 | 47 | 29 | 25 | 48 |
| Left | 17 | 21 | 27 | 23 | 29 | 34 | 27 | 30 | 26 | 37 |
| Esophagus and Stomach | 32 | 55 | 39 | 46 | 22 | 27 | 25 | 40 | 49 | 15 |

* As percent of prepared dose

** As a percentage of the total delivered dose

Conclusion: Use of the BAN resulted in a larger portion of the drug being deposited into the lungs. Since GI distribution was similar between the nebulizers, it appeared that more drug was vented to the surrounding atmosphere with the continuous system.

References: ¹ Hussain S, Zaldonis D, Kusne S et al. Variation in antifungal prophylaxis strategies in lung transplantation. *Transpl Infect Dis* 2006;213-8. ² Drummer JS. A survey of fungal management in lung transplantation. *Journal of Heart and Lung Transplantation* 2004;23:1376-81.

2. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS - THE RELATIONSHIP BETWEEN NEBULIZED DROPLET SIZE AND THE PARTICLE SIZE OF THE SUSPENSION.** JP Mitchell, MW Nagel, and AD Archer. Presented at Drug Delivery to the Lungs-IX, London, UK, 1998, in *J. Aerosol Med.* 12(3), 208, (1999).

A new air entrainment small volume nebulizer (AE-SVN) has been compared with two other SVNs (Neb-U-Mist™ and Misty-Neb™) for the delivery of a suspension of 0.25 µg/ml budesonide. Each SVN was operated at 8 l/min with compressed oxygen (50 psig). The total mass of budesonide was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The time-averaged delivery rate over the period of nebulization ((mean ± 1 S.D.) µg budesonide/min) from the AE-SVN (102 ± 9) was greater than with the Misty-Neb™ (49 ± 2), or Neb-U-Mist™ (25 ± 6). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with the Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). The mass median diameter (MMD) of the droplets from the AE-SVN measured using a laser diffractometer (2.9 ± 0.1 µm), was significantly finer compared with those from the Misty-Neb™ (4.5 ± 0.9 µm) and Neb-U-Mist™ (5.6 ± 0.6 µm) and closest to the size of the micronized budesonide particles in the original suspension. The efficient delivery of medication formulated as micronized powder in aqueous suspension necessitates that the droplets produced upon nebulization are large enough so that single particles are efficiently entrained during atomization, but not so coarse that they cannot leave the nebulizer, extending nebulization time.

3. **EVALUATION OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) WHEN USED WITH OXYGEN AS A DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** JP Mitchell and MW Nagel. Presented at ALA/ATS International Conference, Chicago, 1998.

Purpose: To compare the delivery of saline (0.9% w/v NaCl) by a new AE-SVN (Trudell Medical Int.) with that from two other representative SVNs (UpDraft Neb-U-Mist® (Hudson Oxygen Therapy Sales Co.) and Airlife™ Misty-Neb™ (Baxter Healthcare Corp.)) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 AE-SVNs were tested using a laser diffractometer (Malvern Mastersizer-X) to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. 5, Neb-U-Mist® and a similar number of MistyNeb™ SVNs were also evaluated. **Results:** Total (TM) and respirable ((RM), droplets finer than 4.8 µm diameter) mass output rates and droplet mass median diameter (MMD) were as follows: AE-SVN: TM = 671 ± 26 µg/min, RM = 542 ± 23 µg/min (80.8 ± 1.3% respirable), MMD = 2.88 ± 0.09 µm; Neb-U-Mist™: TM = 266 ± 13 µg/min, RM = 119 ± 16 µg/min (42.1 ± 5.2% respirable), MMD = 5.6 ± 0.6 µm; Misty-Neb™: TM = 336 ± 60 µg/min, RM = 178 ± 43 µg/min (53.1 ± 8.5 % respirable), MMD = 4.5 ± 0.9 µm. **Conclusion:** TM from the new AE-SVN was substantially greater than those from either the Neb-U-Mist® or Misty-Neb™ (1-way ANOVA, $p < 0.001$). The finer MMD produced from the AE-SVN resulted in a significantly greater RM compared with either of the other SVNs ($p < 0.001$).


4. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS: A COMPARATIVE *IN-VITRO* ASSESSMENT.** JP Mitchell, MW Nagel, and AD Archer. Presented at Proc. Ann. Meet. Amer. College of Chest Physicians (ACCP), Toronto, November 1998, in Chest, 114 (4S), 295, (1998); and at the World Asthma Meeting, Barcelona, December 1998, in Eur. Resp. J., 12S29, 7, (1998).

Purpose: To compare the performance of a new air entrainment small volume nebulizer (AE-SVN, Trudell Medical International with other widely used SVNs (LC-Star™ (PARI Respiratory Equipment), Updraft™ Neb-U-Mist™ (Hudson Oxygen Therapy Sales Co.), Circulaire™ (Westmed), Sidestream™ (Medic-Aid), Airlife™ Misty-Neb™ (Baxter Healthcare Corp.)) for the delivery of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). **Methods:** Each SVN (n = 5 devices for each group, 3 replicates per device) was operated with compressed air at 50 psig at a flow rate of 8 l/min. The total mass of budesonide nebulized from 2 x 2 ml ampoules was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The SVN was operated until it spluttered, was then tapped gently to dislodge droplets back to the reservoir. Nebulization was deemed complete 20 seconds later. The mass of budesonide collected was determined by HPLC-UV spectrophotometry. **Results:** The delivery rate ((mean ± 1 S.D) µg budesonide/min) from the AE-SVN (102 ± 9) was significantly greater than with the other groups: (LC-Star™ (91 ± 6), Misty-Neb™ (49 ± 2), Sidestream™ (46 ± 4), Circulaire™ (26 ± 4) and Neb-U-Mist™ (25 ± 6)), (1-way ANOVA, p < 0.02). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with LC-Star™ (229 ± 10 s), Sidestream™ (365 ± 19 s), Circulaire™ (420 ± 84 s), Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). **Conclusions:** The new AE-SVN is highly efficient at entraining the budesonide particles into the liquid droplets at these conditions. **Clinical Implications:** The good delivery rate combined with comparatively short duration of delivery offers the potential for rapid treatment and patient convenience.

Updraft-II Nebu-U-Mist™, Hudson Oxygen Therapy Sales Co.

1. **EFFECT OF NEBULIZER DESIGN ON FINE PARTICLE MASS.** D Hess, J.P Mitchell, D Coppolo, MW Nagel, A.D Archer, R Blacker. Presented at Open Forum, Ann. Meet. of the American Association For Respiratory Care (AARC), Las Vegas, 1999. Published in Resp Care, Oct 1999, 44(10): 1289.

Background: Nebulizer design is known to affect performance. In this study, we compared fine particle mass from nebulizers of four designs. **Methods:** We tested traditional disposable nebulizers (Baxter Misty-Neb™, Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™), nebulizers with collection bags (Westmed Circulaire™), and a Trudell **AeroEclipse*** (with breath actuation disabled). Five of each device with three replicates (n = 15) were tested using an *in-vitro* model of spontaneous breathing. A rigid bar was placed between the two compartments of a test lung (Michigan Instruments TTL). The drive lung was attached to a ventilator (Infrasonics Infant Star™) to simulate spontaneous breathing (tidal volume 0.6 L, rate 10/min, T_i 2 s). A bacterial/viral filter (Trudell MT3000) was placed between the nebulizer and slave lung. Flow was monitored between the test lung and filter (Novamatrix Ventcheck™). Albuterol solution (0.625 mg/mL) was placed into the nebulizers (4 mL), which were powered with air (8 L/min). Filters were replaced at one minute intervals (flow to the nebulizer was discontinued during filter replacement) until sputtering occurred. The filter was washed with methanol (20 mL) and albuterol concentration was measured with HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer. Fine particle mass was calculated as the product of mass % <4.7 µm and total nebulizer output. **Results:** Fine particle mass from the **AeroEclipse*** nebulizer was greater than that from the other nebulizers (P<0.001) (see figure). **Conclusions:** Fine particle mass was affected by nebulizer design. The clinical relevance of this finding awaits further investigation. Further evaluation of the breath-actuated feature of the **AeroEclipse*** is warranted.



2. **THE EFFECT OF SMALL VOLUME NEBULIZER (SVN) DESIGN ON FINE PARTICLE MASS DELIVERY OF A BRONCHODILATOR.** R Blacker, RW Morton, JP Mitchell, MW Nagel and DR Hess. Presented at Drug Delivery to the Lungs-X, London, UK, 1998, J. Aerosol Med., 13(1), 65.

Fine particle mass delivery was compared from six different SVNs, including continuous un-enhanced flow designs (Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™, Medic-Aid Sidestream®), nebulizers with aerosol collection bag (Westmed Circulaire™), and an **AeroEclipse*** with breath actuation disabled (Trudell Medical

International). Five of each type of SVN were tested operating with air (8 l/min , 50 psig), using an *in-vitro* model that simulated spontaneous breathing by an adult (tidal volume 0.6 l, rate 10/min, TI = 2 s). A bacterial/viral filter was placed between the nebulizer and breathing simulator. In each case, salbutamol sulphate (Ventolin®) respirator solution (0.625 mg/ml, 4 ml) was placed into the reservoir of the SVN. The filters were replaced at one-minute intervals until sputtering occurred. The salbutamol collected on the filter was assayed by HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer laser diffractometer. Fine particle mass delivery rates varied significantly from each of the SVN's from more than 110 µg/min (**AeroEclipse***) to ca. 20 µg/min (Circulaire™).

Micromist™, Hudson RCI™

1. **DELIVERY OF ALBUTEROL VIA A NEW BREATH ACTUATED NEBULIZER: COMPARISON WITH CONTINUOUS JET NEBULIZERS** DP Coppola, MW Nagel, CC Doyle, VA Avvakoumova and JP Mitchell. Presented at the American Thoracic Society International Conference 2007, San Francisco, California, USA.

A new breath actuated nebulizer (**AeroEclipse*** II BAN, Monaghan Medical Corp., Plattsburgh, NY) has been developed to deliver medication only when the patient inhales. This study sought to determine the delivery of albuterol (3-ml fill of diluted solution (0.83 mg/ml)) as fine droplets < 4.7 µm aerodynamic diameter, and compare this fine droplet mass (FDM) with equivalent data from 4 widely available continuous jet nebulizers as benchmark devices. Each nebulizer (n=5; 3 replicates/device) was operated with compressed air at 50 psig at ca. 8 L/min to simulate hospital wall outlet conditions. The nebulizer on test was coupled to a breathing simulator set to mimic adult use (tidal volume = 600 ml, rate = 10 breaths/min; duty cycle = 0.33), and the emitted droplets were collected on an electret filter at the mouthpiece. The total mass of albuterol (TM) was assayed subsequently by HPLC-UV spectrophotometry. In a separate study, the droplet size distribution was determined by laser diffractometry so that the fine droplet fraction (FDF) could be obtained. FDM was determined as the product of TM and FDF. FDM (mean SD) from the BAN was 791 84 g, delivered in 8 minutes. Corresponding values (FDM in time from start to sputter) for the VixOne™ (Westmed, Tucson, AZ), MicroMist™ (Hudson RCI, Temecula CA), Misty Max 10™ (Cardinal Health, McGaw Park (IL) and model 8900™ (Salter Labs, Arvin, CA) were 267 11 g in 6 min, 133 8 g in 4 min, 249 10 g in 6 min and 161 10 g in 5 min. Aside from dosage assurance imparted by breath-actuation, the **AeroEclipse*** II BAN delivered substantially more FDM/min than the other devices. The clinician is now able to treat either for extended high dose delivery (potentially eliminating the need for additional therapy), or titrate to a shorter interval based on response.

2. **COMPARISON OF BREATH-ACTUATED JET NEBULIZER (BAN) IN 'CONTINUOUS DELIVERY' MODE WITH OTHER CONTINUOUS DELIVERY NEBULIZERS.** JP Mitchell, KJ Wiersema, CC Doyle and MW Nagel. Respiratory Care 2003;48(11):S1077.

The **AeroEclipse*** BAN (Monaghan Medical Corp., Plattsburgh, N.Y.) has been equipped with an optional blue cap whose purpose is to retain the actuator piston in the position it would occupy during inhalation in breath-actuated mode, so that the nebulizer operates continuously. The present study compared the delivery of a bronchodilator from diluted albuterol sulfate respirator solutions (3-ml of 0.83 and 1-ml of 2.5 mg/ml albuterol in physiologically normal saline (0.9% w/v NaCl)), via this nebulizer, the Micromist® (Hudson RCI, Temecula, CA), Misty-Neb™ (Allegiance Healthcare Corp., McGaw Park, IL) and the LCD™ (PARI Respiratory Equipment, Monterey, CA). Each nebulizer was tested using a breathing simulator set to the following parameters representative of adult use: tidal volume = 600-ml, rate = 10 breaths/min, inspiratory/expiratory ratio 1:2. The total mass of albuterol (TM) delivered to the first sputter was determined by filter collection at the mouthpiece of the nebulizer operated with compressed air supplied at 50 psig at 8 L/min (n = 5 devices/group, 3 replicates/device). The fraction of the aerosol contained in droplets finer than 4.8 µm aerodynamic diameter (FPF) was determined by laser diffractometry in a parallel study, so that the fine droplet mass (FM) could be calculated as the product of TM and FPF. Values of FM (mean ± SD) and time to deliver medication (T_{med}) were as follows:

| Solution(mg/ml) | AeroEclipse* | | LCD™ | | Micromist® | | MistyNeb™ | |
|------------------------|---------------------|-------------|-------------|------------|-------------|---------------|-----------|-----------|
| | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 |
| FM (µg) | 360 ± 22 | 263 ± 26 | 149 ± 16 | 108 ± 4 | 209 ± 12 | 15.4 ± 5.9 | 82 ± 9 | 31 ± 5 |
| T _{med} (min) | 3 | <1 | 2 | <1 | 7 | <1 | 4 | <1 |

The **AeroEclipse*** nebulizer delivered significantly more FM in continuous delivery mode than the other nebulizers when operated in continuous mode with either solution strength (1-way repeated measures ANOVA, p <0.05). T_{med} from the **AeroEclipse*** nebulizer was comparable with the best performing continuous nebulizer (LCD™).

3. **BREATH-ACTUATED NEBULIZER DELIVERS BRONCHO-DILATOR MORE EFFICIENTLY THAN CONVENTIONAL JET NEBULIZER IN A SIMULATION OF AN ADULT TIDAL-BREATHING PATIENT.** MW Nagel and JP Mitchell. Am. J. Resp. Crit. Care Med., 2002;165(8):A189.

Rationale: To compare delivery of albuterol sulfate inhalation solution (2.5 mg/3 ml vial equivalent to 0.083% w/v albuterol, Zenith Goldline Pharmaceuticals, Miami, FL) by conventional and breath-actuated nebulizer (BAN), simulating adult use. **Methods:** Each SVN (n = 5/group, 3 replicates/nebulizer) was operated with 8 l/min air at 50 psig and simulating breathing at tidal volume, I:E ratio and rate of 600-ml, 1:2 and 10/min respectively. Total emitted dose (TED) was determined for 5-**AeroEclipse*** BANs (Monaghan Medical Corp., N.Y., 1.5 ml solution) and 5 Micromist® nebulizers (Hudson RCI, Temecula, CA, 3.0 ml solution) by filter collection, and droplet size distributions were measured in a parallel

study by laser diffractometer. Fine particle dose (FPD) was calculated as the product of TED and the percentage by mass of droplets finer than 4.8 μ m aerodynamic diameter. **Results:** After 3 minutes, the **AeroEclipse*** BAN delivered 282 ± 10 mg FPD (mean \pm SD) and the Micromist[®] delivered 209 ± 12 mg albuterol after 7 minutes. **Conclusion:** Dose delivery and patient compliance are assured by virtue of the breath actuation feature of the **AeroEclipse*** nebulizer and the reduced time to deliver a specific equivalent dose of medication compared with a conventional nebulizer will improve cost effectiveness of treatment.

VixOne™, Westmed™ Corp.

1. **DELIVERY OF ALBUTEROL VIA A NEW BREATH ACTUATED NEBULIZER: COMPARISON WITH CONTINUOUS JET NEBULIZERS** DP Coppola, MW Nagel, CC Doyle, VA Avvakoumova and JP Mitchell. Presented at American Thoracic Society International Conference 2007, San Francisco, California, USA.

A new breath actuated nebulizer (**AeroEclipse*** II BAN, Monaghan Medical Corp., Plattsburgh, NY) has been developed to deliver medication only when the patient inhales. This study sought to determine the delivery of albuterol (3-ml fill of diluted solution (0.83 mg/ml)) as fine droplets < 4.7 μ m aerodynamic diameter, and compare this fine droplet mass (FDM) with equivalent data from 4 widely available continuous jet nebulizers as benchmark devices. Each nebulizer (n=5; 3 replicates/device) was operated with compressed air at 50 psig at ca. 8 L/min to simulate hospital wall outlet conditions. The nebulizer on test was coupled to a breathing simulator set to mimic adult use (tidal volume = 600 ml, rate = 10 breaths/min; duty cycle = 0.33), and the emitted droplets were collected on an electret filter at the mouthpiece. The total mass of albuterol (TM) was assayed subsequently by HPLC-UV spectrophotometry. In a separate study, the droplet size distribution was determined by laser diffractometry so that the fine droplet fraction (FDF) could be obtained. FDM was determined as the product of TM and FDF. FDM (mean SD) from the BAN was 791 84 g, delivered in 8 minutes. Corresponding values (FDM in time from start to sputter) for the VixOne™ (Westmed, Tucson, AZ), MicroMist™ (Hudson RCI, Temecula CA), Misty Max 10™ (Cardinal Health, McGaw Park (IL) and model 8900™ (Salter Labs, Arvin, CA) were 267 11 g in 6 min, 133 8 g in 4 min, 249 10 g in 6 min and 161 10 g in 5 min. Aside from dosage assurance imparted by breath-actuation, the **AeroEclipse*** II BAN delivered substantially more FDM/min than the other devices. The clinician is now able to treat either for extended high dose delivery (potentially eliminating the need for additional therapy), or titrate to a shorter interval based on response.

2. **SIMILAR DELIVERY OF AMPHOTERICIN LIPID COMPLEX IS POSSIBLE AT ONE-HALF DOSE VIA A BREATH-ACTUATED NEBULIZER COMPARED WITH A CONTINUOUSLY OPERATING NEBULIZER.** JP Mitchell, NR MacIntyre, MW Nagel, DP Coppola. Presented at the American Thoracic Society 101st International Congress, May, 2005.

Delivery of aerosolized antibiotics via continuous nebulizers wastes these expensive medications during patient exhalation. Breath-actuated nebulizers (BAN) can minimize waste with significant cost savings in medication, since they only operate when the patient inhales. Furthermore, medication is not emitted into the environment during exhalation. We describe a study in which dose delivery from a BAN (**AeroEclipse***, Monaghan Medical Corp., Plattsburgh, NY) was compared with that from a continuously operating nebulizer (VixOne™, Westmed Corp., Engelwood, CO (VIX)) (n=3/group) for the delivery of amphotericin lipid complex ((AMP) Ablecet, Enzon Pharmaceuticals, Piscataway, NY, 5-mg/ml). Each device was operated with air at 50 psig at 7 L/min (BAN) or 8 L/min (VIX), with the mouthpiece connected to a breathing simulator (Compass, PARI, Germany) set to replicate adult use (500-ml tidal volume, 1:2 inspiratory/expiratory ratio, 20-breaths/min). 5-ml AMP was placed in the BAN and 10-ml in the VIX (5-ml initially, followed by a further 5-ml after 4-min). Each nebulizer was operated for 1-min past first sputter. The mass of AMP collected on a filter at the mouthpiece was determined by HPLC-UV spectrophotometry (3-replicates/nebulizer). Droplet size distributions were determined by laser diffractometer in a separate study. Total emitted mass from the BAN was 7274 123 g, delivered in 10-min, of which 5892 100 g was in fine droplets 4.8 μ m diameter. The VIX delivered a total mass of 5276 557 g in 10-14 min, of which 4326 457 g was contained in fine droplets. The BAN was therefore capable of delivering 36% more medication as fine droplets with only one-half of the dose inserted in the reservoir.

3. **A MECHANICALLY OPERATED BREATH-ACTUATED NEBULIZER ENABLES BOTH IMPROVED CONTROL OF DOSING AND DELIVERY EFFICIENCY.** JP Mitchell, MW Nagel. Presented at Drug Delivery to the Lungs Conference, December 2005.

A mechanically operated, breath-actuated nebulizer (BAN) offers the clinician the prospect of being able to control the rate and duration of medication delivery dosimetrically, providing greater precision when titrating patients to establish an appropriate treatment regimen. We describe an *in vitro* study obtained with two formulations that are representative of formulations available for nebulization (amphotericin-B and ipratropium bromide), in which a BAN (**AeroEclipse***) delivered slightly more medication as fine droplets < 4.8 μ m aerodynamic diameter with approximately one-half of the dose in the reservoir compared with a continuously operating nebulizer (VixOne™). These measurements were made simulating use by an adult (500-ml tidal volume, inspiratory/expiratory ratio 1:2, 20 breaths/minute). Significant cost savings are therefore possible with the BAN with expensive medications, such as antibiotics, if less volume fill is required per treatment.

4. **DELIVERY OF BUDESONIDE INHALATION SOLUTION (BIS) THROUGH AN INFANT UPPER AIRWAY MODEL.** DE Geller, KC Kesser, HM Janssens, HAWM Tiddens. Am J Respir Crit Care Med 2003;167(7):A508.

We investigated variables that may be important in the delivery of BIS to the lungs of infants, a challenging population for aerosol delivery. **METHODS:** The Sophia Anatomical Infant Nose Throat (SAINT) airway model mounted on a breath simulator mimicked the breathing pattern of a 9-mo old infant (RR=30, Vt=100 ml, I:E ratio=1:1.3). Nebulizers were charged with BIS 0.25 mg and run continuously until dry. Drug captured on a filter distal to the SAINT model was the lung dose. Compressor: PARI PRONEB TURBO. Nebulizer/mask systems studied: VIX1/aerosol mask (AM), PediNeb pacifier device (PN) or blow-by (BB); **AeroEclipse*** neb and mask (AE); PARI LC+ and PARI LC*/ PARI Baby mask (PB), Fish mask (FM), and AE masks. The AE neb/mask was also studied with an ill breathing pattern (RR=50, Vt=100, I:E=1:2). **RESULTS:** Lung dose ranged from 2.0 to 7.6% of the neb charge. Lung dose was AE (5.0%) > VIX1 (3.5%), LC+/FM (3.2%), LC*/PB (2.9%), and LC+/PB (2.8%). Also, VIX1/AM (3.5%)>VIX1/PN (2.5%)>VIX1/BB (2.0%). The lung dose of the LC+ and LC* more than doubled (6.8 and 6.3%) when used with the AE mask. Lung dose increased with the ill breath pattern in proportion to increased minute ventilation (7.6%). **CONCLUSION:** 1) The AE system provided higher lung dose than other nebulizers with standard masks. 2) Mask design and fit can substantially impact nebulizer performance. 3) PN performed better than BB, but not as good as a mask. If crying decreases lung dose by 75%, we speculate that the PN and BB (non-crying) may improve lung dose vs mask with a crying infant. 4) An increase in lung dose may occur in ill infants if minute ventilation is elevated. *This abstract is funded by: AstraZeneca*

Airlife™ Misty-Neb™, Baxter™ Healthcare Corp.

1. **PERFORMANCE COMPARISON OF NEBULIZER DESIGNS: CONSTANT-OUTPUT, BREATH-ENHANCED, AND DOSIMETRIC.** JL Rau, A Ari and RD Restrepo. Respiratory Care 2004; 49(2): 174-179.

Introduction: Design differences among pneumatically powered, small-volume nebulizers affect drug disposition (percentage of the dose delivered to the patient, lost to deposition in the equipment, and lost via exhalation to ambient air) and thus affect drug availability and efficacy. **Objective:** Evaluate *in vitro* the dose disposition with 5 nebulizer models, of 3 types (constant-output, breath-enhanced, and dosimetric), using simulated normal, adult breathing. **Methods:** We compared 5 nebulizer models: 2 constant-output (Misty-Neb and SideStream), 1 breath-enhanced (Pari LCD), and 2 dosimetric (Circulaire and **AeroEclipse***). Each nebulizer was filled with a 3-mL unit-dose of albuterol sulfate and powered by oxygen at 8 L/min. The nebulizers were connected to an induction throat, connected to a breathing simulator. We measured (1) inhaled drug (subdivided into mass deposited in the induction throat and mass deposited in the filter at the distal end of the induction throat), (2) exhaled drug (lost to ambient air), (3) drug lost to deposition in the apparatus, and (4) drug left in the unit-dose bottle. The duration of nebulization (until sputter) was measured with a stopwatch. All drug amounts were analyzed via spectrophotometry and expressed as a percentage of the total dose. **Results:** The mean \pm SD inhaled drug percentages were: Misty-Neb $17.2 \pm 0.4\%$, SideStream $15.8 \pm 2.8\%$, Pari LCD $15.2 \pm 4.2\%$, Circulaire $8.7 \pm 1.0\%$, and **AeroEclipse*** $38.7 \pm 1.3\%$. The mean \pm SD percentages of drug lost to ambient air were: Misty-Neb $26.8 \pm 0.7\%$, SideStream $17.3 \pm 0.4\%$, Pari LCD $18.3 \pm 0.8\%$, Circulaire $12.3 \pm 0.8\%$, and **AeroEclipse*** $6.6 \pm 3.3\%$. The mean \pm SD percentages of drug lost to deposition in the apparatus were: Misty-Neb $52.3 \pm 0.6\%$, SideStream $63.4 \pm 3.0\%$, Pari LCD $62.5 \pm 4.0\%$, Circulaire $75.8 \pm 0.5\%$, and **AeroEclipse*** $51.0 \pm 2.1\%$. Duration of nebulization was shortest with the Circulaire and longest with the **AeroEclipse*** ($p < 0.05$ via 1-way analysis of variance). **Conclusions:** The nebulizers we tested differ significantly in overall drug disposition. The dosimetric **AeroEclipse*** provided the largest inhaled drug mass and the lowest loss to ambient air, with the test conditions we used.

2. **COMPARISON OF BREATH-ACTUATED JET NEBULIZER (BAN) IN 'CONTINUOUS DELIVERY' MODE WITH OTHER CONTINUOUS DELIVERY NEBULIZERS.** JP Mitchell, KJ Wiersema, CC Doyle and MW Nagel. Respiratory Care 2003; 48 (11): S1077.

The **AeroEclipse*** BAN (Monaghan Medical Corp., Plattsburgh, N.Y.) has been equipped with an optional blue cap whose purpose is to retain the actuator piston in the position it would occupy during inhalation in breath-actuated mode, so that the nebulizer operates continuously. The present study compared the delivery of a bronchodilator from diluted albuterol sulfate respirator solutions (3-ml of 0.83 and 1-ml of 2.5 mg/ml albuterol in physiologically normal saline (0.9% w/v NaCl)), via this nebulizer, the Micromist® (Hudson RCI, Temecula, CA), Misty-Neb™ (Allegiance Healthcare Corp., McGaw Park, IL) and the LCD™ (PARI Respiratory Equipment, Monterey, CA). Each nebulizer was tested using a breathing simulator set to the following parameters representative of adult use: tidal volume = 600-ml, rate = 10 breaths/min, inspiratory/expiratory ratio 1:2. The total mass of albuterol (TM) delivered to the first sputter was determined by filter collection at the mouthpiece of the nebulizer operated with compressed air supplied at 50 psig at 8 L/min ($n = 5$ devices/group, 3 replicates/device). The fraction of the aerosol contained in droplets finer than $4.8 \mu\text{m}$ aerodynamic diameter (FPF) was determined by laser diffractometry in a parallel study, so that the fine droplet mass (FM) could be calculated as the product of TM and FPF.

Values of FM (mean \pm SD) and time to deliver medication (T_{med}) were as follows:

| Solution(mg/ml) | AeroEclipse* | | LCD™ | | Micromist® | | MistyNeb™ | |
|------------------------|---------------------|-----------------|-----------------|----------------|-----------------|-------------------|---------------|---------------|
| | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 |
| FM (μg) | 360 ± 22 | 263 ± 26 | 149 ± 16 | 108 ± 4 | 209 ± 12 | 15.4 ± 5.9 | 82 ± 9 | 31 ± 5 |
| T_{med} (min) | 3 | <1 | 2 | <1 | 7 | <1 | 4 | <1 |

The **AeroEclipse*** nebulizer delivered significantly more FM in continuous delivery mode than the other nebulizers when operated in continuous mode with either solution strength (1-way repeated measures ANOVA, $p < 0.05$). T_{med} from the **AeroEclipse*** nebulizer was comparable with the best performing continuous nebulizer (LCD™).

3. **CLINICAL AND ECONOMIC IMPACT RESULTING FROM A HOSPITAL-WIDE CONVERSION FROM SMALL VOLUME NEBULIZERS TO THE AEROECLIPSE* BREATH ACTUATED NEBULIZERS.** RS Pikarsky, R Acevedo. Respiratory Care, September 2002; 47(9): 1075.

Purpose: The Respiratory Department converted from the Airlife Misty-Neb (SVN) (Allegiance Healthcare Corporation) to the **AeroEclipse*** Breath Actuated Nebulizer (BAN) (Monaghan Medical Corp. Plattsburgh, N.Y.). This study explores the clinical and economic impact of these interventions. **Methods:** Patients capable of performing aerosol therapy by mouthpiece were converted to BAN. All aerosols treatments, including breakthrough treatments, delivered between 7/1/01 and 2/28/02 were recorded. Fifty-four percent of the nebulizers purchased during this period were BAN. The 5.6-minute timesaving of the BAN over the SVN was from our previous pilot study¹. The FTE average cost (salary/benefits)= \$20.25/hr. The average treatments per patient - day (4.90) were determined from the treatment records. Each patient received an average of 2.5 nebulizers during their admission. **Results & Economic Impact:** The table shows the total number of treatments with estimated timesaving of 0.6 FTE over the 8-month period, or 0.9 FTE annualized (\$37,746). Total number of nebulizers used was estimated at 1,359, 53.8% were BAN. The increase cost of the BAN was \$4,153. Overall savings was \$33,592. **Resource Utilization:** Compared with the same months studied of the prior year, omitted treatments due to "Therapist Unavailable" decreased from 1.79% to 1.40% ($p < .001$). Worked FTEs decreased from 38.05 in 2000 to 33.42 in first quarter 2002. Paid FTEs decreased from 43.76 to 37.59 in the same time periods.

| | | | |
|--|-----------------|----------------------------------|----------------|
| Total treatments | 16,651 | Total treatments | 16,651 |
| % Conversion | 54% | Ave Tx /pt-day | 4.90 |
| # Treatments BAN | 8,958 | Ave pt-days | 3,398 |
| Time saved per treatment (min) | 5.6 | Ave nebs/pt-day | 2.5 |
| Hours saved | 829 | Total nebs used | 1,359 |
| FTE (8 months) | 0.60 | Cost of each Misty Neb | \$0.76 |
| FTE Annualized | 0.90 | Cost of each AeroEclipse* | \$4.55 |
| | | % Conversion | 54% |
| 1.0 FTE RT | \$42,120 | Cost with conversion | \$3,801 |
| 0.9 FTE Savings with AeroEclipse* | \$37,746 | Cost without conversion | \$1,032 |
| Conversion cost difference | -\$4,153 | Conversion cost difference | \$2,769 |
| Total Savings | \$33,592 | Annualized increase | \$4,153 |

Conclusion: Hospital-wide conversion to BAN is cost-effective due to the decrease in administration time. Therapist availability was enhanced, contributing to a significant reduction in omitted treatments. **Clinical Implications:** The conversion to BAN allows the ability to meet our patient care demands and for the reallocation of workforce needs in a manner that is clinically and economically advantageous. This may be one of several strategies to address the problems of Respiratory Care Staff shortages or high seasonal effect in the acute care facility. Quality is enhanced by decreasing the incidence of omitted therapy.

¹ Robert S. Pikarsky, Tracey Farrell, Russ Acevedo, Wendy Fascia and Charles Roman. *The Delivery Time, Efficacy, and Safety Of Beta Agonist Bronchodilator Administration with the Aeroeclipse* Breath - Actuated Nebulizer (BAN)*. CHEST 2001;120(4)218S

4. **THE DELIVERY TIME, EFFICACY, AND SAFETY OF BETA AGONIST BRONCHODILATOR ADMINISTRATION WITH THE AEROECLIPSE* BREATH - ACTUATED NEBULIZER ("BAN").** RS Pikarsky, T Farrell, R Acevedo, W Fascia and C Roman. CHEST 2001; 120(4) 218S.

Purpose: Aerosol delivery consumes the highest level of Respiratory Care resources. This study evaluated the delivery time, efficacy, and safety of rapidly nebulized Albuterol with the use of the **AeroEclipse*** Breath Actuated Nebulizer as compared to both an MDI with **AeroChamber*** VHC (both from Monaghan Medical Corp. Plattsburgh, N.Y.) and the Airlife Misty-Neb Nebulizer (SVN) (Allegiance Healthcare Corporation). **Methods:** A consecutive, non-randomized, mostly COPD population receiving pre & post bronchodilator testing in our Pulmonary Function Lab were studied. Three different Albuterol medication dosages were administered with the BAN: 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline, 1.0 ml (5 mg) of undiluted Albuterol, and 0.75 ml Albuterol (3.75 mg) using an oxygen flow rate of 8 L/min. Two puffs of Albuterol were administered by MDI with **AeroChamber*** VHC. Treatments with the SVN consisted of nebulizing 2.5 mg of Albuterol diluted with 3 ml of Normal Saline Unit Dose (UD) using an oxygen flow rate of 8 L/min. The Sensormedics Vmax 22 Pulmonary Function System was utilized to measure FEV1. A standardized subjective questionnaire to determine side effects was completed.


| Nebulizer (n) | Dose | % Change FEV1 | Time (min) | Tremulousness |
|------------------------------|--------------------|---------------|------------|---------------|
| AeroEclipse* BAN (12) | 0.5 ml + 0.5 ml NS | 8.2% | 2.67* | 0 |
| AeroEclipse* BAN (64) | 1.0 ml undil. | 10.9% | 3.29* | 17 |
| AeroEclipse* BAN (23) | 0.75 ml undil. | 5.6% | 1.30* | 5 |
| MDI (21) | 2 puffs | 8.5% | 2.86** | 1 |
| Misty-Neb (52) | 2.5 mg UD | 9.1% | 8.33 | 2 |

Results: The table shows the Albuterol dosages, mean % change of FEV1 from pre-treatment and 10 minute post treatment, mean administration time and tremulousness. The mean treatment time with all BAN patients was 2.78 minutes as compared to 8.33 minutes with the SVN ($p < .001$) *. The mean treatment time with the MDI was 2.86 minutes as compared to 8.33 minutes with the SVN ($p < .001$) **. The changes in FEV1 were not significant. There was no difference in heart rate, respiratory rate or nausea. Seventeen patients receiving the 1.0 l undiluted Albuterol indicated an increase in tremulousness. **Conclusion:** The rapid administration of Albuterol in the 0.5 ml + 0.5 ml NS and 1.0 ml undiluted doses using the BAN was equally efficacious as the MDI with **AeroChamber*** VHC and SVN UD. The 1.0 ml Albuterol dosage has the highest incidence of tremulousness. The 0.75 ml Albuterol dosage under-performed. Delivering 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline using the BAN offered the best delivery time, efficacy and safety profile of the nebulizer trials. The BAN performance was comparable to the MDI with **AeroChamber*** VHC. **Clinical Implications:** In a health care facility that delivers large volumes of aerosol treatments, the decrease in delivery time could have a significant impact on resource utilization. The results supported changes in the Respiratory Care practice throughout Crouse Hospital. Further studies evaluating additional medication dosing regimens measuring safety, efficacy and resource utilization are needed.

5. **THE CLINICAL EFFICACY OF USING THE AEROECLIPSE* BREATH ACTUATED NEBULIZER ("BAN") IN PULMONARY LAB TESTING AND IMPLICATIONS FOR GENERAL USE.** YM Christensen, CJ Flanagan, SA Ravenscraft. Respiratory Care 2001;46(10):1084.

Purpose: To compare the clinical efficacy and delivery time of nebulization of beta agonist bronchodilator with the use of the **AeroEclipse*** Breath Actuated Nebulizer ("BAN") (Monaghan Medical Corp.) as compared to the AirLife Misty-Neb Nebulizer(SVN) (Allegiance Healthcare Corporation). **Methods:** Adult patients (n=40) presenting with Asthma (50%), COPD (10%) and other pulmonary disorders (40%); receiving pre and post bronchodilator spirometry testing in our Pulmonary Function Lab were included in the study. Each patient received both nebulizers on two separate visits (less than 24 hours apart). Patient received a nebulizer treatment with the BAN (n=40) 2.5mg Albuterol (0.5ml) in 0.5cc saline run to sputter, or the SVN (n=40) 2.5mg Albuterol in 2.5cc saline (3ml unit dose) run to sputter. FVC, FEV1, FEV1% ratio and FEF 25-75% spirometry was conducted using the Medical Graphics 1085DX pre and 5 minutes post treatment with the BAN and 10 minutes post treatment with the SVN. **Results:** The results demonstrated that FVC, FEV1 and FEF 25-75% for patients using the BAN were substantially higher while FEV1% ratio favored the SVN (Table and Chart). Importantly, total nebulization time was reduced from 22 minutes (SVN) to 7 minutes (BAN), and total test time was reduced from 30 minutes (SVN) to 15 minutes (BAN).

| SPIROMETRY RESULTS | | | | |
|-----------------------------|------|------|------------------|-------|
| Absolute % Change by Device | | | % Difference BAN | |
| | SVN | BAN | | BAN |
| FVC | 5.3 | 10.2 | FVC | 91.3 |
| FEV1 | 7.3 | 13.1 | FEV1 | 79.8 |
| FEV1%ratio | 3.0 | 2.3 | FEV1% | -25.1 |
| FEF 25-75% | 29.8 | 57.7 | FEF 25-75% | 93.3 |



Conclusion: The administration of 2.5mg of albuterol with the BAN produced improved results in FVC, FEV1 and FEF 25-75%. Substantially shorter test times delivered by the BAN would allow for more tests and associated revenue. These data support the thesis that the BAN can reduce costs of care by delivering clinically acceptable outcomes in significantly less time.

6. **THE DELIVERY TIME, EFFICACY, AND SAFETY OF BETA AGONIST BRONCHODILATOR ADMINISTRATION WITH THE *AeroEclipse*[®] BREATH ACTUATED NEBULIZER ("BAN") VERSE A CONVENTIONAL T-TYPE SMALL VOLUME NEBULIZER.** RS Pikarsky, T Farrell, R Acevedo, W Fascia and C Roman. Respiratory Care 2001;46(10):1085.

Purpose: Aerosol delivery consumes the highest level of Respiratory Care resources. This study evaluated the delivery time, efficacy, and safety of rapidly nebulized albuterol with the use of a novel breath actuated nebulizer compared to a standard small volume nebulizer. **Methods:** A consecutive, non-randomized, mostly COPD population receiving pre & post bronchodilator testing in our Pulmonary Function Lab were studied. 0.5 ml albuterol (2.5 mg) with 0.5 ml Normal Saline (NS) was administered with the *AeroEclipse*[®] Breath Actuated Nebulizer ("BAN") (Monaghan Medical Corp. Plattsburgh, N.Y.) using an oxygen flow rate of 8 L/min. Treatments with the AirLife[™] brand Misty-Neb[™] small volume nebulizer (SVN) (Allegiance Healthcare Corporation) consisted of nebulizing 2.5 mg of albuterol diluted with 3 ml of Normal Saline Unit Dose (UD) using an oxygen flow rate of 8 L/min. The Sensormedics Vmax 22 Pulmonary Function System was utilized to measure FEV1. A standardized subjective questionnaire to determine side effects was completed. **Results:** The table shows the albuterol dosages, mean % change of FEV1 from pre-treatment and 10 minute post treatment, mean administration time and tremulousness. The mean treatment time with all BAN patients was 2.67 minutes as compared to 8.33 minutes with the SVN (p<.001)*. The changes in FEV1 were not significant. There was no difference in heart rate, respiratory rate or nausea. **Conclusion:** The rapid administration of albuterol in the 0.5 ml + 0.5 ml NS dose using the BAN was equally efficacious as the SVN UD. Delivering 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline using the BAN offered the best delivery time, efficacy and safety profile between the two devices. **Clinical Implications:** In a health care facility that delivers large volumes of aerosol treatments, the decrease in delivery time achieved with the BAN could have a significant impact on resource utilization. The results supported changes in the Respiratory Care practice throughout Crouse Hospital. Further studies evaluating additional medication dosing regimens measuring safety, efficacy and resource utilization are needed.

| Nebulizer (n) | Dose | % Change FEV1 | Time (min) | Tremulousness |
|--|--------------------------------------|---------------|------------|---------------|
| <i>AeroEclipse</i> [®] BAN (12) | 2.5mg (0.5 ml albuterol + 0.5 ml NS) | 8.2% | 2.67* | 0 |
| Misty-Neb [™] (52) | 2.5mg (3 ml unit dose) | 9.1% | 8.33 | 2 |

7. **PREDICTING LUNG DEPOSITION WITH A CASCADE IMPACTOR.** S Sangwan, F Hull, R Condos and GC Smaldone. Journal of Aerosol Medicine, 2001; 14(3):421. *Presented at the 13th International Congress on Aerosols in Medicine, Interlaken, Switzerland, September 17-21, 2001.

Introduction: In recent deposition studies of interferon- β , we failed to predict the deposition pattern from bench studies of aerosols using multistage cascade impaction (MCI). Recent mass balance studies have identified impaction in connecting tubing and effects of breathing on interpretation of cascade data (Gurses BK et al AJRCC 163; 5(A166). 2001). In the present study we related MCI data using our new bench test protocol directly to lung scans in humans. This protocol emphasizes deposition of large particles in connecting tubing and influence of conditions internal to the nebulizer during breathing. **Methods:** Two devices (Misty-Neb and *AeroEclipse*[®] Breath-Actuated Nebulizer ("BAN")) were studied. Mass median aerodynamic diameter (MMAD) and mass balance were measured under standing cloud and ventilation using a piston pump. Deposition images were obtained using gamma camera.

Results:

| Nebulizer & method of assessment | | Respirable Mass [†] (<6 μ m) | Regional Deposition | |
|-------------------------------------|----------------|---|---------------------|---------------------|
| | | | Lung deposition** | Throat deposition** |
| Misty-Neb | Standing Cloud | 46.2% | 32% | 68% |
| | Ventilated | 24.6% | | |
| <i>AeroEclipse</i> [®] BAN | Standing Cloud | 48.3% | 72% | 28% |
| | Ventilated | 71.2% | | |


[†] Calculated by adding T connector deposition to the first stage (>8 μ m) of cascade

** Expressed as Percent of total deposition in the body

Conclusion: Regional deposition (upper airway vs. lung) was predicted by analysis only when effects of both connecting tubing and breathing were considered in the bench protocol.

8. **EFFECT OF NEBULIZER DESIGN ON FINE PARTICLE MASS.** D Hess, JP Mitchell, D Coppolo, MW Nagel, AD Archer, R. Blacker. Presented at Open Forum, Ann. Meet. of the American Association For Respiratory Care (AARC), Las Vegas, 1999. Published in *Resp Care*, Oct 1999, 44(10): 1289.

Background: Nebulizer design is known to affect performance. In this study, we compared fine particle mass from nebulizers of four designs. **Methods:** We tested traditional disposable nebulizers (Baxter Misty-Neb™, Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™), nebulizers with collection bags (Westmed Circulaire™), and a Trudell **AeroEclipse*** (with breath actuation disabled). Five of each device with three replicates (n = 15) were tested using an *in-vitro* model of spontaneous breathing. A rigid bar was placed between the two compartments of a test lung (Michigan Instruments TTL). The drive lung was attached to a ventilator (Infrasonics Infant Star™) to simulate spontaneous breathing (tidal volume 0.6 L, rate 10/min, T_i 2 s). A bacterial/viral filter (Trudell MT3000) was placed between the nebulizer and slave lung. Flow was monitored between the test lung and filter (Novamatrix Ventcheck™). Albuterol solution (0.625 mg/mL) was placed into the nebulizers (4 mL), which were powered with air (8 L/min). Filters were replaced at one minute intervals (flow to the nebulizer was discontinued during filter replacement) until sputtering occurred. The filter was washed with methanol (20 mL) and albuterol concentration was measured with HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer. Fine particle mass was calculated as the product of mass % <4.7 µm and total nebulizer output. **Results:** Fine particle mass from the **AeroEclipse*** nebulizer was greater than that from the other nebulizers (P<0.001) (see figure). **Conclusions:** Fine particle mass was affected by nebulizer design. The clinical relevance of this finding awaits further investigation. Further evaluation of the breath-actuated feature of the **AeroEclipse*** is warranted.



9. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS - THE RELATIONSHIP BETWEEN NEBULIZED DROPLET SIZE AND THE PARTICLE SIZE OF THE SUSPENSION.** JP Mitchell, MW Nagel and AD Archer. Presented at Drug Delivery to the Lungs-IX, London, UK, 1998, in *J. Aerosol Med.* 12(3), 208, (1999).

A new air entrainment small volume nebulizer (AE-SVN) has been compared with two other SVN's (Neb-U-Mist™ and Misty-Neb™) for the delivery of a suspension of 0.25 µg/ml budesonide. Each SVN was operated at 8 l/min with compressed oxygen (50 psig). The total mass of budesonide was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The time-averaged delivery rate over the period of nebulization ((mean ± 1 S.D.) µg budesonide/min) from the AE-SVN (102 ± 9) was greater than with the Misty-Neb™ (49 ± 2), or Neb-U-Mist™ (25 ± 6). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with the Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). The mass median diameter (MMD) of the droplets from the AE-SVN measured using a laser diffractometer (2.9 ± 0.1 µm), was significantly finer compared with those from the Misty-Neb™ (4.5 ± 0.9 µm) and Neb-U-Mist™ (5.6 ± 0.6 µm) and closest to the size of the micronized budesonide particles in the original suspension. The efficient delivery of medication formulated as micronized powder in aqueous suspension necessitates that the droplets produced upon nebulization are large enough so that single particles are efficiently entrained during atomization, but not so coarse that they cannot leave the nebulizer, extending nebulization time.

10. **EVALUATION OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) WHEN USED WITH OXYGEN AS A DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** JP Mitchell and MW Nagel. Presented at ALA/ATS International Conference, Chicago, 1998.

Purpose: To compare the delivery of saline (0.9% w/v NaCl) by a new AE-SVN (Trudell Medical Int.) with that from two other representative SVN's (Updraft Neb-U-Mist® (Hudson Oxygen Therapy Sales Co.) and Airlife™ Misty-Neb™ (Baxter Healthcare Corp.)) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 AE-SVN's were tested using a laser diffractometer (Malvern Mastersizer-X) to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. 5, Neb-U-Mist® and a similar number of MistyNeb™ SVN's were also evaluated. **Results:** Total (TM) and respirable (RM, droplets finer than 4.8 µm diameter) mass output rates and droplet mass median diameter (MMD) were as follows: AE-SVN: TM = 671 ± 26 µg/min, RM = 542 ± 23 µg/min (80.8 ± 1.3% respirable), MMD = 2.88 ± 0.09 µm; Neb-U-Mist™: TM = 266 ± 13 µg/min, RM = 119 ± 16 µg/min (42.1 ± 5.2% respirable), MMD = 5.6 ± 0.6 µm; Misty-Neb™: TM = 336 ± 60 µg/min, RM = 178 ± 43 µg/min (53.1 ± 8.5 % respirable), MMD = 4.5 ± 0.9 µm. **Conclusion:** TM from the new AE-SVN was substantially greater than those from either the Neb-U-Mist® or Misty-Neb™ (1-way ANOVA, p < 0.001). The finer MMD produced from the AE-SVN resulted in a significantly greater RM compared with either of the other SVN's (p < 0.001).

11. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS: A COMPARATIVE *IN-VITRO* ASSESSMENT.** JP Mitchell, MW Nagel and AD Archer. Presented at Proc. Ann. Meet. Amer. College of Chest Physicians (ACCP), Toronto, November 1998, in Chest, 114 (4S), 295, (1998); and at the World Asthma Meeting, Barcelona, December 1998, in Eur. Resp. J., 12S29, 7, (1998).

Purpose: To compare the performance of a new air entrainment small volume nebulizer (AE-SVN, Trudell Medical International with other widely used SVN's (LC-Star™ (PARI Respiratory Equipment), Updraft™ Neb-U-Mist™ (Hudson Oxygen Therapy Sales Co.), Circulaire™ (Westmed), Sidestream™ (Medic-Aid), Airlife™ Misty-Neb™ (Baxter Healthcare Corp.)) for the delivery of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). **Methods:** Each SVN (n = 5 devices for each group, 3 replicates per device) was operated with compressed air at 50 psig at a flow rate of 8 l/min. The total mass of budesonide nebulized from 2 x 2 ml ampoules was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The SVN was operated until it spluttered, was then tapped gently to dislodge droplets back to the reservoir. Nebulization was deemed complete 20 seconds later. The mass of budesonide collected was determined by HPLC-UV spectrophotometry. **Results:** The delivery rate ((mean ± 1 S.D) µg budesonide/min) from the AE-SVN (102 ± 9) was significantly greater than with the other groups: (LC-Star™ (91 ± 6), Misty-Neb™ (49 ± 2), Sidestream™ (46 ± 4), Circulaire™ (26 ± 4) and Neb-U-Mist™ (25 ± 6)), (1-way ANOVA, p < 0.02). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with LC-Star™ (229 ± 10 s), Sidestream™ (365 ± 19 s), Circulaire™ (420 ± 84 s), Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). **Conclusions:** The new AE-SVN is highly efficient at entraining the budesonide particles into the liquid droplets at these conditions. **Clinical Implications:** The good delivery rate combined with comparatively short duration of delivery offers the potential for rapid treatment and patient convenience.

Misty Max 10™, Cardinal Health

1. **DELIVERY OF ALBUTEROL VIA A NEW BREATH ACTUATED NEBULIZER: COMPARISON WITH CONTINUOUS JET NEBULIZERS.** DP Coppola, MW Nagel, CC Doyle, VA Avvakoumova and JP Mitchell. Presented at the American Thoracic Society International Conference 2007, San Francisco, California, USA.

A new breath actuated nebulizer (**AeroEclipse**® II BAN, Monaghan Medical Corp., Plattsburgh, NY) has been developed to deliver medication only when the patient inhales. This study sought to determine the delivery of albuterol (3-ml fill of diluted solution (0.83 mg/ml)) as fine droplets < 4.7 µm aerodynamic diameter, and compare this fine droplet mass (FDM) with equivalent data from 4 widely available continuous jet nebulizers as benchmark devices. Each nebulizer (n=5; 3 replicates/device) was operated with compressed air at 50 psig at ca. 8 L/min to simulate hospital wall outlet conditions. The nebulizer on test was coupled to a breathing simulator set to mimic adult use (tidal volume = 600 ml, rate = 10 breaths/min; duty cycle = 0.33), and the emitted droplets were collected on an electret filter at the mouthpiece. The total mass of albuterol (TM) was assayed subsequently by HPLC-UV spectrophotometry. In a separate study, the droplet size distribution was determined by laser diffractometry so that the fine droplet fraction (FDF) could be obtained. FDM was determined as the product of TM and FDF. FDM (mean SD) from the BAN was 791 84 g, delivered in 8 minutes. Corresponding values (FDM in time from start to sputter) for the VixOne™ (Westmed, Tucson, AZ), MicroMist™ (Hudson RCI, Temecula CA), Misty Max 10™ (Cardinal Health, McGaw Park (IL) and model 8900™ (Salter Labs, Arvin, CA) were 267 11 g in 6 min, 133 8 g in 4 min, 249 10 g in 6 min and 161 10 g in 5 min. Aside from dosage assurance imparted by breath-actuation, the **AeroEclipse**® II BAN delivered substantially more FDM/min than the other devices. The clinician is now able to treat either for extended high dose delivery (potentially eliminating the need for additional therapy), or titrate to a shorter interval based on response.

Pari LC-JET™, PARI™ Respiratory Products Inc.

1. **PERFORMANCE OF A NEW BREATH-ACTUATED SMALL VOLUME NEBULIZER (BA-SVN) WHEN USED WITH OXYGEN AS A DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** AM Verdun, JP Mitchell and MW Nagel. Presented at ALA/ATS International Conference, Chicago, 1998.

Purpose: To compare the delivery of saline (0.9% w/v NaCl) by a new BA-SVN (Trudell Medical Int..) with that from two other representative SVN's (LC-JET™ (PARI Respiratory Products Inc., Canada) and reusable Sidestream™ (MedicAid, UK)) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 BA-SVN's were tested using a laser diffractometer (Malvern Mastersizer-X) to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. 5, LC-JET™ and 5, Sidestream™ SVN's were also tested similarly. The BA-SVN was operated with manual over-ride engaged (continuous delivery of aerosol). **Results:** Total (TM) and respirable (RM), droplets finer than 4.8 µm diameter) mass output rates and droplet mass median diameter (MMD) were as follows: BA-SVN: TM = 672 ± 23 µg/min, RM = 545 ± 31 µg (80.9 ± 2.4% respirable), MMD = 2.79 ± 0.15 µm; LC-JET™: TM = 675 ± 69 µg/min, RM = 449 ± 41 µg/min (66.7 ± 1.8% respirable), MMD = 3.39 ± 0.08 µm; Sidestream™: TM = 442 ± 26 µg/min, RM = 358 ± 38 µg/min (80.8 ± 4.2 %respirable), MMD = 2.94 ± 0.03 µm. **Conclusion:** Although TM from the new BA-SVN was comparable with that from the LC-JET™ (Mann-Whitney rank sum test, p = 0.84), the finer MMAD produced from the BA-SVN resulted in a significantly greater RM (p < 0.001). Both TM and RM from the BA-SVN were greater than those from the Sidestream™ SVN (p < 0.001).

2. **COMPARISON OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) WITH OTHER SVNS WHEN USED WITH OXYGEN AS DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** JP Mitchell and MW Nagel. Presented at Ann. Meet. Amer. Assoc. of Asthma, Allergy and Immunology (AAAAI), Washington D.C., 1998.

The performance of a prototype novel AE-SVN (Trudell Medical International (n = 5)) with normal saline (0.9% w/v NaCl) operating at $20 \pm 2^\circ\text{C}$, $50 \pm 10\%$ RH, has been evaluated with oxygen (50 psig, 8 l/min) as driving gas to simulate hospital use. Comparison testing was also undertaken with two other representative AE-SVNs, (a) LC-JET[™] (Pari Respiratory Equipment Inc.), without inspiratory valve cap which would otherwise restrict aerosol output, (b) SideStream[™] (MedicAid, UK). A laser diffractometer (Malvern Mastersizer-X) was used to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. Total (T) and respirable ((R), droplets finer than $4.8 \mu\text{m}$ aerodynamic diameter) mass output rates and droplet mass median aerodynamic diameter (MMAD) for the new AE-SVN (5 replicate measurements/device) were: $671 \pm 26 \mu\text{g/min}$ (T), $542 \pm 23 \mu\text{g/min}$ (R) and $2.88 \pm 0.09 \mu\text{m}$ (MMAD). Corresponding data for the LC-JET[™] were: $675 \pm 65 \mu\text{g/min}$ (T), $450 \pm 45 \mu\text{g/min}$ (R) and $3.39 \pm 0.14 \mu\text{m}$ (MMAD), and for the SideStream[™] were: $442 \pm 27 \mu\text{g/min}$ (T), $357 \pm 28 \mu\text{g/min}$ (R) and $2.95 \pm 0.13 \mu\text{m}$ (MMAD). The total aerosol delivery rate from the new AE-SVN matched that of the LC-JET[™] (un-paired t-test, $p = 0.79$) and exceeded that from the SideStream[™] ($p < 0.001$). The finer MMAD of the aerosol provided by the new AE-SVN resulted in a significantly greater respirable mass fraction, increasing the respirable mass delivery rate compared with the other SVNs ($p < 0.001$).

3. **COMPARISON OF A NEW BREATH ACTUATED SMALL VOLUME NEBULIZER (BA-SVN) WITH AN SVN SUPPLIED WITH COMPRESSOR INTENDED FOR HOME CARE USE.** AM Verdun, JP Mitchell and MW Nagel. Presented at Ann. Meet. Amer. Assoc. of Asthma, Allergy and Immunology (AAAAI), Washington D.C., 1998.

The performance of a prototype novel BA-SVN (Trudell Medical International (n = 5 devices)) with normal saline (0.9% w/v NaCl) operating at $20 \pm 2^\circ\text{C}$, $50 \pm 10\%$ RH, has been evaluated with an air compressor widely used in home care (Proneb[™], Pari Respiratory Equipment Inc.). A laser diffractometer (Malvern Mastersizer-X) was used to determine the size distribution of droplets emitted at the mouthpiece (5 replicates per device). The total mass output was determined gravimetrically in a parallel series of tests. The BA-SVN was operated with manual over-ride engaged (continuous delivery of aerosol). Total (T) and respirable ((R), droplets finer than $4.8 \mu\text{m}$ aerodynamic diameter) mass output rates, respirable mass fraction (RM) and droplet mass median aerodynamic diameter (MMAD) were $167 \pm 6 \mu\text{g/min}$ (T), $96 \pm 5 \mu\text{g/min}$ (R), $57.5 \pm 2.1\%$ (RM) and $4.40 \pm 0.11 \mu\text{m}$ (MMAD). In comparison, under similar conditions, a Pari LC-JET[™] SVN with Proneb[™] (n = 5 replicate measurements) provided $211 \pm 3 \mu\text{g/min}$ (T), $65 \pm 4 \mu\text{g/min}$ (R), $30.9 \pm 1.5\%$ (RM) and $6.94 \pm 0.20 \mu\text{m}$ (MMAD). The new BA-SVN provided aerosol having a finer MMAD and greater RM (un-paired t-test, $p < 0.001$ for each variable) which resulted in an improved respirable mass output rate compared with the LC-JET[™] SVN. The BA-SVN also has the advantage that no aerosol is produced to waste during the exhalation portion of each breathing cycle.

Pari LC-Star[™], PARI[™] Respiratory Equipment

1. **DELIVERY OF BUDESONIDE INHALATION SOLUTION (BIS) THROUGH AN INFANT UPPER AIRWAY MODEL.** DE Geller, KC Kesser, HM Janssens, HAWM Tiddens. Am J Respir Crit Care Med 2003; 167(7):A508.

We investigated variables that may be important in the delivery of BIS to the lungs of infants, a challenging population for aerosol delivery. **Methods:** The Sophia Anatomical Infant Nose Throat (SAINT) airway model mounted on a breath simulator mimicked the breathing pattern of a 9-mo old infant ($\text{RR}=30$, $\text{Vt}=100 \text{ ml}$, I:E ratio=1:1.3). Nebulizers were charged with BIS 0.25 mg and run continuously until dry. Drug captured on a filter distal to the SAINT model was the lung dose. Compressor: PARI PRONEB TURBO. Nebulizer/mask systems studied: VIX1/aerosol mask (AM), PediNeb pacifier device (PN) or blow-by (BB); **AeroEclipse[®]** neb and mask (AE); PARI LC+ and PARI LC*/ PARI Baby mask (PB), Fish mask (FM), and AE masks. The AE neb/mask was also studied with an ill breathing pattern ($\text{RR}=50$, $\text{Vt}=100$, I:E=1:2). **Results:** Lung dose ranged from 2.0 to 7.6% of the neb charge. Lung dose was AE (5.0%) > VIX1 (3.5%), LC+/FM (3.2%), LC*/PB (2.9%), and LC+/PB (2.8%). Also, VIX1/AM (3.5%)>VIX1/PN (2.5%)>VIX1/BB (2.0%). The lung dose of the LC+ and LC* more than doubled (6.8 and 6.3%) when used with the AE mask. Lung dose increased with the ill breath pattern in proportion to increased minute ventilation (7.6%). **CONCLUSION:** 1) The AE system provided higher lung dose than other nebulizers with standard masks. 2) Mask design and fit can substantially impact nebulizer performance. 3) PN performed better than BB, but not as good as a mask. If crying decreases lung dose by 75%, we speculate that the PN and BB (non-crying) may improve lung dose vs mask with a crying infant. 4) An increase in lung dose may occur in ill infants if minute ventilation is elevated. *This abstract is funded by: AstraZeneca.*

2. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS: A COMPARATIVE IN-VITRO ASSESSMENT.** JP Mitchell, MW Nagel and AD Archer. Presented at Proc. Ann. Meet. Amer. College of Chest Physicians (ACCP), Toronto, November 1998, in Chest, 114 (4S), 295, (1998); and at the World Asthma Meeting, Barcelona, December 1998, in Eur. Resp. J., 12S29, 7, (1998).

Purpose: To compare the performance of a new air entrainment small volume nebulizer (AE-SVN, Trudell Medical International) with other widely used SVNs (LC-Star[™] (PARI Respiratory Equipment), Updraft[™] Neb-U-Mist[™] (Hudson Oxygen Therapy Sales Co.), Circulaire[™] (Westmed), Sidestream[™] (Medic-Aid), AirLife[™] Misty-Neb[™] (Baxter Healthcare Corp.)) for the delivery of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). **Methods:** Each SVN

(n = 5 devices for each group, 3 replicates per device) was operated with compressed air at 50 psig at a flow rate of 8 l/min. The total mass of budesonide nebulized from 2 x 2 ml ampoules was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The SVN was operated until it spluttered, was then tapped gently to dislodge droplets back to the reservoir. Nebulization was deemed complete 20 seconds later. The mass of budesonide collected was determined by HPLC-UV spectrophotometry. **Results:** The delivery rate ((mean \pm 1 S.D) μ g budesonide/min) from the AE-SVN (102 \pm 9) was significantly greater than with the other groups: (LC-Star[™] (91 \pm 6), Misty-Neb[™] (49 \pm 2), Sidestream[™] (46 \pm 4), Circulaire[™] (26 \pm 4) and Neb-U-Mist[™] (25 \pm 6)), (1-way ANOVA, p < 0.02). Duration of nebulization was shortest with the AE-SVN (221 \pm 14 s), compared with LC-Star[™] (229 \pm 10 s), Sidestream[™] (365 \pm 19 s), Circulaire[™] (420 \pm 84 s), Misty-Neb[™] (477 \pm 25 s) and Neb-U-Mist[™] (639 \pm 15 s). **Conclusions:** The new AE-SVN is highly efficient at entraining the budesonide particles into the liquid droplets at these conditions. **Clinical Implications:** The good delivery rate combined with comparatively short duration of delivery offers the potential for rapid treatment and patient convenience.

Pari LC Plus[™], PARI[™] Respiratory Equipment

1. **DELIVERY OF BUDESONIDE INHALATION SOLUTION (BIS) THROUGH AN INFANT UPPER AIRWAY MODEL.** DE Geller, KC Kesser, HM Janssens, HAWM Tiddens. Am J Respir Crit Care Med 2003; 167(7):A508.

We investigated variables that may be important in the delivery of BIS to the lungs of infants, a challenging population for aerosol delivery. **Methods:** The Sophia Anatomical Infant Nose Throat (SAINT) airway model mounted on a breath simulator mimicked the breathing pattern of a 9-mo old infant (RR=30, Vt=100 ml, I:E ratio=1:1.3). Nebulizers were charged with BIS 0.25 mg and run continuously until dry. Drug captured on a filter distal to the SAINT model was the lung dose. Compressor: PARI PRONEB TURBO. Nebulizer/mask systems studied: VIX1/aerosol mask (AM), PediNeb pacifier device (PN) or blow-by (BB); **AeroEclipse**^{*} neb and mask (AE); PARI LC+ and PARI LC*/ PARI Baby mask (PB), Fish mask (FM), and AE masks. The AE neb/mask was also studied with an ill breathing pattern (RR=50, Vt=100, I:E=1:2). **Results:** Lung dose ranged from 2.0 to 7.6% of the neb charge. Lung dose was AE (5.0%) > VIX1 (3.5%), LC+/FM (3.2%), LC*/PB (2.9%), and LC+/PB (2.8%). Also, VIX1/AM (3.5%)>VIX1/PN (2.5%)>VIX1/BB (2.0%). The lung dose of the LC+ and LC* more than doubled (6.8 and 6.3%) when used with the AE mask. Lung dose increased with the ill breath pattern in proportion to increased minute ventilation (7.6%). **Conclusion:** 1) The AE system provided higher lung dose than other nebulizers with standard masks. 2) Mask design and fit can substantially impact nebulizer performance. 3) PN performed better than BB, but not as good as a mask. If crying decreases lung dose by 75%, we speculate that the PN and BB (non-crying) may improve lung dose vs mask with a crying infant. 4) An increase in lung dose may occur in ill infants if minute ventilation is elevated. *This abstract is funded by: AstraZeneca.*

Pari LC-D[™], PARI[™] Respiratory Equipment

1. **PERFORMANCE COMPARISON OF NEBULIZER DESIGNS: CONSTANT-OUTPUT, BREATH-ENHANCED, AND DOSIMETRIC.** JL Rau, A Ari and RD Restrepo. Respiratory Care 2004; 49(2): 174-179.

Introduction: Design differences among pneumatically powered, small-volume nebulizers affect drug disposition (percentage of the dose delivered to the patient, lost to deposition in the equipment, and lost via exhalation to ambient air) and thus affect drug availability and efficacy. **Objective:** Evaluate *in vitro* the dose disposition with 5 nebulizer models, of 3 types (constant-output, breath-enhanced, and dosimetric), using simulated normal, adult breathing. **Methods:** We compared 5 nebulizer models: 2 constant-output (Misty-Neb and SideStream), 1 breath-enhanced (Pari LCD), and 2 dosimetric (Circulaire and **AeroEclipse**^{*}). Each nebulizer was filled with a 3-mL unit-dose of albuterol sulfate and powered by oxygen at 8 L/min. The nebulizers were connected to an induction throat, connected to a breathing simulator. We measured (1) inhaled drug (subdivided into mass deposited in the induction throat and mass deposited in the filter at the distal end of the induction throat), (2) exhaled drug (lost to ambient air), (3) drug lost to deposition in the apparatus, and (4) drug left in the unit-dose bottle. The duration of nebulization (until sputter) was measured with a stopwatch. All drug amounts were analyzed via spectrophotometry and expressed as a percentage of the total dose. **Results:** The mean \pm SD inhaled drug percentages were: Misty-Neb 17.2 \pm 0.4%, SideStream 15.8 \pm 2.8%, Pari LCD 15.2 \pm 4.2%, Circulaire 8.7 \pm 1.0%, and **AeroEclipse**^{*} 38.7 \pm 1.3%. The mean \pm SD percentages of drug lost to ambient air were: Misty-Neb 26.8 \pm 0.7%, SideStream 17.3 \pm 0.4%, Pari LCD 18.3 \pm 0.8%, Circulaire 12.3 \pm 0.8%, and **AeroEclipse**^{*} 6.6 \pm 3.3%. The mean \pm SD percentages of drug lost to deposition in the apparatus were: Misty-Neb 52.3 \pm 0.6%, SideStream 63.4 \pm 3.0%, Pari LCD 62.5 \pm 4.0%, Circulaire 75.8 \pm 0.5%, and **AeroEclipse**^{*} 51.0 \pm 2.1%. Duration of nebulization was shortest with the Circulaire and longest with the **AeroEclipse**^{*} (p < 0.05 via 1-way analysis of variance). **Conclusions:** The nebulizers we tested differ significantly in overall drug disposition. The dosimetric **AeroEclipse**^{*} provided the largest inhaled drug mass and the lowest loss to ambient air, with the test conditions we used.

2. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA A BREATH-ACTUATED SMALL VOLUME NEBULIZER (SVN): A COMPARATIVE *IN-VITRO* ASSESSMENT.** MW Nagel, KJ Wiersema, SL Bates and JP Mitchell. American Journal of Respiratory Care & Critical Care Medicine, 2001; 163(5); A442.

Rationale: To compare the delivery of budesonide suspension in terms of fine particle dose ($< 4.7 \mu\text{m}$ aerodynamic diameter (FPD)) from a breath-actuated (BA) SVN with that from a continuous flow air entrainment (AE) SVN.

Methods: FPD values were determined for 5-**AeroEclipse*** BA SVNs (Monaghan Medical Corp., Plattsburgh, N.Y.) and 5-LC-D™ AE SVNs (PARI Respiratory Equipment, Inc., Monterey, CA), nebulizing 4ml of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). Each SVN was operated with air at 50 psig, 8 l/min until sputtering occurred. Breathing parameters were: tidal volume = 600 ml, I:E = 1:2 rate = 10/min. FPD was determined by cascade impactor at $28.3 \pm 0.5 \text{ l/min}$. **Results:** From the beginning of nebulization until sputtering, the **AeroEclipse*** and the LC-D™ SVNs produced 164 ± 3 and $71 \pm 4 \mu\text{g}$ FPD of budesonide respectively. During the first 5 minutes (after which time the LC-D™ sputtered), values of FPD for the **AeroEclipse*** and the LC-D™ SVNs were 76 ± 4 and $71 \pm 4 \mu\text{g}$ budesonide respectively. **Conclusion:** The **AeroEclipse*** was more efficient than the LCD™ SVN for this suspension formulation [Mann-Whitney rank sum test, $p < 0.001$]. Almost no medication delivery took place from the **AeroEclipse*** SVN during the exhalation portion of the breathing cycle, thereby providing important benefits to both patient and care giver.

3. **COMPARISON OF BREATH-ACTUATED JET NEBULIZER (BAN) IN 'CONTINUOUS DELIVERY' MODE WITH OTHER CONTINUOUS DELIVERY NEBULIZERS.** JP Mitchell, KJ Wiersema, CC Doyle and MW Nagel. Respiratory Care 2003, 48 (11): S1077.

The **AeroEclipse*** BAN (Monaghan Medical Corp., Plattsburgh, N.Y.) has been equipped with an optional blue cap whose purpose is to retain the actuator piston in the position it would occupy during inhalation in breath-actuated mode, so that the nebulizer operates continuously. The present study compared the delivery of a bronchodilator from diluted albuterol sulfate respirator solutions (3-ml of 0.83 and 1-ml of 2.5 mg/ml albuterol in physiologically normal saline (0.9% w/v NaCl)), via this nebulizer, the Micromist® (Hudson RCI, Temecula, CA), Misty-Neb™ (Allegiance Healthcare Corp., McGaw Park, IL) and the LCD™ (PARI Respiratory Equipment, Monterey, CA). Each nebulizer was tested using a breathing simulator set to the following parameters representative of adult use: tidal volume = 600-ml, rate = 10 breaths/min, inspiratory/expiratory ratio 1:2. The total mass of albuterol (TM) delivered to the first sputter was determined by filter collection at the mouthpiece of the nebulizer operated with compressed air supplied at 50 psig at 8 L/min ($n = 5$ devices/group, 3 replicates/device). The fraction of the aerosol contained in droplets finer than $4.8 \mu\text{m}$ aerodynamic diameter (FPF) was determined by laser diffractometry in a parallel study, so that the fine droplet mass (FM) could be calculated as the product of TM and FPF. Values of FM (mean \pm SD) and time to deliver medication (T_{med}) were as follows:

| Solution(mg/ml) | AeroEclipse* | | LCD™ | | Micromist® | | MistyNeb™ | |
|------------------------|---------------------|-----------------|-----------------|----------------|-----------------|-------------------|---------------|---------------|
| | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 | 0.83 | 2.50 |
| FM (μg) | 360 ± 22 | 263 ± 26 | 149 ± 16 | 108 ± 4 | 209 ± 12 | 15.4 ± 5.9 | 82 ± 9 | 31 ± 5 |
| T_{med} (min) | 3 | <1 | 2 | <1 | 7 | <1 | 4 | <1 |

The **AeroEclipse*** nebulizer delivered significantly more FM in continuous delivery mode than the other nebulizers when operated in continuous mode with either solution strength (1-way repeated measures ANOVA, $p < 0.05$). T_{med} from the **AeroEclipse*** nebulizer was comparable with the best performing continuous nebulizer (LCD™).

4. **PERFORMANCE OF A NEW BREATH-ACTUATED SMALL VOLUME NEBULIZER (SVN) FOR THE DELIVERY OF A COMBINATION ANTICHOLINERGIC/BRONCHODILATOR.** MW Nagel, KJ Wiersema, SL Bates and JP Mitchell. American Journal of Respiratory and Critical Care Medicine, 2001, 163 (5): A443.

Purpose: To compare the delivery of ipratropium bromide (IPR) and albuterol sulfate (ALB) as fine droplets ($< 4.8 \mu\text{m}$ diameter (FPD)) and as total emitted dose (ED) from a breath-actuated (BA- SVN) with that from a continuous flow air entrainment (AE-SVN) after 5-minutes of operation. **Methods:** FPD and ED were determined for 5-**AeroEclipse*** BAN (Monaghan Medical Corp., N.Y.) and 5-PARI LCD™ SVNs (PARI Respiratory Equipment, Inc., CA) nebulizing Combivent® (2.5-ml, 0.2 mg/ml IPR and 1.0 mg/ml ALB; Boehringer-Ingelheim (Canada) Inc.). Each SVN was operated with 8 l/min air at 50 psig, simulating breathing at tidal volume, I:E ratio and rate of 750-ml, 1:2 and 10/min respectively. Droplet size distributions were measured by laser diffractometer.

Results: (ED) and (FPD) were as follows:

| | | | |
|-----|-------------------------|-------------------------------|--------------------------------|
| IPR | AeroEclipse* BAN | ED = $102 \pm 7 \mu\text{g}$ | FPD = $82 \pm 6 \mu\text{g}$ |
| IPR | PARI LCD™ SVNs | ED = $55 \pm 7 \mu\text{g}$ | FPD = $45 \pm 5 \mu\text{g}$ |
| ALB | AeroEclipse* BAN | ED = $581 \pm 17 \mu\text{g}$ | FPD = $471 \pm 14 \mu\text{g}$ |
| ALB | PARI LCD™ SVNs | ED = $279 \pm 33 \mu\text{g}$ | FPD = $226 \pm 26 \mu\text{g}$ |

Differences in ED and FPD between SVNs for IPR and ALB components were statistically significant (unpaired t-test for each variable, $p < 0.001$). Mass median aerodynamic diameters were close to $2.8 \mu\text{m}$ for both SVN groups. **Conclusion:** The **AeroEclipse*** BAN is significantly more efficient for the delivery of this combination anticholinergic/bronchodilator than a conventional AE-SVN.

5. **DELIVERY OF A SUSPENSION CORTICOSTEROID FORMULATION BY SMALL VOLUME NEBULIZERS: A COMPARATIVE BENCH STUDY.** JP Mitchell, MW Nagel, KJ Wiersema and SL Bates. Presented at ERS Annual Congress, Berlin, Germany, September 2001.

We report a study of the delivery of 0.25% mg/ml budesonide suspension (Pulmicort[®], Nebuamp[®] (2 x 2-ml), Astra-Zeneca, Canada) by two types of small volume nebulizer (SVN), simulating adult breathing conditions ((tidal volume = 600-ml, duty cycle = 1:2 (2-s inspiration), PIFR = 31 l/min). Each SVN was operated by compressed air (8 l/min at 50 psig). Budesonide mass delivery was determined by filter collection (n = 5 SVNs/group, 3-replicates/device). The **AeroEclipse**[™] BANs (Trudell Medical International, London Canada) delivered 283 ± 32 mg prior to sputtering, and 80 ± 11 mg were lost to the environment. Corresponding data for the LCD[™] SVNs (Pari Respiratory Equipment Inc., Richmond, VA, USA) were 97 ± 7 mg and 305 ± 2 mg respectively. The breath-actuation feature of the **AeroEclipse**[™] SVN minimizes aerosol release to the environment during exhalation, which may cause adverse effects to both patient and health care provider.


| RESULTS | | |
|-------------------------------------|-----------|----------|
| Nebulizer | FILT (µg) | ENV (µg) |
| AeroEclipse [™] BAN | 283 ± 33 | 80 ± 11 |
| LCD [™] | 97 ± 7 | 305 ± 2 |

mean ± SD

6. **EFFECT OF NEBULIZER DESIGN ON FINE PARTICLE MASS.** D Hess, JP Mitchell, D Coppolo, MW Nagel, AD Archer, R Blacker. Presented at Open Forum, Ann. Meet. of the American Association For Respiratory Care (AARC), Las Vegas, 1999. Published in Resp Care, Oct 1999, 44(10): 1289.

Background: Nebulizer design is known to affect performance. In this study, we compared fine particle mass from nebulizers of four designs. **Methods:** We tested traditional disposable nebulizers (Baxter Misty-Neb[™], Hudson Updraft-II Neb-U-Mist[™]), breath-enhanced nebulizers (Pari-LC-D[™]), nebulizers with collection bags (Westmed Circulaire[™]), and a Trudell **AeroEclipse**[™] (with breath actuation disabled). Five of each device with three replicates (n = 15) were tested using an *in-vitro* model of spontaneous breathing. A rigid bar was placed between the two compartments of a test lung (Michigan Instruments TTL). The drive lung was attached to a ventilator (Infrasonics Infant Star[™]) to simulate spontaneous breathing (tidal volume 0.6 L, rate 10/min, T_i 2 s). A bacterial/viral filter (Trudell MT3000) was placed between the nebulizer and slave lung. Flow was monitored between the test lung and filter (Novamatrix Ventcheck[™]). Albuterol solution (0.625 mg/mL) was placed into the nebulizers (4 mL), which were powered with air (8 L/min). Filters were replaced at one minute intervals (flow to the nebulizer was discontinued during filter replacement) until sputtering occurred. The filter was washed with methanol (20 mL) and albuterol concentration was measured with HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer. Fine particle mass was calculated as the product of mass % <4.7 µm and total nebulizer output.

Results: Fine particle mass from the **AeroEclipse**[™] nebulizer was greater than that from the other nebulizers (P<0.001) (see figure). **Conclusions:** Fine particle mass was affected by nebulizer design. The clinical relevance of this finding awaits further investigation. Further evaluation of the breath-actuated feature of the **AeroEclipse**[™] is warranted.



7. **THE EFFECT OF SMALL VOLUME NEBULIZER (SVN) DESIGN ON FINE PARTICLE MASS DELIVERY OF A BRONCHODILATOR.** R Blacker, RW Morton, JP Mitchell, MW Nagel and DR Hess. Drug Delivery to the Lungs-X, London, UK, 1998, J. Aerosol Med., 13(1), 65.

Fine particle mass delivery was compared from six different SVNs, including continuous un-enhanced flow designs (Hudson Updraft-II Neb-U-Mist[™]), breath-enhanced nebulizers (Pari-LC-D[™], Medic-Aid Sidestream[®]), nebulizers with

aerosol collection bag (Westmed Circulaire™), and an **AeroEclipse*** with breath actuation disabled (Trudell Medical International). Five of each type of SVN were tested operating with air (8 l/min, 50 psig), using an *in-vitro* model that simulated spontaneous breathing by an adult (tidal volume 0.6 l, rate 10/min, TI = 2 s). A bacterial/viral filter was placed between the nebulizer and breathing simulator. In each case, salbutamol sulphate (Ventolin®) respirator solution (0.625 mg/ml, 4 ml) was placed into the reservoir of the SVN. The filters were replaced at one-minute intervals until sputtering occurred. The salbutamol collected on the filter was assayed by HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer laser diffractometer. Fine particle mass delivery rates varied significantly from each of the SVN's from more than 110 µg/min (**AeroEclipse***) to ca. 20 µg/min (Circulaire™).

Sidestream™, MedicAid™

1. **PERFORMANCE COMPARISON OF NEBULIZER DESIGNS: CONSTANT-OUTPUT, BREATH-ENHANCED, AND DOSIMETRIC.** JL Rau, A Ari and RD Restrepo. *Respiratory Care* 2004; 49(2): 174-179.

Introduction: Design differences among pneumatically powered, small-volume nebulizers affect drug disposition (percentage of the dose delivered to the patient, lost to deposition in the equipment, and lost via exhalation to ambient air) and thus affect drug availability and efficacy. **Objective:** Evaluate *in vitro* the dose disposition with 5 nebulizer models, of 3 types (constant-output, breath-enhanced, and dosimetric), using simulated normal, adult breathing. **Methods:** We compared 5 nebulizer models: 2 constant-output (Misty-Neb and SideStream), 1 breath-enhanced (Pari LCD), and 2 dosimetric (Circulaire and **AeroEclipse***). Each nebulizer was filled with a 3-mL unit-dose of albuterol sulfate and powered by oxygen at 8 L/min. The nebulizers were connected to an induction throat, connected to a breathing simulator. We measured (1) inhaled drug (subdivided into mass deposited in the induction throat and mass deposited in the filter at the distal end of the induction throat), (2) exhaled drug (lost to ambient air), (3) drug lost to deposition in the apparatus, and (4) drug left in the unit-dose bottle. The duration of nebulization (until sputter) was measured with a stopwatch. All drug amounts were analyzed via spectrophotometry and expressed as a percentage of the total dose. **Results:** The mean ± SD inhaled drug percentages were: Misty-Neb 17.2 ± 0.4%, SideStream 15.8 ± 2.8%, Pari LCD 15.2 ± 4.2%, Circulaire 8.7 ± 1.0%, and **AeroEclipse*** 38.7 ± 1.3%. The mean ± SD percentages of drug lost to ambient air were: Misty-Neb 26.8 ± 0.7%, SideStream 17.3 ± 0.4%, Pari LCD 18.3 ± 0.8%, Circulaire 12.3 ± 0.8%, and **AeroEclipse*** 6.6 ± 3.3%. The mean ± SD percentages of drug lost to deposition in the apparatus were: Misty-Neb 52.3 ± 0.6%, SideStream 63.4 ± 3.0%, Pari LCD 62.5 ± 4.0%, Circulaire 75.8 ± 0.5%, and **AeroEclipse*** 51.0 ± 2.1%. Duration of nebulization was shortest with the Circulaire and longest with the **AeroEclipse*** ($p < 0.05$ via 1-way analysis of variance). **Conclusions:** The nebulizers we tested differ significantly in overall drug disposition. The dosimetric **AeroEclipse*** provided the largest inhaled drug mass and the lowest loss to ambient air, with the test conditions we used.

2. **PERFORMANCE OF A NEW BREATH-ACTUATED SMALL VOLUME NEBULIZER (BA-SVN) WHEN USED WITH OXYGEN AS A DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** AM Verdun, JP Mitchell and MW Nagel. Presented at ALA/ATS International Conference, Chicago, 1998.

Purpose: To compare the delivery of saline (0.9% w/v NaCl) by a new BA-SVN (Trudell Medical Int..) with that from two other representative SVN's (LC-JET™ (PARI Respiratory Products Inc., Canada) and reusable Sidestream™ (MedicAid, UK)) using oxygen delivered at 50 psig at 8 l/min to simulate hospital use. **Methods:** 5 BA-SVN's were tested using a laser diffractometer (Malvern Mastersizer-X) to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. 5, LC-JET™ and 5, Sidestream™ SVN's were also tested similarly. The BA-SVN was operated with manual over-ride engaged (continuous delivery of aerosol). **Results:** Total (TM) and respirable (RM), droplets finer than 4.8 µm diameter) mass output rates and droplet mass median diameter (MMD) were as follows: BA-SVN: TM = 672 ± 23 µg/min, RM = 545 ± 31 µg (80.9 ± 2.4% respirable), MMD = 2.79 ± 0.15 µm; LC-JET™: TM = 675 ± 69 µg/min, RM = 449 ± 41 µg/min (66.7 ± 1.8% respirable), MMD = 3.39 ± 0.08 µm; Sidestream™: TM = 442 ± 26 µg/min, RM = 358 ± 38 µg/min (80.8 ± 4.2% respirable), MMD = 2.94 ± 0.03 µm. **Conclusion:** Although TM from the new BA-SVN was comparable with that from the LC-JET™ (Mann-Whitney rank sum test, $p = 0.84$), the finer MMAD produced from the BA-SVN resulted in a significantly greater RM ($p < 0.001$). Both TM and RM from the BA-SVN were greater than those from the Sidestream™ SVN ($p < 0.001$).

3. **COMPARISON OF A NEW AIR ENTRAINMENT SMALL VOLUME NEBULIZER (AE-SVN) WITH OTHER SVN'S WHEN USED WITH OXYGEN AS DRIVING GAS UNDER CONDITIONS OF HOSPITAL USE.** JP Mitchell and MW Nagel. Presented at Ann. Meet. Amer. Assoc. of Asthma, Allergy and Immunology (AAAAI), Washington D.C., 1998.

The performance of a prototype novel AE-SVN (Trudell Medical International (n = 5)) with normal saline (0.9% w/v NaCl) operating at 20 ± 2°C, 50 ± 10% RH, has been evaluated with oxygen (50 psig, 8 l/min) as driving gas to simulate hospital use. Comparison testing was also undertaken with two other representative AE-SVN's, (a) LC-JET™ (Pari Respiratory Equipment Inc.), without inspiratory valve cap which would otherwise restrict aerosol output, (b) SideStream™ (MedicAid, UK). A laser diffractometer (Malvern Mastersizer-X) was used to determine the size distribution of droplets emitted at the mouthpiece. The total mass output was determined gravimetrically in a parallel series of tests. Total (T) and respirable ((R), droplets finer than 4.8 µm aerodynamic diameter) mass output rates and droplet mass median aerodynamic diameter (MMAD) for the new AE-SVN (5 replicate measurements/device) were: 671 ± 26 µg/min (T), 542 ± 23 µg/min (R) and 2.88 ± 0.09 µm (MMAD). Corresponding data for the LC-JET™ were: 675 ± 65 µg/min (T), 450 ± 45 µg/min (R) and 3.39 ± 0.14 µm (MMAD), and for the SideStream™ were: 442 ± 27 µg/min (T), 357 ± 28 µg/min (R) and 2.95 ± 0.13 µm (MMAD). The total aerosol delivery rate from the new AE-SVN matched that of the LC-JET™ (un-paired t-test, $p = 0.79$).

and exceeded that from the SideStream™ ($p < 0.001$). The finer MMAD of the aerosol provided by the new AE-SVN resulted in a significantly greater respirable mass fraction, increasing the respirable mass delivery rate compared with the other SVNs ($p < 0.001$).

4. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS: A COMPARATIVE *IN-VITRO* ASSESSMENT.** JP Mitchell, MW Nagel and AD Archer. Presented at Proc. Ann. Meet. Amer. College of Chest Physicians (ACCP), Toronto, November 1998, in *Chest*, 114 (4S), 295, (1998); and at the World Asthma Meeting, Barcelona, December 1998, in *Eur. Resp. J.*, 12S29, 7, (1998).

Purpose: To compare the performance of a new air entrainment small volume nebulizer (AE-SVN, Trudell Medical International with other widely used SVNs (LC-Star™ (PARI Respiratory Equipment), Updraft™ Neb-U-Mist™ (Hudson Oxygen Therapy Sales Co.), Circulaire™ (Westmed), Sidestream™ (Medic-Aid), Airlife™ Misty-Neb™ (Baxter Healthcare Corp.)) for the delivery of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). **Methods:** Each SVN ($n = 5$ devices for each group, 3 replicates per device) was operated with compressed air at 50 psig at a flow rate of 8 l/min. The total mass of budesonide nebulized from 2 x 2 ml ampoules was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The SVN was operated until it spluttered, was then tapped gently to dislodge droplets back to the reservoir. Nebulization was deemed complete 20 seconds later. The mass of budesonide collected was determined by HPLC-UV spectrophotometry. **Results:** The delivery rate ((mean \pm 1 S.D) μ g budesonide/min) from the AE-SVN (102 ± 9) was significantly greater than with the other groups: (LC-Star™ (91 ± 6), Misty-Neb™ (49 ± 2), Sidestream™ (46 ± 4), Circulaire™ (26 ± 4) and Neb-U-Mist™ (25 ± 6)), (1-way ANOVA, $p < 0.02$). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with LC-Star™ (229 ± 10 s), Sidestream™ (365 ± 19 s), Circulaire™ (420 ± 84 s), Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). **Conclusions:** The new AE-SVN is highly efficient at entraining the budesonide particles into the liquid droplets at these conditions. **Clinical Implications:** The good delivery rate combined with comparatively short duration of delivery offers the potential for rapid treatment and patient convenience.

5. **THE EFFECT OF SMALL VOLUME NEBULIZER (SVN) DESIGN ON FINE PARTICLE MASS DELIVERY OF A BRONCHODILATOR.** R Blacker, RW Morton, JP Mitchell, MW Nagel and DR Hess. Drug Delivery to the Lungs-X, London, UK, 1998, *J. Aerosol Med.*, 13(1), 65.

Fine particle mass delivery was compared from six different SVNs, including continuous un-enhanced flow designs (Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™, Medic-Aid Sidestream®), nebulizers with aerosol collection bag (Westmed Circulaire™), and an **AeroEclipse*** with breath actuation disabled (Trudell Medical International). Five of each type of SVN were tested operating with air (8 l/min, 50 psig), using an *in-vitro* model that simulated spontaneous breathing by an adult (tidal volume 0.6 l, rate 10/min, TI = 2 s). A bacterial/viral filter was placed between the nebulizer and breathing simulator. In each case, salbutamol sulphate (Ventolin®) respirator solution (0.625 mg/ml, 4 ml) was placed into the reservoir of the SVN. The filters were replaced at one-minute intervals until sputtering occurred. The salbutamol collected on the filter was assayed by HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer laser diffractometer. Fine particle mass delivery rates varied significantly from each of the SVNs from more than 110 μ g/min (**AeroEclipse***) to ca. 20 μ g/min (Circulaire™).

Circulaire™, Westmed™ Corp.

1. **PERFORMANCE COMPARISON OF NEBULIZER DESIGNS: CONSTANT-OUTPUT, BREATH-ENHANCED, AND DOSIMETRIC.** JL Rau, A Ari and RD Restrepo. *Respiratory Care* 2004; 49(2): 174-179.

Introduction: Design differences among pneumatically powered, small-volume nebulizers affect drug disposition (percentage of the dose delivered to the patient, lost to deposition in the equipment, and lost via exhalation to ambient air) and thus affect drug availability and efficacy. **Objective:** Evaluate *in vitro* the dose disposition with 5 nebulizer models, of 3 types (constant-output, breath-enhanced, and dosimetric), using simulated normal, adult breathing. **Methods:** We compared 5 nebulizer models: 2 constant-output (Misty-Neb and SideStream), 1 breath-enhanced (Pari LCD), and 2 dosimetric (Circulaire and **AeroEclipse***). Each nebulizer was filled with a 3-mL unit-dose of albuterol sulfate and powered by oxygen at 8 l/min. The nebulizers were connected to an induction throat, connected to a breathing simulator. We measured (1) inhaled drug (subdivided into mass deposited in the induction throat and mass deposited in the filter at the distal end of the induction throat), (2) exhaled drug (lost to ambient air), (3) drug lost to deposition in the apparatus, and (4) drug left in the unit-dose bottle. The duration of nebulization (until sputter) was measured with a stopwatch. All drug amounts were analyzed via spectrophotometry and expressed as a percentage of the total dose. **Results:** The mean \pm SD inhaled drug percentages were: Misty-Neb $17.2 \pm 0.4\%$, SideStream $15.8 \pm 2.8\%$, Pari LCD $15.2 \pm 4.2\%$, Circulaire $8.7 \pm 1.0\%$, and **AeroEclipse*** $38.7 \pm 1.3\%$. The mean \pm SD percentages of drug lost to ambient air were: Misty-Neb $26.8 \pm 0.7\%$, SideStream $17.3 \pm 0.4\%$, Pari LCD $18.3 \pm 0.8\%$, Circulaire $12.3 \pm 0.8\%$, and **AeroEclipse*** $6.6 \pm 3.3\%$. The mean \pm SD percentages of drug lost to deposition in the apparatus were: Misty-Neb $52.3 \pm 0.6\%$, SideStream $63.4 \pm 3.0\%$, Pari LCD $62.5 \pm 4.0\%$, Circulaire $75.8 \pm 0.5\%$, and **AeroEclipse*** $51.0 \pm 2.1\%$. Duration of nebulization was shortest with the Circulaire and longest with the **AeroEclipse*** ($p < 0.05$ via 1-way analysis of variance). **Conclusions:** The nebulizers we tested differ significantly in overall drug disposition. The dosimetric **AeroEclipse*** provided the largest inhaled drug mass and the lowest loss to ambient air, with the test conditions we used.

2. **BREATH-ACTUATED VS RESERVOIR NEBULIZERS FOR UNDILUTED ALBUTEROL.** D Geller, B Kesser. Presented at the 13th International Congress on Aerosols in Medicine, Interlaken, Switzerland, September 17-21, 2001.

Aim: Some Emergency Departments use undiluted albuterol in nebulizers designed to conserve drug during exhalation. We compared the *in vitro* performance of 4 devices to estimate which would be most effective clinically: **AeroEclipse*** Breath-Actuated Nebulizer ("BAN"); Circulaire® (C) and AeroTee™ (AT) which use a 750 ml reservoir bag to conserve drug during exhalation; and Salter HDN™ (S) with a 50 ml tower reservoir. **Method:** We studied 4 units of each nebulizer type in duplicate, using a Pari Proneb Turbo compressor. Nebulizers were filled with undiluted 0.5% albuterol, 1 ml (5 mg) or 2 ml (10 mg). Particle size distributions were measured by laser diffraction (Malvern SprayTec). Drug output (1 minute after "sputter") was captured on a filter between the device mouthpiece and a Pari breath-simulator, which used a recorded waveform from a 9 yr old male. Albuterol was measured by spectrophotometry, and fine particle dose (FPD) (mg of drug < 5 mm in size) was calculated.

Results:

| Neb | MMAD | FPD (1cc) | Minutes | FPD (2cc) | Minutes |
|-----|------|-----------|---------|-----------|---------|
| AE | 3.9 | 0.60 | 3.8 | 2.41 | 11.0 |
| AT | 4.8 | 0.03 | 2.0 | 0.62 | 3.2 |
| C | 2.5 | 0.09 | 2.0 | 0.65 | 3.7 |
| S | 8.5 | 0.08 | 2.0 | 0.57 | 3.7 |

Conclusions: The AE was superior to the reservoir-type nebulizers in fine-particle output for each fill volume. The AT and C had large dead volumes, and the S produced larger particles. These shortcomings were overcome with larger nominal doses. Each nebulizer produced 0.6-mg FPD of albuterol over 3½ minutes, but the AE required only half the starting dose. Albuterol 0.6 mg is a reasonable clinical respirable dose in a child with acute asthma. These findings must be taken into account when designing clinical treatment protocols for acute asthma. **Background:** Many nebulizers are designed to decrease the amount of drug that is lost during exhalation. The Circulaire® (Westmed) and AeroTee™ (Hudson) incorporate a 750 ml bag on the expiratory side of the nebulizer that collects aerosol while the patient exhales, making it available for inhalation on the next breathing cycle. The Salter HDN™ (Salter) has a 50 ml tower that acts as a reservoir. The **AeroEclipse*** BAN (Trudell/ Monaghan) has a spring mechanism that allows generation of aerosol during inhalation only, so no drug waste occurs during exhalation. We recently reported the aerosol characteristics with these devices nebulizing unit-dose albuterol sulfate (2.5 mg/3 ml).¹ Delivery time with unit-dose (0.083%) albuterol can be long, which may increase personnel costs. To maintain lung-dose delivery and minimize the treatment time, some hospitals use drug-conserving nebulizers with small fill-volumes of undiluted (0.5%) albuterol for patients presenting with acute bronchospasm. We measured the particle size distributions and used a child's breathing pattern to compare albuterol output of these 4 drug-conserving nebulizers, using unit-dose albuterol 2.5 mg (3ml), 0.5% albuterol 5 mg (1ml) and 10 mg (2ml) nominal doses. We calculated the fine particle dose and measured the dose of drug remaining within the nebulizer and all attachments to determine the residual dose. For reference, we compared these results to those of a T-piece (Hudson Micromist) nebulizer using unit-dose albuterol to simulate conventional dosing.

Materials and Methods:

- Drug:
 - Albuterol Sulfate 0.083% unit-dose (2.5mg/3ml)
 - Albuterol Sulfate 0.5% (5mg/ml) 1 and 2cc fill volumes
- Nebulizers: Circulaire® (Model 0260), AeroTee™ with Micromist Nebulizer (Model 1002), Salter HDN™ (Model 8960), and **AeroEclipse*** BAN (Figure 1)
- Compressor: PARI PRONEB TURBO
- 4 nebulizers of each type studied in duplicate
- Particle size by laser diffraction (Malvern Insitac)
- Breathing pattern from 9 year old male volunteer, using the PARI breath simulator (RR 19 bpm, Vt 421 cc, Ti 1.3 seconds)
- Definitions:
 - Inspired dose = drug on inspiratory filter
 - Residual dose = drug collected from nebulizer and accessory components after completion of nebulization.
 - Fine particle dose (FPD) = (Inspired dose) x (% of particles <5 µm) Figure 1.
 - Duration = time (minutes) from the beginning of nebulization to 1 minute past the onset of sputter.
- Samples assayed with spectrophotometer at 228 λ.

Results:

| | | AeroEclipse* BAN | AeroTee™ | Circulaire® | Salter HDN™ |
|----------------------|-------------------------|-------------------------|----------|-------------|-------------|
| Particle Sizing | MMD | 3.87 | 4.80 | 2.47 | 8.46 |
| | GSD | 2.3 | 2.0 | 2.1 | 2.0 |
| | % < 5 µm | 61.7% | 52.9% | 83.6% | 30.0% |
| 2.5 mg Unit Dose* | Duration (minutes) | 14.7 | 7.2 | 7.0 | 3.6 |
| | Inspired Dose (mg) | 0.77 | 0.37 | 0.14 | 0.30 |
| | Residual Dose (mg) | 1.5 | 1.8 | 2.1 | 1.9 |
| | Fine Particle Dose (mg) | 0.52 | 0.19 | 0.12 | 0.10 |
| | | | | | |

| | | | | | |
|-------------------------|-------------------------|------|------|------|------|
| 5 mg (1 ml) Dose | Duration (minutes) | 3.8 | 2.0 | 2.0 | 2.0 |
| | Inspired Dose (mg) | 0.97 | 0.06 | 0.11 | 0.28 |
| | Residual Dose (mg) | 3.5 | 4.9 | 4.6 | 4.4 |
| | Fine Particle Dose (mg) | 0.60 | 0.03 | 0.09 | 0.08 |
| 10 mg (2 ml) Dose | Duration (minutes) | 11.0 | 3.2 | 3.7 | 3.7 |
| | Inspired Dose (mg) | 3.9 | 1.2 | 0.8 | 1.9 |
| | Residual Dose (mg) | 5.8 | 8.7 | 8.6 | 6.9 |
| | Fine Particle Dose (mg) | 2.40 | 0.60 | 0.60 | 0.60 |

* Unit-dose data presented at ATS 2001¹

For comparison, the Hudson Micromist conventional T-Piece Nebulizer (with Unit-Dose 2.5 mg Albuterol) produced a fine-particle dose of 0.14 mg in 7.0 minutes.

Discussion:

- **AeroEclipse*** BAN had highest FPD with all nominal doses:
 - FPD was 2.7 to 5.2 times higher with unit-dose; 6.7 to 20 times higher with 5 mg dose; 4 times higher with 10 mg dose.
 - Lowest residual dose
 - Higher fine particle fraction except for Circulaire®
- **Nebulizer Inefficiencies:**
 - AeroTee™ and Circulaire® had high residual doses in part due to valves and collection bags.
 - Salter HDN™ produces larger particles
 - These inefficiencies were partially compensated for by increasing nominal dose to 10mg (2 ml)
- **Duration of Nebulization:**
 - **AeroEclipse*** BAN had longer delivery time because it is breath actuated; no waste during exhalation
 - Using 0.5% albuterol, all nebulizers produced 0.6 mg fine-particle dose in < 4 minutes, but the **AeroEclipse*** BAN only required half the nominal dose to accomplish this.
- **Comparison to Unit-Dose 2.5 mg:**
 - **AeroEclipse*** BAN produced comparable FPD with unit-dose and 5 mg (1 ml) nominal dose, but delivery time was less than a third with undiluted drug.
- **Comparison to conventional nebulizers:**
 - The FPD with the Hudson and unit-dose drug was 0.14 mg, similar to the reservoir-type nebulizers with unit-dose.
 - The higher FPD with **AeroEclipse*** BAN (all doses) and the reservoir nebs (10 mg dose) may result in better and longer lasting bronchodilation than the Hudson with conventional dosing, thus reducing number of treatments, therapist time, and total costs.

Conclusion:

- **AeroEclipse*** BAN was superior to the reservoir-type nebulizers at all nominal doses.
- **AeroEclipse*** BAN has the additional advantage of being a dosimetric device, i.e. it will not operate or waste drug while the patient is coughing or resting. The patient and health care providers get visual feedback of adequate inspiratory effort necessary to actuate the nebulizer.
- Use of undiluted 0.5% albuterol may result in higher lung doses in a shorter amount of time. These results can be used as a guide when developing bronchodilator protocols for the hospital or E.D. setting.

Funded by the Nemours Foundation

¹ Geller D, Kesser B. Am J Respir Crit Care Med 2001; 163:A444. Journal of Aerosol Medicine 14 (3) 2001; 395:P1-41

3. COMPARISON OF DRUG OUTPUT FROM 4 DIFFERENT RESERVOIR TYPE NEBULIZERS. DE Geller, B Kesser. American Journal of Respiratory Care & Critical Care Medicine, 2001; 163(5); A444.

Rationale: Many nebulizers currently being marketed utilize different techniques to conserve drug that would normally be lost during exhalation. The Circulaire™ and Aero Tee™ nebulizers use a 750 cc reservoir bag to accumulate nebulized drug, while the Salter HDN™ uses a 50ml tower to serve as a reservoir. The **AeroEclipse*** nebulizer uses breath actuated nebulization to deliver drug only during inspiration. We evaluated all 4 nebulizers using a recorded pediatric breathing pattern to measure total drug output. We additionally measured the particle size characteristics of each type with the laser diffraction technique. **Methods:** 4 nebulizers of each type were studied in duplicate for sizing and total output characteristics. Each nebulizer was charged with a unit dose of 2.5 mgs albuterol sulfate in 3cc's. Sizing studies were averaged values preformed over 5 minute runs on each nebulizer with a Malvern Spray Tec laser. Drug output was as calculated as the assayed amount of albuterol collected on a filter distal to the mouthpiece of the nebulizer. Simulated breathing was performed through the nebulizer by a Pari breath simulator from waveforms originally recorded from a healthy 9-year-old male.


Results:

| | Inspired dose | %>1 & <5M | Respirable Dose | Residual Dose |
|---------------------|----------------|------------|-----------------|---------------|
| AeroEclipse* | 0.64 ± 0.06 mg | 52.7 ± 2.0 | 0.34 ± 0.03 mg | 1.27 ± 0.09 |
| Aero Tee | 0.31 ± 0.09 mg | 41.2 ± 7 | 0.13 ± 0.04 mg | 1.51 ± 0.11 |
| Circulaire | 0.12 ± 0.03 mg | 61.9 ± 1 | 0.07 ± 0.02 mg | 1.72 ± 0.13 |
| Salter HDN | 0.25 ± 0.05 mg | 24.7 ± 5 | 0.06 ± 0.02 mg | 1.59 ± 0.10 |

Conclusions: The **AeroEclipse*** delivers a greater total dose of drug as well as a greater amount of drug in the fine particle range, most likely to deposit in the lower airways.

4. **EFFECT OF NEBULIZER DESIGN ON FINE PARTICLE MASS.** D Hess, JP Mitchell, D Coppolo, MW Nagel, AD Archer, R Blacker. Presented at Open Forum, Ann. Meet. of the American Association For Respiratory Care (AARC), Las Vegas, 1999. Published in *Resp Care*, Oct 1999, 44(10): 1289.

Background: Nebulizer design is known to affect performance. In this study, we compared fine particle mass from nebulizers of four designs. **Methods:** We tested traditional disposable nebulizers (Baxter Misty-Neb™, Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™), nebulizers with collection bags (Westmed Circulaire™), and a Trudell **AeroEclipse*** (with breath actuation disabled). Five of each device with three replicates (n = 15) were tested using an *in-vitro* model of spontaneous breathing. A rigid bar was placed between the two compartments of a test lung (Michigan Instruments TTL). The drive lung was attached to a ventilator (Infrasonics Infant Star™) to simulate spontaneous breathing (tidal volume 0.6 L, rate 10/min, T_i 2 s). A bacterial/viral filter (Trudell MT3000) was placed between the nebulizer and slave lung. Flow was monitored between the test lung and filter (Novamatrix Ventcheck™). Albuterol solution (0.625 mg/mL) was placed into the nebulizers (4 mL), which were powered with air (8 L/min). Filters were replaced at one minute intervals (flow to the nebulizer was discontinued during filter replacement) until sputtering occurred. The filter was washed with methanol (20 mL) and albuterol concentration was measured with HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer. Fine particle mass was calculated as the product of mass % <4.7 µm and total nebulizer output. **Results:** Fine particle mass from the **AeroEclipse*** nebulizer was greater than that from the other nebulizers (P<0.001) (see figure). **Conclusions:** Fine particle mass was affected by nebulizer design. The clinical relevance of this finding awaits further investigation. Further evaluation of the breath-actuated feature of the **AeroEclipse*** is warranted.



5. **THE DELIVERY OF BUDESONIDE SUSPENSION VIA SMALL VOLUME NEBULIZERS: A COMPARATIVE IN-VITRO ASSESSMENT.** JP Mitchell, MW Nagel and AD Archer. Presented at Proc. Ann. Meet. Amer. College of Chest Physicians (ACCP), Toronto, November 1998, in *Chest*, 114 (4S), 295, (1998); and at the World Asthma Meeting, Barcelona, December 1998, in *Eur. Resp. J.*, 12S29, 7, (1998).

Purpose: To compare the performance of a new air entrainment small volume nebulizer (AE-SVN, Trudell Medical International) with other widely used SVN (LC-Star™ (PARI Respiratory Equipment), Updraft™ Neb-U-Mist™ (Hudson Oxygen Therapy Sales Co.), Circulaire™ (Westmed), Sidestream™ (Medic-Aid), Airlife™ Misty-Neb™ (Baxter Healthcare Corp.)) for the delivery of a suspension formulation (0.25 mg/ml budesonide (Astra Pharma Inc.)). **Methods:** Each SVN (n = 5 devices for each group, 3 replicates per device) was operated with compressed air at 50 psig at a flow rate of 8 l/min. The total mass of budesonide nebulized from 2 x 2 ml ampoules was determined by filter collection at the mouthpiece at a flow rate of 28.3 l/min. The SVN was operated until it spluttered, was then tapped gently to dislodge droplets back to the reservoir. Nebulization was deemed complete 20 seconds later. The mass of budesonide collected was determined by HPLC-UV spectrophotometry. **Results:** The delivery rate ((mean ± 1 S.D) µg budesonide/min) from the AE-SVN (102 ± 9) was significantly greater than with the other groups: (LC-Star™ (91 ± 6), Misty-Neb™ (49 ± 2), Sidestream™ (46 ± 4), Circulaire™ (26 ± 4) and Neb-U-Mist™ (25 ± 6)), (1-way ANOVA, p < 0.02). Duration of nebulization was shortest with the AE-SVN (221 ± 14 s), compared with LC-Star™ (229 ± 10 s), Sidestream™ (365 ± 19 s), Circulaire™ (420 ± 84 s), Misty-Neb™ (477 ± 25 s) and Neb-U-Mist™ (639 ± 15 s). **Conclusions:** The new AE-SVN is highly efficient at entraining the budesonide particles into the liquid droplets at these conditions. **Clinical Implications:** The good delivery rate combined with comparatively short duration of delivery offers the potential for rapid treatment and patient convenience.

6. **THE EFFECT OF SMALL VOLUME NEBULIZER (SVN) DESIGN ON FINE PARTICLE MASS DELIVERY OF A BRONCHODILATOR.** R Blacker, RW Morton, JP Mitchell, MW Nagel and DR Hess. Drug Delivery to the Lungs-X, London, UK, 1998, J. Aerosol Med., 13(1): 65.

Fine particle mass delivery was compared from six different SVNs, including continuous un-enhanced flow designs (Hudson Updraft-II Neb-U-Mist™), breath-enhanced nebulizers (Pari-LC-D™, Medic-Aid Sidestream®), nebulizers with aerosol collection bag (Westmed Circulaire™), and an **AeroEclipse*** with breath actuation disabled (Trudell Medical International). Five of each type of SVN were tested operating with air (8 l/min, 50 psig), using an *in-vitro* model that simulated spontaneous breathing by an adult (tidal volume 0.6 l, rate 10/min, TI = 2 s). A bacterial/viral filter was placed between the nebulizer and breathing simulator. In each case, salbutamol sulphate (Ventolin®) respirator solution (0.625 mg/ml, 4 ml) was placed into the reservoir of the SVN. The filters were replaced at one-minute intervals until sputtering occurred. The salbutamol collected on the filter was assayed by HPLC-UV spectrophotometry. Particle size was measured using a Malvern Mastersizer laser diffractometer. Fine particle mass delivery rates varied significantly from each of the SVNs from more than 110 µg/min (**AeroEclipse***) to ca. 20 µg/min (Circulaire™).

Salter HDN™, Salter Labs™

1. **BREATH-ACTUATED VS RESERVOIR NEBULIZERS FOR UNDILUTED ALBUTEROL.** D Geller, B Kesser. Presented at the 13th International Congress on Aerosols in Medicine, Interlaken, Switzerland, September 17-21, 2001.

Aim: Some Emergency Departments use undiluted albuterol in nebulizers designed to conserve drug during exhalation. We compared the *in vitro* performance of 4 devices to estimate which would be most effective clinically: **AeroEclipse*** Breath-Actuated Nebulizer ("BAN"); Circulaire® (C) and AeroTee™ (AT) which use a 750 ml reservoir bag to conserve drug during exhalation; and Salter HDN™ (S) with a 50 ml tower reservoir. **Method:** We studied 4 units of each nebulizer type in duplicate, using a Pari Proneb Turbo compressor. Nebulizers were filled with undiluted 0.5% albuterol, 1 ml (5 mg) or 2 ml (10 mg). Particle size distributions were measured by laser diffraction (Malvern SprayTec). Drug output (1 minute after "sputter") was captured on a filter between the device mouthpiece and a Pari breath-simulator, which used a recorded waveform from a 9 yr old male. Albuterol was measured by spectrophotometry, and fine particle dose (FPD) (mg of drug < 5 mm in size) was calculated.

Results:

| Neb | MMAD | FPD (1cc) | Minutes | FPD (2cc) | Minutes |
|------------|-------------|------------------|----------------|------------------|----------------|
| AE | 3.9 | 0.60 | 3.8 | 2.41 | 11.0 |
| AT | 4.8 | 0.03 | 2.0 | 0.62 | 3.2 |
| C | 2.5 | 0.09 | 2.0 | 0.65 | 3.7 |
| S | 8.5 | 0.08 | 2.0 | 0.57 | 3.7 |

Conclusions: The AE was superior to the reservoir-type nebulizers in fine-particle output for each fill volume. The AT and C had large dead volumes, and the S produced larger particles. These shortcomings were overcome with larger nominal doses. Each nebulizer produced 0.6-mg FPD of albuterol over 3½ minutes, but the AE required only half the starting dose. Albuterol 0.6 mg is a reasonable clinical respirable dose in a child with acute asthma. These findings must be taken into account when designing clinical treatment protocols for acute asthma. **Background:** Many nebulizers are designed to decrease the amount of drug that is lost during exhalation. The Circulaire® (Westmed) and AeroTee™ (Hudson) incorporate a 750 ml bag on the expiratory side of the nebulizer that collects aerosol while the patient exhales, making it available for inhalation on the next breathing cycle. The Salter HDN™ (Salter) has a 50 ml tower that acts as a reservoir. The **AeroEclipse*** BAN (Trudell/Monaghan) has a spring mechanism that allows generation of aerosol during inhalation only, so no drug waste occurs during exhalation. We recently reported the aerosol characteristics with these devices nebulizing unit-dose albuterol sulfate (2.5 mg/3 ml).¹ Delivery time with unit-dose (0.083%) albuterol can be long, which may increase personnel costs. To maintain lung-dose delivery and minimize the treatment time, some hospitals use drug-conserving nebulizers with small fill-volumes of undiluted (0.5%) albuterol for patients presenting with acute bronchospasm. We measured the particle size distributions and used a child's breathing pattern to compare albuterol output of these 4 drug-conserving nebulizers, using unit-dose albuterol 2.5 mg (3ml), 0.5% albuterol 5 mg (1ml) and 10 mg (2ml) nominal doses. We calculated the fine particle dose and measured the dose of drug remaining within the nebulizer and all attachments to determine the residual dose. For reference, we compared these results to those of a T-piece (Hudson Micromist) nebulizer using unit-dose albuterol to simulate conventional dosing.

Materials and Methods

- Drug:
 - Albuterol Sulfate 0.083% unit-dose (2.5mg/3ml)
 - Albuterol Sulfate 0.5% (5mg/ml) 1 and 2cc fill volumes
- Nebulizers: Circulaire® (Model 0260), AeroTee™ with Micromist Nebulizer (Model 1002), Salter HDN™ (Model 8960), and **AeroEclipse*** BAN (Figure 1)
- Compressor: PARI PRONEB TURBO
- 4 nebulizers of each type studied in duplicate
- Particle size by laser diffraction (Malvern Insittec)
- Breathing pattern from 9 year old male volunteer, using the PARI breath simulator (RR 19 bpm, Vt 421 cc, Ti 1.3 seconds)

- Definitions:
 - Inspired dose = drug on inspiratory filter
 - Residual dose = drug collected from nebulizer and accessory components after completion of nebulization.
 - Fine particle dose (FPD) = (Inspired dose) x (% of particles <5 µm) Figure 1.
 - Duration = time (minutes) from the beginning of nebulization to 1 minute past the onset of sputter.
- Samples assayed with spectrophotometer at 228 λ

Results:

| | | AeroEclipse* BAN | AeroTee™ | Circulaire® | Salter HDN™ |
|-------------------------|-------------------------|-------------------------|-----------------|--------------------|--------------------|
| Particle Sizing | MMD | 3.87 | 4.80 | 2.47 | 8.46 |
| | GSD | 2.3 | 2.0 | 2.1 | 2.0 |
| | % < 5 µm | 61.7% | 52.9% | 83.6% | 30.0% |
| 2.5 mg Unit Dose* | Duration (minutes) | 14.7 | 7.2 | 7.0 | 3.6 |
| | Inspired Dose (mg) | 0.77 | 0.37 | 0.14 | 0.30 |
| | Residual Dose (mg) | 1.5 | 1.8 | 2.1 | 1.9 |
| | Fine Particle Dose (mg) | 0.52 | 0.19 | 0.12 | 0.10 |
| 5 mg (1 ml) Dose | Duration (minutes) | 3.8 | 2.0 | 2.0 | 2.0 |
| | Inspired Dose (mg) | 0.97 | 0.06 | 0.11 | 0.28 |
| | Residual Dose (mg) | 3.5 | 4.9 | 4.6 | 4.4 |
| | Fine Particle Dose (mg) | 0.60 | 0.03 | 0.09 | 0.08 |
| 10 mg (2 ml) Dose | Duration (minutes) | 11.0 | 3.2 | 3.7 | 3.7 |
| | Inspired Dose (mg) | 3.9 | 1.2 | 0.8 | 1.9 |
| | Residual Dose (mg) | 5.8 | 8.7 | 8.6 | 6.9 |
| | Fine Particle Dose (mg) | 2.40 | 0.60 | 0.60 | 0.60 |

* Unit-dose data presented at ATS 2001¹

For comparison, the Hudson Micromist conventional T-Piece Nebulizer (with Unit-Dose 2.5 mg Albuterol) produced a fine-particle dose of 0.14 mg in 7.0 minutes.

Discussion:

- AeroEclipse* BAN** had highest FPD with all nominal doses:
 - FPD was 2.7 to 5.2 times higher with unit-dose; 6.7 to 20 times higher with 5 mg dose; 4 times higher with 10 mg dose.
 - Lowest residual dose
 - Higher fine particle fraction except for Circulaire®
- Nebulizer Inefficiencies:**
 - AeroTee™ and Circulaire® had high residual doses in part due to valves and collection bags.
 - Salter HDN™ produces larger particles
 - These inefficiencies were partially compensated for by increasing nominal dose to 10 mg (2 ml)
- Duration of Nebulization:**
 - AeroEclipse* BAN** had longer delivery time because it is breath actuated; no waste during exhalation
 - Using 0.5% albuterol, all nebulizers produced 0.6 mg fine-particle dose in < 4 minutes, but the **AeroEclipse* BAN** only required half the nominal dose to accomplish this.
- Comparison to Unit-Dose 2.5 mg:**
 - AeroEclipse* BAN** produced comparable FPD with unit-dose and 5 mg (1 ml) nominal dose, but delivery time was less than a third with undiluted drug.
- Comparison to conventional nebulizers:**
 - The FPD with the Hudson and unit-dose drug was 0.14 mg, similar to the reservoir-type nebulizers with unit-dose.
 - The higher FPD with **AeroEclipse* BAN** (all doses) and the reservoir nebs (10 mg dose) may result in better and longer lasting bronchodilation than the Hudson with conventional dosing, thus reducing number of treatments, therapist time, and total costs.

Conclusion

- AeroEclipse* BAN** was superior to the reservoir-type nebulizers at all nominal doses.
- AeroEclipse* BAN** has the additional advantage of being a dosimetric device, i.e. it will not operate or waste drug while the patient is coughing or resting. The patient and health care providers get visual feedback of adequate inspiratory effort necessary to actuate the nebulizer.
- Use of undiluted 0.5% albuterol may result in higher lung doses in a shorter amount of time. These results can be used as a guide when developing bronchodilator protocols for the hospital or E.D. setting.

Funded by the Nemours Foundation

¹ Geller D, Kesser B. Am J Respir Crit Care Med 2001; 163:A444 (Journal of Aerosol Medicine 14 (3) 2001; 395:P1-41)

2. **COMPARISON OF DRUG OUTPUT FROM 4 DIFFERENT RESERVOIR TYPE NEBULIZERS.** DE Geller, B Kesser. American Journal of Respiratory Care & Critical Care Medicine, 2001; 163(5); A444.

Rationale: Many nebulizers currently being marketed utilize different techniques to conserve drug that would normally be lost during exhalation. The Circulaire™ and Aero Tee™ nebulizers use a 750 cc reservoir bag to accumulate nebulized drug, while the Salter HDN™ uses a 50ml tower to serve as a reservoir. The **AeroEclipse*** nebulizer uses breath actuated nebulization to deliver drug only during inspiration. We evaluated all 4 nebulizers using a recorded pediatric breathing pattern to measure total drug output. We additionally measured the particle size characteristics of each type with the laser diffraction technique. **Methods:** 4 nebulizers of each type were studied in duplicate for sizing and total output characteristics. Each nebulizer was charged with a unit dose of 2.5 mgs albuterol sulfate in 3cc's. Sizing studies were averaged values performed over 5 minute runs on each nebulizer with a Malvern Spray Tec laser. Drug output was as calculated as the assayed amount of albuterol collected on a filter distal to the mouthpiece of the nebulizer. Simulated breathing was performed through the nebulizer by a Pari breath simulator from waveforms originally recorded from a healthy 9-year-old male.

Results:

| | Inspired dose | %>1 & <5M | Respirable Dose | Residual Dose |
|---------------------|----------------|------------|-----------------|---------------|
| AeroEclipse* | 0.64 ± 0.06 mg | 52.7 ± 2.0 | 0.34 ± 0.03 mg | 1.27 ± 0.09 |
| Aero Tee | 0.31 ± 0.09 mg | 41.2 ± 7 | 0.13 ± 0.04 mg | 1.51 ± 0.11 |
| Circulaire | 0.12 ± 0.03 mg | 61.9 ± 1 | 0.07 ± 0.02 mg | 1.72 ± 0.13 |
| Salter HDN | 0.25 ± 0.05 mg | 24.7 ± 5 | 0.06 ± 0.02 mg | 1.59 ± 0.10 |

Conclusions: The **AeroEclipse*** delivers a greater total dose of drug as well as a greater amount of drug in the fine particle range, most likely to deposit in the lower airways.

Salter Model 8900™, Salter Labs™

1. **DELIVERY OF ALBUTEROL VIA A NEW BREATH ACTUATED NEBULIZER: COMPARISON WITH CONTINUOUS JET NEBULIZERS.** DP Coppolo, MW Nagel, CC Doyle, VA Avvakoumova and JP Mitchell. Presented at the American Thoracic Society International Conference 2007, San Francisco, California, USA.

A new breath actuated nebulizer (**AeroEclipse*** II BAN, Monaghan Medical Corp., Plattsburgh, NY) has been developed to deliver medication only when the patient inhales. This study sought to determine the delivery of albuterol (3-ml fill of diluted solution (0.83 mg/ml)) as fine droplets < 4.7 μm aerodynamic diameter, and compare this fine droplet mass (FDM) with equivalent data from 4 widely available continuous jet nebulizers as benchmark devices. Each nebulizer (n=5; 3 replicates/device) was operated with compressed air at 50 psig at ca. 8 L/min to simulate hospital wall outlet conditions. The nebulizer on test was coupled to a breathing simulator set to mimic adult use (tidal volume = 600 ml, rate = 10 breaths/min; duty cycle = 0.33), and the emitted droplets were collected on an electret filter at the mouthpiece. The total mass of albuterol (TM) was assayed subsequently by HPLC-UV spectrophotometry. In a separate study, the droplet size distribution was determined by laser diffraction so that the fine droplet fraction (FDF) could be obtained. FDM was determined as the product of TM and FDF. FDM (mean SD) from the BAN was 791.84 g, delivered in 8 minutes. Corresponding values (FDM in time from start to sputter) for the VixOne™ (Westmed, Tucson, AZ), MicroMist™ (Hudson RCI, Temecula CA), Misty Max 10™ (Cardinal Health, McGaw Park (IL) and model 8900™ (Salter Labs, Arvin, CA) were 267.11 g in 6 min, 133.8 g in 4 min, 249.10 g in 6 min and 161.10 g in 5 min. Aside from dosage assurance imparted by breath-actuation, the **AeroEclipse*** II BAN delivered substantially more FDM/min than the other devices. The clinician is now able to treat either for extended high dose delivery (potentially eliminating the need for additional therapy), or titrate to a shorter interval based on response.

PediNeb™, Westmed™ Corp.

1. **DELIVERY OF BUDESONIDE INHALATION SOLUTION (BIS) THROUGH AN INFANT UPPER AIRWAY MODEL.** DE Geller, KC Kesser, HM Janssens, HAWM Tiddens. Am J Respir Crit Care Med 2003; 167(7):A508.

We investigated variables that may be important in the delivery of BIS to the lungs of infants, a challenging population for aerosol delivery. **Methods:** The Sophia Anatomical Infant Nose Throat (SAINT) airway model mounted on a breath simulator mimicked the breathing pattern of a 9-mo old infant (RR=30, Vt=100 ml, I:E ratio=1:1.3). Nebulizers were charged with BIS 0.25 mg and run continuously until dry. Drug captured on a filter distal to the SAINT model was the lung dose. Compressor: PARI PRONEB TURBO. Nebulizer/mask systems studied: VIX1/aerosol mask (AM), PediNeb pacifier device (PN) or blow-by (BB); **AeroEclipse*** neb and mask (AE); PARI LC+ and PARI LC*/ PARI Baby mask (PB), Fish mask (FM), and AE masks. The AE neb/mask was also studied with an ill breathing pattern (RR=50, Vt=100, I:E=1:2). **Results:** Lung dose ranged from 2.0 to 7.6% of the neb charge. Lung dose was AE (5.0%) > VIX1 (3.5%), LC+/FM (3.2%), LC*/PB (2.9%), and LC+/PB (2.8%). Also, VIX1/AM (3.5%)>VIX1/PN (2.5%)>VIX1/BB (2.0%). The lung dose of the LC+ and LC* more than doubled (6.8 and 6.3%) when used with the AE mask. Lung dose increased with the ill breath pattern in proportion to increased minute ventilation (7.6%). **Conclusion:** 1) The AE system provided higher lung dose

than other nebulizers with standard masks. 2) Mask design and fit can substantially impact nebulizer performance. 3) PN performed better than BB, but not as good as a mask. If crying decreases lung dose by 75%, we speculate that the PN and BB (non-crying) may improve lung dose vs mask with a crying infant. 4) An increase in lung dose may occur in ill infants if minute ventilation is elevated. *This abstract is funded by: AstraZeneca.*

AeroChamber* Valved Holding Chamber and MDI

1. **THE DELIVERY TIME, EFFICACY, AND SAFETY OF BETA AGONIST BRONCHODILATOR ADMINISTRATION WITH THE AEROECLIPSE* BREATH - ACTUATED NEBULIZER ("BAN").** RS Pikarsky, T Farrell, R Acevedo, W Fascia and C Roman. CHEST 2001; 120(4): 218S.

Purpose: Aerosol delivery consumes the highest level of Respiratory Care resources. This study evaluated the delivery time, efficacy, and safety of rapidly nebulized Albuterol with the use of the **AeroEclipse*** Breath Actuated Nebulizer as compared to both an MDI with **AeroChamber*** VHC (both from Monaghan Medical Corp. Plattsburgh, N.Y.) and the Airlife Misty-Neb Nebulizer (SVN) (Allegiance Healthcare Corporation). **Methods:** A consecutive, non-randomized, mostly COPD population receiving pre & post bronchodilator testing in our Pulmonary Function Lab were studied. Three different Albuterol medication dosages were administered with the BAN: 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline, 1.0 ml (5 mg) of undiluted Albuterol, and 0.75 ml Albuterol (3.75 mg) using an oxygen flow rate of 8 L/min. Two puffs of Albuterol were administered by MDI with **AeroChamber*** VHC. Treatments with the SVN consisted of nebulizing 2.5 mg of Albuterol diluted with 3 ml of Normal Saline Unit Dose (UD) using an oxygen flow rate of 8 L/min. The Sensormedics Vmax 22 Pulmonary Function System was utilized to measure FEV1. A standardized subjective questionnaire to determine side effects was completed.

| Nebulizer (n) | Dose | % Change FEV1 | Time (min) | Tremulousness |
|---------------------------------|-----------------------|---------------|------------|---------------|
| AeroEclipse* BAN (12) | 0.5 ml + 0.5 ml NS | 8.2% | 2.67* | 0 |
| AeroEclipse* BAN (64) | 1.0 ml undil. | 10.9% | 3.29* | 17 |
| AeroEclipse* BAN (23) | 0.75 ml undil. | 5.6% | 1.30* | 5 |
| MDI (21) | 2 puffs | 8.5% | 2.86** | 1 |
| Misty-Neb (52) | 2.5 mg UD | 9.1% | 8.33 | 2 |

Results: The table shows the Albuterol dosages, mean % change of FEV1 from pre-treatment and 10 minute post treatment, mean administration time and tremulousness. The mean treatment time with all BAN patients was 2.78 minutes as compared to 8.33 minutes with the SVN (p<.001) *. The mean treatment time with the MDI was 2.86 minutes as compared to 8.33 minutes with the SVN (p<.001) **. The changes in FEV1 were not significant. There was no difference in heart rate, respiratory rate or nausea. Seventeen patients receiving the 1.0 l undiluted Albuterol indicated an increase in tremulousness. **Conclusion:** The rapid administration of Albuterol in the 0.5 ml + 0.5 ml NS and 1.0 ml undiluted doses using the BAN was equally efficacious as the MDI with **AeroChamber*** VHC and SVN UD. The 1.0 ml Albuterol dosage has the highest incidence of tremulousness. The 0.75 ml Albuterol dosage under-performed. Delivering 0.5 ml Albuterol (2.5 mg) with 0.5 ml Normal Saline using the BAN offered the best delivery time, efficacy and safety profile of the nebulizer trials. The BAN performance was comparable to the MDI with **AeroChamber*** VHC. **Clinical Implications:** In a health care facility that delivers large volumes of aerosol treatments, the decrease in delivery time could have a significant impact on resource utilization. The results supported changes in the Respiratory Care practice throughout Crouse Hospital. Further studies evaluating additional medication dosing regimens measuring safety, efficacy and resource utilization are needed.


Large Volume Nebulizers

1. **RAPID DELIVERY OF BRONCHODILATOR MEDICATION IS POSSIBLE USING A BREATH-ACTUATED SMALL VOLUME NEBULIZER AS AN ALTERNATIVE TO EXTENDED.** J Mitchell, D Coppolo, C Doyle, MW Nagel, KJ Wiersema. Presented at the American Association for Respiratory Care (AARC) Open Forum December 2007, Orlando, Florida.

Background: Inhaled beta-2 adrenergic agonist bronchodilators are often given to patients with severe reversible airways disease by continuous nebulization in extended treatments. However data are limited as to whether or not shorter, but higher concentration delivery is as an effective treatment modality. The development of a new breath-actuated nebulizer (**AeroEclipse*** II, Monaghan Medical Corp., Plattsburgh, NY (AEII BAN)) provided an opportunity to compare the two treatment methods in a laboratory study before undertaking a clinical comparison. We investigated the delivery of diluted generic respirator solution albuterol by a widely used continuous jet nebulizer (MiniHeart® Hi-Flo, Westmed Corp., Tucson, AZ (CONT)) with that from the AEII BAN. **Method:** The continuous nebulizers (n=5) were operated with 8 L/min air supplied at 50 psig with a 20-ml fill (albuterol concentration of 0.5 mg/mL). A similar number of AEII BANs were operated with ca. 8.0 L/min air at 50 psi with a 1-ml fill (albuterol concentration of 5 mg/mL). Aerosol from both nebulizers was sampled onto electret filters using a breathing simulator mimicking small child use (250-ml tidal volume,

inspiratory/expiratory ratio 1:2, rate 12 cycles/min) until onset of sputtering. Assay for albuterol was undertaken by UV spectrophotometry. In a parallel study, droplet size distributions were determined by laser diffractometry, so that the fine droplet fraction (mass % < 4.7 μm diameter) likely to penetrate to the airways of the lungs (FDF) could be determined.

Results: Values of FDF for the AEII BAN and were 78.4% and 62.0% respectively. The AEII BAN delivered $758 \pm 36 \mu\text{g}$ as fine droplets after 4-min (delivery rate of $190 \pm 9 \mu\text{g}/\text{min}$), compared to $180 \pm 76 \mu\text{g}$ in the same period by (delivery rate of $45 \pm 19 \mu\text{g}/\text{min}$). Conclusions: The faster delivery rate from the AEII BAN/high albuterol concentration modality (un-paired t-test, $p < 0.001$) may offer an important clinical alternative to CONT/low concentration treatment modality.



2. **A BREATH-ACTUATED SMALL VOLUME NEBULIZER (BAN) OFFERS A RAPID ALTERNATIVE TREATMENT MODALITY FOR THE DELIVERY OF BRONCHODILATORS FOR ASTHMATIC PATIENTS IN A SEVERE EXACERBATION.** DP Coppolo, CC Doyle, JP Mitchell, MW Nagel, KJ Wiersema. Presented at the American Association for Respiratory Care (AARC) Open Forum, 2006.

Large volume continuous nebulizers (LVNs) are often used for the delivery of beta-2 adrenergic agonist bronchodilators in the emergency department to treat severe, reversible airways disease, in particular asthma 1. Treatment time, however, can be lengthy for delivery of the typical LVN fill volume from 20- to 120-ml. Quick delivery of a bronchodilator with an efficient nebulizer may help relieve symptoms from bronchospasm in a shorter period of time. We report a study in which the delivery of diluted generic respirator solution albuterol by LVN (Hope, B&B Medical Technologies Inc., Loomis, CA) was compared with that from a small volume breath-actuated nebulizer (BAN) (**AeroEclipse***, Monaghan Medical Corp., Plattsburgh, NY). The LVNs (n=5) were operated with 10 L/min air supplied at 50 psig with a 20-ml fill (albuterol concentration of 0.167 mg/ml). A similar number of BANs were operated with 8.0 L/min air at 50 psi with a 3-ml fill (albuterol concentration of 0.833 mg/ml). The aerosol from the LVNs was sampled continuously until onset of sputtering at 12 L/min via a Dreschel filter/bottle where the albuterol was captured quantitatively. Aerosol from the BANs was sampled onto electret filters using a breathing simulator (600-ml tidal volume, inspiratory/expiratory ratio 1:2, rate 10 cycles/min) until onset of sputtering, so that operation of the breath actuation mechanism was effected. Assay for albuterol was undertaken by UV spectrophotometry. In a parallel study droplet size distributions were determined by laser diffractometry, so that the fine droplet fraction < 4.8 μm diameter likely to penetrate to the airways of the lungs could be determined. Fine droplet albuterol delivery rates were constant as a function of time for all nebulizers. After 15-min, the LVNs had supplied $127.3 \pm 37.4 \mu\text{g}$ as fine droplets at a rate of $8.5 \pm 2.5 \mu\text{g}/\text{min}$. In contrast, the BANs delivered $810.0 \pm 20.4 \mu\text{g}$ in a 10-min period, equivalent to a rate of $81.0 \pm 2.0 \mu\text{g}/\text{min}$. The significantly higher delivery rate from the BAN group (un-paired t-test, $p < 0.001$) offers an important clinical alternative to the LVN in the emergency department where rapid delivery of a bronchodilator is critical.

Reference: McPeck M, Tandon R, Hughes K, Smaldone GC. Aerosol delivery during continuous nebulization. Chest. 1997;111:1200-1205.

General Information

1. **EUROPEAN PHARMACEUTICAL AEROSOL GROUP (EPAG) NEBULIZER SUB-TEAM: ASSESSMENT OF PROPOSED EUROPEAN PHARMACOPEIAL (PH. EUR.) MONOGRAPH 'PREPARATIONS FOR NEBULIZATION.** E Berg, J Dennis, J Jauernig, M Karlsson, C Kreher, P Lamb, JP Mitchell, S Nichols, D Sandell, M Tservistas. Presented at Drug Delivery to the Lungs, 2006.

Summary – The EPAG Nebulizer Sub-Team was formed to develop industry best practices for the evaluation of nebulizer systems. Its primary objective is to support the development of a new monograph for the European Pharmacopeia concerned with the characterization of preparations for nebulization. Specific tasks are: (1) to establish when it is appropriate to chill the Next Generation Pharmaceutical Impactor (NGI) to avoid bias due to heat-transfer related evaporation; (2) to validate the use of uncoated collection cups for the NGI; (3) to produce a position statement concerning the appropriate use of cascade impaction and laser diffractometry for nebulizer characterization; (4) to establish appropriate breathing patterns for nebulizer mass output testing by breathing simulator. The sub-team is also assisting EPAG develop expert commentary in relation to the development of a proposed international standard (ISO 27427) that will focus on establishing the performance of nebulizing systems during their design verification.

Introduction – Jet and ultrasonic nebulizers continue to be widely used modalities for inhaled aerosol therapy, with new designs such as vibrating mesh and membrane systems being marketed [1], despite the widespread availability of pressurized metered dose inhalers (pMDIs) and dry powder inhalers (DPIs). This is largely because they can be used to deliver almost all therapeutic classes of drugs to the respiratory tract whether the patient is ambulatory or on mechanical ventilation [2]. Nebulizers are typically manufactured for use with a variety of drug products often from different pharmaceutical companies, depending upon the judgment of the prescribing clinician. This contrasts with the situation for pMDIs and DPIs, where the device and drug product are directly linked, and are almost always the responsibility of the pharmaceutical company manufacturing the drug product. As a result, the regulation of nebulizers has traditionally taken place through the devices part of the various agencies, following processes that are separate from those used to regulate the drug products with which they are used. In a departure from this practice, nebulizers are now being included with the other classes of portable inhaler in a joint Health Canada-EMA regulatory guidance on Pharmaceutical Quality of Inhalation and Nasal Products [3, 4].

In terms of nebulizers as drug delivery devices, European-wide guidelines were developed about 5 years ago that established setting uniform standards for their use [5], with performance testing undertaken in accordance with a regional (CEN) standard [6]. At the pharmacopeial compendia level, a monograph on the characterization of preparations for nebulization is in the process of being reviewed by the Inhalanda committee for possible inclusion in the European [7] and US [8] Pharmacopeias. Recognizing the need for harmonization between device- and drug product focused standards where practical, much of the proposed methodology in the draft monograph is based on the procedures described in the CEN standard [6]. However, the advent of the Next Generation Pharmaceutical Impactor (NGI) [9] took place after this standard had been issued. The ability of the NGI to operate at flow rates as low as 15 L/min [10], the flow rate adopted in the CEN standard as representative of adult inhalation, has made it possible to propose this impactor as an apparatus that is suitable for droplet aerodynamic size characterization in the monograph.

The Nebulizer Sub-Team of EPAG was formed in 2005 to reappraise methods that are used for *in vitro* characterization of nebulizers in the context of the above developments. This was deemed both timely and necessary in view of the increased attention being paid to these devices by both the compendia and regulatory agencies, coupled with the development of new types of devices, including breath-actuated and adaptive aerosol delivery-based systems. For instance, methods that are based on constant flow rate sampling are unable to assess properly the function of nebulizers that are either breath-enhanced or breath-actuated. As a further example, optical methods for droplet size characterization, in particular laser diffractometry, are rapid and therefore potentially useful as a tool for assessing quality control of drug product used with a nebulizer. However, without appropriate precautions, such methods are inappropriate for nebulizer designs that allow inherent evaporation of nebulized aerosol, which includes all constant output nebulizers. They are also particularly unsuitable for suspension-based formulations where droplets may include no active drug particles or may contain more than one such particle per droplet. This limitation is not always evident in industry guidance and standards documentation.

Specific Work Being Undertaken Currently by the Nebulizer Sub Team – The sub team started work by addressing four specific work packages that each relate to the proposed pharmacopeial monograph.

1. Assessment of the Need to Chill the NGI to Prevent Heat-Transfer-Related Evaporation – In the late 1990s, Finlay and Stapleton reported that the effect of heat transfer from the impactor to the aerosol droplets being measured can be to bias measurements to finer sizes, when working with the Andersen 8-stage impactor (ACI) [11]. The NGI has significantly greater mass than the ACI, and may therefore be more susceptible to this phenomenon. Attempts to cool the impactor to the temperature of the nebulizer-produced aerosol by immersion in a water bath, though feasible, are awkward and time consuming to perform [12]. Chilling the impactor to a temperature close to +5°C has therefore been proposed as a simpler alternative to water immersion [13]. Although, operating the impactor with air close to saturation is also a practical alternative to minimize evaporative changes [14, 15], there is the risk of condensational growth of droplets and the complication of testing nebulizers in laboratory conditions that do not simulate actual clinical use. Evaporative effects appear also to be dependent upon the nebulizer type, being most apparent with devices that do not entrain air as part of the nebulization process [16]. In summary, there is currently a lack of peer-reviewed experimental data that could be used to develop guidance on when the various techniques are applicable and with which types of nebulizer.

The sub-team has organized a series of experiments to evaluate the effect on NGI-measured droplet size distributions using a cooled impactor (5°C) compared with impactor operated at room ambient temperature (20°C). Included in this experimental are three nebulizer types that represent different categories in terms of droplet formation. These are the MistyMax™ (Cardinal Health, USA), which is a conventional non air entrainment jet nebulizer, the LC-Plus™ (PARI GmbH, Germany), which is an air entrainment jet nebulizer, and the AeroNeb® (Nektar Therapeutics USA), representing newer vibrating mesh/membrane systems. This latter nebulizer is non-air entrainment in design.

Measurements are being made with a generic salbutamol solution formulation, and up to six laboratories are collaborating so that both inter- and intra-nebulizer variability can be assessed. Although data are currently undergoing statistical analysis, preliminary findings are that chilling the NGI may be necessary in determining aerosol size distributions for some nebulizer systems.

2. Assessment of the Need to Coat the Collection Cup Surfaces of the NGI to Mitigate Droplet Bounce – It is well known that impactors are vulnerable to particle bounce and blow-off, biasing particle size distribution data to finer sizes. Various methods have been proposed to mitigate the effect; including coating the collection surfaces with grease or using non-volatile agents that create a tacky surface [17]. Liquid droplet bounce is unlikely but not proven. A study to confirm that coating is not needed was therefore included in the work of the sub-team. This is a two centre experimental evaluation of the behavior of aerosols of a generic salbutamol generated by two different jet nebulizers (Sidestream™ (Respironics Ltd., UK) representing a relatively low output conventional device, LC-plus™ (PARI GmbH, Germany), representing a higher output air entrainment nebulizer). Coating has been undertaken with a thin layer of silicone oil. Although the data are currently undergoing evaluation, initial indications are that coating NGI collection cups is unnecessary, irrespective of nebulizer type, for the collection of aerosol droplets.

3. The Choice of Appropriate Breathing Patterns for Nebulizer Testing – The breathing pattern (tidal volume = 500 mL, inspiratory/expiratory ratio 1:1 (sinusoidal), rate = 15 breaths/min) currently specified in the draft monograph for the determination of active substance delivery rate and total active substance delivered is the same as that adopted in the CEN standard [6], which is based on a normal adult at rest. This pattern was adopted in the CEN standard because it is simple to simulate and reproduce within a test laboratory [18], and the group developing the monograph felt it highly desirable to retain harmonization, given the desirability to have comparable tests for both nebulizers as delivery devices and for the drug products that may be used with them. However, formulations have been marketed that are specifically targeted at infants and young children, whose breathing patterns are very different to those of adults [19]. The sub-team is therefore in the process of developing an evidence-based position statement that will recommend appropriate breathing conditions for these classes of patient. In addition, they are examining the feasibility of capturing exhaled medication during nebulizer operation on a breathing simulator with a view to quantifying mass recovery of active substance, where this is feasible. Such a test might also be indicative of fugitive droplet emission from nebulizers, which is a concern in some countries, post the SARS outbreaks that occurred in 2003-4.

4. The Application of Laser Diffractometry and Cascade Impaction to Nebulizer Testing – The appropriate use of optical methods, in particular low angle laser light scattering (laser diffractometry), would be a significant advantage to industry in the context of routine droplet size testing that involves nebulizer-generated aerosols for drug product quality control purposes [20]. It is well known that the drug product is directly associated with droplet size with solution formulations, so that laser diffractometry should provide a close approximation to the aerodynamic size distribution of the drug product itself [21]. At present, laser diffractometers are used for nebulizer-generated aerosol characterization [22, 23], and the technique has also received limited recognition as a measurement tool in the regulatory literature, provided that it is supported by measurements using cascade impaction, where drug substance traceability is achieved [3, 4]. The sub-team is therefore developing an evidence-based position statement that addresses points to consider for the use of laser diffractometry as an adjunct to support cascade impaction for droplet size distribution measurement where appropriately validated.

The Sub-Team as a Source of Expertise on Nebulizer-Related Issues – In addition to the specific tasks already described, the sub-team which is comprised of individuals with experience in both formulation and device aspects is tasked to provide on-going expert comment as needed to help develop an EPAG position in relation to future regulatory, compendial and national/international standards that relate to nebulizers. In this context, the imminent development of a new international standard for nebulizing systems (ISO 27427) through committee ISO-TC121/SC2 is providing the opportunity to develop a consensus input into the process at the public comment stages via participating national standards bodies of countries which have EPAG members.

References

- Dennis JH. New developments in nebulizer technology. 2004. In: Eds. J Boe, R O'Driscoll, JH Dennis. Practical Handbook of Nebulizer Therapy. London, UK. Martin Dunitz an imprint of Taylor & Francis Group. 41-60.
- Hess DR. 2000. Nebulizers: Principles and performance. *Respir. Care* 45(6):609-622.
- EMA. 2006. Guideline on the Pharmaceutical Quality of Inhalation and Nasal Products. EMEA/CHMP/QWP/49313/2005-final. Available at: www.emea.europa.eu/pdfs/human/qwp/4931305en.pdf
- Health Canada, 2006. Guideline on the Pharmaceutical Quality of Inhalation and Nasal Products. File Number 06-106624-547. Available at: www.hc-sc.gc.ca/dhp-mps/prodpharma/applic-demande/guide-ld/chem/inhalationnas_e.html
- Boe J, Dennis JH, O'Driscoll BR. 2001. European Respiratory Society guidelines on the use of nebulizers. *Eur. Respir. J.* 18:228-242.
- European Committee for standardization (CEN). 2001. Respiratory therapy equipment – Part 1: Nebulizing systems and their components. EN 13544-1.
- European Directorate for the Quality of Medicines (EDQM). 2006 Preparations for Nebulisation: Characterisation (2.9.44). *Pharm. Europa*, 18.2:280-282.
- Truman K, Nichols S, Mitchell J, Vanneste C, Tserivas M, Dennis J. 2006. Preparations for Nebulization: Characterization. *Pharm. Forum*, 32(4):1348-1352.
- Marple VA, Roberts DL, Romay FJ, Miller NC, Truman KG, Van Oort M, Olsson B, Holroyd MJ, Mitchell JP, Hochrainer D. 2003. Next Generation Pharmaceutical Impactor (A new impactor for pharmaceutical inhaler testing) – Part I: Design. *J. Aerosol Med.*, 16(3):283-299.
- Marple VA, Olson BA, Santhanakrishnan K and Mitchell JP. 2004. Next Generation Pharmaceutical Impactor (A new impactor for pharmaceutical inhaler testing) – Part III: Extension of archival calibration to 15 L/min. *J. Aerosol Med.*, 17(4):335-343.
- Finlay WH, Stapleton KW. 1999. Undersizing of droplets from a vented nebulizer caused by aerosol heating during transit through an Andersen impactor. *J Aerosol Sci.* 30:105-109.
- Jauernig J, Ohl S, Luber M, Keller M. 2003. Differences in results obtained with the Next Generation Impactor (NGI) for Pulmicort® suspension and according to the CEN-standard EN-13544-1 for a NaF-solution. *Drug Delivery to the Lungs-14*. London, UK. The Aerosol Society. 41-44.
- Berg E, Asking L. 2004. Nebulizer droplet size distribution - refrigerated NGI at 15 L/min. In: Eds. RN Dalby, PR Byron, J Peart, JD Suman, SJ Farr. *Respiratory Drug Delivery IX*. River Grove IL, USA. Davis Horwood International Publishing LLC. 361-363.
- Jauernig J, Ohl S, Knoch M, Keller M. 2004. Effects of the test set-up, formulation, and nebulizer type on aerodynamic droplet characteristics. In: Eds. RN Dalby, PR Byron, J Peart, JD Suman, SJ Farr. *Respiratory Drug Delivery IX*. River Grove IL, USA. Davis Horwood International Publishing LLC. 609-612.
- Jauernig J, Hug M, Knoch M, Keller M. 2002. Contribution to the aerodynamic particle size assessment of nebulizer drugs using the New Generation Impactor (NGI). *Drug Delivery to the Lungs-13*. London, UK. The Aerosol Society. 44-47.
- Dennis JH. 2004. Theory and science of nebulizer use. In: Eds. J. Boe, R O'Driscoll, JH Dennis. Practical Handbook of Nebulizer Therapy. London, UK. Martin Dunitz an imprint of Taylor & Francis Group. 3-17.
- Mitchell JP, Nagel MW. 2003. Cascade Impactors for the Size Characterization of Aerosols from Medical Inhalers: Their Uses and Limitations. *J. Aerosol Med.* 16(4):341-377.
- Dennis JH, Pieron CA. 2004. Quality control and standards in nebulizer performance and use. In: Eds. J. Boe, R O'Driscoll, JH Dennis. Practical Handbook of Nebulizer Therapy. London, UK. Martin Dunitz an imprint of Taylor & Francis Group. 19-40.
- Stocks J. and Hislop AA. 2002. Structure and function of the respiratory system: developmental aspects and their relevance to aerosol therapy. In: Eds. H. Bisgaard, C. O'Callaghan and G.C. Smaldone. *Drug Delivery to the Lung*. New York, USA, Marcel Dekker Inc. 47-104.
- Nerbrink, O, Dahlback M, and Hansson H-C. 1994. Why do medical nebulizers differ in their output and particle size characteristics? *J. Aerosol Med.*, 3:259-276.
- Clark AR. 1995. The use of laser diffraction for the evaluation of the aerosol clouds generated by medical nebulizers. *Int. J. Pharm.* 115:69-78.
- Kwong WTJ, Ho SL, Coates AL. 2000. Comparison of nebulized particle size distribution with Malvern laser diffraction analyzer versus Andersen cascade impactor and low-flow Marple personal cascade impactor. *J. Aerosol. Med.* 13:304-314.
- Vecellio None L, Grimbard D, Becquemin MH, et al. 2001. Validation of a laser diffraction method as a substitute for cascade impaction in the European project for a nebulizer standard. *J Aerosol Med.* 14:107-114.

Manufactured by:

Trudell Medical International

725 Third Street
London, Ontario,
Canada N5V 5G4

Tel: 1 519 455 7060

Fax: 1 519 455 7858

Email: customerservice@trudellmed.com
www.trudellmed.com

Authorized European Representative:

Trudell Medical International Europe Limited

BioCity Nottingham,
Pennyfoot Street
Nottingham, NG1 1GF, UK

Email: customerservice@trudellmed.com
www.trudellmed.com

Canadian Hospital Distributor:

Trudell Medical Marketing Limited

758 Third Street, London, Ontario
Canada N5V 5J7
Tel: 1-800-265-5494 (In Canada)
Email: tmml@tmml.com
www.tmml.com

UPDATED: October 2008