Instructions for Use
Mode d'emploi
Bedienungsanleitung
Instrucciones de Uso
Istruzioni per l'uso

in-check[™]
DIAL
G16

CE 0120


HS CLEMENT CLARKE
INTERNATIONAL

Tradition and Innovation


How to use the *In-Check DIAL G16*

en
fr
de
es
it


1a


1b


2


3


4

5

6

7


e.g.

* Clinically
Effective
Flow Rate
30-90L/Min

8

Instructions for Use

Caution

Consult 'Instructions for Use'.

This product complies with the essential requirements of the medical devices directive (93/42/EEC). Compliance has been verified by UK notified body per 0120 (SGS United Kingdom Ltd).

Introduction

The **In-Check DIAL G16** is an inhalation airflow training meter that can help educate and assess patients who use inhaler devices.

Inhaler devices are designed to deliver medication to the respiratory tract, and the speed of inhalation through them (the inspiratory flow) can have a significant effect on the quantity of drug delivered and the clinical efficacy of the product.

The **In-Check DIAL G16** is designed to simulate the “internal resistance” of common inhaler devices, and measure inspiratory flow. These measurements enable the healthcare professional to encourage patients to modify their inspiratory technique (by inhaling with more, or less effort), in order to achieve a flow rate consistent with clinical efficacy. The coloured ‘flow’ icons show the clinically effective flow ranges for each different inhaler device. These ‘flow’ icons do not imply any comparison between devices.

Patients that cannot achieve the suggested inspiratory flow for their inhaler may not gain maximum benefit from their prescribed medication, and healthcare professionals may wish to take this factor into account when selecting with the patient, the device that is the most suitable.

Inspiratory Flow and Clinically Effective Flow Range

The inspiratory flow through an inhaler is one of the factors that will influence the clinical effect of the drug delivery from that device. The most effective delivery occurs when the patient achieves a flow within the clinically effective flow range. Flow rates outside this range, may result in a diminished deposition and clinical efficacy.

Inhaler devices

Drug delivery from the various types of inhaler devices is produced by different methods.

Inhaler devices are designed to deliver drug particles of a certain size to the small airways during inhalation. Particles of this size (generally agreed to be approximately between 1 and 5 microns) are known as the “respirable dose”. The particles are either in aerosol (in a suspension or a solution) or dry powder form.

Pressurised Metered Dose Inhalers (pMDIs)

With most pMDIs, the aerosol is delivered under pressure at high speed (often over 90 kilometres per hour). The inhalation should be timed with actuation of the device and should be slow and steady. Inhaling too fast may cause a greater proportion of the aerosol to impact at the back of the throat and be subsequently swallowed, thus reducing the beneficial clinical effect and increasing the potential for local and systemic side effects.

With the majority of pMDIs, manual depression of the drug canister actuates the drug delivery. However, with breath-actuated metered dose inhalers (e.g. Autohaler™), the aerosol is released by mechanical actuation, triggered by the patient breathing in through the device at any flow rate above a minimum level. At inspiratory flows below this, the patient will not receive any of the medication, because a dose will not be released.

pMDIs with Holding Chamber/Spacer

It is recognised that the optimum inhalation technique for using a pMDI with a holding chamber/spacer is a slow inhalation (30 to 60 l/min). As the resistances of the values of most chamber/spacer devices are low, the **In-Check DIAL G16** can be set to “conventional pMDI”  to provide an approximate resistance for inspiratory flow measurements to be made.

Dry Powder Inhalers (DPIs)

Drug delivery from DPIs is triggered by inhaling through the device. Turbulent airflows cause a metered quantity of powder to be drawn into the airflow, and follow a specific pathway within the inhaler, through internal features that create a resistance that is designed to break up the medication into particles of a respirable size.

As the internal design of each DPI is different, and the formulation of medication is not identical, the resistance the patient encounters when inhaling is different from device to device.

The effort required for the patient to achieve a given inspiratory flow will increase as the internal resistance of the device increases. For the same patient effort, the higher the resistance the lower the resulting inspiratory flow.

As the great majority of asthma medication is delivered via inhalers, correct inhaler technique is an important factor in the management of this disease. Patients require their medication for both short-term relief and long-term prevention, and the delivery of these drugs to the lungs is affected by inspiratory flow. Assessing inspiratory flow is an important aspect of device education and review. Determining whether or not a patient can obtain the appropriate inspiratory flow rate for a device can be an important factor in improving inhaled treatment.

In-Check DIAL G16

In-Check DIAL G16 is a low-range inspiratory flow meter (15 to 120 l/min) that has a selectable resistance from high to low, shown by the coloured ‘flow’ icons on the accompanying laminated card, calibrated to enable the measurement of airflow as if the patient was using certain different inhalers. (see laminated card).

The “conventional pMDI” position can be used to approximate the resistance when a conventional pMDI is attached to a holding chamber or spacer device.

In-Check DIAL G16 provides a measurement of inspiratory flow through five classes of device resistance. However, the resistance is not precise for all inhalers, e.g. the resistance of the Genuair is slightly lower than the resistance of the Turbohaler Mk III, but they are both medium resistance inhalers. To mimic the precise resistance of a particular inhaler, specific adaptors are available to use with the DIAL set to ‘conventional pMDI’.

Should new types of inhaler become available, then **In-Check DIAL G16** can be used to measure the new inspiratory flow, by knowing the resistance of the new inhaler.

IMPORTANT

As with any inhalation device, it is important to check for loose foreign objects before the device is used. The transparent material used in the construction of the **In-Check DIAL G16** enables the user to make a visual check before inhalation. Patients should be prevented from exhaling through the device prior to use.

To reset the In-Check DIAL G16

Hold the instrument vertically with the mouthpiece uppermost, so that the rounded end of the meter can be tapped against the other hand or a horizontal surface, such as a table.

A hard tap will dislodge the magnetic resetting weight, which will return the red cursor to a start position. When this has happened, the meter must turn through 180 degrees to return the magnetic weight to its resting position.

 **Do Not** try to reset the **In-Check DIAL G16** as if it were a mercury thermometer – this action causes serious damage to the piston and pointer.

How to use the In-Check DIAL G16 - see diagrams p1

1. Reset the **In-Check DIAL G16**. See above.
2. Align the dial selector with the desired coloured icon – an audible “click” should be heard.
3. Attach a clean mouthpiece. Disposable one-way inspiratory mouthpieces are preferred.
4. Ask the patient to exhale fully.
5. Seal lips around the mouthpiece.
According to the inhaler chosen, instruct the patient to inhale in the manner recommended by the manufacturer Training Device.
6. Record the inspiratory flow from the position of the red cursor against the scale. Reset, and repeat two more times, ensuring correct technique each time.
7. Compare values achieved with target flows for that device. To operate an inhaler device correctly, the patient should be able to achieve a flow rate within the clinically effective range. *e.g Accuhaler: Clinically Effective Flow Rate 30-90L/Min
8. If after repeated training the patient is not able to achieve these values, then the healthcare professional may wish to assess the patient's ability to use an alternative type of inhaler.

Performance Accuracy

Accuracy +/- 10% or 10 l/min (whichever is the greater) and repeatability of +/- 5 l/min.

Cleaning your In-Check DIAL G16

Where local infection control guidelines exist, these should be respected.

Immerse **In-Check DIAL G16** in warm (but not hot) mild detergent solution for 2-3 minutes (maximum 5 minutes). Agitate the meter to ensure thorough cleaning.

Rinse in warm water and shake to remove any excess water. It is important to rinse thoroughly to prevent salt spots appearing on the inside of the body and the spindle.

To shake excess water from the **In-Check DIAL G16**, hold only at the end furthest away from the **DIAL** selector.


Allow to dry thoroughly before using again.

The expected life of the **In-Check DIAL G16**, in normal use, is two years.

 To maintain hygiene the **In-Check DIAL G16** should be used with disposable inspiratory one-way mouthpieces. If use of this product is governed by local infection control guidance, that guidance should be respected in the absence of local guidance.

Disposal Note: The **In-Check DIAL G16** can be disposed of as normal household waste.

Interpreting information from the In-Check DIAL G16

In-Check DIAL G16 was modified in January 2016, in response to customer feedback, to provide information about a wide range of new inhaler devices and the inspiratory flow rates that are associated with clinical efficacy. The new scale is shown above.

The information provided is for guidance only. It does not imply that any particular product will be clinically effective. There are other aspects than inspiratory flow that contribute towards clinical efficacy. The information is based on published clinical studies. In some cases studies may have been performed on a particular formulation in a specific device and may be suggestive that other formulations will behave similarly.

DPIs

Laboratory studies demonstrate that DPIs have a degree of flow-dependence, meaning that improved delivery is achieved at higher flow rates. This does not always equate to improved clinical efficacy. A pale arrow represents the intention that flow rates of a higher level may be beneficial for DPIs.

pMDIs

The pMDIs have been represented in a single block from 20-60 L/Min because the clinical studies support this. The pale arrow is pointed in the opposite direction for pMDIs, indicating that a slower rather than higher flow rate is beneficial in this type of inhaler.


The devices are individually represented on the laminated card.

Adaptors (specific device mimics)

A blue adaptor (example) is included in your pack. It has a code that corresponds to a particular device. Use in accordance with instructions.

These adaptors listed below can be used with the **In-Check DIAL G16** :

Turbohaler Mk II	- T	- 3109781
HandiHaler	- H	- 3109859
Diskhaler 4 Blister	- D4	- 3109788
Clickhaler	- C	- 3109786
Genuair	- G	- 3109861
Breezhaler	- B/F	- 3109783
Easyhaler	- EH	- 3109787
Accuhaler	- A/A/D	- 3109822

Set DIAL to the pMDI to use with device specific adaptor.

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In-Check DIAL G16	Ref. 3109300
In-Check Inhaler Assessment Kit	Ref. 3109760
In-Check Peak Inspiratory Flow Meter (Oral)	Ref. 3109750
Disposable one-way mouthpieces (Inspiratory, 5 x 20)	Ref. 3122061
Disposable one-way mouthpieces (Inspiratory, bulk, 200)	Ref. 3122063
Disposable Cardboard Mouthpieces (Adult, 5 x 20)	Ref. 3125030
Disposable Cardboard Mouthpieces (Adult, 500)	Ref. 3122003
Disposable Cardboard Mouthpieces (Child, 500)	Ref. 3123030
Universal Plastic Mouthpieces (Sterilisable, 15's)	Ref. 3103095
Adult to Paediatric Mouthpiece Adaptor	Ref. 3105027

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Bibliography for In-Check G16

(Please note that references for Clinical Efficacy of different inhaler systems are provided on the laminated card.)

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Capstick, T. G., & Clifton, I. J. (2014). Inhaler technique and training in people with chronic obstructive pulmonary disease and asthma. *Expert review of respiratory medicine*.

Banno M, Iбата H, Niimi T, Sato S, Matsushita R. The Usefulness of Inspiratory Flow Rate during Inhalation Corticosteroid Therapy in Asthma. *Respiration* 2009;78:387-392.

Janssens W, Vanden Brande P, Hardeman E, et al. Inspiratory flow rates at different levels of resistance in elderly COPD patients. *Eur Respir J*. 2008;31:78-83.

Zhang QL, Zheng JP, Pan WH, He H, Chen GL, An JY, Yuan LF, Luo DF. Peak inspiratory flow generated through different analogue dry powder inhalers in Shenzhen healthy preschool children. *Zhonghua Er Ke Za Zhi*. 2008 Feb;46(2): 98-103.

Raid A.M. R.A.M. Al-Showair, Walid Y. W.Y. Tarsin, Khaled H. K.H. Assi, Stanley B. S.B. Pearson, Henry Chrystyn. Can all patients with COPD use the correct inhalation flow with all inhalers and does training help? *Respiratory Medicine - November 2007 Vol. 101, Issue 11, Pages 2395-2401.*

Walid Y. Tarsin, Stanley B. Pearson, Khaled H. Assi, Henry Chrystyn. Emitted dose estimates from Seretide(R) Diskus(R) and Symbicort(R) Turbuhaler(R) following inhalation by severe asthmatics. *International Journal of Pharmaceutics*, Volume 316, Issues 1-2, 19 June 2006, Pages 131-137.

K.N. Desager, J. Geldhof, R. Claes, W. De Backer. Measurement of Inspiratory Flow Through Three Different Dry Powder Inhalation Devices Using In-Check™ in Children with Asthma. *Pediatric Asthma, Allergy & Immunology*. March 2006, Vol. 19, No. 1: 6-13.

Andrea S. Melani; Letizia S. Bracci; Marcello Rossi. Reduced Peak Inspiratory Effort through the Diskus and the Turbuhaler due to Mishandling is Common in Clinical Practice. *Clin Drug Invest.* 2005;25(8):543-549.

Amirav I, Newhouse MT, Mansour Y. Measurement of peak inspiratory flow with In-Check DIAL device to simulate low-resistance (Diskus) and high-resistance (Turbuhaler) dry powder inhalers in children with asthma. *Pediatr Pulmonol*. 2005 May;39(5):447-51.

Y. Adachi, Y.S. Adachi, T. Itazawa, J. Yamamoto, G. Murakami, T. Miyawaki. Measurement of peak inspiratory flow rates with an In-Check™ meter to identify preschool children's ability to use dry powder inhalers; Diskus and Turbuhaler. *The Journal of Allergy and Clinical Immunology - February 2004 (Vol. 113, Issue 2, Supplement, Page S114).*

Please Note

All trademarks and product names are the property of thier respective owners.

Handihaler is a registered trademark of Boehringer Ingelheim Pharma GmbH & Co. KG.

Easyhaler is a registered trademark of Orion Corporation.

NEXThaler is a registered trademark of Chiesi Farmaceutici S.P.A.

Twisthaler is a registered trademark of MSD International Holdings GmbH.

Turbuhaler and Turbohaler are registered trademarks of AstraZeneca AB.

Clickhaler is a trademark of Tianjin Kinnovata Pharmaceutical Company Limetid.

SpiroMax and Respiclick are registered trademarks of IVAX International B.V.

Genuair, Novolizer and Pressair are registered trademarks of AstraZeneca AB.

Ellipta is a registered trademark of GlaxoSmithKline group of companies.

Diskhaler is a registered trademark of GlaxoSmithKline group of companies.

Accuhaler and Discus are registered trademarks of the GlaxoSmithKline group of companies.

Breezhaler and Aerolizer are registered trademarks of Novartis AG.

Clement Clarke International Limited
Edinburgh Way
Harlow
Essex
CM20 2TT UK

Tele: +44 (0)1279 414969
Fax: +44 (0)1279 456300
email: resp@clement-clarke.com
www.clement-clarke.com