
AirFeel IMT

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Medical device for respiratory gymnastics

These operating instructions are provided for the device model P0124EM.

INTENDED USE

Medical device for training the inspiratory muscles. Prevention and/or therapy practices must be prescribed by professional medical personnel.

DIRECTIONS FOR USE

AirFeel IMT is intended for training inspiratory muscles of the respiratory system.

This product is not intended to diagnose, treat, cure, or prevent any disease.

Any changes to the device settings should be made accordingly to what has been prescribed by your doctor/physiotherapist.

Notes for medical personnel: The respiratory gymnastics workload should be based on the patient's respiratory strength, determined by measuring the maximum respiratory pressure (MIP) with a manometer. MIP is the maximum negative pressure that the patient can generate during the inhalation phase.

An initial pressure load of 30% of the measured MIP value is recommended, although it is possible that the initial load for some patients may need to be lower. The pressure load can be increased as the patient's respiratory capacity improves.

CONTRAINDICATIONS

AirFeel IMT is not suitable in the following situations:

- asthmatic patients with low perception of symptoms and/or who suffer from frequent and severe exacerbations;
- people with severe pneumothorax or who cough up large amounts of blood (massive hemoptysis).
- Patients whom have suffered from a ruptured eardrum or any other ear disease
- Patients with exceptionally high left ventricular end-diastolic volume and pressure
- Patients with worsening signs and symptoms of heart failure after IMT session
- Patients suffering of pulmonary hypertension
- Patients showing large bubbles on chest x-ray
- Patients with marked osteoporosis with history of rib fractures
- Patients with desaturation during or after IMT (<94%)

The following side effects may occur: Dizziness

Use of the device is not recommended immediately after lung surgery and in subjects with a history of pneumothorax who require blood pressure monitoring.

If during therapy you experience severe tiredness, shortness of breath or a significant increase in heart rate, stop treatment and contact your doctor.

Inform your doctor or therapist if you notice any side effects not listed in these instructions for use.

INTENDED USERS

The devices are intended to be used by professional medical personnel. The device can be used directly by the patient.

GROUP OF TARGET PATIENTS

Children, men and women of all ages. Children under 16 years of age may only use the device under the supervision of an adult.

OPERATING ENVIRONMENT

This device can be used in healthcare facilities such as hospitals, clinics, etc., or even at home.

WARNINGS

- Use the device only as intended. This medical device is not intended to be a life-saving device. Any other use is considered improper and may be dangerous. The manufacturer is not responsible for any improper use.
- This product must be used only as described in these instructions for use.
- The time required to pass from storage to operating conditions is approximately 2 hours.
- Please keep this manual carefully for future reference.
- When the device is used in combination with another medical device, the instructions for use of the latter must also be followed.
- The materials used in the device are biocompatible materials and comply with mandatory regulations, however it is not possible to completely exclude possible allergic reactions.
- The device does not administer drugs or substances.
- At the end of each use, after cleaning as per the HYGIENIC PREPARATION section, store the device in a dry place, protected from dust.
- It is recommended that you do not use the device if you suffer from cold, sinusitis or respiratory tract infection. In case of doubt, please consult your doctor.
- In case of cold, some users may experience mild ear discomfort during therapy, caused by inadequate equalization of pressure between the mouth and ears. If symptoms persist, consult your doctor.

Choking Hazard:

- Some components of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.

Risk of ineffectiveness of therapy:

- Use the device in the correct position (see instructions for use paragraph)
- Use original Flaem spare parts. Manufacturer declines any responsibility if non-original spare parts or accessories are used.

Risk of infection:

- Single patient use is recommended unless all HYGIENIC PREPARATION, including Sterilization, has been performed.

- Before and after each use, follow the hygienic preparation operations. Make sure that the device is not stored near other accessories or devices for different therapies (e.g. infusions).

FAILURE CASE STUDIES AND THEIR RESOLUTION

Problem	Cause	Remedy
The breathing effort is too heavy or too light	Incorrect respiratory load setting	Check the correct position of the indicator according to the medical prescription. If it is correct, consult your doctor/physiotherapist
The indicator is crooked	Incorrect assembly of the indicator on the internal screw	Disassemble the device and assemble it correctly.
When inhaling, no effort is felt.	Missing inhalation valve or sealing gasket	Check for their presence and, if necessary, reassemble the device with the missing parts.

If, after having checked the conditions described above, the appliance still does not function correctly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest you.

Find the list of all the Service Centers on the page <http://www.flaemnuova.it/it/info/associazione>

DISPOSAL

All components of the product can be disposed of with household waste unless this is prohibited by the disposal regulations in force in the respective member countries.

Packaging materials:

Product box

Product packaging bag

COMMUNICATION OF SERIOUS EVENTS

Please report any serious incidents involving this device to the manufacturer and the competent authority of the Member State in which you reside.

An event is considered serious if it causes or may cause, directly or indirectly, death or an unexpected and serious deterioration in the state of health of a person.

SYMBOLS PRESENT ON THE DEVICE OR PACKAGING

Medical Device

Manufacturer

Production date

CE marking ref. regulation
2017/745 EU and subsequent
updates.

Free from phthalates and
bisphenol

Model number

Attention

See instructions for use

Lot code

Unique identifier of the
device

Atmospheric pressure limits

Temperature limits

Humidity limits

TECHNICAL FEATURES

Model: P0124EM

Weight: 0.05 Kg approx.

MEDICAL DEVICES THAT CAN BE USED IN COMBINATION WITH THE DEVICE

The device can be used in combination with the following medical device: AirFeel respiratory gymnastics device, model P0920EM (not supplied) for IMT therapy and PEP / O-PEP therapy simultaneously.

Combination therapy must be prescribed by medical personnel.

ENVIRONMENTAL CONDITIONS

Operating conditions:

Room temperature

Between +10°C and +40°C

Relative humidity of the air

Between 10% and 95%

Atmospheric pressure

Between 69 KPa and 106 KPa

Storage and transportation conditions:

Room temperature

Between -25°C and +70°C

Relative humidity of the air

Between 10% and 95%

Atmospheric pressure

Between 69 KPa and 106 KPa

DURATION

Model: P0124EM

The expected average life is 5 years.

Warning: Do not remove the valve and the snap ring mounted on the valve holder (4). Do not remove the snap ring mounted on the main body (1).

Disinfection

Sanitize the components by choosing one of the methods described below:

method A: Sanitize the components under hot drinking water (about 40°C) with a mild dish detergent (non-abrasive).

method B: Sanitize the components in the dishwasher on a hot cycle (70°C).

method C: Sanitize the components by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly with hot drinking water (about 40°C).

Disinfection

After sanitizing the components, disinfect them by choosing one of the methods described below. Each method can be practiced for a limited number of times.

method A: Obtain an electrolytic chlor-oxidizing disinfectant (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution:

- Fill a container of a suitable size to contain all the individual components to be disinfected with a solution based on drinking water and disinfectant, respecting the proportions indicated on the disinfectant packaging.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Recover the disinfected components and rinse them thoroughly with warm drinking water.
- Dispose of the solution according to the instructions provided by the disinfectant manufacturer.

method B: Disinfect the components by boiling in water for 10 minutes; use demineralized or distilled water to avoid limescale deposits.

method C: Disinfect the components with a hot steam bottle vaporizer (not microwave). Perform the process by faithfully following the instructions of the vaporizer. For disinfection to be effective, choose a sterilizer with an operating cycle of at least 6 minutes. If you also want to perform Sterilization, skip to the paragraph Sterilization.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a jet of hot air (for example a hairdryer).

Sterilization

Equipment: Steam sterilizer with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in a sterile barrier system or packaging compliant with EN 11607. Insert the packaged components into the steam sterilizer. Perform the sterilization cycle following the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes.

Conservation: Store sterilized components as per the system instructions for use or in selected sterile barrier packaging.

The sterilization procedure has been validated in accordance with ISO 17665-1

At the end of each use, store the appliance, complete with accessories, in a dry place, protected from dust.

HYGIENIC PREPARATION IN THE HOME ENVIRONMENT

Disinfection

method A

✓

Method B

✓

method C

✓

Disinfection

method A

✓ MAX 300 TIMES

Method B

✓ MAX 300 TIMES

method C

✓ MAX 300 TIMES

HYGIENIC PREPARATION IN CLINICAL OR HOSPITAL ENVIRONMENT

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EN> The warranty terms here specified are valid only in Italy for Italian residents. In all other countries, the warranty will be provided by the local dealer that sold you the unit, in accordance with the applicable laws.

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All rights reserved
FLAEM NEW SpA