Cenni accademici sull'Aerosol Terapia / Academic short accounts outlines on Aerosol Therapy

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Aerosol therapy is a method that is occupying an increasingly important role in the treatment of numerous pathologies of the nose, throat, bronchial tubes and lungs as it allows drugs to be directly administered into the respiratory system.

An aerosol is a suspension of solid particles and/or drops of liquid that penetrate into the airways through the mouth and nose, coming into contact with the larynx, pharynx and trachea, and eventually reaching the lungs: here they will come into contact with the bronchial mucosa and air sacs, following precise laws of kinetics such as inertia, gravity and diffusion.

How and where these particles will settle depends on their size.

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




Over 50 years of specific experience have allowed Flaem to realize this new and innovative patented method of performing aerosol treatments targeted to the patient's airways to be treated.

The peculiar characteristic of RF9 nebulizer is the Particle Size Selector System, which allows you to select the desired particle's size.

The 4 mode options, in fact, can dispense different particle's sizes: **Without Selector** larger than 7 µm and high speed for **upper airways**, larger than 5 µm (with **Selector in Position 1**) for the treatment of **tracheobronchial region**, particles between 3 and 4 µm (with **Selector in Position 2**) for the **pulmonary region** and between 2 and 3 µm (with **Selector in Position 3**) for deposition in the **alveoli of peripheral region**. The **valve system** optimizes the delivery of the nebulized medication to the patient and limits its dispersion in the environment.

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Valve System: • Optimize Delivery
• Reduce Dispersion

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A slower flow of medication combined with the synchronised opening and closing of the valve (as shown in pictures A and B) limits the dispersion of medication into the environment.

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The simple extra-rotation of the knob allows the selector group to be removed allowing an easy and effective cleaning in all the parts of the nebulizer.

Gamma accessori completa Complete range of accessories

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
Made with non-allergic materials, phthalates & BPA free, tested in compliance with European standard EN 10993-1 "Biological evaluation of medical devices" and with 93/42/EEC European Medical Device Directive.

The RF9 nebulizer can be used with the mouthpiece, nasal prongs, the Soft Touch adult / child masks or neonatal silicone masks. The nebulizer and the accessories should be replaced every 6 months when frequently used in order to always ensure the maximum therapeutic effect.

The RF9 nebulizer and accessories can be cleaned according to the following procedures: sanitization, disinfection, sterilization.

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The new Soft Touch Masks feature an adaptive edge made of a soft biocompatible material ensuring an excellent adherence to the face and a Dispersion Limiting Vent. These innovations, which make them stand out from the rest, allow increased sedimentation of medication in the lungs and reduced dispersion into the air.

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18733AQ © Copyright by FLAEM NUOVA S.p.A. 2019 - Materiale divulgativo riservato alle categorie professionali specializzate. FLAEM NUOVA si riserva di apportare modifiche tecniche e funzionali al prodotto della presente scheda senza alcun preavviso. FLAEM NUOVA reserves the right to make technical and functional alterations to this product without prior notice.