Mr. CAT

The funny nebulizer

PHTHALATES & BPA-FREE

MANUALE ISTRUZIONI D'USO

INSTRUCTIONS FOR USE MANUAL

MODE D'EMPLOI

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MANUAL DE INSTRUCCIONES

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AEROSOL THERAPY NEBULISER

Congratulations on your purchase and thank you for choosing our unit.

Our goal is to fully satisfy customers by offering them cutting-edge systems for the treatment of respiratory tract ailments. **Carefully read this instruction manual and store it for future reference. Only use the accessories as described in this manual. This is a medical device intended to nebulise and administer medications prescribed or recommended by a doctor.** Visit our Internet site www.flaem.it to view the whole range of Flaem products.

SUPPLIED ACCESSORIES

- 1 Child bi-material face mask
- 2 Mouthpiece
- 3 Kitten Clip
- 4 RF6 Basic² Nebuliser
 - 4.1 - Upper part
 - 4.2 - Complete nozzle
 - 4.3 - Lower part

WARNINGS

- This is a medical device (compliant with Dir. 93/42/EEC) and may only be used with medication prescribed or recommended by a doctor. Before using the accessories, consult the instruction manual of the device. It is important that the patient reads and understands the information on how to use the device and the nebuliser. Contact your retailer or nearest service centre for any questions.
- When used very frequently, we recommend that the nebuliser be replaced every 6 months (or earlier if it becomes clogged), in order to guarantee the maximum therapeutic effect.
- Children and disabled persons should always use the nebuliser under strict supervision of an adult who has read this manual.
- The nebuliser contains small parts that may be swallowed by children; keep therefore the nebuliser out of the reach of children.

INSTRUCTIONS FOR USE


Before each use, thoroughly wash your hands and clean the nebuliser and the accessories, as described in the paragraph "SANITIZATION, STERILIZATION AND DISINFECTION". We recommend you use the nebuliser and the accessories for personal use only to prevent any risk of infection.

1. Open the nebuliser by rotating the upper part (4.1) anti-clockwise.
2. Pour the medication prescribed by the doctor into the lower part (4.3). Close the nebuliser by rotating the upper part (4.1) clockwise.
3. Assemble the accessories as indicated in the 'Connection diagram'. Make sure you are seating comfortably and hold the nebuliser in your hands. Put the mouthpiece in your mouth, or instead use the nasal inhaler, or the face mask.
4. Start the device (see its instruction manual). Inhale gently and deeply; after inspiration we recommend that you hold your breath for a moment in order to allow inhaled aerosol droplets to deposit. Then, exhale slowly.
5. When the treatment is finished, switch the device off.
6. It may happen that after a treatment session visible humidity deposits form in the tube: disconnect the tube from the nebuliser and let it dry using the compressor air flow, as this prevents mould from growing in the tube.


HOW TO USE THE "RF6 Basic?" DUAL SPEED NEBULISER

Due to the geometry inside the nebuliser it delivers a particle size providing effective treatment of the entire respiratory tract.

The RF6 Basic? nebuliser is especially suitable for paediatric use because the non-spill system prevents the medicine from leaking out even if it accidentally tips over (tested to work with up to 5 ml of product, which is the average dose of the most common therapeutic applications). The nebuliser can deliver the medicine in 2 ways:

"standard mode";
for excellent delivery
and minimum
medication waste.

**"high speed
mode";** for faster
administration time.

SANITIZATION, STERILIZATION AND DISINFECTION

Before undertaking any cleaning operation, switch off and unplug the device.

Before and after each use, sanitize the nebuliser and the accessories using one of the methods described below. A personal use of the nebulizer and its accessories is recommended in order to prevent any risk of infection. Open the nebuliser by turning the upper part (4.1) anticlockwise. Detach the nozzle from the upper part (4.1) and remove it by pressing as indicated by the arrows in the 'Connection diagram', section 4.2.

SANITIZING

method A: Sanitize the accessories 1-2-3-4.1-4.2-4.3 using warm tap water (40°C approx.) with mild (non-abrasive) washing-up liquid or wash them in a dishwasher, running a hot cycle.

method B: Sanitize the accessories 1-2-3-4.1-4.2-4.3 by soaking them in a solution of 60% water and 40% white vinegar. Rinse thoroughly with warm tap water (40°C approx.).

method C: Sanitize the accessories 2-3-4.1-4.2-4.3 by boiling them in water for 20 minutes (except for the connection tube) preferably using de-mineralized or distilled water to prevent lime scale deposits.

After sanitizing the accessories shake them vigorously and lay them on a paper towel or dry them with a jet of warm air (e.g. a hair dryer).

STERILIZATION

The accessories subject to sterilization are 2-3-4.1-4.2-4.3.

The sterilization procedure described in this paragraph is validated according to ISO 17665-1 and is effective on all components only if the instructions are followed precisely and if the components have been previously sanitized.

Equipment: Steam sterilizer fractionated vacuum and overpressure compliant with standard EN 13060.

Execution: Pack each component to be treated in a sterile barrier system or packaging in compliance with standard EN 11607. Place the packaged components in the steam sterilizer. Run the sterilization cycle according to the equipment instructions for use at a temperature of 134°C for 10 minutes.

Storage: Store the sterilized components as stated in the instructions for use of the chosen sterile barrier system or packaging.

DISINFECTION -The accessories that can be disinfected are (1-2-3-4.1-4.2-4.3)

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).
- Dispose of the solution as recommended by the disinfectant manufacturer.

TECHNICAL SPECIFICATIONS**RF6 Basic² Nebuliser**

Operating conditions:

Temperature: min 10°C; max 40°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

In compliance with 93/42/EEC Directive

Minimum fill volume: 2 ml

Maximum fill volume: 8 ml

Storage conditions:

Temperature: min -25°C; max 70°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

APPLIED PARTS

The applied parts of type B are:
patient accessories (1,2)

FLAEM compressor range	Operating pressure	MMAD (µm) ⁽²⁾		BREATHABLE FRACTION < 5 µm (FPF) % ⁽²⁾		DELIVERY ml/min ⁽¹⁾	
		standard mode	high speed mode	standard mode	high speed mode	standard mode	high speed mode
F400	0,65	3,8	4,00	62	63	0,25	0,32
F700	0,85	3,6	3,6	66	67	0,32	0,42

⁽¹⁾ Data shown in accordance with Harmonized International Standard EN 12407-2. ⁽²⁾ In vitro characterization carried out at the Biomed Research Center Co. G. on behalf of IGV Research (G. Products GmbH – Germany) in compliance with the new European Standard for aerosol therapy (in its Standard EN 13244-1). Further details are available upon request.

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FLAEM NUOVA S.p.A. - Via Colli Storici, 221