SINCE 1966
MADE IN ITALY

Plus RF7

by **FLAEM**®

Manuale istruzioni d'uso del nebulizzatore "RF7 Dual Speed Plus" con selettore di velocità e sistema valvolare
Instruction manual of the nebulizer "RF7 Dual Speed Plus" with speed selector and valve system
Mode d'emploi du nébuliseur "RF7 Dual Speed Plus" avec sélecteur de vitesse et système à valve
Gebruiksaanwijzing voor de verstuiver "RF7 Dual Speed Plus" met snelheidsregelaar en ventielsysteem
Bedieningsanleitung des Zerstäubers "RF7 Dual Speed Plus" mit Geschwindigkeitswahlschalter und Ventilsystem

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AEROSOL THERAPY NEBULIZER

Congratulations on your purchase and thank you for choosing our unit. Our goal is to fully satisfy customers by offering them cutting-edge systems for the treatment of respiratory tract ailments.

Carefully read this instruction manual and store it for future reference. Only use the accessories as described in this manual. This is a medical device intended to nebulize and administer medications prescribed or recommended by a doctor. Visit our Internet site www.flaem.it to view the whole range of Flaem products.

ENGLISH

ACCESSORIES THAT CAN BE USED

WITH THIS NEBULIZER

- 1 Pediatric nasal prong
- 2 Adult nasal prong
- 3 Pediatric face mask
- 4 Adult face mask
- 5 Mouthpiece with valve
- 6 Non-invasive nasal prong
- 7 Baby mask set (see its instruction manual)
- 8 Nebulization manual control
- 9 Connector
10. 1 or 2 m Connection hose

NEBULIZER

- 11 RF7 Dual Speed Plus Nebulizer
- 11.1 - Lower part
- 11.2 - Nozzle
- 11.3 - Upper part
- 11.4 - Speed selector with valve system

WARNINGS

- This is a medical device (compliant with Dir. 93/42/EEC) and may only be used with medication prescribed or recommended by a doctor. Before using the accessories, consult the instruction manual of the device. It is important that the patient reads and understands the information on how to use the device and the nebulizer. Contact your retailer or nearest service center for any questions.
- When used very frequently, we recommend that the nebuliser be replaced every 6 months (or earlier if it becomes clogged), in order to guarantee the maximum therapeutic effect.
- Children and disabled persons should always use the nebulizer under strict supervision of an adult who has read this manual.
- The nebulizer contains small parts that may be swallowed by children; keep therefore the nebulizer out of the reach of children

OPERATING INSTRUCTIONS

Before each utilisation, wash your hands accurately and clean the device as described in the paragraph "CLEANING, SANITISATION, DISINFECTION AND STERILISATION".

The nebuliser cup and accessories are for personal use only to avoid potential risks of transmitting infection.

During application it is advisable to suitably protect yourself from any dripping.

1. Insert the nozzle (11.2) in the upper part (11.3) pressing as shown by the 2 arrows in the "Connection diagram" in point 11. Insert the Speed selector with valve (11.4) in the upper part (11.3) as shown in the "Connection diagram" in point 11.
2. Pour the medication prescribed by the doctor into the lower part (11.1). Close the nebuliser by turning the upper part (11.3) clockwise.
3. Connect the accessories as shown in the "Assembly diagram". Sit down in a comfortable position holding the nebulizer in your hand, place the mouthpiece over your mouth or use the nasal prong, the face mask, or the Baby mask set.
4. Start the device (see its instruction manual). Inhale gently and deeply; after inspiration we recommend to hold your breath for a moment in order to allow inhaled aerosol droplets to deposit. Then, exhale slowly.
5. Upon completing application, switch off the device and unplug it.



ATTENTION: If after the therapy session an evident deposit of moisture forms within the pipe (10), detach the pipe from the nebuliser and dry it with the very ventilation from the compressor; this operation prevents possible blooms of mould inside the pipe.

HOW TO USE THE "RF7 DUAL SPEED PLUS" NEBULIZER WITH SPEED SELECTOR AND VALVE


The nebulizer is professional fast and suitable for the administration of all types of drugs, including the most expensive ones, also in patients with chronic diseases. The geometry of the internal ducts of the RF7 Dual Speed Plus nebulizer ensures the ideal particle size providing effective treatment right down to the lower airways.

To make inhalation therapy faster, press with your finger on the MAX marking on the speed selector with valve system (11.4). This configuration is also suited for the use of nasal prongs 1, 2.



To make inhalation therapy more effective press with your finger on the speed selector with valve system (11.4) opposite the MAX marking. In this way the valve systems in the nebulizer, in the mouthpiece and in the mask will ensure optimal medication delivery to the lower respiratory tract, limiting dispersion of medication into the environment.

SOFTTOUCH FACE MASKS


The new SoftTouch Masks feature an adaptive edge made of a soft biocompatible material ensuring an excellent adherence to the face and a Dispersion Limiting Vent. These innovations, which make them stand out from the rest, allow increased sedimentation of medication in the lungs and reduced dispersion into the air.


USE OF NEBULIZATION MANUAL CONTROL

To achieve continuous nebulisation action you should not use the manual nebulisation control (8), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (8) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited,

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

CLEANING, SANITISATION, DISINFECTION, STERILISATION

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

DEVICE AND TUBING EXTERIOR

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

ACCESSORIES

Open the nebuliser by turning the upper part (11.3) anticlockwise, remove the nozzle (11.2) from the upper part (11.3) by pressing as shown by the 2 arrows in the "Connection diagram" in point 11.

Then proceed according to the following instructions.

CLEANING AT HOME - SANITISATION AND DISINFECTION

SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **1-2-3-4-5-6-8-9-11.1-11.2-11.3-11.4** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **1-2-3-4-5-6-8-9-11.1-11.2-11.3-11.4** in the dishwasher with a hot cycle.

method C: Sanitise accessories **1-2-3-4-5-6-8-9-11.1-11.2-11.3-11.4** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40°C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

method A: Accessories **1-2-3-4-5-6-8-9-11.1-11.2-11.3-11.4** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

method B: Sanitise the accessories **1-2-5-6-8-9-11.1-11.2-11.3-11.4** by boiling them in water for 20 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: Sanitise the accessories **1-2-5-6-8-9-11.1-11.2-11.3-11.4** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

CLEANING IN A CLINICAL OR HOSPITAL SETTING - DISINFECTION AND STERILISATION

Before disinfection or sterilisation, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **1-2-3-4-5-6-8-9-11.1-11.2-11.3-11.4** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **1-2-3-4-5-6-8-9-11.1-11.2-11.3-11.4** in the dishwasher with

a hot cycle.

DISINFECTION

Accessories **1-2-3-4-5-6-8-9-11.1-11.2-11.3-11.4**. can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

If you want to also perform the STERILISATION, jump to the STERILISATION paragraph.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

STERILISATION

Accessories **1-2-5-6-9-11.1-11.2-11.3-11.4** can be sterilised.

Device: Fractionated vacuum overpressure steam steriliser in accordance with EN 13060.

Implementation: Wrap every single part to be treated with a sterile barrier system or packaging in accordance with Norm EN 11607. Place the packed components in the steam steriliser. Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134°C and a time of 10 minutes first.

Storage: Store the sterilised parts as per the instructions for use of either the sterile barrier system or packaging.

The sterilisation procedure is validated in its conformity to ISO 17665-1.

For the **CLEANING, SANITISATION, DISINFECTION, STERILISATION** of accessories 7 follow the instructions in the instruction manual.

TECHNICAL SPECIFICATIONS

RF7 Dual Speed Plus Nebulizer

Operating conditions:

Temperature: min 10°C; max 40°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

In compliance with 93/42/EEC Directive

Minimum fill volume: 2 ml

Maximum fill volume: 8 ml

Storage conditions:

Temperature: min -25°C; max 70°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

APPLIED PARTS

The applied parts of type BF are:

patient accessories (1,2,3,4,5,6,7)

FLAEM compressor range	Pressure	⁽²⁾ MMAD (µm)		⁽²⁾ BREATHABLE FRACTION < 5 µm (FPF) %		⁽¹⁾ DELIVERY ml/min	
		Min With valve	Max	Min With valve	Max	Min With valve	Max
Portable	0.75 bar	2.95	2.44	74.7	80.7	0.18	0.42
F700/F1000	1.00 bar	2.72	2.53	77.8	79.6	0.23	0.53
F2000	1.40 bar	2.21	2.38	84.2	81.5	0.29	0.65

⁽¹⁾Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5.

In vitro testing⁽¹⁾: certified by TÜV Rheinland LGA Products GmbH – Germany in compliance with the European Standard EN 13544-1 ANNEX CC for nebulizing systems. More details are available on request.

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