RF9

— 50 —
years of italian excellence

**PHTHALATES
& BPA FREE**

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AEROSOL THERAPY NEBULIZER

Congratulations on your purchase and thank you for choosing our unit. Our goal is to fully satisfy customers by offering them cutting-edge systems for the treatment of respiratory tract ailments.

Carefully read this instruction manual and store it for future reference. Only use the accessories as described in this manual. This is a medical device intended to nebulize and administer medications prescribed or recommended by a doctor. Visit our

Internet site www.flaem.it to view the whole range of Flaem products.

ACCESSORIES THAT CAN BE USED

WITH THIS NEBULIZER

- 1 - Pediatric nasal prong
- 2 - Adult nasal prong
- 3 - Pediatric face mask
- 4 - Adult face mask
- 5 - Mouthpiece with valve
 - 5.1 - Exhalation valve
- 6 - Non-invasive nasal prong
- 7 - Baby mask set (see its instruction manual)
- 8 - Nebulization manual control
- 9 - Connector
- 10 - 1 or 2 m Connection hose

NEBULIZER

- 11 - RF9 Nebulizer
 - 11.1 - Complete particle size regulator
 - 11.1a - Regulator
 - 11.1b - Valve
 - 11.1c - Selector
 - 11.2 - Upper part
 - 11.3 - Nozzle
 - 11.4 - Lower part



WARNINGS

- This is a medical device (compliant with Dir. 93/42/EEC) and may only be used with medication prescribed or recommended by a doctor. Before using the accessories, consult the instruction manual of the device. It is important that the patient reads and understands the information on how to use the device and the nebulizer. Contact your retailer or nearest service center for any questions.
- The expected medical life of the accessories is 1 year. It is, however, advisable to replace the nebuliser cup every 6 months in the event of intense use (or earlier if the cup is obstructed) to always guarantee maximum therapeutic efficacy.
- Children and disabled persons should always use the nebulizer under strict supervision of an adult who has read this manual.
- The nebulizer contains small parts that may be swallowed by children; keep therefore the nebulizer out of the reach of children.


ASSEMBLING THE NEBULISER

Follow the assembly instructions below.


1. Assemble the nozzle (11.3) on the upper part (11.2).


2. Place the valve (11.1b) with the flat side down onto the particle size regulator (11.1a) by aligning the hole with the position ribs, then screw the selector (11.1c) onto the threaded pin of the regulator all the way down. Rotate the assembled regulator to make sure that it turns freely and that the selector moves up and down between positions.

3. Insert and press the previously assembled particle size regulator (11.1) all the way down, keeping it in a horizontal position, in the upper part (11.2), and matching the partitions of the selector (11.1c) with the groove in the upper part (11.2).





4. Insert the assembled upper part (11.1 + 11.2 + 11.13) into the lower part (11.4) and hook it by turning it clockwise.


HOW TO USE THE RF9 NEBULISER WITH PARTICLE SIZE REGULATOR AND VALVE SYSTEM

It is professional, quick, recommended for administering all types of medicinal products, even in patients with chronic diseases. The valve system optimises the dispensing of the nebulised medicinal product to the patient and limits its dispersion in the environment. Moreover, thanks to the particle size regulator, the size of nebulised particles can be selected based on the respiratory system area where the medication is to be delivered. Indeed, with the 4 available modes, the bowl can dispense particles of different sizes:


RESPIRATORY SYSTEM	UPPER AIRWAYS	LOWER AIRWAYS		
	Oropharyngeal tract	Tracheobronchial tract	Pulmonary tract	Peripheral tract
CONFIGURATION BOWL	Without Selector	Selector in position 1	Selector in position 2	Selector in position 3
				
FLAEM compressor range - F400				
MMAD (µm)	7.71 ⁽³⁾	5.12 ⁽¹⁾	3.67 ⁽¹⁾	2.49 ⁽¹⁾
Output rate (ml/min)	0.53 ⁽²⁾	0.36 ⁽²⁾	0.30 ⁽²⁾	0.23 ⁽²⁾
Breathable fraction < 5 µm (FPF)	26.7 ⁽³⁾	48.8 ⁽¹⁾	63.4 ⁽¹⁾	83.1 ⁽¹⁾
FLAEM compressor range - F700 / F1000				
MMAD (µm)	6.81 ⁽³⁾	4.56 ⁽³⁾	3.94 ⁽¹⁾	2.70 ⁽¹⁾
Output rate (ml/min)	0.66 ⁽²⁾	0.45 ⁽³⁾	0.41 ⁽²⁾	0.27 ⁽²⁾
Breathable fraction < 5 µm (FPF)	31.3 ⁽³⁾	54.5 ⁽¹⁾	62.0 ⁽¹⁾	80.0 ⁽¹⁾
FLAEM compressor range - F2000				
MMAD (µm)	5.80 ⁽³⁾	4.67 ⁽³⁾	3.62 ⁽¹⁾	2.56 ⁽¹⁾
Output rate (ml/min)	0.75 ⁽²⁾	0.54 ⁽³⁾	0.51 ⁽²⁾	0.38 ⁽²⁾
Breathable fraction < 5 µm (FPF)	40.3 ⁽³⁾	54.4 ⁽¹⁾	66.3 ⁽¹⁾	82.2 ⁽¹⁾
(1) In vitro characterisation certified by: TÜV Rheinland LGA Products GmbH – Germany in compliance with the new European Standard for Aerosol therapy apparatuses EN 13544-1. (2) Data measured according to Flaem I29-P07.5 internal procedure Dispensing speed values may vary according to the patient's respiratory capacity. More details are available on request. (3) Data calculated on the basis of the values obtained with Malvern Spraytec, and compared with the values of the tests performed at the University of Parma.				

Select the desired position by means of the particle size regulator 11.1.

By turning the regulator clockwise, the selector (11.1c) moves and closes the opening, while it opens if the regulator is turned anti-clockwise.

To use the bowl without selector, firmly grasp the bowl body and use your other hand to rotate the particle size regulator (11.1) clockwise until it is ejected from the upper part (11.2).

INSTRUCTIONS FOR USE


Before each use, clean hands thoroughly and clean the device as described in the section on "CLEANING SANITISATION DISINFECTION". During use, it is advisable to protect yourself from any dripping. It is recommended that each person use their own nebulizer cup and accessories to prevent risk of infection due to contamination. This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.

1. Pour the medication prescribed by the physician into the lower part (11.4). Insert the upper part(11.1+11.2+11.3) of the bowl into the lower part (11.4), then close the bowl by turning the upper part (11.2) clockwise. Close the nebuliser by turning the upper part (11.2) clockwise.
2. Connect the accessories as shown in the "Assembly diagram". Sit down in a comfortable position holding the nebulizer in your hand, place the mouthpiece over your mouth or use the nasal prong, the face mask, or the Baby mask set.
3. Set the bowl selector to the desired particle size as described under "How to use the RF9 bowl".



4. Start the device (see its instruction manual). Inhale gently and deeply; after inspiration we recommend to hold your breath for a moment in order to allow inhaled aerosol droplets to deposit. Then, exhale slowly.
5. Upon completing application, switch off the device and unplug it.

ATTENTION: If after the therapy session an evident deposit of moisture forms within the pipe (10), detach the pipe from the nebuliser and dry it with the very ventilation from the compressor; this operation prevents possible blooms of mould inside the pipe.

SOFTTOUCH FACE MASKS


The new SoftTouch Masks feature an adaptive edge made of a soft biocompatible material ensuring an excellent adherence to the face and a Dispersion Limiting Vent. These innovations, which make them stand out from the rest, allow increased sedimentation of medication in the lungs and reduced dispersion into the air.


USE OF NEBULIZATION MANUAL CONTROL

To achieve continuous nebulisation action you should not use the manual nebulisation control (8), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (8) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited,

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

CLEANING SANITISATION DISINFECTION

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

DEVICE AND TUBING EXTERIOR

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

ACCESSORIES

To disassemble the bowl, firmly grasp the bowl body and then use your other hand to rotate the particle size regulator (11.1) clockwise until it is ejected from the upper part (11.2) as a result of unscrewing. The components can then be disassembled (regulator 11.1a, valve 11.1b, selector 11.1c and upper part 11.2) to be cleaned. The nozzle (11.3) is interference fitted in the upper part (11.2) and can be extracted by simply pulling it outwards. Then proceed according to the following instructions.

SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: sanitise accessories **1-2-3-4-5-6-8-9-11.1a-11.1b-11.1c-11.2-11.3-11.4** under potable hot water (approximately 40 °C) with a gentle, non abrasive dish detergent.

method B: sanitise accessories **1-2-3-4-5-6-8-9-11.1a-11.1b-11.1c-11.2-11.3-11.4** in the dishwasher with a hot cycle.

method C: sanitise accessories **1-2-3-4-5-6-8-9-11.1a-11.1b-11.1c-11.2-11.3-11.4** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40 °C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

method A: Accessories **1-2-3-4-5-6-8-9-11.1a-11.1b-11.1c-11.2-11.3-11.4** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

method B: sanitise the accessories **1-2-3-4-5-6-8-9-11.1a-11.1b-11.1c-11.2-11.3-11.4** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: sanitise the accessories **1-2-3-4-5-6-8-9-11.1a-11.1b-11.1c-11.2-11.3-11.4** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

At the end of each use store the device complete with accessories in a dry place away from dust.

For the **CLEANING SANITISATION DISINFECTION** of accessories 7 follow the instructions in the instruction manual.

TECHNICAL SPECIFICATIONS

RF9 Nebulizer

Operating conditions:

Temperature: min 10°C; max 40°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

In compliance with 93/42/EEC Directive

Minimum fill volume: 2 ml

Maximum fill volume: 8 ml

Storage conditions:

Temperature: min -25°C; max 70°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

APPLIED PARTS

The applied parts of type BF are:

patient accessories (1, 2, 3, 4, 5, 6, 7)

SYMBOLS

 EC Marking medical ref. Dir 93/42 EEC and subsequent updates	 Check the instructions for use	 In compliance with: European Standard EN 10993-1 "Biological Evaluation of medical devices" and European Directive 93/42/EEC "Medical Devices." Phthalate-free. In compliance with: Reg. (EC) no. 1907/2006
 Minimum and maximum room temperature	 Manufacturer	
 Minimum and maximum atmospheric pressure	 Minimum and maximum air moisture	

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.