KONICA MINOLTA

PULSE OXIMETER

PULSOX[®]-1

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1. Names of Parts

- A**
- 1 Clip** Squeeze the clip with fingers to open the finger insertion part.
 - 2 Display** Indicates the oxygen saturation, pulse rate, pulse level, error messages, read-out hold status and/or battery level alarm.
 - 3 Backlight sensor** Detects the brightness of the surroundings and if required, turns on the backlight.
 - 4 Display OFF button** Pressing this button during the readout-hold period turns off the LCD and power.
 - 5 Finger insertion part**
 - 6 Finger holder**
 - 7 Protective cover** Protects the instrument body against shock and shuts out ambient light.
- B**
- 8 Battery cover**
 - 9 Rating nameplate**
 - 10 AAA-size alkaline battery (LR03)**
 - 11 Serial number**
 - 12 Strap eyelet** Attach the neck strap to this hook.
- C**
- 1 Hold indicator**
 - 2 Pulse rate (P.R.)**
 - 3 Battery indication** The battery mark comes on when the battery is low.
 - 4 Pulse level meter**
 - 5 Oxygen saturation value (SpO₂)**

Calculation Method for Displayed Values

Pulse level meter

The pulse level is expressed as the ratio of the constant amount of transmitted light to the variable amount (not normalized), as shown in the following formula. It is indicated in eight levels (0 to 9%; full scale for 9% or higher).

$$\text{Pulse level (\%)} = \frac{\text{Transmitted light amount that varies due to the pulse (AC)}}{\text{Transmitted light amount that does not vary due to the pulse (DC)}} \times 100$$

Oxygen saturation (SpO₂)

Updates the moving average of the oxygen saturation over the last 12 seconds at one-second intervals.

Pulse rate (P.R.)

Updates the moving average of the last 8 pulse rates at one-second intervals.

2. Operating Method

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1 Squeeze the clip with fingers to open the finger insertion part.

2 Insert a fingertip until the tip is in contact with the guide in the finger holder; then release the clip.

- The power will turn on and measurement starts.

Tips for accurate measurement

Failure to properly insert a fingertip into the finger holder may result in measurement errors. Follow the instructions in the figures below.

- Position of the fingertip (E)
- Do not insert a finger into the finger holder upside down. (F)

During measurement, softly hold the instrument to avoid unstable measurement. See the following figures. (Do not apply excessive pressure on the instrument nor hold it with excessive force).

- Do not swing the instrument. (G)
- Hold the instrument softly. (H)
- Put the hand and instrument on a table. (I)
- Grip the sides of the instrument lightly by fingers. (J)
- Hold the sides of the instrument lightly by the adjacent two fingers when viewing the readouts for yourself in such a position as shown in the figure below. (K)
- Avoid placing the instrument at a position that is much higher or lower than the height of the heart during measurement. Doing so may result in measurement errors. (L)

Cold fingers will provide very weak pulse signals due to low pulse level and low perfusion.

This may cause the following screens to appear.

In such cases, reinsert the finger into the instrument, rub or warm the finger, or use another finger so that at least two segments of the pulse level meter blink.

Possible causes of measurement errors

Possible causes of measurement errors include the following:

The measurement practice is not correct:

- The finger is not inserted into the finger holder properly.
- Motion of the patient or instrument is excessive.
- Excessive pressure or a grasping force is applied to the instrument.

Environmental conditions are not correct:

- Surrounding light (light from incandescent, fluorescent or infrared lamps, or direct sunlight) is too bright.
- The instrument suffers electromagnetic interference from neighboring electronic equipment (medical devices or electrical apparatus).
- A wireless device such as an FM transmitter is activated in the immediate vicinity of the instrument.
- Blood circulation to the fingertip is insufficient due to pressure on the hand or finger or peripheral circulatory failure, or the fingertip suffers blood stasis.
- A blood pressure cuff is put around the arm on the same side as the finger in question.

The patient is problematic:

- The blood of the patient contains a lot of abnormal hemoglobin such as carbon monoxide hemoglobin (HbCO), methemoglobin or sulfhemoglobin.
- The blood of the patient contains coloring matters such as indocyanine green or methylene blue.
- The patient suffers severe anemia.
- The nail of the finger in question has nail polish or other treatment.

Ending the measurement

Open the finger holder and pull the finger out of it.

- The readouts on the LCD are held for 1 minute; then the power of the instrument turns off automatically.

To clear the readouts being held;

- Press the Display OFF button located above the finger insertion part and hold it in for at least 1 second. The readouts are cleared and the power turns off.

3. Error Messages

The instrument displays error messages when an error occurs. Check and respond to them correctly.

Error Messages	Cause	Solution
 Flashing	Insufficient light for measurement.	<ol style="list-style-type: none">1 Check that the finger holder is properly attached to the patient.2 If the part being measured is thick, insufficient light will be transmitted. Use another, thinner finger.3 Clean the light-source and sensor sections.4 If this message persists, the instrument has likely failed. *1
 Flashing	Pulse signal required for measurement is not being received.	<ol style="list-style-type: none">1 Properly insert a fingertip into the finger holder.2 This is likely caused by insufficient blood circulation to the fingertip. Rub or warm the fingertip or use another finger.
 Flashing	The pulse is too strong for correct measurement.	<ol style="list-style-type: none">1 Properly insert a fingertip into the finger holder.2 This is likely caused by excessive blood circulation to the fingertip. Take a rest for a while and retry again or use another finger.3 If this message persists, the instrument has likely failed. *1

*1 Contact the nearest authorized service facility.

Error Messages	Cause	Solution
 <p>Readout and "A" alternately flashing</p>	<p>The measurement results are unreliable due to motion of the patient during measurement. The stated accuracy of measurement is invalid.</p>	<ol style="list-style-type: none"> 1 Properly insert a fingertip into the finger holder. Avoid moving the finger in question. 2 Warm the fingertip and/or shut out surrounding light or electromagnetic noise.
 <p>Flashing</p>	<p>The pulse rate exceeds the upper limit (230 bpm) of the measurement range.</p>	<p>Measurement cannot be performed when the upper limit of the measurement range is exceeded.</p>
 <p>Flashing</p>	<p>The pulse rate is below the lower limit (30 bpm) of the measurement range.</p>	<p>Measurement cannot be performed when the lower limit of the measurement range is not reached.</p>
 <p>Oxygen saturation readout flashing</p>	<p>The oxygen saturation (SpO₂) is below the lower limit (70%) of the guaranteed accuracy range. The stated accuracy of oxygen saturation measurements is invalid.</p>	<p>This message indicates that the oxygen saturation is out of the range in which the stated accuracy of measurements is valid.</p>

Error Messages	Cause	Solution
 Pulse rate readout flashing	The pulse is weak and pulse signals required for measurement are not obtained. The stated accuracy of pulse rate measurements is invalid.	<ol style="list-style-type: none"> 1 Properly insert a fingertip into the finger holder. 2 This is likely caused by insufficient blood circulation to the fingertip. Rub or warm the fingertip or use another finger.
 Lit	The battery is low.	The battery will run out in less than 6 hours. Replace the battery as early as possible.
 Flashing for about 2 sec. then power off	The battery runs out.	Replace the battery.
  Error indication	The measuring circuitry of the instrument has failed. Measurement is disabled.	If this message persists after the battery is replaced, the instrument has likely failed. *1

*1 Contact the nearest authorized service facility.

4. Safety Precautions

To ensure correct use of this instrument, read the following points carefully and adhere to them. After you have read this manual, keep it in a safe place where it can be referred to anytime a question arises.

Warning

(Failure to adhere to the following points may result in death or serious injury.)

- This device is designed to measure the oxygen saturation level (SpO₂) and the pulse rate. It should not be used as a warning device to detect cessation of breathing or irregular breathing during sleep.
- This instrument has no alarm function. Do not use this instrument if an alarm function is necessary.
- Do not use the instrument in places where flammable or combustible gases including anesthetic gases, high-purity oxygen, or gasoline vapor are present. Doing so may cause a fire.
- Do not disassemble or modify the instrument. Doing so may cause a fire or electric shock.
- The instrument should not be operated if it is damaged, or smoke or odd smells occur. Doing so may result in a fire. In such situations, turn OFF the power immediately, remove the battery, and contact the nearest authorized service facility.
- Do not put the battery on a fire, short-circuit it, heat it or disassemble it. Doing so may cause explosion or heat generation, resulting in fire or injury.
- Do not use the instrument in conjunction with hyperbaric oxygen equipment. Doing so may cause explosion or fire.
- Do not use the instrument in conjunction with an MRI system. Doing so may cause partial heating due to induced electromotive force, resulting in burn.

Caution

(Failure to adhere to the following points may result in injury or damage to the instrument or other property.)

- Do not apply the instrument to infants.
- Do not use batteries other than those specified by KONICA MINOLTA SENSING. When installing a battery in the instrument, make sure that it is correctly oriented according to the (+ ). Failure to adhere to these instructions may cause explosion of the battery or leakage of electrolyte, resulting in fire, injury or corrosion.
- If the alkaline liquid or other substances inside batteries comes in contact with eyes, skin, or clothing, immediately wash the affected area and see a physician for treatment.
- Do not operate the instrument for long periods of time with the finger holder attached to a patient. Low-temperature burn, redness or rash may result. If you feel pain or itchiness, stop use of the instrument immediately and consult a doctor. A doctor should also be consulted before using the instrument on infants, those with peripheral blood circulation problems or those with sensitive skin.

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© 2009 KONICA MINOLTA SENSING, INC.