

Inspiratory Flow: In-Check Dial



The In-Check Dial measures Inspiratory flow. It can be used to determine whether a patient has the correct inhaler technique.

- Quickly identifies whether the patient is using inhaler correctly when their result and the optimum result is compared
- Aids choice of inhaler for new patients
- Reassures patients and carers that inhaler is being used correctly
- Prevents patient misusing inhaler and drug wastage
- Can be set to measure Inspiratory flow with the resistance of a variety of commonly used inhalers
- Individually calibrated to ensure accuracy



www.clement-clarke.com



**CLEMENT CLARKE
INTERNATIONAL**

Engineered for Life

TOTAL RESPIRATORY SOLUTIONS

Purchase this on our E-Shop at www.clement-clarke.com/shop



Technical Specifications

Measurement Range: – 15 –120 l/min

Accuracy: – 10% or +/- 10 l/min

Repeatability: – 5 l/min

Part Codes

In-Check Dial without Information Scale: 3109900

In-Check Dial with Information Scale: 3109910

In Check Dial 5 T&S: 3109969

In-Check Inhaler Assessment Kit: 3109760

Consumables

One Way Inspiratory Mouthpieces

Allows In-Check Dial to be used as a multi-patient device. One Way valve helps to prevent cross infection



Part Codes

Adult Disposable Inspiratory One Way Valve Mouthpieces

Box of 100 (5x20): 3122061

Adult Disposable Inspiratory One Way Valve Mouthpieces

Box of 200: 3122063

Adult Disposable Inspiratory One Way Valve Mouthpieces

Box of 500: 3122065

Paediatric Mouthpiece Adaptor

(Pack of 1): 3122049

Paediatric Mouthpiece Adaptor

(Pack of 20): 3122120



Two Way Mouthpieces

Cost effective mouthpieces that allow the In-Check Dial to be used as a multi-patient device. Cleaning is recommended between each patient when using these mouthpieces.

Part Codes

Adult Disposable Two Way Mouthpieces

Box of 100 (5x20): 3125030

Box of 500: 3122003

Paediatric Disposable Two Way Mouthpieces

Box of 40: 3125005

Box of 500: 3123030

In-Check Paediatric Adapter for use with Paediatric

Disposable Two Way Mouthpieces (Pack of 1): 3105027

Clement Clarke are well established in offering accurate, precise and dependable products for the diagnosis, monitoring and treatment of respiratory disease

Clement Clarke continues to develop broad product ranges to strengthen our portfolio in **Asthma Management, COPD, Smoking Cessation, Cystic Fibrosis and Allergic Rhinitis** in providing the total respiratory solution.

www.clement-clarke.com

DISTRIBUTED BY: