

## Inspiratory Flow: In-Check Nasal



The In-Check Nasal is a portable Inspiratory flow meter that measures nasal passage obstruction.

- Measures Inspiratory flow through the nose
- Accurate assessment of nasal congestion
- Assists in determining the effectiveness of treatments
- Can be used during assessment of Allergic Rhinitis
- Lightweight and easy to use
- Portable
- In-Check Nasal comes with medium masks. Small masks can also be purchased



[www.clement-clarke.com](http://www.clement-clarke.com)



**CLEMENT CLARKE  
INTERNATIONAL**

Engineered for Life

TOTAL RESPIRATORY SOLUTIONS

Purchase this on our E-Shop at [www.clement-clarke.com/shop](http://www.clement-clarke.com/shop)



**Technical Specifications**

Measurement Range: 15 –120 l/min

Accuracy: 10% or +/- 10 l/min

Repeatability: 5 l/min

**Part Codes**

In-Check with Medium Mask: 3109751

**Consumables****In-Check Masks**

In-Check Nasals masks are available in two sizes.

- Small
- Medium

**Part Codes**

In-Check Nasal Small Masks (Pack of 1): 3109756

In-Check Nasal Medium Mask (Pack of 1): 3109753



[www.clement-clarke.com](http://www.clement-clarke.com)

**Clement Clarke** are well established in offering accurate, precise, dependable products for the diagnosis, monitoring and treatment of respiratory disease

Clement Clarke continues to develop broad product ranges to strengthen our portfolio in **Asthma Management, COPD, Smoking Cessation, Cystic Fibrosis and Allergic Rhinitis** in providing the total respiratory solution.

DISTRIBUTED BY: